



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

November 2, 2016

YTY Industry (Manjung) Sdn. Bhd.,
Punitha Samy
Deputy Manager-DC/RA
Lot 1422-1424
Batu 10 Lekir
Sitiawan, Perak 32020 MY

Re: K162095

Trade/Device Name: Non-sterile, Powder Free Nitrile Examination Gloves, Tested For Use
With Chemotherapy Drugs (Original Blue, Cobalt Blue)

Regulation Number: 21 CFR 880.6250

Regulation Name: Patient Examination Glove

Regulatory Class: I

Product Code: LZA, LZC

Dated: September 19, 2016

Received: September 23, 2016

Dear Punitha Samy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in

the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

Michael J. Ryan -S

for Tina Kiang, Ph.D.

Acting Director

Division of Anesthesiology,

General Hospital, Respiratory,

Infection Control, and Dental Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K162095

Device Name

Non-Sterile, Powder Free Nitrile Examination Gloves Tested For Use With Chemotherapy Drugs - Original Blue

Indications for Use (Describe)

A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. This glove is also tested for use against Chemotherapy Drugs. The Chemotherapy Drugs and its permeation time is listed as below.

Tested Chemotherapy Drug And Concentration

Average Breakthrough Detection Time

Carboplatin, 10 mg/ml	> 240 minutes
Carmustine (3.3 mg/ml)	1.84 minutes
Cisplatin (BCNU), 1.0 mg/ml	> 240 minutes
Cyclophosphamide (Cytosan), 20.0 mg/ml	> 240 minutes
Dacarbazine (DTIC), 10.0 mg/ml	> 240 minutes
Doxorubicin Hydrochloride, 2.0 mg/ml	> 240 minutes
Etoposide (Toposar), 20.0 mg/ml	> 240 minutes
Fluorouracil, 50.0 mg/ml	> 240 minutes
Ifosfamide, 50.0 mg/ml	> 240 minutes
Methotrexate, 25 mg/ml	> 240 minutes
Mitomycin C, 0.5 mg/ml	> 240 minutes
Mitoxantrone, 2 mg/ml	> 240 minutes
Paclitaxel (Taxol), 6.0 mg/ml	> 240 minutes
ThioTEPA (10.0 mg/ml)	0.76 minutes
Vincristine Sulfate, 1.0 mg/ml	> 240 minutes

The following chemotherapy drugs and concentration have extremely low permeation time.

Carmustine (3.3 mg/ml): 1.84 minutes

ThioTEPA (10.0 mg/ml): 0.76 minutes

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Indications for Use

510(k) Number (if known)

K162095

Device Name

Non-Sterile, Powder Free Nitrile Examination Gloves Tested For Use With Chemotherapy Drugs - Cobalt Blue

Indications for Use (Describe)

A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. This glove is also tested for use against Chemotherapy Drugs. The Chemotherapy Drugs and its permeation time is listed as below.

Tested Chemotherapy Drug And Concentration

Average Breakthrough Detection Time

Carboplatin, 10 mg/ml	> 240 minutes
Carmustine (3.3 mg/ml)	1.82 minutes
Cisplatin (BCNU), 1.0 mg/ml	> 240 minutes
Cyclophosphamide (Cytosan), 20.0 mg/ml	> 240 minutes
Dacarbazine (DTIC), 10.0 mg/ml	> 240 minutes
Doxorubicin Hydrochloride, 2.0 mg/ml	> 240 minutes
Etoposide (Toposar), 20.0 mg/ml	> 240 minutes
Fluorouracil, 50.0 mg/ml	> 240 minutes
Ifosfamide , 50.0 mg/ml	> 240 minutes
Methotrexate, 25 mg/ml	> 240 minutes
Mitomycin C, 0.5 mg/ml	> 240 minutes
Mitoxantrone, 2 mg/ml	> 240 minutes
Paclitaxel (Taxol), 6.0 mg/ml	> 240 minutes
ThioTEPA (10.0 mg/ml)	0.93 minutes
Vincristine Sulfate, 1.0 mg/ml	> 240 minutes

The following chemotherapy drugs and concentration have extremely low permeation time.

Carmustine (3.3 mg/ml): 1.82 minutes

ThioTEPA (10.0 mg/ml): 0.93 minutes

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510 (K) SUMMARY SHEETS

1.0**510 (K) SUMMARY****2.0** Submitter

YTY INDUSTRY (MANJUNG) SDN. BHD.,
 Lot 1422-1424, Batu 10 Lekir
 32020 Sitiawan, Perak
 Malaysia
 Tel 605-6908533
 Fax 605-6908533
 Name of Contact Person Punitha Samy (Ms)
 E-mail: punitha.samy@ytygroup.com.my
 Date Summary Prepared November 2, 2016

3.0 Name of Device

Trade Name: Non-Sterile, Powder Free Blue Nitrile Patient Examination Gloves,
 Tested for use with Chemotherapy Drugs (Original Blue, Cobalt Blue)
 Common Name: Nitrile Examination Gloves
 Classification Name: Patient Examination Gloves, Powder Free
 Device Classification: I

Regulation No. & Classification Name: Patient Examination Gloves Specialty (21 CFR
 880.6250 product code LZC)
 Patient Examination Gloves (21 CFR 880.6250 product
 code LZA)

Panel: General Hospital

4.0 Identification of The Legally Marketed Devices that equivalency is claimed:

Device Name	Dermagrip Powder Free Blue Nitrile Patient Examination Gloves, Non- Sterile, Tested for use with Chemotherapy Drugs	Nitrile Powder Free Examination Gloves Tested for Use with Chemotherapy Drugs - Violet Blue (VBLU)
Predicate 510(K) number:	K141982	K151997
Device Classification:	1	1
Product Code:	LZC	LZC

5.0 Description of The Device

Non-Sterile, Powder Free Blue Nitrile Patient Examination Gloves, Tested for use with Chemotherapy Drugs (Original Blue, Cobalt Blue) meets all the current specifications listed under the ASTM Specification D6319-10, Standard Specification for Nitrile Examination Gloves for Medical Application and D6978-05 Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs and FDA 21 CFR 880.6250. The principle operation and mechanism of this device is to prevent contamination between patient and examiner and this principle is achieved through testing of barrier, physical properties and other testing stated in the performance data. This device is for over-the counter single use.

6.0 The Intended Use of Glove

A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. This glove is also tested for use against Chemotherapy Drugs. The Chemotherapy Drugs and its permeation time is listed as below.

Characteristic	Concentration	K162095 Original Blue	K162095 Cobalt Blue
		Average Breakthrough Detection Time (minutes)	
Carboplatin	10mg/ml	> 240	> 240
*Carmustine	3.3mg/ml	1.84	1.82
Cisplatin (BCNU)	1.0mg/ml	> 240	> 240
Cyclophosphamide (Cytosan)	20.0mg/ml	> 240	> 240
Dacarbazine (DTIC)	10.0mg/ml	> 240	> 240
Doxorubicin Hydrochloride	2.0mg/ml	> 240	> 240
Etoposide (Toposar)	20.0mg/ml	> 240	> 240
Fluorouracil	50.0mg/ml	> 240	> 240
Ifosfamide	50.0mg/ml	> 240	> 240
Methotrexate	25mg/ml	> 240	> 240
Mitomycin C	0.5mg/ml	> 240	> 240
Mitoxantrone	2 mg/ml	> 240	> 240
Paclitaxel (Taxol)	6.0mg/ml	> 240	> 240
Thiotepa	10.0mg/ml	0.76	0.93
Vincristine Sulfate	1.0mg/ml	> 240	> 240
Warning Statement		*WARNING: Please note that the following drugs have extremely low permeation times: Carmustine: 1.84 minutes and ThioTepa: 0.76 minutes. NOT TO BE USED WITH CARMUSTINE OR THIOTEPa	* WARNING: Please note that the following drugs have extremely low permeation times: Carmustine: 1.82 minutes and ThioTepa: 0.93 minutes. NOT TO BE USED WITH CARMUSTINE OR THIOTEPa

7.0 Summary of the Technological Characteristic of the Device compared to the Predicate Device for substantial equivalent discussion

There is no difference in technology characteristic compared to the predicate device. Gloves are made from nitrile latex compound. Non-Sterile, Powder Free Blue Nitrile Patient Examination Gloves, Tested for use with Chemotherapy Drugs (Original Blue, Cobalt Blue) have the below technological characteristic compared to ASTM or Equivalent standards.

This glove is compliant with the ASTM standards and similar to the predicate device in the market, which has a glove thickness less than 0.10mm and length less than 270mm. The chemotherapy claim is similar to Predicate.

Table 1 – Predicate 1

Characteristic		Standards	Device Performance	
			Predicate	K162095 – Original Blue
510(K) Number			K141982	K162095
Dimension		ASTM D6319-10	Min 240mm	Min 240mm
Physical Properties		ASTM D6319-10	Meets	Meets
Thickness	Finger	ASTM D6319-10	0.07 – 0.10mm	0.09 – 0.14mm
	Palm		0.07 – 0.09mm	0.07 – 0.10mm
	Cuff		0.06 – 0.08mm	0.06 - 0.09mm
Powder-free		ASTM D6124-06 (≤ 2mg/glove)	Meets	Meets
Bio-compatibility		Primary skin irritation ISO 10993-10 (2010)	Not a primary skin irritant under the conditions of the study	Under the condition of the study the device is non-irritant
		Dermal Sensitization ISO 10993-10 (2010)	Not a contact sensitizer under the conditions of the study	Under the condition of the study the device is non-sensitizer

Table 2 – Predicate 1

Characteristic	Tested according to ASTM D6978-05	Device Performance	
		Predicate-K141982	K162095– Original Blue
		Average Breakthrough Detection Time (minutes)	
	Concentration	Predicate	Current
Carboplatin	10mg/ml	-	> 240
*Carmustine	3.3mg/ml	15	1.84
Cisplatin	1.0mg/ml	> 240	> 240
Cyclophosphamide (Cytosan)	20.0mg/ml	> 240	> 240
Dacarbazine (DTIC)	10.0mg/ml	> 240	> 240
Doxorubicin Hydrochloride	2.0mg/ml	> 240	> 240
Etoposide (Toposar)	20.0mg/ml	> 240	> 240
Fluorouracil	50.0mg/ml	> 240	> 240
Ifosfamide (Ifex)	50.0mg/ml	> 240	> 240
Methotrexate	25mg/ml	-	> 240
Mitomycin C	0.5mg/ml	-	> 240
Mitoxantrone	2 mg/ml	> 240	> 240
Paclitaxel (Taxol)	6.0mg/ml	> 240	> 240
Thiotepa	10.0mg/ml	2.0	0.76
Vincristine Sulfate	1.0mg/ml	> 240	> 240
Warning Statement		*WARNING: Please note that the following drugs have extremely low permeation times: Carmustine (BCNU) : 15 minutes and Thiotepa : 2 minutes	*WARNING: Please note that the following drugs have extremely low permeation times: Carmustine: 1.84 minutes and ThioTepa: 0.76 minutes. NOT TO BE USED WITH CARMUSTINE OR THIOTEPA

Table 3 – Predicate 2

Characteristic		Standards	Device Performance	
			Predicate	K162095 – Cobalt Blue
510(K) Number			K151997	K162095
Dimension		ASTM D6319-10	Min 240mm	Min 240mm
Physical Properties		ASTM D6319-10	Meets	Meets
Thickness	Finger Palm Cuff	ASTM D6319-10	0.10min 0.07 – 0.09mm -	0.10min 0.07 – 0.09mm 0.05 - 0.08mm
Powder-free		ASTM D6124-06 ($\leq$ 2mg/glove)	Meets	Meets
Bio-compatibility		Primary skin irritation ISO 10993-10 (2010)	Not a primary skin irritant under the conditions of the study	Under the condition of the study the device is non-irritant
		Dermal Sensitization ISO 10993-10 (2010)	Not a contact sensitizer under the conditions of the study	Under the condition of the study the device is non-sensitizer

Table 4 – Predicate 2

Characteristic	Tested according to ASTM D6978-05	Device Performance	
		Predicate – K151997	K162095 – Cobalt Blue
		Average Breakthrough Detection Time (minutes)	
	Concentration	Predicate	Current
Carboplatin	10mg/ml	-	> 240
*Carmustine	3.3mg/ml	10.2	1.82
Cisplatin	1.0mg/ml	> 240	> 240
Cyclophosphamide (Cytoxan)	20.0mg/ml	> 240	> 240
Dacarbazine (DTIC)	10.0mg/ml	> 240	> 240
Doxorubicin Hydrochloride	2.0mg/ml	> 240	> 240
Etoposide (Toposar)	20.0mg/ml	> 240	> 240
Fluorouracil	50.0mg/ml	> 240	> 240
Ifosfamide (Ifex)	50.0mg/ml	-	> 240
Methotrexate	25mg/ml	> 240	> 240
Mitomycin C	0.5mg/ml	> 240	> 240
Mitoxantrone	2 mg/ml	-	> 240
Paclitaxel (Taxol)	6.0mg/ml	> 240	> 240
Thiotepa	10.0mg/ml	30.2	0.93
Vincristine Sulfate	1.0mg/ml	> 240	> 240
Warning Statement		* Please note that Carmustine and Thiotepa have extremely low permeation times of 10.2 minutes and 30.2 minutes	* WARNING: Please note that the following drugs have extremely low permeation times: Carmustine: 1.82 minutes and ThioTepa: 0.93 minutes. NOT TO BE USED WITH CARMUSTINE OR THIOTEPA

Table 5

Characteristics	Device Performance				Medical Glove Manual (1661) and Standards
	Predicate K141982	K162095 Original Blue	Predicate K151997	K162095 Cobalt Blue	
Water tightness (1000ml)	Passes	Passes	Passes	Passes	ASTM D5151-06
Color	Blue	Original Blue	Violet Blue	Cobalt Blue	
Material	Nitrile	Nitrile	Nitrile	Nitrile	ASTM D6319-10 (Reapproved 2015)
Size	Extra Small Small Medium Large Extra Large	Extra Small Small Medium Large Extra Large	Extra Small Small Medium Large Extra Large	Extra Small Small Medium Large Extra Large	Labeling
Texture	Finger Textured	Finger Textured	Finger Textured	Finger Textured	
Single Use	Single Use	Single Use	Single Use	Single Use	Labeling
Bio-compatibility	Passes	Passes	Passes	Passes	Primary Skin Irritation -
	Not primary skin irritant under the conditions of the study.	Under the condition of the study the device is non-irritant	Under the condition of the study the device is non-irritant	Under the condition of the study the device is non-irritant	ISO 10993-10: 2010(E) Test for Irritation and Skin Sensitization & Consumer Product Safety Commission II, Part 1500
	Passes	Passes	Passes	Passes	Dermal Sensitization –
	Not a contact sensitizer under the conditions of the study.	Under the condition of the study the device in non-sensitizer	Under the condition of the study the device in non-sensitizer	Under the condition of the study the device in non-sensitizer	ISO 10993-10(E) & Consumer Product Safety Commission, Title 16, Chapter II, Part 1500.3 (c) (4)

8.0 Substantial Equivalent Based on Assessment of Non-Clinical Performance Data

The non-clinical tests performance data support a determination of substantial equivalence is the same as and in compliance to ASTM requirements.

9.0 Conclusion

The conclusions drawn from the intended uses, technological characteristics and non-clinical tests demonstrate that the subject device is substantially equivalent to the predicate device K141982 and K151997.