



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

January 11, 2017

Tianjin UWell Medical Device Manufacturing Co. Ltd.  
Mr. Tao Fan  
General Manager  
A02, Plant B, No. 278, Hangkong Rd, Tianjin Free Trade Zone  
Tianjin, China 300308 CN

Re: K162387  
Trade/Device Name: U-IGNITE Bladeless Trocar  
Regulation Number: 21 CFR 876.1500  
Regulation Name: Endoscope and accessories  
Regulatory Class: Class II  
Product Code: GCJ  
Dated: August 10, 2016  
Received: August 25, 2016

Dear Mr. Fan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and Part 809), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Post Market Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

**Jennifer R. Stevenson -A**

For Binita S. Ashar, M.D., M.B.A., F.A.C.S.  
Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K162387

Device Name

U-IGNITE Bladeless Trocar

Indications for Use (Describe)

The U-IGNITE Bladeless Trocar has applications in endoscopic procedures to provide a port of entry for endoscopic instruments.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

### FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

## 510(k) Summary K162387

**Submitter:** Tianjin UWell Medical Device  
Manufacturing Co., Ltd.  
A02, Plant B, No. 278, Hangkong Road, Tianjin Free Trade Zone(Air Port  
Industrial Park), Tianjin, 300308

**Contact:** Mr. Fan, Tao      General Manager  
Phone: (86)18616806984  
Email: [fant01@uwellmed.com](mailto:fant01@uwellmed.com)

**Date Prepared:** Jan., 4, 2017

**Device Trade Name:** U-IGNITE Bladeless Trocar

**Device Common Name:** Disposable Surgical Trocar / Cannula

**Classification Name:** Laparoscope, General & Plastic Surgery

**Class:** II

**Regulation Number:** 876.1500

**Panel:** General & Plastic Surgery

**Product Code:** GCJ

### Primary Predicate Device:

| Device   | Company                      | Product Code | 510(k) Number |
|--|------------------------------|--------------|---------------|
| ENDOPATH XCEL<br>Bladeless Trocar with<br>OPTI VIEW Technology | Ethicon Endo-Surgery,<br>LLC | GCJ          | K122511       |

### Reference Device:

| Device                                       | Company   | Product Code | 510(k) Number |
|--|---|--------------|---------------|
| Unimicro Trocar Kit,<br>models: Auto-Locking | Unimicro Medical Systems<br>(ShenZhen) Company, | GCJ          | K141594       |

**3. Indications for Use:**

The U-IGNITE Bladeless Trocar has applications in endoscopic procedures to provide a port of entry for endoscopic instruments.

**3. Product Description:**

It is a rigid hand-held surgical instrument with a rounded conical point used to manually puncture a bodily entry point to assist in the location and positioning of gently part the surrounding soft tissue. The trocar assembly punctures the muscle and tissue layers in the body. The trocar assembly's obturator is removed leaving the cannula to provide a port through which laparoscopes and other laparoscopic instruments can be introduced. The cannula sealing system shall be able to maintain pneumoperitoneum.

The U-IGNITE Bladeless Trocar is a sterile single patient use instrument consisting of a radiolucent sleeve and obturator in various sizes diameter. The obturator contains a clear, taped optical element. The obturator accommodate an appropriate sized 0 degree endoscope and provide visibility of individual tissue layers during insertion. Basic type do not have the optical element for its obturator.

In addition, the U-IGNITE Bladeless Trocar with Endo-assistant part is a design enhancement that reduces the incidence of trocar-induced endoscope lens smudging during endoscope insertion. Endoscope smudging occurs when bodily fluids and debris smear across the endoscope lens during a laparoscopic procedure. Trocar-induced smudging occurs when these bodily fluids and debris are deposited within the trocar's seal system when an endoscope or instrument has been exchanged through the trocar.

The trocar sleeves for the 12 mm devices contain two seals, an outer integrated removable self-adjusting seal that accommodates instruments ranging from 5mm to 12 mm in diameter where indicated and an internal seal. Together, these two seals minimize gas leakage when instruments are inserted or withdraw through the trocar. The 5 mm trocar sleeve does not contain an integrated removable outer seal and accommodates only 5 mm instruments.

A stopcock valve is compatible with standard luer lock fitting and provides attachment for gas insufflation and desufflation. The stopcock is in the closed position when it is parallel to the sleeve.

This is an irradiation sterilized single use disposable device to avoid cross infection.

## Product Specification

| Trocar       |                |             |                              |
|--------------|----------------|-------------|------------------------------|
| Product Code | Allow Size(mm) | Length (mm) | Comments                     |
| IG5ST        | 5              | 75          | Without Endo. Assistant part |
| IG5LT        | 5              | 100         | Without Endo. Assistant part |
| SIG5ST       | 5              | 75          | With Endo. Assistant part    |
| SIG5LT       | 5              | 100         | With Endo. Assistant part    |
| BIG5ST       | 5              | 75          | Basic Type                   |
| BIG5LT       | 5              | 100         | Basic Type                   |
| BIG8LT       | 8              | 100         | Basic Type                   |
| IG10LT       | 10             | 100         | Without Endo. Assistant part |
| IG10XT       | 10             | 150         | Without Endo. Assistant part |
| SIG10LT      | 10             | 100         | With Endo. Assistant part    |
| SIG10XT      | 10             | 150         | With Endo. Assistant part    |
| IG12ST       | 12             | 75          | Without Endo. Assistant part |
| IG12LT       | 12             | 100         | Without Endo. Assistant part |
| SIG12ST      | 12             | 75          | With Endo. Assistant part    |
| SIG12LT      | 12             | 100         | With Endo. Assistant part    |
| BIG12ST      | 12             | 75          | Basic Type                   |
| BIG12LT      | 12             | 100         | Basic Type                   |

Basic type is a simplified version of normal type without Endo-Assistant part.

### 4. Non-Clinical Testing

A series of safety tests were performed to assess the safety and effectiveness of the U-IGNITE Bladeless trocar. The safety tests were conducted in accordance with ISO 10993-1, ISO 10993-5, ISO 10993-10, ISO 10993-12, and ISO 11137-1/-2. The performance testing conducted on subject device and predicate device are listed below:

- Obturator & Cannula Compatibility
- Insertion & Cannula Stability
- Air Leakage as a whole device
- Trocar System Puncture Performance (animal simulation test)
- Trocar sleeve retention force (animal simulation test)
- Air Leakage with obturator withdrawn
- Endoscope Visualization Image Quality
- Endo-Assistant Part Liquid Absorbing Performance

All the test results demonstrate U-IGNITE Bladeless trocar meet the requirements of its pre-defined acceptance criteria and intended uses.

### 5. Clinical Data Evaluation:

Based on this literature review no new risks that have not already been identified in the instructions for use or the previously reviewed literature have been noted.

After search from websites cited, post manufacturing experience and data, no clinical incident is reported that is the cause of the trocar device itself.

Identified hazards, such as the occurrence of incisional hernias and injuries to internal viscera and vessels are hazards that are ordinarily associated with operative procedures in the anatomical regions examined regardless of method.

These post-manufacturing experience/ risks are considered & evaluated periodically in our own risk analysis management system, e.g. the Instruction for Use updating based on the literature review result.

Summary of the referred literature and expert opinions that shows that the method is safe and that the device fulfils its intended purpose and clinical performance.

## 6. Substantial Equivalence Determination

The indications for use and technology characteristics of the proposed U-IGNITE Bladeless Trocar, and the substantial equivalence to the predicate devices have been demonstrated via data collected in design verifications. The results of these tests provide reasonable assurance that the proposed device has been designed and tested to assure conformance to the requirements for its indications for use. And the differences between the devices cited in this section do not raise any new issues of safety or effectiveness.

The comparison to predicate devices as below table.

**Comparison to Predicate Device and Reference Device**

| Item                              | Proposed Device<br>(U-IGNITE Bladeless Trocar)  | Primary Predicate Device<br>(ENDOPATH XCEL Bladeless Trocar with OPTI VIEW Technology)  | Reference Device<br>(Unimicro Trocar Kit, models: Auto-Locking Trocar, Hasson)                           |
|-----------------------------------|---|---|--|
| <b>Models</b>                     | Bladeless without Endo-Assistant Part (BIG & IG models)                                     | Bladeless   | Bladeless  |
|                                   | Bladeless with Endo-Assistant Part (SIG model)  | Bladeless with OPTIVIEW™ Technology   | /  |
|                                   | /   | /   | Auto-Locking   |
|                                   | /   | Dilating Tip  | /  |
|                                   | /   | Blunt Tip (Hasson)  | Hasson   |
| <b>Comprised Elements</b>         | Cannula sleeve, Obturator   | Cannula sleeve, Obturator   | Cannula sleeve, Obturator  |
| <b>Dimension</b>                  | Diameter:   | Diameter:   | Diameter:  |
|                                   | 5-12 mm   | 5-12 mm   | 5-12 mm  |
|                                   | Length:   | Length:   | Length:  |
|                                   | 75-150 mm   | 75-150 mm   | 70-120 mm  |
| <b>Classification</b>             | 21 CFR 876.1500   | 21 CFR 876.1500   | 21 CFR 876.1500  |
| <b>Classification and Code</b>    | Class II, GCJ   | Class II, GCJ   | Class II, GCJ  |
| <b>Device Classification Name</b> | Laparoscope, General & Plastic Surgery  | Laparoscope, General & Plastic Surgery  | Laparoscope, General & Plastic Surgery   |
| <b>Indications for Use</b>        | applications in endoscopic procedures to provide a port of entry for endoscopic instruments | Applications in abdominal, thoracic, and gynecologic minimally invasive surgical procedures to establish a path of entry for endoscopic Instruments. The trocar may be used with or without visualization for primary and secondary insertions. | Application in a variety of endoscopic procedures to provide a port of entry for endoscopic instruments. |
| <b>Sterilization</b>              | Gamma sterilization   | Gamma sterilization   | EO sterilization   |

|                                 |  |   |  |
|---------------------------------|--|---|--|
| <b>Safety standards</b>         | ISO 10993-1<br>ISO 10993-5<br>ISO 10993-10<br>ISO 10993-12<br>ISO 11137-1  | ISO 10993-1<br>ISO 10993-5<br>ISO 10993-10<br>ISO 10993-12<br>ISO 11137-1 | ISO 10993-1<br>ISO 10993-5<br>ISO 10993-10<br>ISO 10993-7<br>ISO 10993-12<br>ISO 11135 |
| <b>Performance Standards</b>    | Not Applicable   | Not Applicable  | Not Applicable   |
| <b>Performance Testing Item</b> | Obturator & Cannula Compatibility  | Obturator & Cannula Compatibility   | Obturator & Cannula Compatibility  |
|                                 | Insertion & Cannula Stability  | Insertion & Cannula Stability   | Insertion & Cannula Stability  |
|                                 | Trocar System Puncture Performance   | Trocar System Puncture Performance  | Trocar System Puncture Performance   |
|                                 | Trocar sleeve retention force  | Trocar sleeve retention force   | Trocar sleeve retention force  |
|                                 | Air Leakage as a whole device  | Air Leakage as a whole device   | Air Leakage as a whole device  |
|                                 | Air Leakage with obturator withdrawn   | Air Leakage with obturator withdrawn                                      | Air Leakage with obturator withdrawn   |
|                                 | Endoscope Visualization Quality  | Endoscope Visualization Quality   | Not Applicable   |
|                                 | Endo-Assistant Part Liquid Absorbing Performance   | Endo-Assistant Part Liquid Absorbing Performance                          | Not Applicable   |
| <b>Comparison Summary</b>       | <p>The subject device has the same or simpler intended use comparing to predicate device and reference device.</p> <p>The subject device has same design, construction, and operation principle comparing to predicate device, K122511 and reference device K141594.</p> <p>For both subject device and predicate &amp; reference device, all materials and final assembled product used can comply with the biocompatibility international standard ISO 10993-1, with the evidence from Third Party report. And the test results show the materials used and final product safe biologically.</p> <p>The technology and performance specification of our product (type SIG series and IG series), Bladeless Trocar with Endo-assistant part, are same in Essential Components, physical and general safety feature comparing to the primary predicate device, K122511 from ENDOPATH XCEL Bladeless Trocar with OPTIVIEW Technology, and Bladeless XCEL.</p> <p>Compare to the reference device specified in K141594, our product (type IG, and type BIG series) and the reference device are same in physical and general safety features based on bench test results.</p> <p>The subject device utilizes the same sterilization method, Gamma irradiation sterilization, as what used in primary predicate device, K122511. This sterilization method is safer than EO sterilization method used in K141594.</p> <p>With the same / smaller intended use scope &amp; same principle of operation, and complying with the similar bench tests, the proposed device, Tianjin UWell Medical Device Manufacturing Co. Ltd., performs in same manner to the predicate device, ENDOPATH XCEL Bladeless Trocar with OPTI VIEW Technology, and reference device,</p> |   |  |

|  |  |
|--|--|
|  | <p>Unimicro Trocar Kit, with respect to the following characteristics:</p> <ul style="list-style-type: none"><li>● Design</li><li>● Operation principle</li><li>● Constructed materials,</li><li>● Performance characteristics</li><li>● Perceived safety of the device</li></ul> <p>Any minor structure, material, dimensional or labeling differences between the predicate &amp; reference device and the subject device do not pose risk to their performance and usage.</p> |
|--|--|

## 7. Conclusion

After analyzing bench tests, safety testing data, it can be concluded that: U-IGNITE Bladeless Trocar is as safe and effective as the predicate device and reference device.