



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

May 12, 2017

Avazzia, Inc.  
Tammy Lahutsky  
Regulatory Affairs  
13140 Coit Road  
Suite 515  
Dallas, Texas 75240  
(214) 575-2820 [voice]  
(214) 575-2824 [fax]  
tlahutsky@avazzia.com

Re: K162392

Trade/Device Name: Avazzia OTC TENS Model Best-AV1™; Med-Best™, Med-Sport™,  
Avazzia Blue™, Avazzia Star™

Regulation Number: 21 CFR 882.5890

Regulation Name: Transcutaneous Electrical Nerve Stimulator For Pain Relief

Regulatory Class: Class II

Product Code: NUH

Dated: April 12, 2017

Received: April 14, 2017

Dear Tammy Lahutsky:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

William J.  
Heetderks -S

Digitally signed by William J. Heetderks -S  
DN: c=US, o=U.S. Government, ou=HHS,  
ou=FDA, ou=People,  
0.9.2342.19200300.100.1.1=0010149848,  
cn=William J. Heetderks -S  
Date: 2017.05.12 09:54:28 -0400

for Carlos L. Peña, PhD, MS  
Director  
Division of Neurological  
and Physical Medicine Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K162392

Device Name

Avazzia OTC TENS Model Best-AV1™: Med-Best™, Med-Sport™, Avazzia Blue™, Avazzia Star™

Indications for Use (Describe)

Transcutaneous Electrical Nerve Stimulation (TENS) for the temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, back of the neck, upper extremities (arm), and lower extremities (leg) due to strain from exercise or normal household work activities.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

### Submitter's Identifications:

#### Manufacturer and Sponsor:

Avazzia, Inc.  
13140 Coit Road., Suite 515  
Dallas, TX 75240 USA  
Tel. 214-575-2820  
Establishment Registration # 3004839404  
Official Contact Person: Tammy Lahutsky  
Date of Summary Preparation: August 18, 2016

### Trade Names of Candidate Devices:

Avazzia OTC TENS Model Best-AV1™: Med-Best™, Med-Sport™, Avazzia Blue™, Avazzia Star™

**Classifications:** TENS, Class II 21 CFR 882.5890 NUH

### Information of the 510(k) Cleared Devices (Predicate Devices):

#### Primary Predicate:

Prospera OTC TENS Electronic Pulse Massager, Models PL009, PL009A, and PL029, K122744  
03/28/2013, TENS Class II 21 CFR 882.5890 NUH, NGX

#### Reference Devices:

Avazzia TENS Model BEST-AV1™, K062641, 4/30/2007, TENS, Class II 21 CFR 882.5890 GZJ  
Avazzia TENS Model Best-AV2™, K123099, 05/23/2014, TENS, Class II, 21 CFR 882.5890, GZJ  
Gemore Technology Co. LTd. Low back pain relief system, models model GM310PP,  
GM320PP, GM321PP, K060222, 04/28/2006, TENS Class II 21 CFR 882.5890 NUH

### Intended Use:

Transcutaneous Electrical Nerve Stimulation (TENS) for the temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, back of the neck, upper extremities (arm), and lower extremities (leg) due to strain from exercise or normal household work activities.

### Description of Candidate Devices:

#### Description of Candidate Devices as might be found in the labeling or promotional material:

Avazzia OTC TENS electronic pulse massage devices are electrically powered device intended for over-the-counter use and used to apply a microcurrent electrical pulses to electrodes on a user's skin to relieve pain. They are easy-to-use, handheld, AA battery-operated portable devices for use in the home or clinic. The candidate devices have pre-set modes.

#### Description of the Candidate Devices – How the devices work:

The devices are portable, battery operated microcurrent electrical stimulation device that applies charge and power to the tissue through electrodes where maximum power delivered to the load is controlled and limited, and an automatic shut off is implemented.

The user can apply the built-in onboard electrodes or place self-adhesive conductive electrodes where indicated and apply stimulation for a period of time.

The user controls the output by selecting the mode and power setting.

**Scientific Concepts that form the basis of the devices:**

Various modes in the candidate device are suited for TENS applications.

**Design Description: Significant physical and performance characteristics of the device including stimulation output specifications & Summary of the Technological Characteristics Compared to the Predicates are shown below:**

Discussion of any different technological characteristics between the devices and the predicate devices:

The candidate devices have the same technological characteristics (i.e., design, material, chemical composition, energy source) as the predicate device as summarized in the table of the technological characteristics.

**Accessories**

- Self-adhesive conductive electrode pads
- Soft tissue electrodes with common name: Y electrodes
- Lead wires

**Nonclinical tests**

Verification and validation tests as well as certificates and test reports contained in this submission demonstrate that the submitted models are substantially equivalent to the safety and effectiveness as that of the cleared devices.

The software verification is conducted according to the Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices – Guidance for Industry and FDA Staff – May 11, 2005.

The candidates devices and accessories meet general safety requirements IEC 60601-1 and requirements EMC: IEC 60601-1-2.

**Safety and effectiveness**

The Avazzia device technological specifications are the same; therefore, they do not pose new questions regarding safety and effectiveness.

Indications for use for TENS OTC has been established, therefore, indications for use do not pose new questions regarding safety and effectiveness.

**Basis for a determination of substantial equivalence:**

- *Do the candidate Devices have the same indications for use? Yes*
- *Do the candidate Devices have the same Technological Characteristics? Yes*

- Do the candidate devices pose any new questions regarding safety of effectiveness? No
- Do the candidate devices have accepted scientific methods for evaluating safety and effectiveness?  
Yes
- Do the candidate devices demonstrate no diminishment of performance? Yes

### Conclusions:

The Avazzia OTC TENS Model BEST – AV1™ devices have intended uses and technological characteristics that are substantially equivalent to predicate devices.

Summary of Technical Characteristics		
	Candidates	Predicate
<b>Device Name and Model</b>	BEST-AV1, Med-Best, Med-Sport, Avazzia Blue, Avazzia Star	Models PL009, PL009A, and PL029
Manufacturer	Avazzia, Inc.	Propera
Power Source	Two 1.5 V AA batteries	Four 1.5 V AAA batteries
Method of Line Current Isolation	Battery Supply N/A	Battery Supply N/A
Patient Leakage Current		
- Normal Condition	13.7	1.8
- Single Fault Condition	27	3.3
Average DC current through electrodes when device is on but no pulses are being applied	0	0
Number of output modes	2, 4	3
Number of output channels:	1	Synchronous
- Synchronous/Alternating?	n/a	ON/Off Switch or By
- Method of Channel Isolation	n/a	Software
Regulated Current or Regulated Voltage	Voltage control	Voltage control
Software control	Yes	Yes
Automatic Overload Trip?	No	No
Automatic No-Load Trip?	Yes	No
Automatic Shut Off?	Yes	Yes
User Override Control?	Yes	Yes
Indicator Display		
- On/Off Status?	Yes	Yes
- Low Battery?	Yes	No
- Voltage/Current intensity?	Yes	No
Timer Range	60	5, 10, or 15
Compliance with Voluntary Standards?	Yes	Yes
Compliance with 21 CFR 898	Yes	Yes
Weight	7	8
Dimensions	4.6 x 2.5 x 1.3 (117mm X 64mm X 33mm)	2.2 x 7.8 x 0.9
Housing materials and construction	Enclosure: ABS	Enclosure: ABS
Wavform	Biphasic	monophasic
Shape	positive square wave followed by a damped sinusoidal waveform of variable duration depending on damping and body loading	regular

Summary of Technical Characteristics			
		Candidates	Predicate
<b>Device Name and Model</b>		BEST-AV1, Med-Best, Med-Sport, Avazzia Blue, Avazzia Star	Models PL009, PL009A, and PL029
Max output voltage (+/- 20%)			
- at 500 $\Omega$	V	-42	12.8
- at 2,000 $\Omega$	V	-122	51
- at 10,000 $\Omega$	V	-348	368
- Max output current (+/- 20%)			
- at 500 $\Omega$	$\mu$ A	3080	15,000
- at 2,000 $\Omega$	$\mu$ A	451	3,200
- at 10,000 $\Omega$	$\mu$ A	535	600
Duration of primary (depolarizing phase)	$\mu$ Sec	506	0.040
Pulse Duration	$\mu$ Sec	1100	120 to 6800
Frequency	Hz	20 to 185	1 to 100
Net Charge per pulse at 500 $\Omega$	$\mu$ C	4.0	18,000
Max Charge per pulse at 500 $\Omega$	$\mu$ C	10.12	23
Max current density at 500 $\Omega$	mA/cm <sup>2</sup>	1.6	1.4
Max average power density at 500 $\Omega$	mW/cm <sup>2</sup>	2.445	230
Burst mode			
(a) pulse per burst		Up to 8	1
(b) burst per second		Up to 77	0-25
(c) burst duration (sec)		Up to 6mSec	1
(d) Duty Cycle		46%	5
On Time	Sec	0.46	Potentiometer switch
Off time	Sec	0.54	Potentiometer switch
Additional Features		N/A	N/A