



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

February 8, 2017

Shenzhen Dongdixin Technology Co., Ltd.
Truman Shen
R.A. Specialist
No. 3 Building
Xilibaimang Xusheng
Industrial Estate Nanshan
Shenzhen, 518108 CN

Re: K162479

Trade/Device Name: Smart Pain Reliever, Model LT5019

Regulation Number: 21 CFR 882.5890

Regulation Name: Transcutaneous Electrical Nerve Stimulator For Pain Relief

Regulatory Class: Class II

Product Code: NUH, NGX, NYN

Dated: August 15, 2016

Received: January 9, 2017

Dear Truman Shen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR

Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,


Michael J. Hoffmann -S

for Carlos L. Peña, PhD, MS
Director
Division of Neurological
and Physical Medicine Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K162479

Device Name

Smart Pain Reliever

Model: LT5019

Indications for Use (Describe)

TENS: The device is designed to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household work activities. And to be used for the symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis.

EMS: The device is designed to be used for stimulate healthy muscles in order to improve and facilitate muscle performance.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) SUMMARY

as required by section 21 CFR 807.92

Smart Pain Reliever

Date of Submission: 08/15/2016

Submitter's Name: Shenzhen Dongdixin Technology Co., Ltd

Address: No. 3 Building, Xilibaimang Xusheng Industrial
Estate Nanshan, Shenzhen, CHINA 518108

E-mail: shenc@dundex.com

Tel: +86(755) 27652471

FAX: +86(755) 27652674

Contact: Truman Shen

1. Proposed Device:

Device classification Name: Stimulator, Nerve, Transcutaneous, Over-The-Counter

Regulation Description: Transcutaneous electrical nerve stimulator for pain relief

Regulation Medical Specialty: Neurology

Review Panel: Neurology

Regulation Number: 882.5890

Product Code: NUH

Device Class: II

Device classification Name: Stimulator, Muscle, Powered, For Muscle Conditioning

Regulation Description: Powered muscle stimulator.

Regulation Medical Specialty: Physical Medicine

Review Panel: Physical Medicine

Regulation Number: 890.5850

Product Code: NGX

Device Class: II

Device classification Name: Stimulator, Electrical, Transcutaneous, For Arthritis

Regulation Description: Transcutaneous electrical nerve stimulator for pain relief.

Regulation Medical Specialty: Neurology

Review Panel: Neurology

Regulation Number: 882.5890

Product Code: NYN

Device Class: II

2. Predicate Device:

No.	Legally Marketed Device:	510(k) Number:	Manufacturer:
1	PulseRelief	K151035	Philips Consumer Lifestyle
2	SmartRelief	K131159	Chattem, Inc.
3	Compex® Sport Plus “MASSAGE” program	K083140	Compex S.A.

3. Description of Proposed Device:

LT5019 is designed to be used at home, by adults of all genders.

Explanation of how the device functions:

The Smart Pain Reliever includes TENS, EMS (including MASSAGE Program). Transcutaneous Electrical Nerve Stimulation (TENS) is a noninvasive, drug free method of controlling pain. EMS works by sending electronic pulses to the muscle needing treatment; this causes the muscle to exercise passively.

Scientific concepts that form the basis for the device:

TENS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. TENS does not cure any physiological problem; it only helps control the pain, this activates the underlying sensory nerves. EMS is a product derived from the square waveform (ladder - shaped), through the square wave pattern it is able to work directly on muscle motor neurons, contraction helps to condition the muscle in order to facilitate performance.

For Transcutaneous Electrical Nerve Stimulation (TENS) Self-adhesive electrodes are placed on the skin close to the area of pain. The user can choose 11 pre-set TENS programs with different pulse settings. In each program, the intensity of the pulse can be adjusted.

For Electrical Muscle Stimulation (EMS), the electrodes are placed near the muscle to be stimulated. The user can choose 8 pre-set EMS programs with different pulse settings. In each program, the intensity of the pulse can be adjusted.

The Smart Pain Reliever device only have one massage program, and this massage stimulation program is a stimulation program that is used with the EMS mode of stimulation to facilitate recovery from muscle fatigue and to help recover muscle strength after training sessions and competitions.

The functions of the LT5019 device:

The Smart Pain Reliever device is controlled by means of an APP on a mobile device (phone or table). The communication is done via Bluetooth Low Energy. The APP operates on IOS and Android platforms (IOS 8.0 or greater, Android 5.0 or greater).

Table 1: Smart Pain Reliever device controls and indicators

Control/Indicator	User interface	Function
On/off Button	Press button (when device is off)	Turn on the device
	Press button (when device is on)	Turn off the device (at any time, before, during or after treatment)
LED indicator	off	Device off
	Green flashing	It means that the device is in standby mode
	Orange continuous	It means in treating mode or the battery is being charging
	Orange flashing	It means that the device is in paused mode
	Green continuous	It means it has connected with APP or the battery is full when charging

4. Proposed Device Intended for Use Statement:

Device Name: Smart Pain Reliever, Model: LT5019

Indications for Use:

- ◆ TENS: The device is designed to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities(arm) and lower extremities (leg) due to strain from exercise or normal household work activities. And to be used for the symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis.
- ◆ EMS: The device is designed to be used for stimulate healthy muscles in order to improve and facilitate muscle performance.

5. Technological Characteristics and Substantial Equivalence

Both the LT5019 and Predicate Device utilize the application of electrical current through electrodes placed on the skin for pain control, or electrical muscle stimulation (EMS) the elicitation of muscle contraction using electric impulses. The impulses are generated by the device and delivered through electrodes on the skin in direct proximity to the (painful) muscles to be stimulated.

Basic technological characteristics, new device vs. Predicate device

Table 2: Substantial Equivalence Comparison Table

		New device		Predicate device		Remark
1	510K#	K162479	K151035	K131159	K083140	--
2	Device Name and Model	Smart Pain Reliever Mode: LT5019	PulseRelief	SmartRelief	Compex® Sport Plus “MASSAGE” program	--
3	Manufacturer	Dongdixin Technology Co., Ltd	Philips Consumer Lifestyle	Chattem, Inc.	Compex S.A.	--
4	Intended for use	TENS: The device is designed to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities(arm) and lower extremities (leg) due to strain from exercise or normal household work activities. And to be used for the symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis. EMS: The device is designed to be used for stimulate	The OTC TENS/EMS stimulator PulseRelief is designed to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household work activities. It should be applied to normal, healthy, dry and clean skin of adult patients, and is to be used for stimulate healthy muscles in	To be used for temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities, It is also intended for symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis	Compex® Sport Plus is intended to stimulate healthy muscles in order to improve or facilitate muscle performance.	Similar

		healthy muscles in order to improve and facilitate muscle performance.	order to improve and facilitate muscle performance.				
5	Power Source	3.7V battery supply	3.7V battery supply	3V battery supply	NIMH rechargeable battery (4.8V \geq 1200mA/h) (4cells pack of 1.2V)	SE Note 1	
	-Method of Line current isolation	N/A	N/A	N/A	N/A	Same	
	- Patient Leakage Current (μ A) -Normal condition -Single fault condition	0 μ A 2.0 μ A	<10 μ A <50 μ A	N/A	N/A N/A	SE Note 1	
6	Average DC current through electrodes when device is on but no pulses are being applied (μ A)	0	0	0	0	Same	
7	Number of Output Modes	2 (TENS/EMS(including MASSAGE program))	2 (TENS/EMS)	1 (TENS)	9 modes/programs	SE Note 2	
8	Number of Output Channels:	Synchro nous or Alternati -ng?	N/A	N/A	N/A	Synchronous	SE Note 2
		Method of Channel Isolation	N/A	N/A	N/A	N/A	Same
9	Regulated Current	Current control	Current control	Current control	Current control	Same	

	or Regulated Voltage?						
10	Software/Firmware /Micro processor Control?	Yes	Yes	Yes	Yes	Same	
11	Automatic Overload Trip	Yes	Yes	No	No	SE Note 2	
12	Automatic No Load contact Trip	Yes	Yes	No	Yes	SE Note 2	
13	Automatic Shut off	Yes	Yes	Yes	Yes	Same	
14	User Override Control?	Yes Power on/off button on the device, and power on/off in the APP software	Yes Power on/off button on the device, and power on/off in the APP software	No	Power on/off button	SE Note 2	
15	Indicator Display:	On/Off Status?	Yes	Yes	Yes	Yes	Same
		Low Battery?	Yes	Yes	No	Yes	SE Note 2
		Voltage/Current Level?	Yes	Yes	Active/inactive output	Yes	SE Note 2
16	Timer Range (minutes)	1 ~ 59 minutes and continuous	1 ~ 59 minutes and continuous	30 minutes	3-55.5min	SE Note 2	
17	Weight (grams.)	28grams (main device)	62 grams (main device)	20 grams	350g	SE Note 2	
18	Dimensions (cm.) H*W * L	13x41x67mm(main device)	2 units, each 54x54x14 (excl electrodes)	64x38x13[mm]	L142*W99*H36	SE Note 2	
19	Housing Materials & Construction	PC + ABS	PC + ABS	Molded Plastic	/	SE Note 3	

Table 3: Output Specification TENS mode

		New device	Predicate device		Remark
1	510K#	K162479	K151035	K131159	--
2	Device Name	Smart Pain Reliever Mode: LT5019	PulseRelief	SmartRelief	--
3	Waveform	Biphasic	Biphasic	Asymmetrical Biphasic	Similar Note 4
4	Shape	Rectangular	Rectangular	Rectangular	Same
5	Max Output Voltage (V) $\pm 20\%$				Similar Note 5
6	500 Ω	31.2	31	30	
7	2k Ω	69.6	69	70	
8	10k Ω	69.6	70	70	
9	Pulse Duration (μ sec)	150~250us	60~350us	30-220 μ s at 50% of peak amplitude	Similar Note 5
10	Frequency (Hz)	2~100Hz	1~100Hz	1~100Hz	Similar Note 5
11	Maximum Phase Charge (μ C) 500 Ω	30	1.6~6.8	13.2	Similar Note 5
12	Maximum Current Density 500 Ω	0.32mA/cm ²	0.002~0.045mA/cm ²	0.1mA/cm ²	Similar Note 5
13	Maximum Average Current (average absolute value), mA, 500 Ω	1.2	0.06~1.36	2.06	Similar Note 5
14	Maximum Average Power Density, (mW/cm ²),500 Ω	1.9	0.24~1.69	1.5	Similar Note 5

Table 4: Output Specification EMS mode

		New device	Predicate device	Remark
1	510K#	K162479	K151035	--
2	Device Name or Program Name	Smart Pain Reliever Mode: LT5019	PulseRelief	--
3	Waveform	Biphasic	Biphasic	Same
4	Shape	Rectangular	Rectangular	Same
5	Max Output Voltage (V) $\pm 20\%$			Similar Note 5
6	500 Ω	31.2	31	
7	2k Ω	69.6	69	
8	10k Ω	69.6	70	
9	Pulse Duration (usec)	200~370us	150~350us	Similar Note 5
10	Frequency (Hz)	3~75Hz	40~65Hz	Similar Note 5
11	Maximum Phase Charge (μ C) 500 Ω	24	4.7~10.9	Similar Note 5
12	Maximum Current Density	0.248mA/cm ²	0.019~0.037mA/cm ²	Similar

	500Ω			Note 5
13	Maximum Average Current (average absolute value), mA, 500Ω	0.72	0.47~0.93	Similar Note 5
14	Maximum Average Power Density, (mW/cm ²),500Ω	1.2	0.62~1.15	Similar Note 5

Table 5: Output Specification MASSAGE mode

		New device	Predicate device	Remark
1	510K#	To be assigned	K083140	--
2	Device Name or Program Name	Smart Pain Reliever Mode: LT5019	Compex® Sport Plus “MASSAGE” program	--
3	Waveform	Biphasic	Biphasic	Same
4	Shape	Rectangular	Rectangular	Same
5	Max Output Voltage (V) ±20%			Similar Note 5
6	500Ω	31.2	60	
7	2kΩ	69.6	136	
8	10kΩ	69.6	137	
9	Pulse Duration (usec)	250us	200~400us	Similar Note 5
10	Frequency (Hz)	10~80Hz	2~100Hz	Similar Note 5
11	Maximum Phase Charge (uC) 500Ω	30	48	Similar Note 5
12	Maximum Current Density 500Ω	0.32mA/cm ²	Not available	Similar Note 5
13	Maximum Average Power Density, (mW/cm ²),500Ω	1.9	<12mW/cm ² (microwatt) @500Ω	Similar Note 5

Comparison in Detail(s):

Note1:

Although the “Power Source(s)” and “Patient Leakage Current” of subject device are a little different from the predicate devices, they all comply with IEC 60601-1 requirements. So the differences will not raise any safety issue.

Note2:

Although the “Number of Output Modes” “Output Intensity Level”, “Method of Channel Isolation”, “Timer Range”, “Weight” and “Dimensions” of subject device are different from the predicate devices, they are all comply with IEC 60601-1 and IEC 60601-2-10 requirements. So the differences of the function specifications will not raise any safety or effectiveness issue. **Note**

3:

Although the “Housing Materials & Construction” of subject device is different from the predicate devices, they are all comply with ISO 10993-1, ISO 10993-5 and ISO 10993-10 requirements. So the differences of the Housing Materials & Construction will not raise any biological hazard.

Note4:

Although the “Waveform” of subject device is different from the predicate devices, the

waveform is the key factors of the Net Charge (per pulse), and all the Net Charge (per pulse) are zero. The “Biphasic” waveform won’t create Net Charge (per pulse), and not cause charge accumulation on the skin and cause burns. So the differences of the “Waveform” will not raise any safety or effectiveness issue.

Note5:

Although the “Maximum Output Voltage”, “Pulse Duration”, “Maximum pulse frequency”, “Maximum Phase Charge”, “Maximum Average Current”, “Maximum Current Density”, “Maximum Average Power Density of subject device”, are a little different from the predicate devices, they all comply with the FDA guidance requirement for Transcutaneous Electrical Nerve Stimulator for Pain Relief and FDA guidance requirement for Powered Muscle Stimulator for Muscle Conditioning. So the differences of function specification will not raise any safety or effectiveness issue.

6. Performance Data:

The following performance data are provided in support of the substantial equivalence determination:

6.1 Biocompatibility testing

The biocompatibility evaluation for the LT5019 was conducted in accordance with the International Standard ISO 10993-1 “Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing Within a Risk Management Process”. As dictated by the application and duration of contact with the intact skin, the shell and electrode of testing included the following tests:

- Cytotoxicity
- Sensitization
- Irritation

6.2 Electrical safety and electromagnetic compatibility (EMC)

Electrical safety and EMC testing were conducted on the LT5019. The system complies with the IEC 60601-1, IEC 60601-1-6, IEC 60601-1-11 and IEC 60601-2-10 standards for safety and the IEC 60601-1-2 standard for EMC, and IEC 62133 for battery safety.

6.3 Software Verification and Validation Testing

Software verification and validation testing were conducted and documentation was provided as recommended by FDA’s Guidance for Industry and FDA Staff, “Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices.” The software for this device was considered as a “Moderate” level of concern. IEC 62304 was followed.

6.4 Usability Testing

Usability testing according to IEC 62366, Applying Human Factors and Usability Engineering to Medical Device, was conducted.

7. Conclusions

The intended use and basic technological characteristics of the LT5019 device are equivalent with those of the referenced Predicate device K151035, K083140 and K131159.

The LT5019 device complies with the requirement of IEC60601-1, IEC 60601-2-10, IEC60601-1-2 and IEC60601-1-11. The bench testing and safety report documentation supplied in this submission demonstrates that any differences in their technological characteristics do not raise any new questions of safety or effectiveness.

Moreover, the bench testing include the testing performed to characterize the waveform and stimulus parameters, and the performance testing proved that the ability of the LT5019 device to deliver a stimulus that meets design specifications.

In the other words, those engineering difference do not affect the intended use or alter the fundamental scientific technology of the device.

Thus, the Smart Pain Reliever device (K162479) is substantially equivalent to the predicate devices (K151035, K131159 and K083140).