



April 28, 2017

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

IOB Medical, Inc.  
% Joe Shia  
Regulatory Consultant  
LSI International, Inc.  
504 East Diamond Ave., Suite F  
Gaithersburg, Maryland 20877

Re: K162679

Trade/Device Name: IOB Temperature Management System  
Regulation Number: 21 CFR 870.5900  
Regulation Name: Thermal Regulating System  
Regulatory Class: Class II  
Product Code: DWJ  
Dated: March 10, 2017  
Received: March 20, 2017

Dear Joe Shia:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in

the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

 **Fernando  
Aguel-S**

for Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K162679

Device Name

IOB Temperature Management System

Indications for Use (Describe)

The IOB Temperature Management system is indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) SUMMARY

1. Date: April 28, 2017
2. Submitter: IOB Medical Inc  
504E Diamond Ave., Suite I  
Gaithersburg, MD 20877
3. Contact person: Joe Shia  
LSI International Inc.  
504 East Diamond Ave., Suite F  
Gaithersburg, MD 20877  
Telephone: 240-505-7880  
Fax: 301-916-6213  
Email: shiajl@yahoo.com
4. Device Name: IOB Temperature Management System
5. Classification:  
**Class: Class II**

Product Code	CFR #	Product Name
DWJ	870.5900	Thermal Regulating System

6. Predicate Devices:  
K021473  
Bair Hugger Temperature Management System
7. Intended Use  
The IOB Temperature Management system is indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated.
8. Device Description  
The IOB Temperature Management System consists of the IOB warming unit and the IOB warming blankets. The IOB warming unit draws ambient- temperature air through a 0.2 micron air filter. The filtered air is warmed to a selected temperature. The warmed air enters the IOB Warming Blanket through the hose and is distributed through delivery channels. Perforations on the patient side of the air delivery channels in the warming blanket gently disperse the warmed air over and around the patient. The warming unit has three temperature settings of 32°C, 38°C, and 43°C. These temperature settings are servo-controlled by a thermistor placed at the end of hose where the hose connects to the blanket. The unit can also deliver ambient-temperature air. The temperature indicated on the control panel is the temperature of air being delivered to the blanket at the end of the hose. A control thermistor in the warming unit adjusts the power applied to the heater to maintain the selected temperature. This enables the system to maintain the selected temperature under variations in ambient temperature. Besides, the warming unit has high and low air flow for choose.

The IOB Warming Blankets in this submission are the following:

- Torso IOB-001
- Lower Body IOB-002
- Upper Body IOB-003
- Full Body IOB-004
- Pediatric Underbody IOB-005
- Adult Underbody IOB-006
- Pediatric Full Body IOB-007
- Full Body Surgical IOB-008
- Large Pediatric Underbody IOB-009
- Spinal Underbody IOB-010
- Lithotomy Underbody IOB-011
- Pediatric Lower Body IOB-012
- Pediatric Long IOB-014
- Cath Lab IOB-015
- Surgical Access IOB-016
- Chest Access IOB-017
- Multi-Access IOB-018
- Dual Port Torso IOB-019
- Cardiac Access IOB-020
- XL Upper Body IOB-021
- Outpatient Care IOB-022
- Cardiac IOB-023

These blankets are single-use and disposable. Each blanket consists of two layers of non-woven polypropylene fabric coated with a layer of polyethylene. The layers are bonded together to form a distribution network of air delivery channels. The warm air is distributed around the blanket through the delivery channels and exits the blanket through a specially designed series of perforations in the patient side of the blanket.

9. Substantial Equivalence Information

A summary comparison of features of the IOB Temperature Management system and the predicate devices is provided in following Tables.

**Table 1: Comparison between IOB Warming Blankets and predicate devices.**

Parameters	Predicate devices K021473	Proposed devices
	Bair Hugger blankets (models 540,525,522,300,555,635,310,610,550,575,585,537,530,560,570,305,315,542,645,523,110,630)	IOB Warming Blankets (IOB-001,IOB-002,IOB-003,IOB-004,IOB-005,IOB-006,IOB-007,IOB-008,IOB-009,IOB-010,IOB-011, IOB-012,IOB-014,IOB-15,IOB-16,IOB-17, IOB-18, IOB-019,IOB-020,IOB-021, IOB-022,IOB-023)

INDICATIONS FOR USE	The Bair Hugger temperature management system is indicated for hyper- or hypothermic patients or normothermic patients for whom induced hyper- or hypothermia or localized temperature therapy is clinically indicated. In addition, the Bair Hugger temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to become too warm or too cold. The Bair Hugger temperature management system can be used with adult and pediatric patients.	The IOB Temperature Management system is indicated for hypothermic patients or normothermic patients for whom induced hyperthermia or localized increase in temperature is clinically indicated.
MATERIAL DESIGN	Consists of two layers of non- woven polypropylene fabric bonded to a fusion layer of polyethylene. The layers are bonded together to form a distribution network of air delivery channels The warm air is distributed around the patient's body through the delivery channels and exits the blanket through a specially designed series of perforations in the patient side layer of the blanket. The distribution of air is designed to minimize temperature differences of delivered air at different blanket locations.	SAME
Shelf Life	3 years	SAME
Sterility	Non-sterile except 645 and 630	Non-sterile except IOB-14, IOB-16, IOB-020 and IOB-023

Blanket Dimensions (approximate)	540	107×91cm	IOB-001	142×120cm
	525	152×91cm	IOB-002	142×120cm
	522	188×61cm	IOB-003	195×80cm
	300	213×91cm	IOB-004	210×120cm
	635	221×91cm	IOB-006	215×100cm
	610	183×91cm	IOB-008	210×120cm
	575	188×91cm	IOB-010	215×100cm
	585	188×91cm	IOB-011	170×100cm
	560	23×335cm	IOB-015	17x180cm
	570	213×91cm	IOB-016	210×120cm
	305	183×91cm	IOB-017	180×120cm
	315	213×91cm	IOB-018	210×120cm
	542	107×91cm	IOB-019	142×120cm
	523	213×91cm	IOB-021	215×80cm
	110	213×91cm	IOB-022	210×120cm
	310	152×91cm	IOB-007	170×100cm
	555	91×84cm	IOB-005	91×80cm
	550	152×81cm	IOB-009	160×80cm
	537	89×61cm	IOB-012	142×100cm
	530	188×21cm	IOB-014	110×17cm
645	23x274cm	IOB-020	17x150cm	
630	152x91cm	IOB-023	142x120cm	
Air Velocity	28-30cfm		Same	
Air Temperature	Hose End Temperature 43°C +/- 3.0°C 38°C +/- 3.0°C 32°C +/- 3.0°C Ambient		Hose End Temperature 43°C +/- 3.0°C 38°C +/- 3.0°C 32°C +/- 3.0°C Ambient	
System Power	110-120VAC,60Hz,10A 220-240VAC,50Hz,4.5A		110-120V, 60 Hz, 12A 230V, 50/60 Hz, 8A	
Heater Power	850 Watts		1000 W	
Dimensions	33X 25 X 28cm		28 X 22X 22cm	
Weight	5.2kg		4.5kg	
EMI/EMC Compliant	IEC60601-1, IEC 60601-1-2		IEC60601-1, IEC 60601-1-2	
Forced air Over Temperature	Auto-shuts heater off at 47°C +/- 2°C		Auto-shuts heater off at 47°C +/- 2°C .	

Hose with Secure Locking	Yes	Yes
Air Filter	Replaceable 0.2 micron	Replaceable 0.2 micron
Temperature Display	Front panel LCD display	Front panel digital display

10. Safety and Performance Characteristics

1. Temperature uniformity tests show equivalence between the IOB warming blankets and the predicate.
2. Accelerated stability tests show three years shelf-life of the IOB warming blankets.
3. Electrical safety/EMC testing, and software validation show that the device meets its design specifications, performs as intended.
4. Biocompatibility testing (cytotoxicity, irritation and sensitivity) according to ISO 10993 for a limited contact device was demonstrated to be suitable for the intended use of the product.
5. Clinical Studies  
Not applicable

11. Conclusion

Based on the information presented in this 510K premarket notification, the IOB Temperature Management System is substantially equivalent to the predicate.