



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

April 20, 2017

Jiangxi Hongda Medical Equipment Group, Ltd.
% Diana Hong
General Manager
Mid-Link Consulting Co., Ltd.
P.O. Box 120-119
Shanghai, 200120
CHINA

Re: K163160
Trade/Device Name: Sterile Single-use Infusion Set
Regulation Number: 21 CFR 880.5440
Regulation Name: Intravascular Administration Set
Regulatory Class: Class II
Product Code: FPA
Dated: March 10, 2017
Received: March 21, 2017

Dear Diana Hong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

 Tina Kiang -S

Tina Kiang, Ph.D.
Acting Director
Division of Anesthesiology,
General Hospital, Respiratory,
Infection Control, and Dental Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K163160

Device Name

Sterile Single-use Infusion Set

Indications for Use (Describe)

The device is intended to administer fluids from a container to a patient's vascular system through a needle or catheter inserted into the vein.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Exhibit #2 510(k) Summary

This 510(k) Summary is being submitted in accordance with requirements of Title 21, CFR Section 807.92.

The assigned 510(k) Number: K163160

1. Date of Preparation: 04/17/2017
2. Sponsor Identification

JiangXi HongDa Medical Equipment Group Ltd.

No. 39 South Shengli Road, Jinxian County, Nanchang City,
Jiangxi Province, 331700, China

Establishment Registration Number: Not yet registered.

Contact Person: Xiaohua Yang
Position: Vice-general Manager
Tel: +86-791-85628918
Fax: +86-791-85628881
Email: yxh2821006@126.com

3. Designated Submission Correspondent

Ms. Diana Hong (Primary Contact Person)
Ms. Ying Xu (Alternative Contact Person)

Mid-Link Consulting Co., Ltd

P.O. Box 120-119, Shanghai, 200120, China

Tel: +86-21-22815850,
Fax: 240-238-7587
Email: info@mid-link.net

4. Identification of Proposed Device

Trade Name: Sterile Single-use Infusion Set

Common Name: Disposable Infusion Set

Regulatory Information

Classification Name: Set, Administration, Intravascular;

Classification: II;

Product Code: FPA;

Regulation Number: 21CFR 880.5440

Review Panel: General Hospital;

Indication for Use:

The device is intended to administer fluids from a container to a patient's vascular system through a needle or catheter inserted into the vein.

Device Description

The proposed device is intended to administer fluids from a container to a patient's vascular system through a needle or catheter inserted into the vein under the action of gravity. The proposed device is provided with/without injection needle. It is comprised of protective cap of spike, spike, air vent, air filter, drip chamber, fluid filter, flexible tube, injection site, roller clamp, luer lock connector, protector cap of luer lock connector. The devices are provided sterile and single use.

5. Identification of Predicate Device

Premark Notification

510(k) Number: K112204

Trade Name: KDL Disposable Infusion Set

Classification Name: Set, Administration, Intravascular

Manufacturer: Shanghai Kindly Enterprise Development Group Co., Ltd

6. Non-Clinical Test Conclusion

Non clinical tests were conducted to verify that the proposed device met all design specifications as was Substantially Equivalent (SE) to the predicate device. The test provided in this submission include:

Physical, Mechanical and Chemical Tests performed on the proposed device

Particulate contamination	Clause 6.1 of ISO 8536-4:2010 AMD1 2013
Leakage	Clause 6.2 of ISO 8536-4:2010 AMD1 2013
Tensile strength	Clause 6.3 of ISO 8536-4:2010 AMD1 2013
Closure-piercing device	Clause 6.4 of ISO 8536-4:2010 AMD1 2013
Air-inlet device	Clause 6.5 of ISO 8536-4:2010 AMD1 2013
Tubing	Clause 6.6 of ISO 8536-4:2010 AMD1 2013
Fluid filter	Clause 6.7 of ISO 8536-4:2010 AMD1 2013
Drip chamber and drip tube	Clause 6.8 of ISO 8536-4:2010 AMD1 2013
Flow regulator	Clause 6.9 of ISO 8536-4:2010 AMD1 2013
Flow rate of infusion fluid	Clause 6.10 of ISO 8536-4:2010 AMD1 2013
Injection site	Clause 6.11 of ISO 8536-4:2010 AMD1 2013
Male conical fitting	Clause 6.12 of ISO 8536-4:2010 AMD1 2013
Protective caps	Clause 6.13 of ISO 8536-4:2010 AMD1 2013
Reducing (oxidizable) matter	Clause 7.1 of ISO 8536-4:2010 AMD1 2013
Metal ions	Clause 7.2 of ISO 8536-4:2010 AMD1 2013
Titration acidity or alkalinity	Clause 7.3 of ISO 8536-4:2010 AMD1 2013
Residue on evaporation	Clause 7.4 of ISO 8536-4:2010 AMD1 2013
UV absorption of extract solution	Clause 7.5 of ISO 8536-4:2010 AMD1 2013
Materials	Clause 5 of ISO 8536-14:2016
Tubing specifications	Clause 6.1 of ISO 8536-14:2016
Operating temperature	Clause 6.2 of ISO 8536-14:2016
Construction	Clause 6.3 of ISO 8536-14:2016
Flow rate	Clause 6.4 of ISO 8536-14:2016
Materials	Clause 3 of ISO 9626:1991/AMD-1:2001
Surface finish	Clause 4 of ISO 9626:1991/AMD-1:2001
Cleanliness	Clause 5 of ISO 9626:1991/AMD-1:2001
Limits for acidity and alkalinity	Clause 6 of ISO 9626:1991/AMD-1:2001
Size designation	Clause 7 of ISO 9626:1991/AMD-1:2001
Dimensions	Clause 8 of ISO 9626:1991/AMD-1:2001
Stiffness	Clause 9 of ISO 9626:1991/AMD-1:2001
Resistance to breakage	Clause 10 of ISO 9626:1991/AMD-1:2001
Resistance to corrosion	Clause 11 of ISO 9626:1991/AMD-1:2001
Cleanliness	Clause 4 of ISO 7864:1993
Limits for acidity or alkalinity	Clause 5 of ISO 7864:1993
Limits for extractable metals	Clause 6 of ISO 7864:1993

Size designation	Clause 7 of ISO 7864:1993
Colour coding	Clause 8 of ISO 7864:1993
Needle hub	Clause 9 of ISO 7864:1993
Sheath	Clause 10 of ISO 7864:1993
Needle tube	Clause 11 of ISO 7864:1993
Needle point	Clause 12 of ISO 7864:1993
Performance	Clause 13 of ISO 7864:1993

Gauge	Clause 4.1 of ISO 594-1:1986
Liquid leakage	Clause 4.2.1 of ISO 594-2:1998
Air leakage	Clause 4.2.2 of ISO 594-2:1998
Separation force	Clause 4.3 of ISO 594-2:1998
Unscrewing torque	Clause 4.4 of ISO 594-2:1998
Ease of assembly	Clause 4.5 of ISO 594-2:1998
Resistance to overriding	Clause 4.6 of ISO 594-2:1998
Stress cracking	Clause 4.7 of ISO 594-2:1998

Sterile Barrier Packaging Testing performed on the proposed device:

Seal strength	ASTM F88/F88-09
Internal pressure	ASTM F1140/F1140M-13
Dye Penetration	ASTM F 1929-12

Sterilization and Shelf Life Testing performed on the proposed device:

EO residue	ISO 10993-7:2008
ECH residue	ISO 10993-7:2008
Bacteria Endotoxin Limit	USP 36-NF 31 <85>
Shelf Life Evaluation	Physical, Mechanical, Chemical, Package and Sterility Tests were performed on accelerated aging samples to verify the claimed shelf life of the device

Biocompatibility Testing:

The patient-contact materials of blood collection sets are identified and biocompatibility testing is performed according to ISO 10993 standards.

Transportation Testing

Transportation test is performed on the final product to verify its package integrity during transportation.

7. Clinical Test Conclusion

No clinical study is included in this submission.

8. Substantially Equivalent (SE) Comparison

Table 1 Comparison of Technology Characteristics

Item	Proposed Device		Predicate Device K112204	
Product Code	FPA		Same	
Regulation Number	21 CFR 880.5440		Same	
Indication for Use	The device is intended to administer fluids from a container to a patient's vascular system through a needle or catheter inserted into the vein.		KDL Disposable Infusion Set is intended to administer fluids from a container to a patient's vascular system through a needle or catheter inserted into the vein.	
Configuration and material	Configuration	Material	Configuration	Material
	Protector Cap of Spike	HDPE	Protective Cap of Closure-piercing Device	Unknown
	Spike	ABS	Closure-piercing Device	
	Air Vent	PVC	Air Filter	
	Drip Chamber	PVC	Drip Chamber	
	Fluid Filter	ABS	Fluid Filter	
	Flexible Tube	PVC	Flexible Tube	
	Roller Clamp	HDPE	Roller Clamp	
	Latex Tube	Polyisoprene	Check Valve	
	Y Injection Site	ABS	Y Injection Set	
	Luer Lock Connector	ABS	Clamp	
	Protector Cap of Luer Lock Connector	HDPE	Injection Site	
	Needle Hub	PP	Luer Lock Connector	
	Needle Tube	Stainless Steel	Infusion Needle	
Tubing Diameter	3.9mm		Unknown	
Filter Characteristics	15µm		15µm	

Needle Gauge	21G	Unknown
Sterile	EO sterilized	Same
	10^{-6}	Same
Single Use	Single Use	Same
Infusion Set Performance		
Particulate contamination	Contamination index limit is less than 90	Comply with ISO 8536-4
Leakage	No leakage	
Tensile strength	Withstand a static tensile force of not less than 15N for 15s	
Closure piercing device	Comply the dimension of ISO 8536	
Air-inlet device	Flow of fluid is reduced less than 20% of that from a freely ventilated container	
Tubing	Length is not less than 1500mm	
Fluid filter	Retention of latex particles is not less than 80%	
Drip chamber and drip tube	20 drops of distilled water delivered by the drip are equivalence to a volume of $(1\pm 0.1)\text{ml}(1\pm 0.1)\text{g}$	
Flow rate	Deliver not less than 1000ml	
Injection site	No leakage	
Construction	The clamps can resist the flow of fluid and air at an applied pressure of 50 kPa when closed.	
Needle Performance		
Cleanliness	Free from particles and extraneous matter	Comply with ISO 7864 and ISO 9626
Needle length	38mm	
Bond of hub and needle tube	Not broken by the minimum force	
Stiffness	Deflection is less than 0.50mm	
Dimensions of needle tubing	0.8mm	
Resistance to breakage	The tubing is not break	

Resistance to corrosion	No evidence of corrosion	
Biocompatibility		
Cytotoxicity	No Cytotoxicity	Conform with ISO 10993 requirements
Intracutaneous Reactivity	No Intracutaneous Reactivity	
Skin Sensitization	No Skin Sensitization	
Acute Systemic Toxicity	No Acute Systemic Toxicity	
Hemolysis	No Hemolysis	
Pyrogen	No pyrogen	

9. Substantially Equivalent (SE) Conclusion

Based on the performance testing, comparison and analysis above, the proposed devices are determined to be Substantially Equivalent (SE) to the predicate devices.