



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

March 24, 2017

Brightway Holdings Sdn Bhd
G. Baskaran
Group Managing Director
Lot 1559, Jalan Istimewa, Batu Belah
Klang, 42100 MY

Re: K163267

Trade/Device Name: Brightway Nitrile Examination Glove, powder free (blue)
Regulation Number: 21 CFR 880.6250
Regulation Name: Patient Examination Glove
Regulatory Class: Class I
Product Code: LZA
Dated: November 7, 2016
Received: November 21, 2016

Dear G. Baskaran:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,


Michael J. Ryan -S

for Tina Kiang, Ph.D.

Acting Director

Division of Anesthesiology,

General Hospital, Respiratory,

Infection Control, and Dental Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K163267

Device Name

BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE (BLUE)

Indications for Use (Describe)

BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE (BLUE) is a disposable device intended for medical purpose to be worn on the examiners hand or finger to prevent contamination between the patient and the examiner.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) SUMMARY

K163267

BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE, [BLUE]

1. Submitter :

Company Name : **BRIGHTWAY HOLDINGS SDN. BHD.**

Street Address : Lot 1559, Jalan Istimewa,
Batu Belah, 42100 Klang
Selangor Darul Ehsan.

Country : Malaysia

Phone No. : 603-3343 1007 & 603-3343 1094.

Fax No. : 603-3341 4800

E-mail Address : brightway@brightway919.com

Contact Person : Mr. G. Baskaran (Group Managing Director)
baskar@brightway919.com
Mr. Felix Darrel (Group Marketing Manager)
felix.marketing@brightway919.com

2. Preparation Date : 18th March 2017

3. Name of the Device :

Device trade or proprietary name: BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE, [BLUE]

Device Classification Name: Polymer Patient Examination Glove
(21 CFR 88.6250)

Device common or usual name : NITRILE EXAMINATION GLOVE

FDA Device Class : Class 1

Product Code : LZA

4. Identification of the Device :

Class I patient Examination Nitrile gloves, Powder Free, LZA, which meets all the requirement of ASTM D 6319-10 and FDA 21 CFR 880.6250.

Predicate Device:

Legally Marketed Devices to which Substantial Equivalence is claimed ;
K143477- Nitrile Patient Examination Gloves, Powderfree, Blue Color

5. Device Description :

The subject device in this 510(k) Notification is a Blue Examination Glove.

The subject device is a patient examination glove made from a Nitrile compound, Blue in colour, powder free and non sterile (as per 21 CFR 880.6250, class I).

The principle operation of the medical device to provide single use barrier protection for the wearer and the device that meets all the requirement specifications for Barrier Protection, tensile properties as defined in ASTM D6319-10 ; Standard specification for Nitrile Examination Gloves.

This device is manufactured in facilities compliant to ISO 9001:2008 certified in Manufacture of Non Sterile Natural (Latex) and Synthetic Latex (Nitrile) Examination, Surgical and Industrial Gloves & Nitrile Sheath

The device is manufactured to comply with ISO 13485:2003 / EN ISO 13485:2012; Manufacture of Non Sterile Natural (Latex) and Synthetic Latex (Nitrile) Examination & Sterile Surgical Gloves

6. Intended use of the Device

BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE, [BLUE] is a disposable device intended for medical purpose, to be worn on the examiners hand or finger to prevent contamination between patient and examiner.

7. Specification For Nitrile gloves

7.1. Dimension and Thickness of Gloves

Certificate No : 6H22 06

Date : 22/08/2016

Brand : 9.5" Textured Non Sterile Blue Nitrile Exam Gloves.

Prd Description : Brightway Brand Textured Non Sterile Nitrile Examination Gloves, Powder Free (Blue)

Standard Specification for Nitrile Examination Gloves for Medical Application
Test Method: ASTM D6319-10

Dimension Measurement

Test Method: ASTM D6319-10

Sample tested: 20 pieces per Batch.

Random Sampling based on ISO2859-1:1999; S2 AQL 2.5, Ac=1 Rej=2

Batch No: 6H10 03

Sample No.	Size	Weight (g)	Length (mm)	Width (mm)	Thickness (mm)		
					Finger Tip	Palm	Cuff
1	XS	4.5	240	75	0.16	0.13	0.085
2	XS	4.6	241	76	0.155	0.125	0.09
3	XS	4.6	241	76	0.16	0.13	0.085
4	XS	4.5	240	78	0.165	0.12	0.095
5	S	5.3	242	87	0.155	0.125	0.085
6	S	5.2	239	87	0.155	0.13	0.09
7	S	5.4	242	87	0.165	0.13	0.085
8	S	5.3	242	86	0.155	0.125	0.095
9	M	6.0	240	97	0.16	0.13	0.085
10	M	6.1	239	96	0.155	0.125	0.09
11	M	5.9	241	96	0.16	0.125	0.085
12	M	6.1	240	97	0.165	0.13	0.095
13	L	7.2	240	115	0.155	0.12	0.085
14	L	7.1	241	116	0.16	0.125	0.09
15	L	7.0	241	117	0.165	0.13	0.08
16	L	7.2	240	115	0.155	0.125	0.095
17	XL	7.3	240	122	0.165	0.125	0.085
18	XL	7.4	239	121	0.155	0.12	0.085
19	XL	7.2	242	120	0.16	0.13	0.09
20	XL	7.3	241	122	0.155	0.12	0.095
<i>Specification For XS</i>		$4.6 \pm 0.3g$	<i>Min - 230 Tar - 242 Max - 258</i>	<i>Min - 60 Tar - 70 Max - 80</i>	<i>Min - 0.12 Tar - 0.15 Max - 0.19</i>	<i>Min - 0.08 Tar - 0.12 Max - 0.16</i>	<i>Min - 0.06 Tar - 0.09 Max - 0.13</i>
<i>Specification For Small</i>		$5.2 \pm 0.3g$	<i>Min - 230 Tar - 242 Max - 258</i>	<i>Min - 70 Tar - 80 Max - 90</i>	<i>Min - 0.12 Tar - 0.15 Max - 0.19</i>	<i>Min - 0.08 Tar - 0.12 Max - 0.16</i>	<i>Min - 0.06 Tar - 0.09 Max - 0.13</i>
<i>Specification For Medium</i>		$6.1 \pm 0.3g$	<i>Min - 230 Tar - 242 Max - 258</i>	<i>Min - 85 Tar - 95 Max - 105</i>	<i>Min - 0.12 Tar - 0.15 Max - 0.19</i>	<i>Min - 0.08 Tar - 0.12 Max - 0.16</i>	<i>Min - 0.06 Tar - 0.09 Max - 0.13</i>
<i>Specification For Large</i>		$7.0 \pm 0.3g$	<i>Min - 230 Tar - 242 Max - 258</i>	<i>Min - 100 Tar - 110 Max - 120</i>	<i>Min - 0.12 Tar - 0.15 Max - 0.19</i>	<i>Min - 0.08 Tar - 0.12 Max - 0.16</i>	<i>Min - 0.06 Tar - 0.09 Max - 0.13</i>
<i>Specification For XL</i>		$7.2 \pm 0.3g$	<i>Min - 230 Tar - 242 Max - 258</i>	<i>Min - 110 Tar - 120 Max - 130</i>	<i>Min - 0.12 Tar - 0.15 Max - 0.19</i>	<i>Min - 0.08 Tar - 0.12 Max - 0.16</i>	<i>Min - 0.06 Tar - 0.09 Max - 0.13</i>
<i>Number of non-conformance</i>			<i>0</i>	<i>0</i>	<i>0</i>	<i>0</i>	<i>0</i>
<i>Status</i>			<i>Pass</i>	<i>Pass</i>	<i>Pass</i>	<i>Pass</i>	<i>Pass</i>

7.2. Gloves Physical Properties and Freedom From Holes

Water-tight Test

Test method: ASTM D5151- 06

Sample tested: 200 pieces per Batch.

Random Sampling based on ISO2859-1:1999; G2 AQL 1.0, Ac=5 Rej=6

Location of holes	Size : X - Small	Size : Small	Size : Medium	Size : Large	Size : X - Large
Cuff	0	0	0	0	0
Palm	0	1	0	0	0
Finger-tip	0	0	0	0	0
Total defective Found	1				
Acceptance number	5				
Status	Pass				

BRIGHTWAY HOLDINGS SDN BHD Testing Services Physical Properties Analysis Report

Certificate No : 6H22 168 Date : 22/08/2016
 Customer : N/A
 Brand : 9.5" Textured Non Sterile Blue Nitrile Exam Glove.
 PO No : N/A
 PI No : N/A
 Batch No : 6H10 03
 Product Description : Brightway Brand Textured Non Sterile Nitrile Examination Glove, Powder Free [Blue]
 Sample Receiving Date : 14.08.2016
 Sampling Plan : Single Normal S2 AQL: 2.5 Acc/Rej : 1/2 .
 Sizes Comprising Of : XS, S, M, L & XL
 Test Method : ASTM D 6319 - 10
 Test Conducted on : 21/08/2016
 Aging : 70 +/- 2 deg C for 168 Hrs .

Test Environment Condition		Room Temperature (°C) : 23 RH (%) : 48			Room Temperature (°C) : 23 RH (%) : 48	
Sample No	Size	Before Aging			After Aging	
		Modulus At 500% (MPa)	Tensile Strength (MPa)	Ultimate Elongation (%)	Tensile Strength (MPa)	Ultimate Elongation (%)
1	XS	16.5	23.546	575.0	21.155	495.5
2	XS	18.2	24.852	546.2	22.459	481.2
3	XS	17.4	24.865	537.8	21.398	477.2
4	XS	19.3	23.654	562.3	23.546	499.3
5	S	16.5	24.784	543.7	21.654	487.7
6	S	18.4	23.123	569.5	22.478	454.0
7	S	17.7	23.321	544.4	21.398	451.2
8	S	18.6	24.512	562.3	23.485	486.8
9	M	17.4	25.845	543.7	21.145	495.6
10	M	19.4	24.956	569.5	22.409	485.5
11	M	16.9	23.698	544.4	21.545	479.6
12	M	17.6	24.654	575.4	23.348	496.5
13	L	18.2	25.415	568.7	21.145	492.9
14	L	17.8	23.154	522.2	22.365	457.3
15	L	16.3	24.856	519.4	21.398	489.3
16	L	18.7	23.156	530.5	23.348	487.7
17	XL	18.5	25.852	517.0	21.145	496.9
18	XL	17.8	23.465	596.0	22.409	498.5
19	XL	16.9	25.865	543.1	21.478	488.5
20	XL	17.8	23.456	541.8	23.369	489.7
Average Result		17.8	24.351	550.6	22.134	484.5
Cust. Requirement		N/A	14 (Min)	500 (Min)	14 (Min)	400 (Min)

Gloves meet all the specification listed in ASTM D 6319-10

8. Substantial Equivalence Based on Assessment of Non-Clinical Performance Data

Testing was performed per ASTM D6319-10 , ASTM D5151-06, ASTM D6124- 06,ISO 10993-10:2010 and 16 CFR Part 1500.41. The glove meet standards requirement referenced in section 6.0 above.

Biocompatibility tests indicated that under the conditions of the studies, the gloves were non-sensitizing, non-irritating, and did not reveal any systemic toxicity

9. Summary of Technical Characteristics and Substantial Equivalence Comparison Table

The predicate device in scope is as follows:

- 1) K143477- Nitrile Patient Examination Gloves, Powderfree, Blue Color

BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE, [BLUE] is substantially equivalent in safety and effectiveness to K143477- Nitrile Patient Examination Gloves, Powderfree, Blue Color.

The subject device and predicate device use a similar Nitrile barrier film to achieve a device for the intended use. The properties between the subject device and the predicate device are compared in the following table:

		Device Performance		
Characteristics	Standard	Predicate Device: K143477- Nitrile Patient Examination Gloves, Powderfree, Blue Color	Subject Device : K161215 - BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE, [BLUE]	Result of Comparison
Labeling	N/A	There are no special labeling claims and do not claim gloves as hypoallergenic on labels.	There are no special labeling claims and do not claim gloves as hypoallergenic on labels.	Same
Device Materials	N/A	Nitrile Compound	Nitrile Compound	Same

Device Description		The subject device is a patient examination glove made of synthetic nitrile latex compound. It is non-sterile, powderfree and is Blue in color. The device is ambidextrous and can be worn on either the left or right hand. The device meets ASTM 06319-10: Standard specification for Nitrile Examination Gloves for Medical Application. The subject device is substantially equivalent to legally marketed Nitrile examination gloves identified as Product code LZA. The device is for over-the counter single use.	The subject device in this 510(k) Notification is a Blue Examination Glove. The subject device is a patient examination glove made from a Nitrile compound, Blue in colour, powder free and non sterile (as per 21 CFR 880.6250, class I). The principle operation of the medical device to provide single use barrier protection for the wearer and the device that meets all the requirement specifications for Barrier Protection, tensile properties as defined in ASTM D6319-10 ; Standard specification for Nitrile Examination Gloves. This device is manufactured in facilities compliant to ISO 9001:2008 certified in Manufacture of Non Sterile Natural (Latex) and Synthetic Latex (Nitrile) Examination,	
----------------------------------	--	---	---	--

			Surgical and Industrial Gloves & Nitrile Sheath The device is manufactured to comply with ISO 13485:2003 / EN ISO 13485:2012; Manufacture of Non Sterile Natural (Latex) and Synthetic Latex (Nitrile) Examination & Sterile Surgical Gloves	
Indications for use	N/A	The Nitrile Patient Examination gloves, Powderfree, Blue color, is a disposable device intended for medical purposes that is worn on the examiners' hand or finger to prevent contamination between patient and examiner	BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE, [BLUE] is a disposable device intended for medical purpose, to be worn on the examiners hand or finger to prevent contamination between patient and examiner.	Same
Color	N/A	Blue	Blue	Same
Device Tolerances and Specifications & Performance Data				
Tensile strength : before and after ageing	ASTM 6319-10	> 14Mpa ; PASS	> 14Mpa ; PASS	Same
Ultimate Elongation :	ASTM 6319-10	> 500 % ; PASS	> 500 % ; PASS	Same

Before and after ageing				
Freedom from Pinholes	ASTM 6319-10 (FDA 1000 ml water leak test)	Pass	Pass	Same
Dimensions:				
Length	ASTM 6319-10	> 230 mm (240-400mm)	>230 mm (230-258mm)	Same
Width	ASTM 6319-10	70±10 mm to 120±10 mm (sizes XS to XL)	70±10 mm to 120±10 mm (sizes XS to XL)	Same
Thickness	ASTM 6319-10	> 0.05 mm (palm & finger)	> 0.05 mm (Palm :0.08mm-0.16mm Finger : 0.12mm-0.19mm)	Same
Residual Powder :	ASTM 6319-10, ASTM D6124	Less than 2 mg per glove ; PASS	Less than 2 mg per glove ; PASS	Same
Biocompatibility				
Primary Skin Irritation test	ISO 10993-10	Under conditions of the study, not an irritant	Under conditions of the study, not an irritant	Same
Dermal sensitization assay	ISO 10993-10	Under conditions of the study, not a sensitizer	Under conditions of the study, not a sensitizer	Same
Systemic Toxicity	ISO 10993-11	N/A	Under conditions of study, did not reveal any systemic toxicity	Predicated Device did not undertake the testing. Subject Device met requirements under ISO 109930-11

10.0 Conclusion

Based on intended uses, technological characteristics and Non-Clinical performance data, the subject device; BRIGHTWAY BRAND NITRILE EXAMINATION GLOVES, POWDER FREE, [BLUE] is substantially equivalent to the predicate device K143477- Nitrile Patient Examination Gloves, Powder Free, Blue Color

The subject device meets the requirements of ASTM D 6319- 10 standards as well as applicable 21 CFR references, and meets FDA recognized standards for physical properties requirements, pinhole requirements, biocompatibility requirements which are as shown and discussed above.

Based on the complete list of non-clinical tests, the subject device herein mentioned, is as safe, as effective, and performs as well as the legally marketed predicate devices.