



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

Interson Corporation  
% David Asarnow, Ph.D.  
Quality/Regulatory Manager  
7150 Koll Center Parkway  
PLEASANTON CA 94566

April 13, 2017

Re: K163443

Trade/Device Name: Interson USB Ultrasound System  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulatory Class: II  
Product Code: IYN, IYO, ITX  
Dated: March 17, 2017  
Received: March 17, 2017

Dear Dr. Asarnow:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



For

Robert Ochs, Ph.D.  
Director  
Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K163443

Device Name

Interson USB Ultrasound System

Indications for Use (Describe)

The Interson USB Ultrasound System is intended for diagnostic ultrasound imaging in B, color Doppler, or Combined (B + Color) modes. It is indicated for diagnostic ultrasound imaging and fluid flow analysis in the following applications:

Fetal/Obstetric

Abdominal

Pediatric

Small Organ

Musculo-skeletal (conventional)

Musculo-skeletal (superficial)

Urology

Gynecology

Pelvic Floor

Neuro-muscular

Peripheral Vessel

The system is intended for use by healthcare professionals.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**SECTION 1.3 - INDICATIONS FOR USE**

The Interson USB Ultrasound System is intended for diagnostic ultrasound imaging in B, color Doppler, or Combined (B + Color) modes. It is indicated for diagnostic ultrasound imaging and fluid flow analysis in the following applications:

Fetal/Obstetric  
Abdominal  
Pediatric  
Small Organ  
Musculo-skeletal (conventional)  
Musculo-skeletal (superficial)  
Urology  
Gynecology  
Pelvic Floor  
Neuro-muscular  
Peripheral Vessel

The system is intended for use by healthcare professionals.

The Indications for Use statement on Form 3881 is provided in Tab 4, followed by the clinical applications and modes of operation applicable to each application for each transducer and for the system.

Diagnostic Ultrasound Indications for Use

Interson USB Ultrasound System

**Device Name: GP-C01 Transducer**

Mode: B, color Doppler, B + color Doppler

Intended Use: Diagnostic ultrasound imaging and fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 only)	Specific (Tracks 1&3)	B	M	PWD	CWD	Color Doppler (CD)	Combined (B+CD)	Other* (Specify)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal/Obstetric	N				N	N	
	Abdominal	N				N	N	
	Intra-Operative (Specify)							
	Intra-Operative Neurological							
	Laparoscopic							
	Pediatric (excluding transcranial & neonatal)							
	Small Organ (Scrotum, prostate, lymph nodes, thyroid, breast)	N					N	N
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-card.)							
	Muscular-Skeletal (Conventional)	N					N	N
	Muscular-Skeletal (Superficial)	N					N	N
	Intravascular							
	Other (Urology)	N					N	N
Other (Gynecology)	N					N	N	
Other (Pelvic Floor)	N					N	N	
Other (Neuromuscular)	N					N	N	
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
	Other (Specify)							
Peripheral Vessel	Peripheral vessel							
	Other (Specify)							

N = New Indication, P=Previously cleared by FDA; E=added under this appendix

\*Examples of other modes may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler and Color Velocity Imaging

Diagnostic Ultrasound Indications for Use

Interson USB Ultrasound System

**Device Name: SP-L01 Transducer**

Modes: B, color Doppler, B + color Doppler

Intended Use: Diagnostic ultrasound imaging and fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation							
General (Track 1 only)	Specific (Tracks 1&3)	B	M	PWD	CWD	Color Doppler (CD)	Combined (B+CD)	Other* (Specify)	
Ophthalmic	Ophthalmic								
Fetal Imaging & Other	Fetal/Obstetric								
	Abdominal	N				N	N		
	Intra-Operative (Specify)								
	Intra-Operative Neurological								
	Laparoscopic								
	Pediatric (excluding transcranial & neonatal)	N				N	N		
	Small Organ (Scrotum, prostate, lymph nodes, thyroid, breast)	N				N	N		
	Neonatal Cephalic								
	Adult Cephalic								
	Trans-rectal								
	Trans-vaginal								
	Trans-urethral								
	Trans-esoph. (non-card.)								
	Muscular-Skeletal (Conventional)	N					N	N	
	Muscular-Skeletal (Superficial)	N					N	N	
	Intravascular								
	Other (Urology)	N					N	N	
	Other (Gynecology)	N					N	N	
Other (Pelvic Floor)	N					N	N		
Other (Neuromuscular)	N					N	N		
Cardiac	Cardiac Adult								
	Cardiac Pediatric								
	Intravascular (Cardiac)								
	Trans-esoph. (Cardiac)								
	Intra-cardiac								
	Other (Specify)								
Peripheral Vessel	Peripheral vessel	N				N	N		
	Other (Specify)								

N = New Indication, P=Previously cleared by FDA; E=added under this appendix

\*Examples of other modes may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler and Color Velocity Imaging

Diagnostic Ultrasound Indications for Use

**Device Name: Interson USB Ultrasound System**

Modes: B, color Doppler, B + color Doppler

Intended Use: Diagnostic ultrasound imaging and fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation							
General (Track 1 only)	Specific (Tracks 1&3)	B	M	PWD	CWD	Color Doppler (CD)	Combined (B+CD)	Other* (Specify)	
Ophthalmic	Ophthalmic								
Fetal Imaging & Other	Fetal/Obstetric	N				N	N		
	Abdominal	N				N	N		
	Intra-Operative (Specify)								
	Intra-Operative Neurological								
	Laparoscopic								
	Pediatric (excluding transcranial & neonatal)	N				N	N		
	Small Organ (Scrotum, prostate, lymph nodes, thyroid, breast)	N				N	N		
	Neonatal Cephalic								
	Adult Cephalic								
	Trans-rectal								
	Trans-vaginal								
	Trans-urethral								
	Trans-esoph. (non-card.)								
	Muscular-Skeletal (Conventional)	N					N	N	
	Muscular-Skeletal (Superficial)	N					N	N	
	Intravascular								
	Other (Urology)	N					N	N	
Other (Gynecology)	N					N	N		
Other (Pelvic Floor)	N					N	N		
Other (Neuromuscular)	N					N	N		
Cardiac	Cardiac Adult								
	Cardiac Pediatric								
	Intravascular (Cardiac)								
	Trans-esoph. (Cardiac)								
	Intra-cardiac								
Other (Specify)									
Peripheral Vessel	Peripheral vessel	N				N	N		
	Other (Specify)								

N = New Indication, P=Previously cleared by FDA; E=added under this appendix

\*Examples of other modes may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler and Color Velocity Imaging

**510(k) Summary K163443**

Provided in accordance with 21CFR 807.92 (c).

**Submitter Information- 21 CFR 807.92 (a)(1)**

Date of submission: 11/28/2016

Submitter information: Interson Corporation  
7150 Koll Center Parkway  
Pleasanton, CA 94566

Establishment Registration No: 2939830

Contact person: David Asarnow  
Quality/Regulatory Manager  
(925) 462-4948 (tel)  
(925) 462-4833 (fax)  
dasarnow@interson.com

**Name of Device and Classification – 21 CFR 807.92 (a)(2)**

Device trade name: USB Ultrasound System  
Common name: Diagnostic ultrasound system and transducers  
Classification: Class II

<u>21 CFR Section</u>	<u>Classification Name</u>	<u>Product Code</u>
892.1550	Ultrasonic pulsed doppler imaging system	IYN
892.1560	Ultrasonic pulsed echo imaging system	IYO
892.1570	Diagnostic ultrasonic transducer	ITX

**Track**

Track 3

**Predicate Device – 21 CFR 807.92 (a)(3)**

Philips Lumify Diagnostic Ultrasound System	K152899	10/30/2015
Interson USB Ultrasound Probe System	K070907	4/17/2007

**Device Description – 21 CFR 807.92 (a)(4):**

The Interson USB Ultrasound Imaging System is a self-contained portable, solid-state, multiple-mode and multiple-application ultrasound imaging system. The system comprises a series of handheld probes containing an ultrasound generator/receiver, analog to digital converter, microcontroller, control logic, USB 2.0 interface, and



software-based controls offering a full complement of conventional operating modes, parameter controls, and recording functions. The selection of transducers offered with the system permit a wide range of clinical applications including Fetal/Obstetric, Abdominal, Pediatric, Small Organ, Musculo-skeletal, Urology, Gynecology, Pelvic Floor, Neuro-muscular, and Peripheral Vessel.

The initial operational settings for each transducer are preprogrammed in the system. User-customized parameter settings for each transducer may be set by the operator and stored for recall as needed via the system control panel. Customization includes transmit power, images controls selection, and Time Gain Compensation (TGC). Controls are also provided to select display format and to utilize the Cine function.

The Interson USB Ultrasound transducers operate in B-Mode, color Doppler, and combined mode (B + color Doppler), providing high resolution and high penetration performance. The product family uses curved and flat linear array transducers and includes a General Purpose (GP) probe and a Small Parts (SP) probe.

The SeeMore software application used in the Interson USB Ultrasound System supports static image acquisition and a cine function capable of storing up to 32 to 1024 sequential images. Management of patient history is enabled by the image or cine-storage function.

The SeeMore control software runs on standard Microsoft Windows platforms.

#### **Intended Use/Indications for Use – 21 CFR 807.92 (a)(5)**

The Interson Ultrasound Imaging System is intended for diagnostic ultrasound imaging in B, color Doppler, or Combined (B + Color) modes. It is indicated for diagnostic ultrasound imaging and fluid flow analysis in the following applications:

- Fetal/Obstetric
- Abdominal
- Pediatric
- Small Organ
- Musculo-skeletal
- Urology
- Gynecology
- Pelvic Floor
- Neuro-muscular
- Peripheral Vessel

The system is intended for use by healthcare professionals.

**Product Models**

Model	Common Applications
Transducer Model GP- C01	Fetal/Obstetric, Abdominal, Urology, Musculo-skeletal, Small Organs, Neuromuscular
Transducer Model SP-L01	Pediatric, Small Organs, Peripheral Vessel

All models are used with SeeMore control software.

**Summary of technological characteristics of the device compared to the predicate device- 21 CFR 807.92 (a)(6)**

Device Features	Subject Device: Interson Ultrasound System	Predicate 1: Philips Lumify K152899	Predicate 2: Interson Ultrasound System K070907
Intended Use	Diagnostic ultrasound imaging in B, color Doppler and Combined (B + Color) modes.	Diagnostic ultrasound imaging in B, color Doppler, and Combined (B + Color) modes.	Diagnostic ultrasound imaging in B mode (all transducers), A mode (ophthalmic)
Indications for Use	Indicated for diagnostic ultrasound imaging and fluid flow analysis in specified applications	Indicated for diagnostic ultrasound imaging and fluid flow analysis in specified applications	Indicated for diagnostic ultrasound imaging in specified applications
Array Geometry	Curved and linear	Curved and linear	Non-Array
Mechanics	Solid State	Solid State	Mechanical
Software platform	Commercial off-the-shelf operating system (Windows)	Commercial off-the-shelf operating system (Android)	Commercial off-the-shelf operating system (Windows)
Software control	Standalone	Standalone	Standalone
Measurement function	2D measurement and area measurement	2D measurement tool	2D measurement tool
Wireless networking	Not supported	Supported	Not supported
Connector	USB	USB	USB

## **Determination of Substantial Equivalence**

### **Non-clinical Performance Data**

Non-clinical testing relied on in this premarket notification submission for a determination of substantial equivalence include tests which show compliance with the following standards:

1. Recognition Number 12-105: NEMA UD 2-2004 (R2009), Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment Revision 3. (Radiology)
2. Recognition Number 19-4: AAMI ANSI ES60601-1:2005/(R)2012 and A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012, (Consolidated Text) Medical electrical equipment - Part 1: General requirements for basic safety and essential performance (IEC 60601-1:2005, MOD)
3. Recognition Number 12-293: IEC 60601-2-37 (2015), Amendment 2, Medical electrical equipment - Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment. (Radiology)
4. Recognition Number 19-8: IEC 60601-1-2 Edition 4.0 2014-02, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests
5. Recognition Number 5-89: IEC 60601-1-6 Edition 3.1 2013-10, Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
6. Recognition Number 5-40: ISO 14971 Second edition 2007-03-01, Medical devices - Application of risk management to medical devices
7. Recognition Number 2-173: AAMI ANSI ISO 10993-10:2010/(R)2014, Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
8. Recognition Number 2-153: AAMI ANSI ISO 10993-5:2009/(R)2014, Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity

### **Summary of Clinical Tests**

The Interson USB Ultrasound System product family introduces no new modes, features, or technologies relative to the predicate devices (Philips Lumify K152899 and Interson K070907) that require clinical testing. The clinical safety and effectiveness of ultrasound systems with these characteristics are well accepted for both predicate and subject devices.

### **514 Performance Standards**

There are no Sec. 514 performance standards for this device.

Prescription Status

This is a prescription device. The prescription device statement appears in the labeling.

Sterilization Site(s)

Not applicable. No components are supplied sterile.

Conclusions

Interson Corporation concludes that the subject device, the Interson Ultrasound USB Imaging System, has been shown to be substantially equivalent to the predicate devices identified above.