



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

May 25, 2017

Stryker Corporation
David Colao
Senior Manager Regulatory Affairs
4100 E. Milham Ave
Kalamazoo, Michigan 49001

Re: K170242

Trade/Device Name: MultiGen™ 2 RF Generator System
Regulation Number: 21 CFR 882.4400
Regulation Name: Radiofrequency Lesion Generator
Regulatory Class: Class II
Product Code: GXD, GXI
Dated: April 25, 2017
Received: April 26, 2017

Dear Mr. Colao:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-

related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

Michael J. Hoffmann -S

for Carlos L. Peña, PhD, MS
Director
Division of Neurological
and Physical Medicine Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K170242

Device Name

Stryker MultiGen™ 2 RF Generator System

Indications for Use (Describe)

Stryker MultiGen™ 2 RF Generator

The Stryker MultiGen™ 2 RF Generator, in combination with the Stryker MultiGen 2 Electrodes and RF Cannulae, are intended for coagulation of soft tissues in orthopedic, spinal, and neurosurgical applications. Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.

Venom MultiGen 2 Electrodes

The Stryker Venom MultiGen 2 Electrodes, in combination with the Stryker RF Cannulae and Stryker MultiGen™ 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. These products are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerve, and nerve roots for the relief of pain.

Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.

Standard MultiGen 2 Standard Electrodes

The Stryker MultiGen 2 Electrodes, in combination with the Stryker RF Cannulae and Stryker MultiGen™ 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. These products are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerve, and nerve roots for the relief of pain.

Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.

Venom and Standard Cannulae (RF Cannulae)

The Stryker RF Cannulae, in combination with Stryker Electrodes and the Stryker RF Generator, Stryker MultiGen, or Stryker MultiGen™ 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. They are also used for selective denervation and destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerve, and nerve roots for the relief of pain.

Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

SECTION
510(k) Summary

Table of Contents

510(k) Summary	3
Submitter.....	3
Device	3
Predicate Devices	3
Device Description	3
Indications for Use	4
Comparison of Technological Characteristics with the Predicate Device.....	5
Performance Data.....	14
Conclusions	15

List of Figures and Tables

Table 1: Generator Comparison Table	8
Table 2: Electrode and Cannula Comparison Table	13

510(k) Summary

Submitter

Stryker Instruments
4100 E Milham Ave
Kalamazoo, MI 49001
(269) 389-4285

Contact Person: David Colao
Date Prepared: May 23, 2017

Device

Name of Device: MultiGen 2 RF Generator system
Common or Usual Name: RF Generator
Classification Regulation: 882.4400 / 882.4725
Panel: Neurology
Product Code: GXD / GXI

Predicate Devices

The following devices are predicates for the system:
Stryker RF Multi-lesion Generator (K071482)
Stryker RF Electrodes and Cannulae (K032406)
Stryker Venom Electrodes and Cannulae (K123178)

The following reference device was used in this submission:
Stryker Interventional Pain RF Generator (K032601)

Device Description

The MultiGen™ 2 RF Generator system is a bipolar and monopolar, high frequency electrosurgical system. The MultiGen 2 RF Generator will be used in conjunction with the MultiGen 2 Splitter Cable, MultiGen 2 Electrodes, RF Cannulae, and other accessories to produce lesions by the direct application of radiofrequency currents. The generator applies temperature-controlled, radiofrequency (RF) energy into the electrode. During lesion creation, targeted nerve tissue is exposed to RF energy using an active electrode inserted into an insulated cannula with an uninsulated tip.

The application of RF energy causes a thermal reaction at the targeted nerve tissue site to create a lesion.

Indications for Use

Stryker MultiGen™ 2 RF Generator

The Stryker MultiGen™ 2 RF Generator, in combination with the Stryker MultiGen 2 Electrodes and RF Cannulae, are intended for coagulation of soft tissues in orthopedic, spinal, and neurosurgical applications. Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.

Venom MultiGen 2 Electrodes

The Stryker Venom MultiGen 2 Electrodes, in combination with the Stryker RF Cannulae and Stryker MultiGen™ 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. These products are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerve, and nerve roots for the relief of pain.

Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.

Standard MultiGen 2 Standard Electrodes

The Stryker MultiGen 2 Electrodes, in combination with the Stryker RF Cannulae and Stryker MultiGen™ 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. These products are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerve, and nerve roots for the relief of pain.

Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.

Venom and Standard Cannulae (RF Cannulae)

The Stryker RF Cannulae, in combination with Stryker Electrodes and the Stryker RF Generator, Stryker MultiGen, or Stryker MultiGen™ 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. They are also used for selective denervation and destruction

procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerve, and nerve roots for the relief of pain.

Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.

Comparison of Technological Characteristics with the Predicate Device

The predicate device is a previous generation of Stryker radiofrequency generator, the Stryker RF Multi-Lesion Generator (MultiGen). The intended uses of the subject and predicate devices are similar. The systems are intended for coagulation of soft tissues in orthopedic, spinal, and neurosurgical applications using radiofrequency energy. Minor modifications were made to the indication for use statement to provide clarification and consistency, but do not change the meaning of the indication for use statement.

The technological characteristics of the subject and predicate devices are equivalent. Both the subject device and the primary predicate device use same modes of operation, types of control, energy delivery and type of energy output. None of the changes alter the operating principle, modes of operation, temperature range. Both systems use the same type of accessories (a RF generator, connecting cables, electrodes, cannulae, and neutral electrode). The user interaction with the devices are similar. As demonstrated by the performance testing, the subject device has similar performance specifications as the predicate device.

Based on the comparison of intended use and technological characteristics, the device is similar to the predicate device. The hardware and software verification and validation testing demonstrate that the subject device meets its performance specifications and will perform as intended in the specified use conditions and that any differences between the subject device and predicate device do not raise new questions of safety and effectiveness. Therefore, the subject device can be found substantially equivalent to the predicate device.

FEATURE	<i>MultiGen™ 2 RF Generator</i> SUBJECT DEVICE	<i>Stryker RF Multi-Lesion Generator</i> PREDICATE DEVICE	Equivalency
Generator			
510(k) Clearance	Subject Device of Current Submission	K071482	N/A.
Indication for Use	The Stryker MultiGen™ 2 RF Generator, in combination with the Stryker MultiGen 2 Electrodes and Cannulae are intended for coagulation of soft tissues in orthopedic, spinal, and neurosurgical applications. Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.	The Stryker Interventional Pain RF MultiGen, in combination with the Stryker RF Electrodes and Cannulae are intended for coagulation of soft tissues in orthopedic, spinal, and neurosurgical applications. Examples include, but are not limited to: Facet Denervation, Percutaneous Chordotomy/Dorsal Root Entry Zone (DREZ) Lesion, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy. The Stryker Interventional Pain RF MultiGen in combination with a Smith & Nephew SPINECATH™ & Acutherm™ catheters are intended for coagulation and decompression of disc material to treat symptomatic patients with annular disruption of contained herniated discs.	Equivalent. Indications within scope of predicate indications.
Modes of Operation	Sensory and Motor Stimulation. Pulsed and Continuous Lesioning	Sensory and Motor Stimulation. Pulsed and Continuous Lesioning	Identical.
Type of Control	Manual and Automatic	Manual and Automatic	Identical.
Temperature measurement technology	K-type thermocouple in electrode Thermocouple interface electronics in RF Generators	K-type thermocouple in electrode Thermocouple interface electronics in RF Generators	Identical.
Temperature measurement location	Distal Tip of the Nitinol Electrode	Distal Tip of the Nitinol Electrode	Identical.
Temperature Range	37 °C - 95 °C	37 °C - 97 °C	Equivalent. Within range of predicate.
Temperature Accuracy	37 °C - 95 °C +2°C	37 °C - 95 °C +2°C	Identical
Lesion Time	0-999 seconds	0-999 seconds	Identical.

FEATURE	<i>MultiGen™ 2 RF Generator</i> SUBJECT DEVICE	<i>Stryker RF Multi-Lesion Generator</i> PREDICATE DEVICE	Equivalency
Stimulation Frequency	Sensory = 50 Hz Motor = 2 Hz	Sensory: 10, 20, 30, 50, 75, 100, 150, 200 (Hz) Motor: 1, 2, 3, 4, 5, 7, 10, 20 (Hz)	Equivalent. Reduce complexity. Within range of predicate.
Stimulation Pulse Width	1 ms	.1,.2,.5,.75,1,2,3,5 (ms)	Similar. Reduce complexity. Within range of predicate.
Stimulation Amplitude Voltage Regulation Mode	0-10 V peak; accuracy + 10%, current limited to 25 mA	0-10 V peak; accuracy + 10%, current limited to 40 mA	Equivalent. Reduction in voltage.
Impedance Range	0-2000 ohms	0-2000 ohms	Identical.
Type of User Interface	Touch screen graphical user interface	Touch screen graphical user interface	Identical.
Type of Probe Recognition	Automatic	Automatic	Identical.
Number of Amp	4	1	Equivalent. Simplifies control of energy.
Energy delivery during multi-channel RF treatment	Sequential non simultaneous energy delivery	Sequential non simultaneous energy delivery	Identical.
Output Energy	0-100 watts	0-50 watts	Equivalent. The maximum current the patient is exposed to remains the same (700 mA).
Output Waveform(s)	500kHz sinusoid	Malis Dual/Wave Waveform	Equivalent. Sinusoid waveform requires lower voltage to achieve equivalent energy transfer.
Monopolar/Bipolar	Monopolar and Bipolar	Monopolar and Bipolar	Identical.
Ability to perform multiple lesions simultaneously	Yes	Yes	Identical.

FEATURE	<i>MultiGen™ 2 RF Generator</i> SUBJECT DEVICE	<i>Stryker RF Multi-Lesion Generator</i> PREDICATE DEVICE	Equivalency
Number of electrode connections	4	4	Identical.
Printer Output	No	Yes	Equivalent. Feature is not frequently used.
Accessories	• Splitter Box • 3M Neutral Electrode (Ground Pad) • (Neutral Electrode) Ground Pad Cable • Stryker MultiGen 2 (Standard) Electrodes • Stryker RF (Standard) Cannulae • Venom MultiGen 2 Electrodes • Venom Cannulae	• MultiGen Cable • Connector Cable • Hand Controller • 3M Neutral Electrode (Ground Pad) • (Neutral Electrode) Ground Pad Cable • Stryker RF (Standard) Electrodes • Stryker RF (Standard) Cannulae • Venom Electrodes • Venom Cannulae	Similar. The Splitter Cable replaces the need for two MultiGen Cables or other connecting cables.

TABLE 1: GENERATOR COMPARISON TABLE

FEATURE	<i>Standard and Venom MultiGen 2 Electrodes</i> SUBJECT DEVICE	<i>Standard Electrodes</i> PREDICATE DEVICE	<i>Venom Electrodes</i> PREDICATE DEVICE	Equivalency
Intended Use				
Electrodes				
510(k) Clearance	Subject Device	K032406	K123178	N/A.
Indication for Use (Electrodes)	Standard MultiGen 2 Electrodes The Stryker MultiGen 2 Electrodes, in combination with the Stryker RF Cannulae and Stryker MultiGen 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. These products are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerve, and nerve roots for the relief of pain. Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy. Venom MultiGen 2 Electrodes The Stryker Venom MultiGen 2 Electrodes, in combination with the Stryker RF Cannulae and Stryker MultiGen 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. These products are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and	The Stryker RF Electrodes and Cannulae, in combination with the Stryker RF Generator, are intended for coagulation of soft tissues in orthopedic, arthroscopic, spinal, and neurosurgical applications. They are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the spinal cord, peripheral nerves, and nerve roots for the relief of pain. Examples include, but are not limited to, Facette Denervation, Percutaneous Chordotomy/Dorsal Root Entry Zone (DREZ) Lesion, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.	The Stryker RF electrodes and cannulae, in combination with the Stryker RF Generator/Multigen, are intended for coagulation of soft tissues in orthopedic, spinal, and neurosurgical applications. These products are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerves, and nerve roots for the relief of pain. Examples include, but are not limited to, Facette Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.	Equivalent. Indications within scope of predicate indications.

FEATURE	<i>Standard and Venom MultiGen 2 Electrodes</i> SUBJECT DEVICE	<i>Standard Electrodes</i> PREDICATE DEVICE	<i>Venom Electrodes</i> PREDICATE DEVICE	Equivalency
	cervical regions of the peripheral nerve, and nerve roots for the relief of pain. Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.			
RF Cannulae				
Indication for Use (Cannulae)	Standard RF Cannulae / Venom RF Cannulae The Stryker RF Cannulae, in combination with Stryker Electrodes and the Stryker RF Generator, Stryker MultiGen, or Stryker MultiGen 2 RF Generator, are intended for coagulation of soft tissue in orthopedic, spinal, and neurosurgical applications. They are also used for selective denervation and destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerve, and nerve roots for the relief of pain. Examples include, but are not limited to: Facet Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.	The Stryker RF Electrodes and Cannulae, in combination with the Stryker RF Generator, are intended for coagulation of soft tissues in orthopedic, arthroscopic, spinal, and neurosurgical applications. They are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the spinal cord, peripheral nerves, and nerve roots for the relief of pain. Examples include, but are not limited to, Facette Denervation, Percutaneous Chordotomy/Dorsal Root Entry Zone (DREZ) Lesion, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.	The Stryker RF electrodes and cannulae, in combination with the Stryker RF Generator/Multigen, are intended for coagulation of soft tissues in orthopedic, spinal, and neurosurgical applications. These products are also used for selective denervation and tissue destruction procedures which may be performed on the lumbar, thoracic, and cervical regions of the peripheral nerves, and nerve roots for the relief of pain. Examples include, but are not limited to, Facette Denervation, Trigeminal Neuralgia, Peripheral Neuralgia and Rhizotomy.	Equivalent. Indications within scope of predicate indications.

FEATURE	<i>Standard and Venom MultiGen 2 Electrodes</i> SUBJECT DEVICE	<i>Standard Electrodes</i> PREDICATE DEVICE	<i>Venom Electrodes</i> PREDICATE DEVICE	Equivalency
Technological Indications				
Electrodes				
Electrode Lengths	Standard 50mm, 100mm, 150mm, 200mm Venom 100mm, 150mm	50mm, 100mm, 150mm,	100mm, 150mm	Equivalent. Additional length.
Nitinol Electrode Radial Dimensions	Standard Inside Diameter: 0.25mm Outside Diameter: 0.40mm Venom Inside Diameter: 0.25mm Outside Diameter: 0.50mm	Inside Diameter: 0.25mm Outside Diameter: 0.40mm	Inside Diameter: 0.25mm Outside Diameter: 0.50mm	Identical.
Design	Standard 27 gauge (0.4mm OD) Nitinol electrode PEEK electrode hub 4 conductor silicone cable (approximately 5 feet long) K-type thermocouple ODU cable connector assembly (diameter: 8.9mm, length 34.9 mm) Venom 25 gauge (0.5mm OD) Nitinol electrode PEEK electrode hub 4 conductor silicone cable (approximately 5 feet long) K-type thermocouple ODU cable connector assembly (diameter: 8.9mm, length 34.9 mm)	27 gauge (0.4mm OD) Nitinol electrode PEEK electrode hub 4 conductor silicone cable (approximately 4 feet long) K-type thermocouple ODU cable connector assembly (diameter: 13.9mm, length 47.5 mm)	25 gauge (0.5mm OD) Nitinol electrode PEEK electrode hub 4 conductor silicone cable (approximately 4 feet long) K-type thermocouple (diameter: 13.9mm, length 47.5 mm)	Equivalent.

FEATURE	<i>Standard and Venom MultiGen 2 Electrodes</i> SUBJECT DEVICE	<i>Standard Electrodes</i> PREDICATE DEVICE	<i>Venom Electrodes</i> PREDICATE DEVICE	Equivalency
Temperature Measurement Technology	K-type thermocouple in electrode Thermocouple interface electronics in RF Generators	K-type thermocouple in electrode Thermocouple interface electronics in RF Generators	K-type thermocouple in electrode Thermocouple interface electronics in RF Generators	Identical.
Temperature Measurement Location	Distal Tip of the Electrode	Distal Tip of the Electrode	Distal Tip of the Electrode	Identical.
Sterility	Steam or Vaporized Hydrogen Peroxide	Steam	Steam	Equivalent.
Materials	Nitinol Silicone Rubber PEEK(Polyether ether ketone) Stainless Steel	Nitinol Silicone Rubber PEEK(Polyether ether ketone) Stainless Steel	Nitinol Silicone Rubber PEEK(Polyether ether ketone) Stainless Steel	Identical.
Cannulae				
Sizes	<u>Lengths</u> Standard 50mm, 100mm, 150mm, 200mm Venom 100mm, 150mm <u>Diameters</u> Standard 18G, 20G, 22G Venom 18G, 20G, <u>Active Tip</u> Standard 5.0mm, 7.5mm, 10.0mm, 15.0mm Venom 10mm	<u>Lengths</u> 50mm, 100mm, 150mm <u>Diameters</u> 20G, 22G <u>Active Tip</u> 2.5mm, 4.0mm, 5.0mm, 7.5mm, 10.0mm, 15.0mm	<u>Lengths</u> 100mm, 150mm <u>Diameters</u> 18G, 20G <u>Active Tip</u> 10mm	Equivalent. Longer Length. Gages and Active Tips within cleared range.

FEATURE	<i>Standard and Venom MultiGen 2 Electrodes</i> SUBJECT DEVICE	<i>Standard Electrodes</i> PREDICATE DEVICE	<i>Venom Electrodes</i> PREDICATE DEVICE	Equivalency
Design	Standard The Standard RF Cannula consists of 18, 20 or 22 gauge stainless steel tubing cut to length. A bevel tip is created via an electro-chemical grinding process. The cannula is insulated with siliconised polyester. A polycarbonate cannula hub with inert black ink marking is attached to the proximal end. Venom The Venom™ Cannula consists of 18 or 20 gauge stainless steel tubing cut to length. A bevel tip and a side port are created at the distal end via an electro-chemical grinding process. The cannula is insulated with siliconised polyester. A polycarbonate cannula hub with inert black ink marking is attached to the proximal end.	The Standard RF Cannula consists of 20 or 22 gauge stainless steel tubing cut to length. A bevel tip is created via an electro-chemical grinding process. The cannula is insulated with siliconised polyester. Additionally the cannulae include a small amount of liquid silicone. A polycarbonate cannula hub with inert black ink marking is attached to the proximal end.	The Venom™ Cannula consists of 18 or 20 gauge stainless steel tubing cut to length. A bevel tip and a side port are created at the distal end via an electro-chemical grinding process. The cannula is insulated with siliconised polyester. A polycarbonate cannula hub with inert black ink marking is attached to the proximal end.	Equivalent.
Sterility	EtO	EtO	EtO	Identical.
Materials	Stainless Steel Siliconised Polyester Polycarbonate plastic Inert black ink Polyethylene	Stainless Steel Siliconised Polyester Liquid Silicone Polycarbonate plastic Inert black ink Polyethylene	Stainless Steel Siliconised Polyester Polycarbonate plastic Inert black ink Polyethylene	Equivalent.

TABLE 2: ELECTRODE AND CANNULA COMPARISON TABLE

Performance Data

The following performance data were provided in support of the substantial equivalence decision:

Electrical Safety and Electromagnetic Compatibility (EMC) Testing

Electrical safety and EMC testing were conducted on the subject device in accordance with the following standards:

- AAMI / ANSI ES60601-1:2005/(R)2012 and C1:2009/(R)2012 and, A2:2010/(R)2012 (IEC 60601-1:2005, mod).
- IEC 60601-2-2 Edition 5.0 2009-02 [including: technical corrigendum 1 (2014)]
- IEC 60601-1-2 Edition 3: 2007-03 and
- IEC 60601-1-2 Edition 4.0 2014-02

Software Verification and Validation Testing

Software verification and validation testing were conducted and documentation was provided as recommended by FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices."

Performance Testing – Bench

The following design verification activities have been performed to ensure the correct functionality of the system as it has been specified:

- Substantial equivalence testing verifying the radiofrequency ablations created by the subject device are equivalent to the radiofrequency ablations created by the predicate device.
- Multi-Lesion testing verifying that the MultiGen 2 RF Generator system functions properly with all possible set ups and procedural variations.
- MultiGen 2 Splitter Cable (Splitter Cable) life testing.
- MultiGen 2 Electrode NOVRAM testing showing that the NOVRAM meets specifications.
- MultiGen 2 RF Generator current limitation testing.

Animal Study

No animal studies were performed to support substantial equivalence.

Clinical Studies

No clinical studies were performed to support substantial equivalence.

Conclusions

Based on the comparison of intended use and technological characteristics, the device is similar to the predicate device. The hardware and software verification and validation testing demonstrate that the subject device meets its performance specifications and will perform as intended in the specified use conditions and that any differences between the subject device and predicate device do not raise new questions of safety and effectiveness. Therefore, the subject device can be found substantially equivalent to the predicate device.