



November 22, 2017

Intrinity Global Limited
% Elly Xu
Consultant
Shenzhen Joyantech Consulting Co., Ltd
NO. 55 Shizhou middle road, Nanshan District
Shenzhen, 518000 China

Re: K170662
Trade/Device Name: Non Contact Infrared Forehead Thermometer
Regulation Number: 21 CFR 880.2910
Regulation Name: Clinical Electronic Thermometer
Regulatory Class: Class II
Product Code: FLL
Dated: October 15, 2017
Received: October 26, 2017

Dear Elly Xu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

 Tina Kiang
-S

Tina Kiang, Ph.D.
Acting Director
Division of Anesthesiology,
General Hospital, Respiratory,
Infection Control, and Dental Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K170662

Device Name

Non Contact Infrared Forehead Thermometer

Indications for Use (Describe)

Non Contact Infrared Forehead Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in clinic as reference. It is intended for measuring human body temperature of all ranges of people by detecting infrared heat from the forehead.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

VOL_05: 001 510(k) Summary

K170662

5.1 Administrative Information**Date of Summary prepared**

Nov, 14, 2017

Manufacturer information

Intrinity Global Limited

Unit 1707, Floor 17, Westley Square, 48 Hoi Yuen Road, Kwun Tong, Hong Kong

Contact person: Alex Kwok

Phone: 852-35231874

Fax: 852-37484128

E-mail: alex.kwok@intrinity.com.hk

Submission Correspondent

卓远天成

Shenzhen Joyantech Consulting Co., Ltd.

Address: Room 1122, International Mayors Communication Centre, NO. 55 Shizhou middle road , Nanshan District, Shenzhen

Contact person: Ms. Elly Xu; Mr. Field Fu

E-Mail: elly@cefd.com; cefd13485@163.com**Establishment registration number****5.2 Device Information****Type of 510(k) submission:**

Traditional

Trade Name:

Non Contact Infrared Forehead Thermometer

Model:

TVT-200, TVT-200 PLUS

Classification name:

thermometer, electronic, clinical

Review Panel:

General Hospital

Product Code:

FLL

Device Class:

II

Regulation Number: | 880.2910

5.3 Predicate Device Information

Sponsor: | KAZ USA, Inc.
Device: | No Touch + Forehead Thermometer
510(K) Number: | K134043

5.4 Device Description

The thermometer measures the temperatures of people by detecting the energy.

The external probe plays an important role in the measuring process.

As soon as the thermometer is placed near the body and the radiation sensor is activated, the measurement will be taken instantly by detection of the infrared heat.

The thermometer includes two models: TVT-200, TVT-200 PLUS. The TVT-200 model includes four colors: pink, grey, orange and purple. While TVT-200 PLUS has only orange color.

The Non Contact Infrared Forehead Thermometer is composed by an IR sensor, a Human temperature measurement button, an object temperature button, Battery compartment, Buzzer, °C /°F button, a LCD and an enclosure.

The functions of TVT-200 and TVT-200 PLUS are the same, except the Back Light. The detailed functions see below table:

Models	Functions							
	Memory	°C/°F	Scan Mode	LCD	Buzzer	Fever alert	Power off/ Auto Off	Back Light
TVT-200	•	•	•	•	•	•	•	•
TVT-200 PLUS	•	•	•	•	•	•	•	○

5.5 Intended Use/ Indications for Use

Non Contact Infrared Forehead Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in

clinic as reference. It is intended for measuring human body temperature of all ranges of people by detecting infrared heat from the forehead.

5.6 Technological characteristics of the subject device compared to the predicate device

Items	Predicate Device (K134043), KAZ	Subject Device	Remarks
Indications for use	The No Touch + Forehead Thermometer (Model NTF3000US) is a non-sterile, reusable clinical thermometer intended for the intermittent determination of human body temperature in a touch and no touch on the centre of the forehead as the measurement site on people of all ages.	Non Contact Infrared Forehead Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in clinic as reference. It is intended for measuring human body temperature of all ranges of people by detecting infrared heat from the forehead.	Equivalent: intended use remains unchanged
Measurement method	Infrared radiation detection	Infrared radiation detection	Same
Measurement mode	Forehead measure mode	Forehead measurement mode	Same
Measuring range	Forehead temperature mode: 34.4°C -42.2°C (93.9°F - 108°F)	Forehead measurement mode: 32° to 43°C (89.6°F to 109.4°F)	Similar (Note01)
Display resolution	0.1°C (0.1°F)	0.1°C/0.1°F	Same
C/F switchable	Yes	Yes	Same
Measuring accuracy	Forehead temperature mode: ±0.2°C (0.4°F)	Forehead measurement mode: ±0.2°C (0.4°F)	Same
Display	LCD display	LCD display	Same
Measurement distance	0-5cm	1cm	Similar (Note02)
Memory	Not available.	16 sets.	(Note03)

Items	Predicate Device (K134043), KAZ	Subject Device	Remarks
Power source	Two AA batteries	One 1.5V AAA alkaline battery	Similar (Note04)
Low battery indication	Yes	Yes	Same
Waterproof	No, IP20	No, IP22	Similar (Note05)
Operating condition	15 °C – 40 °C (59 °F – 104 °F), 15–95% non condensing	15°C~ 40°C (59°F~ 104°F); ≤95% RH	Similar (Note06)
Patient contact materials	metals and resins	ABS with colorants (pink, grey, orange and purple), Glass and Metal	Different (Note07)
Cleaning/ disinfection	The thermometer casing and the measuring probe are cleaned and disinfected by alcohol (70% Isopropyl).	The thermometer casing and the measuring probe are cleaned and disinfected by 70% alcohol.	Different (Note08)
Biocompatibility	Comply with ISO 10993-5:2009, ISO 10993-10:2010	Comply with ISO 10993-5:2009, ISO 10993-10:2010	Same
Electric Safety and EMC	IEC 60601-1 3rd edition: 2005, IEC 60601-1-2: 2007, IEC 60601-1-11: 2010	IEC 60601-1: 2005+CORR.1 (2006)+ CORR.2 (2007), IEC 60601-1-2: 2014, IEC 60601-1-11: 2010, ISO 80601-2-56: 2009.	Similar
Performance	ASTM E1965-98 (2009)	ASTM E1965-98 (2009), ISO 80601-2-56: 2009.	Similar

Note 01:

The measuring range of subject device meet the minimum rated output range of clinical thermometer requirement, from 35°C to 42°C, which is stated in standard ISO 80601-2-56. The differences between the predicate device and subject device will not affect the safety and effectiveness of the subject device.

Note 02:

The difference in measurement distance will not affect the safety and effectiveness of the subject device. And the use method are described detailed in IFU.

Note 03:

The memory capacity will not affect the safety and effectiveness of the subject device.

Note 04:

The main difference between the batteries of the Predicate and the subject device is the size, which will not affect the safety and effectiveness.

Furthermore, the electric safety of the battery has been validated with the subject device according to the IEC 60601-1. The differences between the predicate device and subject device will not affect the safety and effectiveness of the subject device.

Note 05:

The difference in waterproof will not affect the safety and effectiveness of the subject device, which has passed IEC 60601-1 and IEC 60601-1-11 safety test. The differences between the predicate device and subject device will not affect the safety and effectiveness of the subject device.

Note 06:

The operating condition of subject device has passed the safety test, and the Instructions for Use in VOL_13_Labeling provides the operating condition, so the difference between the operating conditions of subject device and predicate device will not affect the safety and effectiveness of subject device.

Note 07:

The patient contact materials including colorants (pink, grey, orange and purple) have passed the ISO 10993-1:2009, ISO 10993-5:2009, ISO 10993-10:2010 Biocompatibility Test. The differences between the predicate device and subject device will not affect the safety and effectiveness of the subject device.

Note 08:

The cleaning has been validated according to the standards ASTM E2314 - 03(2014). The differences between the predicate device and subject device will not affect the safety and effectiveness of the subject device.

The subject device and the predicate device have the same intended use and similar technological characteristics, they both use infrared radiation detection method to detect human body forehead temperature. Their design is compact, small and light-weight. They are same in measuring accuracy, and similar in measuring range. Thus, the subject device is substantially equivalent to the predicate devices.

5.7 Brief discussion of nonclinical tests

Nonclinical tests of the Infrared Thermometer are listed as below table:

Tests	Test Standards	Results
Electric Safety	IEC 60601-1:2005+A1:2012, Medical electrical equipment - Part 1: General requirements for basic safety and essential performance; IEC 60601-1-11: 2010, Medical electrical equipment –Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment	Pass
EMC	IEC 60601-1-2:2014, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic compatibility - Requirements and tests;	Pass
Bench performance	ISO 80601-2-56 First Edition 2009-10-01, Medical electrical equipment - Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement	Pass

Tests	Test Standards	Results
	ASTM E1965-98 (Reapproved 2009): Standard Specification for Infrared Thermometers for Intermittent Determination of Patient Temperature	
Clinical accuracy	ASTM E1965-98 (Reapproved 2009): Standard Specification for Infrared Thermometers for Intermittent Determination of Patient Temperature	Pass
Biological Evaluation	AAMI ANSI ISO 10993-1:2009/(R)2013, Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Pass
<i>In vitro</i> Cytotoxicity	ISO 10993-5:2009, Biological evaluation of medical devices - Part 5: Tests for <i>in vitro</i> cytotoxicity	Pass
irritation and skin sensitization	ISO 10993-10: 2010, Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization	Pass

5.8 Brief discussion of clinical tests

The clinical performance test protocol and data analysis is conducted as the requirement of ASTM E1965-98 (2009). The test report showed the clinical performance of the subject device complied with the requirement of ASTM E1965-98 (2009).

5.9 Conclusions

Based on the above information, we conclude that the subject device, Non Contact Infrared Forehead Thermometer, is substantially equivalent to the predicate device.