



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

June 13, 2017

Cynosure, Inc  
Ms. Amy Tannenbaum  
Regulatory Affairs Specialist  
5 Carlisle Road  
Westford, Massachusetts 01886

Re: K171111

Trade/Device Name: Sculpsure  
Regulation Number: 21 CFR 878.5400  
Regulation Name: Low level laser system for aesthetic use  
Regulatory Class: Class II  
Product Code: PKT  
Dated: April 13, 2017  
Received: April 14, 2017

Dear Ms. Tannenbaum:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-

related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

**Jennifer R. Stevenson -S3**

For Binita S. Ashar, M.D., M.B.A., F.A.C.S.  
Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K171111

Device Name

SculpSure™

Indications for Use (Describe)

The Cynosure SculpSure™ is intended for non-invasive lipolysis of the abdomen, flanks, back, and thighs in individuals with a Body Mass Index (BMI) of 30 or less. The device is intended to affect the appearance of visible fat bulges in the abdomen, flanks, back, and thighs.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

Section 5	510(k) Summary
-----------	----------------

**510(k) Summary for SculpSure Laser System**

A summary of 510(k) safety and effectiveness information in accordance with the requirements of 21 CFR 807.92.

<b>807.92(a)(1) Submitter Information</b>	
Applicant	Cynosure, Inc.
Address	5 Carlisle Road Westford, MA 01886
Phone Number	(781)-993-2454
Fax Number	(978) 256-6556
Establishment Registration Number	1222993
Contact Person	Amy Tannenbaum
Preparation Date	April 13, 2017
<b>807.92(a)(2) Name of Device</b>	
Trade or Proprietary Name	SculpSure
Common or Usual Name	Laser
Classification Name	Laser for disruption of adipocyte cells for aesthetic use
Classification Panel	General & Plastic Surgery
Regulation	21 CFR 878.5400
Regulatory Class	II
Product Code(s)	PKT
<b>807.92 (a)(3) Legally marketed device(s) to which equivalence is claimed</b>	
Predicate Devices	Cynosure SculpSure K160470 Zeltiq Coolsculpting K162050 The predicate devices have not been subject to a design-related recall.
<b>807.92(a)(4) Device Description</b>	
	The Cynosure SculpSure is a diode laser system. The main components of SculpSure are a console and four applicators that deliver the laser energy to the patient. Electrically efficient semiconductors generate optical radiation (1060 nm) which is used to deliver laser energy to subcutaneous tissue layers.
<b>807.92(a)(5) Intended Use of the Device</b>	
	The SculpSure laser system is intended for non-invasive lipolysis of the abdomen, flanks, back, and thighs, in individuals with a

Section 5	510(k) Summary
-----------	----------------

	Body Mass Index (BMI) of 30 or less. The device is intended to affect the appearance of visible fat bulges in the abdomen, flanks, back, and thighs.
--	--

**807.92(a)(6) Summary of the Technological Characteristics of the Device Compared to the Predicate**

	<b>Proposed SculpSure Laser System</b>	<b>Cynosure SculpSure Laser System</b>	<b>Zeltiq Coolsculpting Device</b>
510(k) Number	Pending	K160470	K162050
Manufacturer	Cynosure, Inc.	Cynosure, Inc.	Zeltiq Aesthetics, Inc.
Lipolysis Method	Heat-assisted	Heat-assisted	Cold-assisted
Device Type	Diode Laser	Diode Laser	N/A
Wavelength	1060 ±20 nm (infrared)	1060 ±20 nm (infrared)	N/A
Spot Size	4 x 6 cm <sup>2</sup> on each of the Applicator heads	4 x 6 cm <sup>2</sup> on each of the Applicator heads	4.5 x 7 cm <sup>2</sup> (2x) approximate
Pulse Width (laser ON time)	CW	CW	N/A
Power Density	Up to 1.4 W/cm <sup>2</sup>	Up to 1.4 W/cm <sup>2</sup>	N/A
Attachment to patient	Belt	Belt	Belt
Voltage	200-240V~, Single Phase	200-240V~, Single Phase	N/A
Current	20A	20A	N/A

**807.92(b)(1) Non-clinical tests submitted**

There have been minor changes to the software from the previously cleared SculpSure in K160470. These changes have been addressed and additional information is available in Section 16. There have been no changes to the biocompatibility or electrical safety testing from the previously cleared SculpSure device. Additional information is available in Section 15 and Section 17, respectively.

**807.92(b)(2) Clinical tests submitted**

Prospective, controlled studies were conducted at 3 study centers to evaluate the safety and efficacy of Cynosure SculpSure laser for the non-invasive fat reduction of the inner thigh, outer thigh and back. The study was conducted in the same manner as the previous abdomen and flanks studies. The treatment techniques used and subsequent follow up schedule for the study were kept consistent between the studies. A total of 168 subjects, making up 214 treatment areas, were enrolled in the study. This is broken down into 55 subjects/55 treatment areas for back, 52 subjects/52 treatment areas for outer thigh, and 61 subjects/107 treatment areas for inner thighs

Each subject received 1 treatment with the SculpSure. Treatment effectiveness was assessed through blind evaluation of pre and post treatment (12-week) photographs; percentage change in the adipose tissue thickness from baseline to the 12-week follow up measured through ultrasound imaging; and patient satisfaction measured at the 12-week post treatment follow up visit through a 6-point Likert scale. All subjects that were treated in the study were included in the safety analysis and adverse events were assessed at all visits. On an average, blind evaluators were able to identify 86% of the post treatment photographs (back 91%, outer thigh 87% and inner thigh 83%); ultrasound imaging at 12-week post follow up visit showed 8.6% normalized fat reduction (back 10.6%, outer thigh 7.2%, and inner thigh 8.0%). All treatments had a p-value of <0.01 indicating that results were statistically significant. The patient satisfaction survey at 12 week post follow up visit showed 83% satisfaction rate (back 81%, outer thigh 88% and inner thigh 80%). All 168 subjects were included in the safety analysis and adverse events reporting. All events were transient and majority was mild in nature. There were no serious adverse events. The study concluded that the SculpSure is safe and effective for non-invasive fat reduction for the back and thighs.

**807.92(b)(3) Conclusions drawn from clinical and non-clinical tests submitted**

Cynosure believes that SculpSure is substantially equivalent to its predicate devices with same intended use and similar technological characteristics. The non-clinical data for SculpSure supports the safety of the device and the biocompatibility, hardware and software verification and validation demonstrate that the SculpSure should perform as intended in the specified use conditions. Additionally, the successful clinical performance of SculpSure as documented in the clinical study demonstrate that the SculpSure has a safety and effectiveness profile that is similar to its predicate devices with a successful treatment of the back and thighs.