



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

May 25, 2017

Medacta International SA  
% Elizabeth Rose  
Manager, Regulatory Affairs  
Mapi USA, Inc.  
2343 Alexandria Drive, Suite 100  
Lexington, Kentucky 40504

Re: K171170

Trade/Device Name: M.U.S.T. Pedicle Screw System Extension - Straight Connectors and  
Additional Screws

Regulation Number: 21 CFR 888.3070

Regulation Name: Thoracolumbosacral Pedicle Screw System

Regulatory Class: Class III

Product Code: NKB, KWQ, KWP

Dated: April 20, 2017

Received: April 21, 2017

Dear Ms. Rose:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply

with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

Vincent J. Devlin -S

for  
Mark N. Melkerson  
Director  
Division of Orthopedic Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K171170

Device Name

M.U.S.T. Pedicle Screw System Extension - Straight Connectors and Additional Screws

Indications for Use (Describe)

The M.U.S.T. Pedicle Screw System is intended for posterior non-cervical pedicle fixation (T1-S2/ilium) and nonpedicle fixation, or anterolateral fixation (T8-L5). These devices are indicated as an adjunct to fusion for all of the following indications: degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor; pseudoarthrosis and failed previous fusion in skeletally mature patients.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

### I. Submitter

Medacta International SA  
Strada Regina  
6874 Castel San Pietro (CH)  
Switzerland  
Phone (+41) 91 696 60 60  
Fax (+41) 91 696 60 66

Contact Person: Stefano Baj, Regulatory Affairs Manager  
Consultant: Elizabeth Rose, Regulatory Affairs Manager at Mapi USA, Inc.  
Date Prepared: April 20, 2017  
Date Revised: May 23, 2017

### II. Device

Device Proprietary Name: M.U.S.T. Pedicle Screw System Extension - Straight  
Connectors and Additional Screws  
Common or Usual Name: Pedicle Screw Spinal System  
Classification Name: Thoracolumbosacral Pedicle Screw Spinal System  
Primary Product Code: NKB  
Secondary Product Code: KWQ, KWP  
Regulation Number: 21 CFR 888.3070, 21 CFR 888.3060, 21 CFR 888.3050  
Device Classification 2

### III. Predicate Device

Substantial equivalence is claimed to the following devices:

Primary Predicate:

- M.U.S.T. Pedicle Screws System, K121115, Medacta International SA

Additional Predicates:

- M.U.S.T. Extension, K132878, Medacta International SA (also known as M.U.S.T. Pedicle Screws System)
- M.U.S.T. Pedicle Screws System, K141988, Medacta International SA
- M.U.S.T. Pedicle Screws System, K153664, Medacta International SA

#### IV. Device Description

The M.U.S.T. Extension - Straight Connectors and Additional Screws are intended to be used as part of the M.U.S.T. Pedicle Screws System (K121118, K132878, K141988, K153664) for the stabilization and the fusion of the lumbar and thoracic spine. The M.U.S.T. Pedicle Screws System includes: cannulated or non-cannulated poly-axial pedicle screws (K121115 and K132878), cannulated or non-cannulated mono-axial pedicle screws (K132878), set screws (K121115), straight and pre-bent rods (K121115), cross connectors (K132878), enhanced screws and rods designed for percutaneous surgery (K141988), and cannulated or non-cannulated reduction screws (K153664).

The M.U.S.T. Additional Pedicle Screws (standard and enhanced pedicle) are made of Ti6Al4V ELI (ISO 5832-3, ASTM F136-13) and CoCrMo (ISO 5832-12, ASTM F1537-11), while the M.U.S.T. Additional Reduction Screws are made of Ti6Al4V ELI (ISO 5832-3, ASTM F136-13). The screw shaft is color anodized to simplify the identification of the screw diameter. The pedicle screw has a dual lead thread to simplify the screw insertion and reduce the number of turns. The threads are designed with a cylindrical diameter.

The M.U.S.T. Straight Connectors and Additional Screws introduces the following additional components:

| Component                                                          | Diameter | Length     | Material                     |
|--------------------------------------------------------------------|----------|------------|------------------------------|
| Pedicle screw<br>(Sterile and non-sterile)                         | 4 mm     | 20 - 50 mm | Ti6Al4V ELI<br>and<br>CoCrMo |
| Enhanced Solid<br>Pedicle screw<br>(Sterile and non-sterile)       | 4 mm     | 20 - 50 mm | Ti6Al4V ELI<br>and<br>CoCrMo |
| Enhanced Cannulated<br>Pedicle screw<br>(Sterile and non-sterile)  | 4.5 mm   | 20 - 50 mm | Ti6Al4V ELI<br>and<br>CoCrMo |
| Reduction Pedicle<br>screw (Sterile and<br>non-sterile)            | 4 mm     | 20 - 50 mm | Ti6Al4V ELI                  |
| Reduction Cannulated<br>Pedicle screw<br>(Sterile and non-sterile) | 4.5 mm   | 20 - 50 mm | Ti6Al4V ELI                  |
| Straight Cross<br>Connector                                        | N/A      | 20–100 mm  | Ti6Al4V ELI<br>and<br>CoCrMo |

## **V. Indications for Use**

The M.U.S.T. Pedicle Screws System is intended for posterior non-cervical pedicle fixation (T1-S2/ilium) and non-pedicle fixation, or anterolateral fixation (T8-L5). These devices are indicated as an adjunct to fusion for all of the following indications: degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor; pseudoarthrosis and failed previous fusion in skeletally mature patients.

## **VI. Comparison of Technological Characteristics**

The M.U.S.T. Extension - Straight Connectors and Additional Screws and the predicate devices share the following characteristics:

- indications for use;
- design – screws only;
- materials – screws only;
- packaging;
- sterile.

The M.U.S.T. Extension - Straight Connectors and Additional Screws are technologically different from the predicate devices as follows:

- new diameters added;
- addition of a straight connector;
- materials - straight connector.

The M.U.S.T. Extension - Straight Connectors and Additional Screws components are manufactured from Titanium-6Aluminum-4Vanadium Extra Low Interstitial Alloy (ISO 5832-3:1996 Implants For Surgery – Metallic Materials – Part 3: Wrought Titanium 6-Aluminum 4-Vanadium Alloy and ASTM F136-13 Standard Specification For Wrought Titanium-6 Aluminum-4 Vanadium ELI (Extra Low Interstitial) Alloy For Surgical Implant Applications) and Cobalt-Chromium-Molybdenum alloy (Co-Cr-Mo) (ISO 5832-12 Second Edition 2007-05-01: Implants For Surgery -- Metallic Materials -- Part 12: Wrought Cobalt-Chromium-Molybdenum Alloy [Including: Technical Corrigendum 1 (2008)] and ASTM F1537-11 Standard Specification for Wrought Cobalt-28Chromium-6-Molybdenum Alloys for Surgical Implants (UNS R31537, UNS R31538, and UNS R31539).

Biocompatibility testing was conducted on the predicate devices for the same material and testing supports the biological safety of the M.U.S.T. Extension - Straight Connectors and Additional Screws. Additional biocompatibility testing was deemed unnecessary because the materials and manufacturing process are identical to the predicate devices described below.

A comparison of the subject and predicate devices are provided in the tables below:

## Technological comparison

| <b>Parameters</b> | <b>M.U.S.T. Extension - Pedicle Screw<br/>(Subject Device)</b> | <b>M.U.S.T. Pedicle Screw<br/>K121115 and K124034<br/>(Predicate Device)</b> |
|-------------------|----------------------------------------------------------------|------------------------------------------------------------------------------|
| Design/Types      | Posterior fixation devices to be used with rods and set screws | Identical                                                                    |
| Material          | Ti6Al4V ELI<br>CoCrMo                                          | Identical                                                                    |
| Diameter          | 4.0 mm                                                         | 4.5 – 10 mm                                                                  |
| Length            | 20 – 50 mm                                                     | 20 – 100 mm                                                                  |
| Device usage      | Single Use                                                     | Identical                                                                    |
| Shelf Life        | 5 years                                                        | Identical                                                                    |
| Biocompatibility  | Implant with permanent >30 day                                 | Identical                                                                    |
| Sterilization     | Sterile: Gamma<br>Non-sterile                                  | Identical                                                                    |
| Packaging         | Individual packaging                                           | Identical                                                                    |

| <b>Parameters</b> | <b>M.U.S.T. Extension - Solid Enhanced Pedicle Screw<br/>(Subject Device)</b> | <b>M.U.S.T. Solid Enhanced Pedicle Screw<br/>K141988<br/>(Predicate Device)</b> |
|-------------------|-------------------------------------------------------------------------------|---------------------------------------------------------------------------------|
| Design/Types      | Posterior fixation devices to be used with rods and set screws                | Identical                                                                       |
| Material          | Ti6Al4V ELI<br>CoCrMo                                                         | Identical                                                                       |
| Diameter          | 4.0 mm                                                                        | 4.5 – 7 mm                                                                      |
| Length            | 20 – 50 mm                                                                    | 20 – 90 mm                                                                      |
| Device usage      | Single Use                                                                    | Identical                                                                       |
| Shelf Life        | 5 years                                                                       | Identical                                                                       |
| Biocompatibility  | Implant with permanent >30 day                                                | Identical                                                                       |
| Sterilization     | Sterile: Gamma<br>Non-sterile                                                 | Identical                                                                       |
| Packaging         | Individual packaging                                                          | Identical                                                                       |

| <b>Parameters</b> | <b>M.U.S.T. Extension - Cannulated Enhanced Pedicle Screw<br/>(Subject Device)</b> | <b>M.U.S.T. Cannulated Enhanced Pedicle Screw<br/>K141988<br/>(Predicate Device)</b> |
|-------------------|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| Design/Types      | Posterior fixation devices to be used with rods and set screws                     | Identical                                                                            |
| Material          | Ti6Al4V ELI<br>CoCrMo                                                              | Identical                                                                            |

| <b>Parameters</b> | <b>M.U.S.T. Extension - Cannulated Enhanced Pedicle Screw<br/>(Subject Device)</b> | <b>M.U.S.T. Cannulated Enhanced Pedicle Screw<br/>K141988<br/>(Predicate Device)</b> |
|-------------------|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| Diameter          | 4.5 mm                                                                             | 5 – 7 mm                                                                             |
| Length            | 20 – 50 mm                                                                         | 20 – 90 mm                                                                           |
| Device usage      | Single Use                                                                         | Identical                                                                            |
| Shelf Life        | 5 years                                                                            | Identical                                                                            |
| Biocompatibility  | Implant with permanent >30 day                                                     | Identical                                                                            |
| Sterilization     | Sterile: Gamma<br>Non-sterile                                                      | Identical                                                                            |
| Packaging         | Individual packaging                                                               | Identical                                                                            |

| <b>Parameters</b> | <b>M.U.S.T. Extension - Solid Reduction Pedicle Screw<br/>(Subject Device)</b> | <b>M.U.S.T. Solid Reduction Pedicle Screw<br/>K153664<br/>(Predicate Device)</b> |
|-------------------|--------------------------------------------------------------------------------|----------------------------------------------------------------------------------|
| Design/Types      | Posterior fixation devices to be used with rods and set screws                 | Identical                                                                        |
| Material          | Ti6Al4V ELI                                                                    | Identical                                                                        |
| Diameter          | 4.0 mm                                                                         | 4.5 – 7 mm                                                                       |
| Length            | 20 – 50 mm                                                                     | 20 – 90 mm                                                                       |
| Device usage      | Single Use                                                                     | Identical                                                                        |
| Shelf Life        | 5 years                                                                        | Identical                                                                        |
| Biocompatibility  | Implant with permanent >30 day                                                 | Identical                                                                        |
| Sterilization     | Sterile: Gamma<br>Non-sterile                                                  | Identical                                                                        |
| Packaging         | Individual packaging                                                           | Identical                                                                        |

| <b>Parameters</b> | <b>M.U.S.T. Extension - Cannulated Reduction Pedicle Screw<br/>(Subject Device)</b> | <b>M.U.S.T. Cannulated Reduction Pedicle Screw<br/>K153664<br/>(Predicate Device)</b> |
|-------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Design/Types      | Posterior fixation devices to be used with rods and set screws                      | Identical                                                                             |
| Material          | Ti6Al4V ELI                                                                         | Identical                                                                             |
| Diameter          | 4.5 mm                                                                              | 5 – 7 mm                                                                              |
| Length            | 20 – 50 mm                                                                          | 25 – 90 mm                                                                            |
| Device usage      | Single Use                                                                          | Identical                                                                             |
| Shelf Life        | 5 years                                                                             | Identical                                                                             |
| Biocompatibility  | Implant with permanent >30 day                                                      | Identical                                                                             |

| Parameters    | M.U.S.T. Extension - Cannulated Reduction Pedicle Screw<br>(Subject Device) | M.U.S.T. Cannulated Reduction Pedicle Screw<br>K153664<br>(Predicate Device) |
|---------------|-----------------------------------------------------------------------------|------------------------------------------------------------------------------|
| Sterilization | Sterile: Gamma<br>Non-sterile                                               | Identical                                                                    |
| Packaging     | Individual packaging                                                        | Identical                                                                    |

| Parameters                    | M.U.S.T. Extension – Straight Connectors<br>(Subject Device)   | M.U.S.T. Cross Connectors<br>K132878<br>(Predicate Device) |
|-------------------------------|----------------------------------------------------------------|------------------------------------------------------------|
| Design/Types                  | Posterior fixation devices to be used with rods and set screws | Identical                                                  |
| Material                      | Ti6Al4V ELI<br>CoCrMo                                          | Ti6Al4V ELI                                                |
| Length                        | 20 – 50 mm                                                     | 20 – 90 mm                                                 |
| Device usage                  | Single Use                                                     | Identical                                                  |
| Shelf Life                    | 5 years                                                        | Identical                                                  |
| Biocompatibility              | Implant with permanent >30 day                                 | Identical                                                  |
| Sterilization                 | Sterile: Gamma<br>Non-sterile                                  | Identical                                                  |
| Packaging                     | Individual packaging                                           | Identical                                                  |
| Additional technical features | Length is not adjustable                                       | Length is adjustable from 35 – 98 mm                       |

### *Discussion*

As seen above, the differences between the subject and predicate devices are the addition of new screw diameters (4.0 mm for solid screws and 4.5 mm for cannulated screws) to the range of sizes and the addition of a straight connector. These technological differences do not raise new questions of safety or effectiveness. A comparison evaluation shows there are no new risks associated with the subject device design.

## **VII. Performance Data**

The addition of the subject devices was evaluated by risk analysis to identify any new risks associated with the extension of the pedicle screw sizes and the inclusion of the straight connectors. Based on the risk analysis, design verification was conducted to written protocols with pre-defined acceptance criteria. The protocols and predefined acceptance criteria were based on the standards, FDA guidance, and comparison to the predicate device system.

The M.U.S.T. Extension - Straight Connectors and Additional Screws was tested per ASTM F1717-12 Standard Test Methods for Spinal Implant Constructs in a Vertebrectomy Model using the worst-case device for each of the following tests:

1. Dynamic Compression - Bending
2. Static Compression – Bending
3. Static Torsion

The testing met all acceptance criteria and verifies that the performance of the M.U.S.T. Extension - Straight Connectors and Additional Screws is substantially equivalent to the predicate devices.

### **VIII. Conclusion**

The information provided above supports that the M.U.S.T. Extension - Straight Connectors and Additional Screws are as safe and effective as the predicate devices. Therefore, it is concluded that the M.U.S.T. Extension - Straight Connectors and Additional Screws is substantially equivalent to the predicate devices.