



July 7, 2017

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

ECO Medi Glove Sdn. Bhd.
Suresh Kumar
Official Correspondent
Lot 23826, Jalan Tembaga Kuning
Kamunting Raya Industrial Estate
Kamunting Perak, 34600 Malaysia

Re: K171339

Trade/Device Name: EMG Blue Nitrile Examination Gloves Powder free with tested for use with chemotherapy Drugs

Regulation Number: 21 CFR 880.6250

Regulation Name: Patient Examination Glove

Regulatory Class: Class I

Product Code: LZA, LZC

Dated: May 15, 2017

Received: May 19, 2017

Dear Suresh Kumar:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply

with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

Mark S. Fellman -S

for

Lori A. Wiggins, MPT, CLT
Acting Director
Division of Anesthesiology,
General Hospital, Respiratory,
Infection Control, and Dental Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K171339

Device Name
EMG Blue Nitrile Examination Glove Powder Free with tested for use with Chemotherapy Drugs

Indications for Use (Describe)

A powder-free patient examination glove is a disposable device intended for use medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner and use for use with chemotherapy drugs.

In addition these gloves were tested for use with chemotherapy drugs in accordance ASTM D6978-05 standards- Practice for assessment of Medical Glove to Permeation by chemotherapy drugs

| Chemotherapy Drugs and Concentration | Minimum Breakthrough Detection Time(Min), $\mu\text{g}/\text{cm}^2/\text{minutes}$ |
|--------------------------------------|--|
| 1) Carmustine (BCNU)(3.3 mg/ml) | 24.0 Minutes |
| 2)Cyclophosphamide (20 mg/ml) | > 240 Minutes |
| 3)Doxorubicin Hydrochloride (2mg/ml) | > 240 Minutes |
| 4)Etoposide (20mg/ml) | > 240 Minutes |
| 5)Fluorouracil (50mg/ml) | > 240 Minutes |
| 6)Methorexate (25mg/ml) | > 240 Minutes |
| 7)Paclitaxel (6mg/ml) | > 240 Minutes |
| 8)Thiotepa (10mg/ml) | 56.9 Minutes |
| 9)Carboplatin (10mg/ml) | > 240 Minutes |
| 10)Cisplatin (1 mg/ml) | > 240 Minutes |
| 11)Dacarabazine (1 mg/ml) | > 240 Minutes |
| 12)Ifosfamide (5 mg/ml) | > 240 Minutes |
| 13)Mitomycin (0.5mg/ml) | > 240 Minutes |
| 14)Mitoxantrone (2 mg/ml) | > 240 Minutes |
| 15)Vincristine Sulfate (1 mg/ml) | > 240 Minutes |

The Maximum testing is 240 minutes. Please note that the following drugs have extremely low Permeation time,
Carmustine (BCNU)(3.3 mg/ml) - Minimum breakthrough detection time : 24.0 minutes
Thiotepa (10mg/ml) - Minimum breakthrough detection time : 56.9 minutes

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510K No : K171339

ECO Medi Glove Sdn. Bhd. (815262-D)
(formerly known as Sinetimed Consumables Sdn. Bhd.)

Lot 23826, Jalan Tembaga Kuning, Kamunting Raya Industrial Estate, 34600 Taiping, Perak Darul Ridzuan. **MALAYSIA.**
TEL +60-5-891 2777 FAX +60-5-891 2999

510(K) Summary
EMG Blue Nitrile Examination Gloves Powder Free with
Tested for use with Chemotherapy Drugs

1.0 Submitter :

Company Name : ECO MEDI GLOVE SDN. BHD.

Company Address : Lot 23826, Jalan Tembaga Kuning
Kamunting Raya Industrial Estate,
34600, Kamunting Perak
Malaysia.

Contact Person : Mr Suresh Kumar

Telephone No : 603-60283033

Email : suresh@ecomediglove.com.my

2.0 Preparation Date : 4th July 2017

3.0 Name of the Device

Trade Name / Proprietary Name : EMG Blue Nitrile Examination Gloves
Powder Free with tested for use with
Chemotherapy drugs.

Device Name : Nitrile Patient Examination gloves.

Device Classification Name : Patient Examination gloves (21 CFR 880.6250).

Device Class : Class I.

Product Code : Nitrile-LZA and LZC.

ECO Medi Glove Sdn. Bhd. (815262-D)
(formerly known as Sinetimed Consumables Sdn. Bhd.)

Lot 23826, Jalan Tembaga Kuning, Kamunting Raya Industrial Estate, 34600 Taiping, Perak Darul Ridzuan. **MALAYSIA.**
TEL +60-5-891 2777 FAX +60-5-891 2999

4.0) Identification of The Legally Marketed Device :

EMG Blue Nitrile Examination Gloves Powder Free with tested for use with Chemotherapy drugs,LZC, which meets all the requirement of ASTM D 6319-10 and FDA 21 CFR 880.6250.It is equivalent to K161187, EMG Blue Nitrile Examination Glove Powder Free with tested for use with Chemotherapy Drugs.

5.0 Device Description

The subject device in this 510(k) Notification is EMG Blue Nitrile Examination Gloves Powder Free with tested for use with Chemotherapy drugs. The subject device is a patient examination glove made from nitrile latex compound, Blue color, powder free and non-sterile (Per 21 CFR 880.6250, class I). The device meets all the specifications in ASTM D6319-10, Standard specification for Nitrile Examination Gloves. Additionally, the gloves have been tested for biocompatibility and permeability to chemotherapy drugs.

The Blue Nitrile Medical Examination Gloves ,Powder Free, is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner (product Code LZA) and is used with chemotherapy drugs (Product code LZC). The subject device is substantially equivalent to the legally marketed Nitrile Medical Examination Gloves (product Code LZA and LZC).

6.1 Specification for Nitrile gloves:

6.2 Dimension and Thickness of Gloves

| Dimension | Size S | Size M | Size L | Size XL |
|---------------------------------------|---------------|---------------|---------------|----------------|
| Overall Length (mm) (Minimum) | 230 | 230 | 230 | 230 |
| Width (\pm 5mm) | 85 | 95 | 105 | 115 |
| Thickness at Palm (mm) (Minimum) | 0.05 | 0.05 | 0.05 | 0.05 |
| Thickness at Finger Tip (mm)(Minimum) | 0.05 | 0.05 | 0.05 | 0.05 |

ECO Medi Glove Sdn. Bhd. (815262-D)
(formerly known as Sinetimed Consumables Sdn. Bhd.)

Lot 23826, Jalan Tembaga Kuning, Kamunting Raya Industrial Estate, 34600 Taiping, Perak Darul Ridzuan, MALAYSIA.
TEL +60-5-891 2777 FAX +60-5-891 2999

6.1.2 Gloves Physical Properties and Holes

| Measurement | Before Ageing | After Aging at 70°C for 168 hrs @ 100°C for 22 hrs |
|-------------------------|---------------------------------|--|
| Tensile Strength (MPa) | 14min | 14 Min |
| Ultimate Elongation (%) | 500min | 400min |
| Pin-hole Level | AQL 2.5 Inspection Level G-1 | AQL 2.5 Inspection Level G-1 |

Gloves meet all the specification listed in ASTM D 6319-10

| Characteristics | Acceptance Criteria | EMG Blue Nitrile Examination Gloves Powder Free with tested for use with chemotherapy drugs. K171339 | EMG Blue Nitrile Examination Gloves, Powder Free with tested for use with chemotherapy drugs, K161187 | Comparison |
|-----------------|---|---|---|------------|
| Product Code | LZA and LZC | LZA and LZC | LZA and LZC | Same |
| Intended use | A powder free patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. The device is for over-the-counter use. Intended use. | A powder free patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. The device is for over-the-counter use. | A powder free patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. The device is for over-the-counter use. | Same |
| Material use | Nitrile latex compound | Nitrile latex compound | Nitrile latex compound | Same |
| Color | Blue | Blue | Blue | Same |
| Sterility | Non sterile | Non sterile | Non sterile | Same |
| Single used | Single used | Single used | Single used | Same |

ECO Medi Glove Sdn. Bhd. (815262-D)
(formerly known as Sinetimed Consumables Sdn. Bhd.)

Lot 23826, Jalan Tembaga Kuning, Kamunting Raya Industrial Estate, 34600 Taiping, Perak Darul Ridzuan, MALAYSIA.

TEL +60-5-891 2777 FAX +60-5-891 2999

| | | | | |
|---|---|---|---|------|
| Dimensions | Overall Length (mm) Min 230 mm Width (\pm 5mm) Size S = 85mm Size M = 95mm Size L = 105mm Size XL = 115mm Thickness at Palm (mm) Min; 0.05 mm Thickness at Finger Tip (mm) Min 0.05 mm | Size : Medium Palm Width: 98 mm Length : 247 mm Finger Thickness : 0.08mm Palm Thickness : 0.05mm Meets ASTM D6319-10 | Size : Medium Palm Width: 99 mm Length : 247 mm Finger Thickness : 0.11mm Palm Thickness : 0.07mm Meets ASTM D6319-10 | Same |
| Physical properties | Before Ageing Tensile Strength (MPa) = Min 14 Mpa Ultimate Elongation (%) = Min 500 % After Aging at 70°C for 168 hrs @ 100°C for 22 hrs Tensile Strength (MPa) = Min 14 Mpa Ultimate Elongation (%) = Min 400 % | Tensile strength (Mpa) : Before aging : 25.45Mpa After Accelerated aging: 26.33Mpa Ultimate Elongation(%): Before aging:600% After aging : 550% Meets ASTM D6319-10 | Tensile strength (Mpa) : Before aging : 25.45Mpa After Accelerated aging: 26.33Mpa Ultimate Elongation(%): Before aging:600% After aging : 550% Meets ASTM D6319-10 | Same |
| Freedom from pinholes | AQL 2.5 Inspection Level G-1 | AQL 2.5 Meets ASTM D5151-06 | AQL 2.5 Meets ASTM D5151-06 | Same |
| Residual Powder | \leq 2.0 mg/glove | Average 0.17 mg/glove Meets ASTM D6124-06 | Average 0.19 mg/glove Meets ASTM D6124-06 | Same |
| Biological Evaluation on Medical Device -Primary Skin Irritation Test | ISO 10993-10:2010 | Under the conditions of this study, the test article was a non-irritant. | Under the conditions of this study, the test article was a non-irritant. | Same |
| Biological Evaluation on Medical Device- Dermal Sensitization Assay | ISO 10993-10:2010 | Under the conditions of this study, the test article was a non-sensitizer. | Under the conditions of this study, the test article was a non- sensitizer. | Same |

ECO Medi Glove Sdn. Bhd. (815262-D)
(formerly known as Sinetimed Consumables Sdn. Bhd.)

Lot 23826, Jalan Tembaga Kuning, Kamunting Raya Industrial Estate, 34600 Taiping, Perak Darul Ridzuan, MALAYSIA.

TEL +60-5-891 2777 FAX +60-5-891 2999

| | | | | |
|--|---|--|---|---|
| <p>Resistance against Chemotherapy Drugs</p> | <p>Standards Practice for Assessment of resistance of Medical Glove to Permeation by Chemotherapy drugs ASTM D6978-05(2013)</p> | <p>1) Carmustine (3.3mg/ml or 3000ppm), Minimum Breakthrough detection time : 24 minutes</p> <p>2) Cyclophosphamide (20mg/ml or 20,000ppm), Minimum Breakthrough detection time: >240 minutes</p> <p>3) Cisplatin (1.0mg/ml or 1,000ppm), Minimum Breakthrough detection time: >240 minutes</p> <p>4)Doxorubicin Hydrochloride (2.0mg/ml or 2000ppm), Minimum Breakthrough detection time: >240 minutes</p> <p>5) Etoposide (20mg/ml or 20,000ppm), Minimum Breakthrough detection time: >240 minutes</p> <p>6) Flourouracil (50mg/ml or 50,000), Minimum Breakthrough detection time: >240 minutes</p> | <p>1) Carmustine (3.3mg/ml or 3000ppm), Breakthrough : 1.3 minutes.</p> <p>2) Cyclophosphamide (20mg/ml or 20,000ppm), Breakthrough time : >240 minutes.</p> <p>3) Cytarabine (100mg/ml or 100,000ppm), Breakthrough time : >240 minutes.</p> <p>4)Doxorubicin Hydrochloride (2.0mg/ml or 2000ppm), Breakthrough time : >240 minutes.</p> <p>5) Etoposide (20mg/ml or 20,000ppm), Breakthrough time : >240 minutes.</p> <p>6) Flourouracil (50mg/ml or 50,000), Breakthrough time : >240 minutes.</p> <p>7) Methorexate (25mg/ml or 25,000ppm), Breakthrough time : >240 minutes.</p> <p>8) Paclitaxel (6mg/ml or 6,000ppm), Breakthrough time : >240 minutes.</p> | <p>Subject device was tested with 7 mandatory chemotherapy drugs with additional 8 chemotherapy drugs. Predicate device was tested with 7 mandatory chemotherapy drugs with additional of 2 chemotherapy drugs.</p> |
|--|---|--|---|---|

ECO Medi Glove Sdn. Bhd. (815262-D)
(formerly known as Sinetimed Consumables Sdn. Bhd.)

Lot 23826, Jalan Tembaga Kuning, Kamunting Raya Industrial Estate, 34600 Taiping, Perak Darul Ridzuan, **MALAYSIA.**

TEL +60-5-891 2777 FAX +60-5-891 2999

| | | | | |
|--|--|--|---|--|
| | | <p>7) Methorexate (25mg/ml or 25,000ppm), Minimum Breakthrough detection time: > 240 minutes</p> <p>8) Paclitaxel (6mg/ml or 6,000ppm), Minimum Breakthrough detection time: >240minutes</p> <p>9) Thiotepa (10mg/ml or 10,000ppm), Minimum Breakthrough detection time: 56.9 minutes</p> <p>10) Ifosfamide (50mg/ml) Minimum Breakthrough detection time: >240 minutes</p> <p>11) Mitoxantrone (2mg/ml) , Minimum Breakthrough detection time: >240 minutes</p> <p>12) Vincristine Sulfate (1mg/ml) , Minimum Breakthrough detection time: >240 minutes</p> <p>13) Carboplatin (10.0mg/ml), Minimum Breakthrough detection time: >240 minutes</p> | <p>9) Thiotepa (10mg/ml or 10,000ppm), Breakthrough time: 67.8 minutes.</p> | |
|--|--|--|---|--|

ECO Medi Glove Sdn. Bhd. (815262-D)
(formerly known as Sinetimed Consumables Sdn. Bhd.)

Lot 23826, Jalan Tembaga Kuning, Kamunting Raya Industrial Estate, 34600 Taiping, Perak Darul Ridzuan, MALAYSIA.

TEL +60-5-891 2777 FAX +60-5-891 2999

| | | | | |
|--|--|---|--|--|
| | | <p>14) Dacarbazine (DTIC) (10.0mg/ml) Minimum Breakthrough detection time: >240 minutes</p> <p>15) Mitomicin C (0.5mg/ml) Minimum Breakthrough detection time: >240 minutes</p> | | |
|--|--|---|--|--|

7.0 Indications for Use

A powder-free patient examination glove is a disposable device intended for medical Purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner and for use with chemotherapy drugs .

In addition these gloves were tested for use with Chemotherapy drugs in accordance With ASTM D6978-05 standards Practice for assessment of Medical Glove to Permeation by chemotherapy drugs.

| Chemotherapy Drug and concentration | Minimum Breakthrough detection time in Minutes, $\mu\text{g}/\text{cm}^2/\text{minute}$ |
|--|---|
| 1) Carmustine (BCNU) (3.3mg/ml) | 24 minutes |
| 2) Cyclophosphamide (20mg/ml) | > 240 minutes |
| 3) Cytarabine (10mg/ml) | > 240 minutes |
| 4) Doxorubicin Hydrochloride (2 mg/ml) | > 240 minutes |
| 5) Etoposide (20mg/ml) | > 240 minutes |
| 6) Fluorouracil (50mg/ml) | > 240 minutes |
| 7) Methorexate (25mg/ml) | > 240 minutes |
| 8) Paclitaxel (6mg/ml) | > 240 minutes |
| 9) Thiotepa (10mg/ml) | 56.9 minutes |
| 10) Cisplatin (1.0mg/ml) | > 240 minutes |

ECO Medi Glove Sdn. Bhd. (815262-D)
(formerly known as Sinetimed Consumables Sdn. Bhd.)

Lot 23826, Jalan Tembaga Kuning, Kamunting Raya Industrial Estate, 34600 Taiping, Perak Darul Ridzuan, MALAYSIA.
TEL +60-5-891 2777 FAX +60-5-891 2999

| | |
|---|-------------------------|
| 11) Dacarbazine (DTIC) (10.0mg/ml) | > 240 minutes |
| 12) Ifosfamide (50.0mg/ml) | > 240 minutes |
| 13) Mitomycin C (0.5mg/ml) | > 240 minutes |
| 14) Mitoxantrone (2.0mg/ml) | > 240 minutes |
| 15) Vincristine Sulfate (1.0mg/ml) | > 240 minutes |

The maximum testing time is 240 minutes. Please note that the following drugs have extremely low permeation time,

Carmustine (BCNU) (3.3mg/ml) - Minimum Breakthrough detection time 24 Minutes

Thiotepa (10mg/ml) – Minimum Breakthrough detection time 56.9 minutes.

8.0 Summary of the Technological Characteristics of the Device compared to the Predicate Device for substantial equivalent discussion

There are no differences in technological characteristics of the subject device compare with the predicate device.

The gloves are made from nitrile latex compound, Blue color, Powder free and non-sterile. The gloves met all the specifications in ASTM D6319-10 Standard specification for Nitrile Examination Gloves as well Biological Evaluation on medical device. Additionally, the gloves have been tested for permeability to Chemotherapy drugs.

9.0 Conclusion

Based on intended uses, technological characteristics and non – clinical performance Data, the EMG Blue Nitrile Examination Gloves Powder Free with tested for use with Chemotherapy Drugs is as safe, as effective, and performs as well as the legally marketed predicate device, K161187.