



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

Arterys, Inc.  
% Ms. Golnaz Moeini  
Director of Quality and Regulatory Affairs  
51 Federal Street, Suite 305  
SAN FRANCISCO CA 94107

July 18, 2017

Re: K171544

Trade/Device Name: Arterys Viewer  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Picture archiving and communications system  
Regulatory Class: II  
Product Code: LLZ  
Dated: May 24, 2017  
Received: May 26, 2017

Dear Ms. Moeini:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

A handwritten signature in blue ink that reads "Michael D. O'Hara". The signature is written in a cursive style. To the right of the signature, the word "For" is printed in a small, black, sans-serif font.

Robert Ochs, Ph.D.  
Director  
Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K171544

Device Name

Arterys Viewer

Indications for Use (Describe)

Arterys Viewer is intended to be used as a Digital Imaging and Communications in Medicine (DICOM) and non-DICOM information and data management system. The Arterys Viewer displays, processes, stores, and transfers medical data from original equipment manufacturers (OEMs) that support the DICOM standard, with the exception of mammography. It provides the capability to store images and patient information from OEM equipment, and perform filtering, digital manipulation and quantitative measurements.

The client software is designed to run on standard personal and business computers. The product is intended to be used by trained medical professionals, including but not limited to radiologists, oncologists, and physicians. It is intended to provide image and related information that is interpreted by a trained professional to render findings and/or diagnosis, but it does not directly generate any diagnosis or potential findings.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## Section 5. 510(k) Summary

### 1. General Information

<b>510(k) Sponsor</b>	Arterys Inc.
<b>Address</b>	51 Federal St. Suite 305 San Francisco, CA 94107
<b>Correspondence Person</b>	Golnaz Moeini, RAC Director of Quality and Regulatory Arterys Inc.
<b>Contact Information</b>	Email: quality@arterys.com Phone: 408-504-3187
<b>Date Prepared</b>	May 24, 2017

### 2. Proposed Device

<b>Proprietary Name</b>	<i>Arterys Viewer</i>
<b>Common Name</b>	<i>Viewer</i>
<b>Classification Name</b>	System, Image Processing, Radiological
<b>Regulation Number</b>	21 CFR 892.2050
<b>Product Code</b>	LLZ
<b>Regulatory Class</b>	II

### 3. Predicate Device

#### Primary Predicate

<b>Proprietary Name</b>	<i>HealthMyne PACS</i>
<b>Premarket Notification</b>	K152186
<b>Classification Name</b>	System, Image Processing, Radiological
<b>Regulation Number</b>	21 CFR 892.2050
<b>Product Code</b>	LLZ
<b>Regulatory Class</b>	II

#### Reference Device

<b>Proprietary Name</b>	<i>Arterys Software 2.0</i>
<b>Premarket Notification</b>	K162513
<b>Classification Name</b>	System, Image Processing, Radiological
<b>Regulation Number</b>	21 CFR 892.2050
<b>Product Code</b>	LLZ
<b>Regulatory Class</b>	II

### 4. Device Description

*Arterys Viewer* is intended to be used as a Digital Imaging and Communications in Medicine (DICOM) and non-DICOM information and data management system. The Arterys Viewer displays, processes, stores, and transfers medical data from original equipment manufacturers (OEMs) that support the DICOM standard, with the exception of mammography. It provides

the capability to store images and patient information from OEM equipment, and perform filtering, digital manipulation and quantitative measurements.

The client software is designed to run on standard personal and business computers. The product is intended to be used by trained medical professionals, including but not limited to radiologists, oncologists, and physicians. It is intended to provide image and related information that is interpreted by a trained professional to render findings and/or diagnosis, but it does not directly generate any diagnosis or potential findings.

The following visualization, quantification and data-reporting functionalities are provided by the software:

Visualization:

- 2D image review
- 3D image review by means of MIP, MinIP, Surface or Average
- Multi-planar reconstruction (MPR) views (axial, coronal, sagittal, and oblique)
- Image navigation tools
- Cine play

Quantification:

- Distance and area measurements

Data reporting:

The distance, area, and user-selected annotated images are displayed to the user within the software client web browser. The user has the option to save the data for later use. The user can also send the data to PACS, for review.

## **5. Indications for Use**

*Arterys Viewer* is intended to be used as a Digital Imaging and Communications in Medicine (DICOM) and non-DICOM information and data management system. The Arterys Viewer displays, processes, stores, and transfers medical data from original equipment manufacturers (OEMs) that support the DICOM standard, with the exception of mammography. It provides the capability to store images and patient information from OEM equipment, and perform filtering, digital manipulation and quantitative measurements.

The client software is designed to run on standard personal and business computers. The product is intended to be used by trained medical professionals, including but not limited to radiologists, oncologists, and physicians. It is intended to provide image and related information that is interpreted by a trained professional to render findings and/or diagnosis, but it does not directly generate any diagnosis or potential findings.

## 6. Comparison of Technological Characteristics with the Predicate and Reference Device

<b>Feature/ Function</b>	<b>Proposed Device: Arterys Viewer</b>	<b>Primary Predicate Device: HealthMyne PACS (K152186)</b>	<b>Reference Device: Arterys Software v2.0 (K162513)</b>
<b>Support Mammography</b>	No	No	
<b>Operating System</b>	Client server architecture using Linux server and web browser client (Web-based only)	Client server architecture utilizing Windows and Linux Platforms (Installed-client only)	Client server architecture using Linux server and web browser client (Web-based only)
<b>Image storage/ compression</b>	Support JPEG2000 and compression	Support JPEG2000 and compression	
<b>DICOM Compliant</b>	Yes	Yes	
<b>Worklists</b>	Yes	Yes	
<b>Filter and Search capabilities</b>	Yes	Yes	
<b>Ability to search studies</b>	Yes	Yes. Dynamic (freeform) search and matching. Also customized, indexed search parameters.	
<b>View study-related documents</b>	Yes. Can view reports/documents (except RT Struct)	Yes. Presentation states and RT Struct (regions of interest)	
<b>Priority “stat” studies</b>	No	Yes. Studies with a DICOM priority tag (“Stat” studies) are given priority order (top of the list) in the exam view.	No
<b>Managed pushed studies</b>	Yes. This product supports only pushed studies. No studies originate within the product	Yes. This product support only pushed studies. No studies originate within the product	
<b>View current and prior studies at the same time</b>	No	Yes. Called current study and prior study. Prior studies have a large “PRIOR” label in each viewport.	No

<b>Feature/ Function</b>	<b>Proposed Device: Arterys Viewer</b>	<b>Primary Predicate Device: HealthMyne PACS (K152186)</b>	<b>Reference Device: Arterys Software v2.0 (K162513)</b>
<b>Re-organize series in a study (for viewing)</b>	Yes	Yes. Has a thumbnail view with the ability to drag and drop the thumbnail into a viewport. Cannot save the order.	
<b>Create separate displays</b>	Yes	Yes. Can display a viewport as a single viewport, can select a viewport layout and add series to it.	
<b>Cycle through series</b>	Yes	Yes. Can show the “next” and “previous” sets of series	
<b>Image display modes</b>	Yes. Static and cine	Yes. Static and manual cine.	
<b>Select images</b>	Yes. The images currently being viewed are highlighted, and can be chosen by the user.	Yes. There is an active image indicator. The active image can be chosen or is automatically set based on tool use.	
<b>Delete Images</b>	No	Yes. Only with PACS admin privileges and from the admin console.	No
<b>Sort Images</b>	Yes	Yes. Sorting and grouping are by system-defined rules	
<b>Scrolling through slices</b>	Yes	Yes. Linked series are scrolled together. Can “swipe” on a scroll bar to move through slices quickly. Can lock/unlock scrolling through every image.	
<b>Zoom in/out</b>	Yes. Default settings can zoom interactively	Yes. Default settings can zoom interactively	

<b>Feature/ Function</b>	<b>Proposed Device: Arterys Viewer</b>	<b>Primary Predicate Device: HealthMyne PACS (K152186)</b>	<b>Reference Device: Arterys Software v2.0 (K162513)</b>
<b>Pan an image</b>	Yes	Yes	
<b>Standard viewport layouts</b>	Yes	Yes. Viewport layouts that are independent of any modality or common features of series	
<b>Labels</b>	Yes	Yes. There are labels in the viewport for patient, study, and image information.	
<b>Orientation labels</b>	Yes	Yes	
<b>Cross-reference indicator</b>	Yes	Yes. There is a cutline in linked viewports to indicate intersection	
<b>View DICOM data</b>	Yes. You can view the DICOM information about the patient and study, and the pixel information	Yes. You can view the DICOM information about the patient and study, and the pixel information	
<b>Create MPR images</b>	Yes	No. MPRs from the external source are supported/displayed.	Yes



<b>Feature/ Function</b>	<b>Proposed Device: Arterys Viewer</b>	<b>Primary Predicate Device: HealthMyne PACS (K152186)</b>	<b>Reference Device: Arterys Software v2.0 (K162513)</b>
<b>Window/level determination</b>	Yes Same approach as predicate but we let the user adjust W/L	Determined by a lookup table function (linear) and the W/L values of the image. If the image has a custom lookup table or a fixed W/L, those settings are used instead of allowing changing of W/L. If no W/L, then a histogram is used.	
<b>Window/level access series</b>	Yes	Yes. W/L settings are applied to the active image and any linked images.	
<b>Window/level presets</b>	Yes	Yes. Factory default	
<b>Adjust window/level</b>	Yes. Can interactively adjust the window and level	Yes. Can interactively adjust the window and level	
<b>Annotation</b>	Yes. Display Only	Yes. Display Only	
<b>Measuring tools</b>	Yes. Linear, area, and pixel intensity and location of a point	No. Pixel intensity and location.	Yes. Linear, area, and pixel intensity and location of a point

<b>Feature/ Function</b>	<b>Proposed Device: Arterys Viewer</b>	<b>Primary Predicate Device: HealthMyne PACS (K152186)</b>	<b>Reference Device: Arterys Software v2.0 (K162513)</b>
<b>Detect image/patient issues</b>	No	Yes. Can view patients and studies with errors	No
<b>Custom filters</b>	Yes. Can set filters to affect the studies listed	Yes. Can set filters to affect the studies listed	
<b>Set reading state</b>	Yes	Yes. Can mark a study as read.	
<b>Custom search groups</b>	Yes	Yes. Can set “codes” to index elements for searching (for example, referring physician or sets of exam types_ for faster auto-complete during search	
<b>Display radiation therapy information</b>	No	Yes. Supports DICOM RT Structures	No
<b>Surface Rendering</b>	Yes	No	Yes
<b>Annotation Propagation</b>	Yes	No	Yes

<b>Feature/ Function</b>	<b>Proposed Device: Arterys Viewer</b>	<b>Primary Predicate Device: HealthMyne PACS (K152186)</b>	<b>Reference Device: Arterys Software v2.0 (K162513)</b>
<b>Time Curve Display</b>	Yes	No	Yes
<b>Contour Creation</b>	Yes	No	Yes
<b>Report Creation</b>	Yes	No	Yes

## 7. Substantial Equivalence Summary

*Arterys Viewer* has the same intended use, indications for use and technological characteristics as the HealthMyne PACS predicate device (K152186). The added features to allow for usability enhancements is similar to the features in the reference device, *Arterys Software v2.0*, which was cleared by the FDA under K162513. Any noted minor differences have been explained and do not raise any different questions of safety or effectiveness. The implemented design controls, risk management activities, labeling and performed verification and validation tests demonstrate the safety and efficacy of the proposed device. Based on the comparison information provided above, *Arterys Viewer* is substantially equivalent to the predicate device.

## 8. Performance Data

Safety and performance of *Arterys Viewer* has been evaluated and verified in accordance with software specifications and applicable performance standards through software verification and validation testing. Additionally, the software validation activities were performed in accordance with *IEC 62304:2006/AC: 2008- Medical device software – Software life cycle processes*, in addition to the FDA Guidance documents, “*Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices*” and “*Content of Premarket Submission for Management of Cybersecurity in Medical Devices.*”

## 9. Conclusion

Based on the information submitted in this premarket notification, and based on the indications for use, technological characteristics and performance testing, *Arterys Viewer* raises no new

questions of safety and effectiveness and is substantially equivalent to the predicate device in terms of safety, efficacy and performance.