



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

July 3, 2017

Codman & Shurtleff, Inc.  
Yoon Hee Beatty  
Senior Regulatory Affairs Specialist  
325 Paramount Drive  
Raynham, Massachusetts 02767

Re: K171653  
Trade/Device Name: YOGA Microcatheter  
Regulation Number: 21 CFR 870.1250  
Regulation Name: Percutaneous Catheter  
Regulatory Class: Class II  
Product Code: DQY  
Dated: June 2, 2017  
Received: June 5, 2017

Dear Ms. Yoon Hee Beatty:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-

related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

Carlos L. Pena -S 

Carlos L. Peña, PhD, MS  
Director  
Division of Neurological  
and Physical Medicine Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K171653

Device Name

YOGA Microcatheter

Indications for Use (Describe)

The YOGA Microcatheter is intended for use in peripheral, coronary, and neuro vasculature for the intravascular introduction of interventional/diagnostic devices.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

### I. Submitter

Codman & Shurtleff, Inc.  
325 Paramount Drive  
Raynham, MA 02767

Tel: (305) 265-2919  
Fax: (305) 265-6889

Contact Person: Yoon Hee Beatty  
Date Prepared: June 2, 2017

### II. Device

<b>Table 1. Device</b>	
<b>Device Proprietary Name</b>	YOGA Microcatheter
<b>Common or Usual Name</b>	Catheter, Percutaneous
<b>Classification Name</b>	Catheter, Percutaneous, Class II, 21 CFR 870.1250
<b>Regulatory Classification</b>	II
<b>Product Code</b>	DQY

### III. Predicate Device

The predicate device listed in **Table 2** below is applicable to the device in this submission.

<b>Table 2. Predicate 510(k) Clearance</b>			
<b>510(k) Number</b>	<b>Date Cleared</b>	<b>Device Name</b>	<b>Manufacturer</b>
K162563	01/05/2017	YOGA Microcatheter	Codman & Shurtleff, Inc.

### IV. Device Description

The YOGA Microcatheter is a variable stiffness, single lumen catheter designed to access small, tortuous vasculature. The catheter shaft is composed of a variable pitch stainless steel braid with a PTFE inner liner to facilitate movement of guide wires and other devices. The exterior of the catheter shaft is covered with a polymer material, which encapsulates the stainless steel braid construction. The distal end of the catheter has a radiopaque marker band and has a hydrophilic coating to provide lubricity for navigation of vessels. The proximal end of the catheter has a hub and an ID band is placed at the distal end of the hub over a strain relief. A steam shaping mandrel is provided in the package. The two new YOGA Microcatheters are provided in back up (XB) and extra backup (XXB) configurations, i.e. a stiffer distal end to provide additional support.

### V. Indications for Use

The YOGA Microcatheter is intended for use in peripheral, coronary, and neuro vasculature for the intravascular introduction of interventional/diagnostic devices.

*Continued on next page*

**510(k) Summary, Continued****VI.  
Comparison of  
Technological  
Characteristics  
with the  
Predicate  
Device**

**Table 3** below provides a comparison of technological characteristics of the subject and predicate devices.

<b>Table 3. Comparison of the Subject and Predicate Device</b>		
<b>Characteristics</b>	<b>Predicate Device: YOGA Microcatheter (K162563)</b>	<b>This Submission: YOGA 32XB and YOGA 32XXB Microcatheters</b>
<b>Intended Use</b>	The YOGA Microcatheter is intended for use in peripheral, coronary, and neuro vasculature for the intravascular introduction of interventional/diagnostic devices.	Same as Predicate
<b>Product Code</b>	DQY	Same as Predicate
<b>Classification</b>	21 CFR 870.1250 – Class II	Same as Predicate
<b>Sterilization Method</b>	Ethylene Oxide	Same as Predicate
<b>Sterilization Assurance Level (SAL)</b>	10 <sup>-6</sup>	Same as Predicate
<b>Length</b>	150cm	Same as Predicate
<b>Internal Diameter</b>	0.024", 0.028", 0.032"	Same as Predicate
<b>Outer Diameter</b>	2.4F – 3.4F	Same as Predicate
<b>Polymers</b>	Pebax, Vestamid, Nylon compounds	Same as Predicate
<b>Reinforcement shaft</b>	Stainless Steel/PTFE	Same as Predicate
<b>Marker Band Material</b>	90% Platinum / 10% Iridium	Same as Predicate
<b>Hub</b>	Grilamid	Same as Predicate
<b>Strain Relief</b>	Pebax	Same as Predicate
<b>Packaging</b>	hoop, pouch, carton	Same as Predicate
<b>Shelf Life</b>	3 years	Same as Predicate

**VII.  
Performance  
Data****Performance Testing - Bench**

Appropriate testing was identified based on design, risk analyses and the intended use of the predicate YOGA Microcatheter which was cleared under K162563. The following performance data are being provided in support of the substantial equivalence determination. All testing was conducted using sampling methods as required by Codman & Shurtleff, Inc. Design Control procedures. The bench testing included the following tests:

<b>Table 4: Bench Test Summary</b>		
<b>Test</b>	<b>Test Method Summary</b>	<b>Result</b>
Linear Stiffness Test	Linear Stiffness was measured using Instron tester, holding fixture and ruler	PASS: Samples met the established acceptance criteria
Lateral Stiffness Test	Lateral stiffness was measured using the dynamic three point bend tester and a pin gage	PASS: Samples met the established acceptance criteria
Track Testing	Trackability measured the force to push each device through a representative tortuous anatomical model	PASS: Samples met the established acceptance criteria

*Continued on next page*

---

## 510(k) Summary, Continued

---

**VII.  
Performance  
Data, Continued**

**Performance Testing - Animal**

No *in vivo* testing was required as appropriate verification and validation of the catheter and packaging modifications were achieved based on the similarities of the proposed device to the predicate device, and from results of bench testing.

**Performance Testing - Clinical**

No clinical studies were required as appropriate verification and validation of the catheter and packaging modifications were achieved based on the similarities of the proposed device to the predicate device, and from results of bench testing.

**Sterilization**

The YOGA Microcatheter is sterilized using a validated Ethylene Oxide sterilization process to ensure sterility assurance level (SAL) of  $10^{-6}$  in accordance with ISO 11135-1.

**Shelf-Life Testing**

The YOGA Microcatheter will have a shelf life of 3 years based on the shelf life of currently cleared products.

**Biocompatibility Testing**

Based on the same materials used to manufacture the YOGA Microcatheters (K162563), it was determined that the biocompatibility data from the predicate device can be leveraged. The tests were previously completed in accordance with International Standard ISO 10993-1 "*Biological Evaluation of Medical Devices – Part 1: Evaluation of Testing within a Risk Management Process*", FDA Bluebook Memorandum G95-1, and FDA's Guidance Document entitled "Use of International Standard ISO 10993, Biological Evaluation of Medical Devices Part 1: Evaluation and Testing" issued April 23, 2013.

**Conclusion**

Based upon the intended use, design, materials, function, side-by-side *in-vitro* testing and packaging testing it is concluded that the subject device, YOGA Microcatheter, is substantially equivalent to the predicate device YOGA Microcatheter cleared under 510(k) K162563. Therefore, the modified YOGA Microcatheter does not raise different issues of safety and effectiveness.

---