



October 20, 2017

Grason-Stadler Inc.  
Amy Yanta  
Director of Regulatory Affairs  
10395 West 70th Street  
Eden Prairie, MN 55344

Re: K172403  
Trade/Device Name: GSI Novus  
Regulation Number: 21 CFR 874.1050  
Regulation Name: Audiometer  
Regulatory Class: Class II  
Product Code: EWO, GWJ  
Dated: August 9, 2017  
Received: August 9, 2017

Dear Amy Yanta:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education (DICE) at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education (DICE) at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

  
**Srinivas Nandkumar -S**

for Malvina B. Eydelman, M.D.  
Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K172403

Device Name

GSI Novus

Indications for Use (Describe)

The GSI Novus is intended to be used for the measurement and automated analysis of auditory evoked responses (auditory brainstem responses, ABR) and/or otoacoustic emissions (distortion product, DPOAE and transient evoked, TEOAE). These measures are useful in the screening evaluation, identification, documentation and diagnosis of auditory and hearing related disorders. The auditory evoked response (ABR) measurement is intended for newborns and infants up to 6 months of age. The otoacoustic emissions (DPOAE and/or TEOAE) measurement is intended for use in patients of all ages.

The GSI Novus is intended to be used by a healthcare professional such as an ENT doctor, nurse or audiologist or by a trained technician under the supervision of a professional. The device is intended to be used in a hospital, clinic, or other facility with a suitable quiet testing environment.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**Section 5: 510(k) Summary as required by 21 CFR 807.92.**

Administrative Information

---

Submitter: Grason-Stadler Inc.  
10395 West 70<sup>th</sup> Street  
Eden Prairie, MN 55344  
Tel: 952-947-6097  
Fax: 952-278-4401

Contact Person: Amy Yanta  
Director of Regulatory Affairs  
10393 West 70<sup>th</sup> Street  
Eden Prairie, MN 55344  
952-947-6097  
amy@diagnostic-group.us

Date Summary Prepared: August 7, 2017

Device Identification

---

Trade Name: Novus  
Common Name: otoacoustic emission device/evoked response auditory stimulator and measurements  
Device Classification Name: Audiometer/Stimulator, Auditory, Evoked Response  
Device classification: Class II  
Panel: Ear Nose & Throat / Neurology  
Classification Regulation: 874.1050/882.1900  
Product Code: EWO/GWJ

Predicate Device 1: Titan (with DPOAE440 and ABRIS440), cleared on 05/05/2011 via K103760  
Predicate Device 2: Titan (with TEOAE440), cleared on 06/20/2013 via K130795

Device Description

---

The device is audiometric equipment used for testing of inner ear and auditory brainstem abnormalities.

Novus features a touch-screen display and user-friendly software in a compact hardware design. Novus can be purchased with various licenses allowing you to perform different hearing screening tests.

Novus uses auditory brainstem response (ABR) technology to screen patients for hearing loss. A modified click stimulus, the CE-Chirp<sup>®</sup>, of 35 dB nHL is delivered into the patient's ear while electrodes placed on the patient's head measure EEG activity.

The EEG is processed and analyzed automatically using the Novus's response detection algorithm. When a response is detected, the screening is stopped automatically and a Pass result is assigned to the test ear. When no response is detected after 3 minutes of EEG activity has been processed, a Refer result is assigned.

Auditory brainstem response (ABR) test produces a short acoustic stimulus and measures via transcutaneous electrodes the auditory evoked potentials from the inner ear, the auditory nerve and the brainstem.

Distortion product otoacoustic emissions (DPOAE) technology uses pairs of pure tones presented in sequence to screen patients for cochlear hearing loss. Responses to the stimulus are predictable and therefore can be measured via a sensitive microphone placed in the patient's ear canal.

Transient otoacoustic emissions (TEOAE) technology uses a click stimulus to screen patients for cochlear hearing loss. Responses to the stimulus are predictable and therefore can be measured via a sensitive microphone placed in the patient's ear canal. The response can be divided into frequency bands for assessment.

#### Device Intended Use

---

The GSI Novus is intended to be used for the measurement and automated analysis of auditory evoked responses (auditory brainstem responses, ABR) and/or otoacoustic emissions (distortion product, DPOAE and transient evoked, TEOAE). These measures are useful in the screening evaluation, identification, documentation and diagnosis of auditory and hearing related disorders.

The auditory evoked response (ABR) measurement is intended for newborns and infants up to 6 months of age. The otoacoustic emissions (DPOAE and/or TEOAE) measurement is intended for use in patients of all ages.

The GSI Novus is intended to be used by a healthcare professional such as an ENT doctor, nurse or audiologist or by a trained technician under the supervision of a professional. The device is intended to be used in a hospital, clinic, or other facility with a suitable quiet testing environment.

#### Technological Characteristics

---

The Novus consists of a handheld unit that utilizes a touchscreen display and a rechargeable battery. A simple cradle is included to support charging of the device's battery. The device supports Bluetooth® communication with a label printer for the purpose of printing screening results.

A comparison between the new and predicate devices shows that the technological characteristics and indications for use are equivalent. The device employs similar technology to accomplish the same tasks as the predicates. A detailed table is provided below.

Equivalence Predicate Chart 1:

Description	Titan with TEOAE440 (k130795)	Novus
Type	Audiometer	Audiometer
Regulation Number	21 CFR 874.1050 (otoacoustic emission device)	21 CFR 874.1050 (otoacoustic emission device)
Classification Product Code	EWO	EWO
Regulatory Class	Class II	Class II
Indications for Use	The Titan with TEOAE440 is intended for use in the audiologic evaluation and documentation of ear disorders	The GSI Novus is intended to be used for the measurement and automated analysis of

	<p>using Transient Evoked Otoacoustic Emissions.</p> <p>The target population for Titan with TEOAE440 includes all ages.</p>	<p>auditory evoked responses (auditory brainstem responses, ABR) and/or otoacoustic emissions (distortion product, DPOAE and transient evoked, TEOAE). These measures are useful in the screening evaluation, identification, documentation and diagnosis of auditory and hearing related disorders. The auditory evoked response (ABR) measurement is intended for newborns and infants up to 6 months of age. The otoacoustic emissions (DPOAE and/or TEOAE) measurement is intended for use in patients of all ages. The GSI Novus is intended to be used by a healthcare professional such as an ENT doctor, nurse or audiologist or by a trained technician under the supervision of a professional. The device is intended to be used in a hospital, clinic, or other facility with a suitable quiet testing environment.</p> <p>Intended use differs due to the device encompassing all functions into one device.</p>
Target Population	The devices are suitable for all populations including new-born infants	The device are suitable for all populations including new-born infants
Intended User	The Titan System is to be used by trained personnel only such as audiologists, ENT surgeons, doctors, hearing healthcare professionals or personnel with a similar level of education.	To be used by trained personnel only such as audiologists, ENT surgeons, doctors, hearing healthcare professionals or personnel with similar level of education. (or trained user with supervision of a professional)
TEOAE Stimulus		
Frequency Range	500 to 5500Hz	500 to 5500Hz

Stimuli Type	Non-Linear and Linear Short duration signal (Click) According to IEC 60645-3	Non-Linear and Linear Short duration signal (Click) According to IEC 60645-3
Level	30 to 90 dB peSPL	30 to 90 dB peSPL
Level Step	1 dB SPL	1 dB SPL
Transducer	Dedicated OAE Probe	Dedicated OAE Probe
Probe Detection	Auto detection	Auto detection
Recording		
A/D Resolution	24 bit	24 bit
Artifact Reject System	0 -> +60 dB SPL or off	0 -> +60 dB SPL or off
Automatic test with display of PASS-REFER	Yes	Yes

## Equivalence Predicate Chart 2:

Description	Titan (K103760)		Novus	
	With ABRIS440	With DPOAE440	With ABRIS	With DPOAE
Type	Auditory Brainstem Response – Audiometric equipment	Audiometer	Auditory Brainstem Response – Audiometric equipment	Audiometer
Regulation Number	21 CFR 882.1900 (Evoked response auditory stimulator)	21 CFR 874.1050 (otoacoustic emission device)	21 CFR 882.1900 (Evoked response auditory stimulator)	21 CFR 874.1050 (otoacoustic emission device)
Classification Product Code	GWJ	EWO	GWJ	EWO
Regulatory Class	Class II	Class II	Class II	Class II
Indications for Use	The Titan with ABRIS440 is intended for use in the audiologic evaluation and documentation of ear and nerve disorders using auditory evoked potentials from the inner ear, the auditory nerve and the brainstem.	The Titan with DPOAE440 is intended for use in the audiologic evaluation and documentation of ear disorders using Distortion Product Otoacoustic Emissions.	The GSI Novus is intended to be used for the measurement and automated analysis of auditory evoked responses (auditory brainstem responses, ABR) and/or otoacoustic emissions (distortion product, DPOAE and transient evoked, TEOAE). These measures are useful in the screening evaluation, identification, documentation and diagnosis of auditory	The GSI Novus is intended to be used for the measurement and automated analysis of auditory evoked responses (auditory brainstem responses, ABR) and/or otoacoustic emissions (distortion product, DPOAE and transient evoked, TEOAE). These measures are useful in the screening evaluation, identification, documentation and diagnosis of auditory

			<p>and hearing related disorders. The auditory evoked response (ABR) measurement is intended for newborns and infants up to 6 months of age. The otoacoustic emissions (DPOAE and/or TEOAE) measurement is intended for use in patients of all ages. The GSI Novus is intended to be used by a healthcare professional such as an ENT doctor, nurse or audiologist or by a trained technician under the supervision of a professional. The device is intended to be used in a hospital, clinic, or other facility with a suitable quiet testing environment.</p> <p>Intended use differs due to the device encompassing all functions into one device.</p>	<p>and hearing related disorders. The auditory evoked response (ABR) measurement is intended for newborns and infants up to 6 months of age. The otoacoustic emissions (DPOAE and/or TEOAE) measurement is intended for use in patients of all ages. The GSI Novus is intended to be used by a healthcare professional such as an ENT doctor, nurse or audiologist or by a trained technician under the supervision of a professional. The device is intended to be used in a hospital, clinic, or other facility with a suitable quiet testing environment.</p> <p>Intended use differs due to the device encompassing all functions into one device.</p>
Target Population	Children and newborn	The patient group includes all ages and sexes.	Children and newborn	The patient group includes all ages and sexes.
Anatomical Sites	Examination of Ear and hearing nerves	Examination of Ear	Examination of Ear and hearing nerves	Examination of Ear
Safety Standards	IEC 60601-1	IEC 60601-1	IEC 60601-1	IEC 60601-1
Performance standard	IEC 60645-7	IEC 60645-6	IEC 60645-7	IEC 60645-6
Device Type	Screening device (PASS/REFER) result	Screening and diagnostic	Screening and diagnostic device (PASS/REFER) result	Screening and diagnostic device (PASS/REFER) result
System Configuration	1 -channel ABR system operated through a handheld base unit. The base unit can be operated stand		1 -channel ABR system operated through a handheld base unit. The	



	alone or PC controlled through USB or Bluetooth.		base unit can be operated as stand-alone device.  This allows the freedom of flexibility. The data can be transferred to the PC via USB port for secure transmission.	
Display Information	PASS/REFER status, indicated with value between 0 and 100% where 100% indicates a pass. EEG peak or RMS value, rejection status, residual noise and what transducer(s) are detected.		Noise status, detected transducer, artifact %, PASS/REFER status.  This creates simplicity for user use.	
Stimulus	Click and Chirps	2 pure tones	Click and Chirps	2 pure tones
Electrode quality check	YES		YES	
Impedance Test	Before recording: Electrode impedance is measured if they are above 10kOhm, below 10 kOhm or below 3 kOhm.		Similar impedance test; acceptable impedance <40kOhm.  Below <40kOhm is acceptable value and one test value provides quicker values to reference.	
Binaural screening	YES		YES	
Pre-amplifier channels	1		1	
Stimulus rate	90/s		90/s	
Pre-amplifier Gain	64 dB (fixed) + 64dB (Variable)		64 dB (fixed) + 64dB (Variable)	
Stimulus Level	30,40 and 45dB HL	30 dB SPL to 80 dB SPL	30,40 and 45dB HL	30 dB SPL to 80 dB SPL
Masking	None		None	
Transducers	Titan Probe (mono) Stereo headset: TDH39 and EAR3A		OAE probe, RadioEar IP30 insert earphone  Titan used different accessories but Novus accessories perform same functions.	
				500Hz – 10kHz

### Summary of Non-Clinical Testing

---

Design verification and validation were performed according to current standards for OAE and ABR to assure the device meets its performance specifications. EMC and Safety was performed in compliance with recognized standards IEC 60601-1 series, Medical Electrical Equipment – General requirements for basic safety and essential performance. The product meets the requirements from the international standard for OAE measurements IEC 60645 series. Software verification and validation testing were conducted and documentation was provided as recommended by FDA’s Guidance for Industry and FDA Staff, “Guidance for the Content of Premarket Submissions for Software Contained in medical Devices.” The software for this device was considered as a “moderate” level of concern since a malfunction of, or a latent design flaw in, the Software Device could lead to an erroneous diagnosis or a delay in delivery of appropriate medical care that would likely lead to Minor Injury. The OAE and ABRIS measurements were divided into 3 phases. Phase 1 included when optimization occurred and involved feedback to the operator so that they could adjust such as probe fit, electrode impedance, ambient electrical and acoustic noise etc. Once the pre-test conditions were optimized, phase 2 of data collection proceeded as rapid as possible to allow the maximum quantity of good quality data to be collected in the shortest possible time. Phase 3 proceeded into the data assessment and decision stage and this ran concurrently with Phase 2 once the predetermined minimum amount of data had been collected. Phase 3 then went into the algorithm descriptions for each TEOAE, DPOAE and ABRIS measurements modes. The detailed information about the validation and verification of PASS/REFER algorithms for the OAE and ABR modules is provided in the GSI Novus Manual, e.g., PASS/REFER Criteria, Sensitivity and Specificity etc. No clinical tests were performed, but based on the fulfillment of the international standards for OAE and ABR we believe the device is safe and effective. The auditory impedance testing characteristics and safety systems were compared and found to be comparable.

### Summary of Clinical Testing

---

Not applicable. Not required to establish substantial equivalence.

### Conclusion

---

We have compared the intended use and performance characteristics with the predicate devices. The Novus was tested according to current standards and the differences found between the devices were related to functionality, not in relation to safety and efficiency. The Novus conforms to the current standards. After analyzing bench testing, safety, EMC, and software validation (with risk analysis) testing we conclude that the Novus is found to be substantially equivalent to the predicate devices in technological characteristics and indications for use.