



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

October 5, 2017

Dornoch Medical Systems  
Michael Wolford  
Regulatory Affairs Sr. Specialist  
200 West Parkway  
Riverside, Missouri 64150

Re: K172481

Trade/Device Name: IntelliCart System  
Regulation Number: 21 CFR 878.4780  
Regulation Name: Powered Suction Pump  
Regulatory Class: Class II  
Product Code: JCX  
Dated: August 11, 2017  
Received: August 16, 2017

Dear Michael Wolford:

This letter corrects our substantially equivalent letter of September 14, 2017

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-

related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education (DICE) at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education (DICE) at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

Michael J. Ryan -S

for Tina Kiang, PhD  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K172481

Device Name  
IntelliCart™ System

Indications for Use (Describe)

The IntelliCart™ System is a self-powered suction/vacuum pump intended to collect and dispose of liquid waste within Hospital Operating Rooms, Pathology Labs, Surgical Outpatient Centers, and Doctor's Offices.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

<b>510(k) Summary</b>	<b>Sponsor:</b>	Dornoch Medical Systems, Inc. 200 North West Parkway Riverside, MO 64150 Establishment Registration Number: 1954182
	<b>Contact:</b>	Michael T. Wolford Regulatory Affairs Sr. Specialist Phone: (330) 364-9411
	<b>Date:</b>	8/11/2017
	<b>Trade Name:</b>	<i>IntelliCart™ System</i>
	<b>FDA Product Code/Device Common Name:</b>	JCX – Apparatus, Suction, Ward Use, Portable, Ac-Powered
	<b>Regulation Number and Description; Device Classification:</b>	21 CFR 878.4780 – Powered suction pump; Class II
	<b>Predicate Device</b>	<i>IntelliCart™ System (K162421)</i>
	<b>Device Description</b>	<p>The IntelliCart™ System is a self-powered suction / vacuum pump intended to collect and dispose of liquid waste within hospital operating rooms, pathology labs, surgical outpatient centers, and doctor's offices.</p> <p>The IntelliCart™ System (System) consists of a Fluid Cart (Cart) and Evac Station (Evac).</p> <p>Carts are mobile devices used during surgery to collect patient fluids. Evacs are stationary devices used to process and clean Carts. The System's closed design significantly reduces employee exposure to potentially infectious surgical fluids, while eliminating up to 70% of operating room red bag waste.</p> <p>Carts are equipped with large volume fluid reservoirs, a single regulated vacuum pump, and a height-adjustable IV pole. A single-use manifold (REF 00-5140-201-00) serves as the interface between a Cart's fluid reservoir and sterile patient suction tubing. A Cart model equipped with a bracket to hold a Smoke Evacuator is also available.</p>
	<b>Intended Use:</b>	The <i>IntelliCart™ System</i> is a self-powered suction/vacuum pump intended to collect and dispose of liquid waste within Hospital Operating Rooms, Pathology Labs, Surgical Outpatient Centers, and Doctor's Offices.
	<b>Technological Characteristics</b>	The <i>IntelliCart™ System</i> is a closed self-powered suction device providing various levels of suction, as selected by the user, ranging from 0-600 mmHg

at maximum flow volume, and collection of up to 34L of liquid surgical waste into cart canisters equipped with a viewing window for display of canister contents shown in 100ml increments. The suction cart is portable, equipped with caster type wheels with a profile able to pass through 36 inch minimum wide doorways and allows for connection to a standard 120V/60Hz/15A facility power supply and connection to facility wall suction when necessary. Cart Canisters are integrated with a window tint feature that visually restricts canister contents during transport.

The *IntelliCart™ System* is operated using a back lit LCD touch screen display monitor equipped with lateral adjustment capabilities and provides digital display confirmation of total fluid capacity used as well as a functional interface for suction selection/adjustment and cart cleaning.

Single use manifolds are equipped with hydrophobic filters which turn color after wash cycle is completed, giving visual indication of a manifold which must be discarded.

The Evac Station is equipped with a coupler that connects to the Suction Cart and, using compatible bleach and enzymatic solutions to break down solid waste, empties, cleans and disinfects cart canisters. After emptying and cleaning cart canisters, the Evac Station allows for draining surgical waste into facility drains.

The *IntelliCart™ System* can be equipped with a powered or manual IV pole, supporting up to 12,000mL of bagged fluid and allows mounting of a compatible smoke evacuation device.

**Comparison to Predicate:**

The *IntelliCart™ System* is substantially equivalent to the legally marketed predicate device *IntelliCart™ System* in that these devices are identical in intended use, materials and design. The following tables provide a comparison between the predicate device and modified device showing similarities and differences:

**System Use**

Property	Predicate Device ( IntelliCart™ System (K162421)	Modified Device (IntelliCart™ System)
<b>Intended Use</b>	The IntelliCart™ System is a self-powered suction/vacuum pump intended to collect and dispose of liquid waste within Hospital Operating Rooms, Pathology Labs, Surgical Outpatient Centers, and Doctor's Offices.	Same
<b>Single Use or Reusable</b>	Cart, Evac – Reusable Manifold, Single Use	Same
<b>Target Population</b>	Ward Use, General Surgical Application	Same

**System Operating Principle**

Property	Predicate Device (IntelliCart™ System (K162421)	Modified Device (IntelliCart™ System)
<b>Display</b>	LCD, Backlit Touch	Same
<b>Suction/Exhaust</b>	Closed suction system with variable suction levels based on user selection providing regulated vacuum from 0-600mmHg using AC powered twin piston vacuum pump passing through single use manifolds and into canisters.	Same
<b>Filtration</b>	HEPA Filtration 99.97%	Same
<b>Cleaning</b>	Cart Canister Processing using automated validated cycle	Cart Canister Processing using automated validated cycle to low level of disinfection
<b>Manifold</b>	Single Use, disposable manifold with Hydrophobic Filter	Same
<b>Software</b>	Validated Software	Same
<b>Electrical Requirements</b>	<b>(Cart)</b> 120VAC, 60Hz, 4.75amp <b>(Evac)</b> Same	Same
<b>Portable</b>	Locking casters and push handle allow cart to be easily moved	Same

**System Design**

Property	Predicate Device (IntelliCart™ System (K162421)	Modified Device (IntelliCart™ System)
<b>Materials</b>	No direct or indirect patient contact	Same
<b>Touch Screen</b>	LCD backlit color touch screen display	Same
<b>Total Fluid Capacity</b>	34L	Same

<b>Single Use Manifold</b>	Manifold hydrophobic filter turns blue after use, indicating it must be replaced	Same
<b>Suction Range</b>	0-600mmHg	Same
<b>Footprint</b>	55”H x 23.5”W x 23”D	Same
<b>Safety Alert</b>	Audible and visual display for early warning of overflow protection system and vacuum start up alert	Same
<b>Compatible with Wall Suction</b>	Yes	Same
<b>ConSeal™ Tint</b>	Yes	Same
<b>Canister Exterior Graduations</b>	Graduations marked on canister exterior provide approximate visual indication of contents in 100mL increments.	Same
<b>IV Pole</b>	Carts available with powered or manual IV pole, capable of 12,000mL maximum fluid capacity	Same
<b>External Mounting of compatible smoke evacuation</b>	Yes	Same
<b>HEPA Filtration</b>	99.97%	Same
<b>Cleaning</b>	Multiple user selected Enzymatic and Bleach cleaning cycles from 3-12 minutes	Same
<b>Safety Features</b>		
<b>Feature</b>	<b>Predicate Device (IntelliCart™ System (K162421))</b>	<b>Modified Device (IntelliCart™ System)</b>
<b>Closed System Design</b>	Confine and contain fluids during collection and disposal	Same
<b>Interactive Controls</b>	Graphical User Interface (GUI) guides user through system operation	Same
<b>Bleach Cycle Monitoring</b>	Continuous automatic monitoring of bleach cycle using electronic sensors	Same
<b>Fluid Backflow Protection</b>	Positive air gap maintained between collected fluids and single-use Manifolds	Same
<b>Vacuum System Protection</b>	Hydrophobic filters shut off vacuum supply when wet, preventing liquids from contaminating the Cart vacuum supply	Same

<b>Air Exhaust Protection</b>	Replaceable HEPA rated filter	Same
<b>Overflow Protection</b>	Warning to alert user when reservoir is 3000mL and 1500mL from full capacity – automatic shut off when canisters have reached full capacity	Same
<b>ES 60601-1 Classification</b>	Class 1, Type B applied part (suction tubing – not supplied by Zimmer Biomet)	Same

### Non-Clinical Performance Data as Compared with the Predicate Device

The following performance standards are applicable to the predicate and modified devices:

Standard or Test Type	Predetermined Acceptance Criteria	Results	
		Predicate Device (IntelliCart™ System (K162421))	Modified Device (IntelliCart™ System)
<b>IEC 60601-1</b> Medical electrical equipment - Part 1: General requirements for basic safety and essential performance	Yes	Pass	Pass
<b>IEC 60601-1-2</b> Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and Tests	Yes	Pass	Pass
<b>IEC 60601-1-6</b> Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability	Yes	Pass	Pass
<b>ISO 10079-1</b> Medical suction equipment - Part 1: Electrically powered suction equipment	Yes	Pass	Pass
<b>ISO 10993-1</b> Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Not required – no direct or indirect patient contact	N/A	N/A

### Electrical Safety, Electromagnetic compatibility (EMC)

Medical device safety testing was conducted in accordance with IEC 60601-1 standard for safety, IEC 60601-1-2 standard for EMC. Human factors were also considered in the design of the IntelliCart™ System and testing was performed in accordance with IEC 60601-1-6 and Applying

Human Factors and Usability Engineering to Optimize Medical Device Design, Guidance for Industry and Food and Drug Administration Staff, February 3, 2016. Critical and essential tasks were identified as part of the studies. All identified issues were managed and mitigated to an acceptable level. The IntelliCart™ System passed all electrical safety, and electromagnetic compatibility tests to predetermined acceptance criteria.

### **Software and Display Verification and Validation**

Software verification and validation were performed at a systems and unit level and information supplied in accordance with guidance as recommended by FDA Guidance for Industry “Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices”. Software passed all predetermined acceptance criteria including reservoir overflow protection, cart reservoir processing and low level disinfection using validated cycles, watchdog safety transducer and safety start up warning screen.

### **Mechanical Verification and Validation**

#### **Suction Performance**

The IntelliCart™ System suction performance was designed and evaluated according to ISO 10079-1 Electrically powered suction equipment and as recommended by FDA “Guidance for Industry and FDA Reviewers/Staff: Guidance Document for Powered Suction Pump 510(k)s”. All critical and essential tasks were identified and mitigated through risk analysis to an acceptable level.

The IntelliCart™ System was evaluated to provide adjustable vacuum pressure as measured at the manifold port to predetermined acceptance criteria to demonstrate operation capable of providing adjustable vacuum up to 12 hours of continuous use at 600mmHg.

#### **Safety Features**

The IntelliCart™ System was tested to predetermined acceptance criteria for Closed System Design, Interactive controls including de-bounce feature, Fluid Backflow Protection, Watchdog Safety Transducer, Wash Fluid leak free seals around covers and manifolds, and User Interface Vacuum Start Up Warning. The IntelliCart™ System met all predetermined acceptance criteria with passing results.

The IntelliCart™ fluid path was tested to predetermined acceptance criteria for chemical resistance against compatible enzymatic cleaners, bleach and Bactisure™. The IntelliCart™ met all predetermined acceptance criteria with passing results.

The IntelliCart™ System disposable manifold was tested to predetermined acceptance criteria to provide vacuum seal in the reservoir when installed. The IntelliCart™ met all acceptance criteria with passing results.

### **IV Pole Function and Support**

Manual and Power IV Poles were tested to predetermined acceptance criteria for maximum bagged fluid support and automatic function (power IV pole only). The IntelliCart™ met all acceptance criteria with passing results.

### **Outer Body Shell and Fluid Capacity**

The IntelliCart™ was tested to predetermined acceptance criteria for portability using casters and brakes including profile requirements to pass through standard 36" door thresholds. The IntelliCart™ met all acceptance criteria with passing results.

IntelliCart™ was tested to predetermined acceptance criteria to include directly viewable fluid measurements, easily viewable touchscreen monitor and exterior surfaces which provide chemical resistance to standard hospital wipes. The IntelliCart™ met all acceptance criteria with passing results.

### **Fluid Reservoir Concealment**

ConSeal™ technology providing concealment of canister contents was evaluated to predetermined acceptance criteria for the concealment of collected waste fluid. The IntelliCart™ met all acceptance criteria with passing results.

### **Conclusion**

The subject device is substantially equivalent to the predicate device.