



December 14, 2017

Toshiba Medical Systems Corporation
% Ms. Janine Reyes
Manager, Regulatory Affairs
Toshiba America Medical Systems, Inc.
2441 Michelle Drive
TUSTIN CA 92780

Re: K172863

Trade/Device Name: Infinix-i, INFX-8000V, V7.0
Regulation Number: 21 CFR 892.1650
Regulation Name: Image-intensified fluoroscopic x-ray system
Regulatory Class: II
Product Code: OWB
Dated: September 20, 2017
Received: September 20, 2017

Dear Ms. Reyes:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Robert Ochs, Ph.D.
Director
Division of Radiological Health
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES
Food and Drug Administration

Form Approved: OMB No. 0910-0120
Expiration Date: January 31, 2017
See PRA Statement below.

Reset Form

Indications for Use

510(k) Number (if known)

K172863

Device Name

Infinix-i, INFX-8000V, V7.0

Indications for Use (Describe)

This device is a digital radiography/fluoroscopy system used in a diagnostic interventional angiography configuration. The system is indicated for use in diagnostic and angiographic procedures for blood vessels in the heart, brain, abdomen and lower extremities.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”

510(k) SUMMARY

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of Safe Medical Device Act 1990 and 21 CFR § 807.92

1. CLASSIFICATION and DEVICE NAME

Classification Name	Image-intensified fluoroscopic X-ray system
Primary Product Code	OWB
Regulation Number	21 CFR 892.1650
Regulatory Class	Class II
Trade Proprietary Name/Model Number	Infinix-i, INFX-8000V, V7.0

2. SUBMITTER’S NAME

Toshiba Medical Systems Corporation (TMSC)
1385 Shimoishigami
Otawara-Shi, Tochigi-ken, Japan 324-8550

3. OFFICIAL CORRESPONDENT

Naofumi Watanabe
Senior Manager, Regulatory Affairs and Vigilance Department

4. CONTACT PERSON, U.S. AGENT and ADDRESS

Contact Person

Janine Reyes
Manager, Regulatory Affairs
Toshiba America Medical Systems, Inc.
2441 Michelle Drive, Tustin, CA 92780
Phone: (714) 669-7853
Fax: (714) 730-1310
jfreyes@tams.com

Official Correspondent/U.S. Agent

Paul Biggins
Director, Regulatory Affairs
Toshiba America Medical Systems, Inc.
2441 Michelle Drive, Tustin, CA 92780
Phone: (714) 730-7808
Fax: (714) 730-1310
pbiggins@tams.com

5. MANUFACTURING SITE

Toshiba Medical Systems Corporation (TMSC)
1385 Shimoishigami
Otawara-shi, Tochigi 324-8550, Japan

6. ESTABLISHMENT REGISTRATION

9614698

7. DATE PREPARED

September 20, 2017

8. DEVICE MODEL NAME

INFX-8000V, V7.0 [Infinix CF-i/SP, Infinix CF-i/BP, Infinix VF-i/SP, Infinix VF-i/BP]

9. TRADE NAME

Infinix-i, INFX-8000V, V7.0

10. DEVICE NAME

Interventional Fluoroscopic X-ray System

11. CLASSIFICATION PANEL

Radiology

12. DEVICE CLASSIFICATION

Class II (per 21 CFR 892.1650)

13. PRODUCT CODE / DESCRIPTION

Product Code: OWB – Image-intensified fluoroscopic X-ray system

14. PERFORMANCE STANDARD

This device conforms to applicable Performance Standards for Ionizing Radiation Emitting Products [21 CFR Subchapter J, Federal Diagnostic X-ray Equipment Standard].

15. PREDICATE DEVICE

Infinix, INFX-8000V, V6.35 (K162614)

Product	Marketed by	510(k) Number	Clearance Date
INFX-8000V, V6.35	Toshiba America Medical Systems	K162614	October 17, 2016

PREDICATE DEVICE CLASSIFICATION and DEVICE NAME

Regulation Description	Image-Intensified Fluoroscopic X-ray System
Regulation Number	21 CFR 892.1650
Regulatory Class	Class II
Primary Product Code	OWB

16. REASON FOR SUBMISSION

Modification of a cleared device

17. SUBMISSION TYPE

Traditional 510(k)

18. DEVICE DESCRIPTION

INFX-8000V, V7.0, is an X-ray system that is capable of radiographic and fluoroscopic studies and is used in an interventional setting. The system consists of a C-arm/ Ω -arm which is equipped with an X-ray tube, beam limiter and X-ray receptor, X-ray controller, computers with system and processing software, and a patient radiographic table.

19. INDICATIONS FOR USE

This device is a digital radiography/fluoroscopy system used in a diagnostic and interventional angiography configuration. The system is indicated for use in diagnostic and angiographic procedures for blood vessels in the heart, brain, abdomen and lower extremities.

20. SUMMARY OF CHANGE(S)

This submission is to report image quality improvements, workflow improvements, and DFP (Digital Fluoroscopy Processor) hardware changes.

- **Image Quality Improvements:**
 - **Time Axis SNRF:** Improvement to SNRF processing. Utilizing chronological frame data, image noise is further reduced with minimal lag.
 - **F-SUB (fluoroscopic subtraction) improvement:** artifact reduction during fluoroscopic subtraction.
 - **Real Time Auto Pixel Shift (RAPS):** Real Time Auto Pixel Shift automatically corrects misalignment between the contrast image and mask image during DSA and F-SUB acquisition.
 - **Enhanced Live Zoom image quality**
 - **UNSUB DSA halation reduction:** enables the ability to evaluate positional relationship between vessel and bone in the UNSUB (unsubtracted) display of the DSA image.
 - **16-bit data processing:** maximized FPD output data use.
- **Workflow Improvements:**
 - DFP supports the concurrent operation of up to four streams of display (three streams dynamic display plus export).
 - Reduction of steps and time during F-SUB sequence.
 - Reduction in processing time for both system startup and map saving.
- **DFP (Digital Fluoroscopy Processor):**
 - Hardware changes to enhance both system operability and image quality. Changes include: host system PC, real time controller CPU board, image processing unit and storage.

21. SAFETY

The device is designed and manufactured under the Quality System Regulations as outlined in 21 CFR § 820 and ISO 13485 Standards. This device is in conformance with the applicable parts of the IEC60601-1 standards, its collateral standards and particular standards; IEC 60601-2-43 and IEC60601-2-28. All requirements of the Federal Diagnostic Equipment Standard, as outlined in 21 CFR §1020, that apply to this device, will be met and reported via product report.

LIST OF APPLICABLE STANDARDS

- | | |
|----------------------|-----------------------|
| ● IEC60601-1:2005 | ● IEC60601-2-28: 2010 |
| ● IEC60601-1-2:2007 | ● IEC60601-2-43:2010 |
| ● IEC60602-1-3: 2008 | ● IEC62304:2006 |
| ● IEC60601-1-6: 2010 | ● IEC62366:2007 |

22. TESTING

Risk analysis and verification/validation testing conducted through bench testing demonstrate that the established specifications for the device have been met.

System evaluation of image quality:

Image quality metrics were performed, utilizing phantoms, to assess spatial resolution, low-contrast resolution/dynamic range, fluoroscopic still image resolution/dynamic image resolution/afterimage, artifacts/contrast/dynamic range of DSA, reconstructed image spatial resolution in 3D-DA acquisition, reconstructed image spatial resolution in 3D-DSA acquisition, reconstructed image spatial resolution in 3D-LD acquisition, reconstructed image density resolution in LCI acquisition.

Evaluation of items supported to improve image quality:

Image quality metrics were performed, utilizing phantoms, to assess image quality improvements, DSA halation reduction, fluoroscopic subtraction artifact reduction, Live Zoom interpolation method comparison, SNRF noise reduction improvement comparison, real-time auto pixel shift comparison, 16-bit data LCI image reconstruction improved visualization comparison.

This submission contains test data that demonstrates that the system modifications result in performance that is equal to or better than the predicate system. Testing of the modified system was conducted in accordance with the applicable standards published by the International Electromechanical Commission (IEC) for Medical Devices and XR Systems.

Software Documentation for a Moderate Level of Concern, per the FDA guidance document, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices Document" issued on May 11, 2005, is also included as part of this submission.

Cybersecurity documentation, per the FDA cybersecurity premarket guidance document "Content of Premarket Submissions for Management of Cybersecurity in Medical Devices" issued on October 2, 2014, is also included as part of this submission.

Additionally, the design controls used for this device included risk management and all known risks were mitigated to an acceptable level.

23. SUBSTANTIAL EQUIVALENCE

This device is substantially equivalent to the INFX-8000V, V6.35, (K162614), marketed by Toshiba America Medical Systems. INFX-8000V, V7.0, includes image quality improvements, Time Axis SNRF, F-SUB improvement, Real Time Auto Pixel Shift, enhanced live zoom image quality, UNSUB DSA halation reduction, 16-bit data processing, workflow improvements, and DFP hardware changes (host system PC, real time controller CPU board, image processing unit

and storage). The basic system configuration, method of operation, base software and manufacturing process remain unchanged from the cleared device. There are no new indications for use or intended use of the device.

24. CONCLUSION

The subject device, INFX-8000V, V7.0, is substantially equivalent to the INFX-8000V, V6.35, which was cleared via Pre-Market Notification 510(k), K162614. The INFX-8000V, V7.0, incorporates modifications to the cleared device which include image quality improvements, Time Axis SNRF, F-SUB improvement, Real Time Auto Pixel Shift, enhanced live zoom image quality, UNSUB DSA halation reduction, 16-bit data processing, workflow improvements, and DFP hardware changes (host system PC, real time controller CPU board, image processing unit and storage). The changes to this device do not alter the indications for use or the intended uses associated with the previously cleared device, as described in the labeling. It is the contention of Toshiba that all new safety issues have been addressed in the design of this change and that adequate evidence of this is presented with this submission.