



June 14, 2018

Jiangsu Caina Medical Co., Ltd.
% Diana Hong
General Manager
Mid-link Consulting Co., Ltd
P.O.Box 120-119
Shanghai, 200120, China

Re: K172938

Trade/Device Name: Disposable Sterile Needle
Regulation Number: 21 CFR 880.5570
Regulation Name: Hypodermic Single Lumen Needle
Regulatory Class: Class II
Product Code: FMI
Dated: May 4, 2018
Received: May 8, 2018

Dear Diana Hong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820);

and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Alan M.
Stevens -S

Digitally signed by Alan M.
Stevens -S
DN: c=US, o=U.S. Government,
ou=HHS, ou=FDA, ou=People,
0.9.2342.19200300.100.1.1=130
0189211, cn=Alan M. Stevens -S
Date: 2018.06.14 13:35:06 -04'00'

for Tina Kiang, Ph.D.
Acting Director
Division of Anesthesiology,
General Hospital, Respiratory,
Infection Control, and Dental Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K172938

Device Name

Disposable Sterile Needle

Indications for Use (Describe)

The Disposable Sterile Needle is intended for use with syringes and injection devices for general purpose fluid injection/aspiration.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K172938 510(k) SUMMARY

Submitter: Jiangsu Caina Medical Co., Ltd.
No.23, Huanxi Road, Zhutang Town, Jiangyin, Jiangsu, 214425,
China

Application Correspondent: Ms. Diana Hong
General Manager
Mid-Link Consulting Co., Ltd
P.O. Box 120-119, Shanghai, 200120, China
Tel: +86-21-22815850,
Fax: 240-238-7587
Email: info@mid-link.net

Preparation Date: June 13, 2018

Trade Name: Disposable Sterile Needle

Common or Usual Name: Needle, Hypodermic, Single Lumen

Regulation Name: Hypodermic Single Lumen Needle.

Regulation Number: 21 CFR 880.5570

Product Code: FMI

Device Class: Class II

Primary Predicate Device: K072739; Jierui Syringes and Needles

Note – K072739 is a bundled 510(k) submission containing a sterile hypodermic needle. K172938 is substantially equivalent to the sterile hypodermic needle component of K072739.

Indications for Use:

The Disposable Sterile Needle is intended for use with syringes and injection devices for general purpose fluid injection/aspiration.

Device Description:

The Disposable Sterile Needle consists of three components: needle tube, needle hub, and the needle cap. The Disposable Sterile Needle is supplied sterile for single use and has a shelf life of 5 years. The proposed device is available in various combinations of needle gauge, length, wall type and Edges. The hub color of the needle complies with ISO 6009: 2016 - Hypodermic needles for single use — Colour coding for identification. The product specifications are listed in following table.

Gauge	Hub Color	Product specification		Gauge	Hub Color	Product specification
33G	Black	33Gx6mm RWxLB		22G	Black	22Gx25mm RWxLB
		33Gx6mm RWxSB				22Gx25mm RWxSB
		33Gx6mm TWxLB				22Gx25mm TWxLB
		33Gx6mm TWxSB				22Gx25mm TWxSB
32G	Deep Green	32Gx13mm RWxLB				22Gx25mm ETWxLB
		32Gx13mm RWxSB				22Gx25mm ETWxSB
		32Gx13mm TWxLB				22Gx38mm RWxLB
		32Gx13mm TWxSB				22Gx38mm RWxSB
		32Gx13mm ETWxLB				22Gx38mm TWxLB
		32Gx13mm ETWxSB				22Gx38mm TWxSB
30G	Yellow	30Gx25mm RWxLB				22Gx38mm ETWxLB
		30Gx25mm RWxSB				22Gx38mm ETWxSB
		30Gx25mm TWxLB		21G	Deep Green	21Gx38mm RWxLB
		30Gx25mm TWxSB				21Gx38mm RWxSB
		30Gx25mm ETWxLB				21Gx38mm TWxLB
		30Gx25mm ETWxSB				21Gx38mm TWxSB
25G	Orange	25Gx25mm RWxLB				21Gx38mm ETWxLB
		25Gx25mm RWxSB				21Gx38mm ETWxSB
		25Gx25mm TWxLB		20G	Yellow	20Gx38mm RWxLB
		25Gx25mm TWxSB				20Gx38mm RWxSB

		25Gx38mm RWxLB				20Gx38mm TWxLB
		25Gx38mm RWxSB				20Gx38mm TWxSB
		25Gx38mm TWxLB				20Gx38mm ETWxLB
		25Gx38mm TWxSB				20Gx38mm ETWxSB
Note: RW = Regular Wall; TW = Thin Wall; ETW= Extra Thin Wall LB=Long Bevel; SB=Short Bevel			18G	Pink		18Gx50mm RWxLB
						18Gx50mm RWxSB
						18Gx50mm TWxLB
						18Gx50mm TWxSB

Device Comparison of Technology

A technological comparison table is provided below that compares the subject device and predicate device:

Item		Proposed Device –K172938 Disposable Sterile Needle	Predicate Device - K072739 Jierui Syringes and Needles: Sterile Hypodermic Needle Component Only	Comparison
Product Code		FMI	FMI	Same
Intended Use		The Disposable Sterile Needle is intended for use with syringes and injection devices for general purpose fluid injection/aspiration.	The Sterile Hypodermic Needle for single use is intended for use with syringes and injection devices for general purpose fluid injection/aspiration	Same
Single Use		Yes	Yes	Same
Sterility Condition		EO Sterilization	EO Sterilization	Same
Material	Needle	304 Stainless Steel	Stainless Steel	Different
	Tube			
	Cap	Polypropylene	Unknown	Different
	Hub	Polypropylene	Unknown	Different
Needle Gauge		18G,20G,21G,22G,25G,30G,32 G,33G	16G,18G,19G,20G,21G,22G,23G,24G,25G,26G,27G,29G	Different

Needle Length	6mm, 13mm, 20mm, 25mm, 38mm, 45mm, 50mm	13 mm, 16 mm, 19 mm, 25 mm, 32 mm, 38 mm	Different
Needle Edge/Bevel Angle	Long Bevel: angle 11 ° Short Bevel: angle 17 °	Long bevel Short Bevel	Different
Needle Performance	Complies with ISO 7864: 2016 ISO 9626: 2016	Complies with ISO 7864: 1993 ISO 9626: 1991	Same
Luer Connector Performance	Complies with ISO 594-1:1986 ISO 594-2:1998	Complies with ISO 594-1:1986 ISO 594-2:1998	Same
Needle Hub Color	Complies with ISO 6009: 2016	Complies with ISO 6009: 2016	Same

Substantial Equivalence Discussion

The indications for use and intended use of the hypodermic needle component in the predicate bundled 510(k) K072739, are equivalent to the subject 510(k) device.

The differences between the subject device and the predicate device are the following:

- Needle length
Although the needle length range of the subject device is larger than needle length range of the predicate device, the needle length of the subject device meets 4.10.2 Tolerances requirements of ISO 7864. This difference does not raise different questions of safety and effectiveness
- Needle Gauge
The subject device has 30G, 32G and 33G needle gauges, which are not included in the predicate device. All needle gauges of proposed device meet 5.4 Dimensions requirements of ISO 9626. This difference does not raise different questions of safety and effectiveness
- Edge/Bevel angle
Although the /bevel angle of the predicate device is unknown, the edge/bevel angle of subject device are listed in the product labeling, and the end user can choose the needle per their preference. This difference does not raise different questions of safety and effectiveness

- **Materials**

Although the materials of the subject device are different from those of the predicate device, the subject device has been tested for biocompatibility and meets the biocompatibility requirements of the applicable standards in the ISO 10993 series. This difference does not raise different questions of safety and effectiveness

Based on the aforementioned modifications to the subject device, the subject device does not raise different types of safety and effectiveness questions when compared to the predicate device.

Performance Testing Summary

Non-clinical tests were conducted to verify that the proposed devices met all design specifications and are substantially equivalent (SE) to the hypodermic needle component of the predicate device. The test results demonstrated that the proposed devices comply with the following standards:

- ISO 7864:2016, Sterile Hypodermic Needles for Single Use
- ISO 6009: 2016 - Hypodermic needles for single use — Colour coding for identification
- ISO 9626:2016, Stainless Steel Needle Tubing for the Manufacture of Medical Devices
- ISO 594-1:1986 Conical Fittings with A 6% (Luer) Taper for Syringes, Needles and Certain Other Medical Equipment - Part 1: General Requirements.
- ISO 594-2:1998 Conical Fittings with A 6% (Luer) Taper for Syringes, Needles and Certain Other Medical Equipment - Part 2: Lock Fittings.
- ASTM F1929-15 Standard Test Method for Detecting Seal Leaks in Porous Medical Package by Dye Penetration
- ASTM F88/F88M-15 Standard Test Method for Seal Strength of Flexible Barrier Materials.
- USP <85> Bacterial Endotoxins Test
- ISO 11135-1:2007 – Sterilization of Health Care Products- Ethylene Oxide- Part 1: Requirements for Development, Validation and Routine Control of a Sterilization Process for Medical Devices

The Disposable Sterile Needle is supplied sterile for single use and are sterilized by Ethylene Oxide to achieve a SAL of 10^{-6} and supplied in sterility maintenance package which could maintain the sterility of the device during the shelf life of 5 years. Sterility was validated in accordance with ISO 11135-1:2007. Packaging tests were performed per ASTM F88/F88M-15 and ASTM F1929-15.

The Disposable Sterile Needle has been tested for biocompatibility and meets the biocompatibility requirements of the applicable standards in the ISO 10993 series. The device

was tested to the following biocompatibility endpoints: cytotoxicity, sensitization, irritation sensitivity, acute systemic toxicity, and material mediated pyrogenicity, hemocompatibility.

A risk analysis was conducted in accordance with ISO 14971: 2007 – Medical devices — Application of risk management to medical devices.

In all testing, the pre-determined acceptance criteria were met.

Substantially Equivalence Conclusion

The subject device has the same intended use and technological characteristics as the predicate device. The performance of the device is supported by non-clinical testing and risk management activities. The Disposable Sterile Needle is Substantially Equivalent (SE) to the Sterile Hypodermic Needle component of the bundled 510(k), Jierui Syringes and Needles, cleared under K072739.