



April 4, 2018

DeRoyal Industries, Inc.
Sarah Bennett
Regulatory Affairs Specialist
200 DeBusk Lane
Powell, Tennessee 37849

Re: K173072

Trade/Device Name: Algidex Ag Silver Alginate Wound Dressing
Regulatory Class: Unclassified
Product Code: FRO
Dated: March 5, 2018
Received: March 5, 2018

Dear Sarah Bennett:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jennifer R. Stevenson -S3

For Binita S. Ashar, M.D., M.B.A., F.A.C.S.
Director
Division of Surgical Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K173072

Device Name

Algidex Ag Silver Alginate Wound Dressing

Indications for Use (Describe)

Abrasions and lacerations, dermal wounds, donor sites, first- and second-degree burns, surgical incisions, vascular access sites, pressure ulcers stages I-IV, stasis ulcers, venous ulcers.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



Corporate

200 DeBusk Lane
Powell, TN 37849 USA
865.938.7828
800.251.9864
www.deroyal.com

DeRoyal Industries, Inc.
Traditional 510(k) – Algidex Ag Silver Alginate Wound Dressing
September 28, 2017

510(k) Summary

Date prepared: September 28, 2017

510(k) Owner: DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849
Owner/Operator #1044833

510(k) Contact: Sarah Bennett
Regulatory Affairs Specialist
P: 865-362-6112
F: 865-362-3741
sabennett@deroyal.com

Contract Manufacturer: DeRoyal Industries, Inc.
1595 Highway 33 South
New Tazewell, TN 37825
FDA Establishment # 2320762

Trade Name: Algidex Ag Silver Alginate Wound Dressing

Common Name: Wound Dressing

Classification Name: Dressing, Wound, Drug (Pre-Amendment)

Device Product Code: FRO

Regulatory Class: Unclassified

Classification Panel: General & Plastic Surgery

Predicate Devices: SilverSite [K041268] (Primary Predicate)
Calgitrol Ag Silver Alginate Wound Dressing [K011618]

Device Description

Algidex Ag[®] silver alginate wound dressings are an antimicrobial barrier for up to three (3) days. The dressings utilize a formulation of ionic silver combined in an alginate and maltodextrin matrix.



Corporate

200 DeBusk Lane
 Powell, TN 37849 USA
 865.938.7828
 800.251.9864
 www.deroyal.com

The silver is intended to prevent microbes from passing through the dressing during use. The dressings are offered in two configurations:

1. a paste;
2. a paste applied to a non-adherent polyurethane foam layers in various sizes and shapes;

The device is sterilized by gamma irradiation.

Intended Use

Algidex Ag[®] silver alginate wound dressings are intended to manage abrasions and lacerations, dermal wounds, donor sites, first- and second-degree burns, surgical incisions, vascular access sites, pressure ulcers stage I-IV, stasis ulcers, and venous ulcers.

Summary of Technological Characteristics

| Characteristic | Algidex Ag Silver Alginate Wound Dressing | SilverSite (K041268) (Primary Predicate) | Calgitrol Ag Silver Alginate Wound Dressing (K011618) |
|---------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Indications for Use | Abrasions and lacerations, dermal wounds, donor sites, first- and second-degree burns, surgical incisions, vascular access sites, pressure ulcers stages I-IV, stasis ulcers, venous ulcers | Pressure ulcers (stages I-IV), dermal lesions (or secreting skin injuries), venous ulcers, stasis ulcers, intended to protect vascular access sites, intramuscular sites, and surgical incisions, 1 st and 2 nd degree burns, donor sites | Pressure ulcers (stages I-IV), dermal lesions (or secreting skin injuries), venous ulcers, stasis ulcers, 1 st and 2 nd degree burns, donor sites |
| Prescription Only | Yes | Yes | Yes |
| Design | Formulation of ionic silver combined in an alginate and maltodextrin matrix. Placed on the wound in a paste or on a dressing. | Formulation of ionic silver combined in an alginate and maltodextrin matrix. Placed on the wound on a dressing. | Formulation of ionic silver combined in an alginate and maltodextrin matrix. Placed on the wound in a paste or on a dressing. |
| Materials | Silver alginate paste placed on the wound on its own or on a polyurethane foam | Silver alginate paste placed on the wound on a polyurethane foam | Silver alginate paste placed on the wound on its own or on a polyurethane foam |
| Microbial Barrier Effectiveness | Effective for 3 days | Effective for 3 days | Effective for 3 days |



Corporate

200 DeBusk Lane
Powell, TN 37849 USA
865.938.7828
800.251.9864
www.deroyal.com

| | | | |
|------------------|---------------------------------|---------------------------------|---------------------------------|
| Log Reduction | 4 log reduction | 4 log reduction | 4 log reduction |
| Adherence | Non-adherent | Non-adherent | Non-adherent |
| Moisture Content | 2.46 mL | 2.46 mL | 2.46 mL |
| Silver Content | Contains 141 mg silver/100sq cm | Contains 141 mg silver/100sq cm | Contains 141 mg silver/100sq cm |
| Biocompatibility | ISO 10993-1 compliant | ISO 10993-1 compliant | ISO 10993-1 compliant |

Summary of Performance Tests

To ensure performance, Algidex Ag was tested when it was developed and originally cleared for the market in 2004, under the brand names SilverSite (K041268) and Calgitrol AG (K011618). Algidex Ag contains the exact same silver alginate and maltodextrin matrix as the predicate devices, and therefore, the testing performed and submitted in the Calgitrol AG 510(k)s continues to support the performance of this device.

The following tests have been performed and are contained within this premarket submission: AATCC 100 test method, antimicrobial barrier testing, and absorbency test. The results indicate the inhibition of microorganism growth within the dressing over the seven-day wear-time of the dressing. The antimicrobial barrier was shown for three days. The absorbency information confirms the dressing absorbs fluid. The proposed device, predicate devices, and other silver-containing devices are substantially equivalent.

Conclusion

The results of performance testing demonstrate that Algidex Ag silver alginate wound dressings are substantially equivalent to the predicate. It is the exact same silver alginate and maltodextrin matrix that was cleared for market in K041268 and K011618. Thus, Algidex Ag is substantially equivalent to the predicate devices.