



February 12, 2018

STERIS Corporation  
Anthony Piotrkowski  
Senior Manager, Regulatory Affairs  
5976 Heisley Road  
Mentor, Ohio 44060

Re: K173670

Trade/Device Name: CELERITY 20 Steam Incubator  
Regulation Number: 21 CFR 880.2800  
Regulation Name: Sterilization Process Indicator  
Regulatory Class: Class II  
Product Code: FRC  
Dated: November 29, 2017  
Received: November 30, 2017

Dear Anthony Piotrkowski:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820);

and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

  
Michael J. Ryan -S

for Tina Kiang, Ph.D.  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K173670

Device Name

Celerity 20 Steam Incubator

Indications for Use (Describe)

Use the Celerity 20 Steam Incubator to incubate and automatically read Celerity 20 Steam Biological Indicators at 57 °C for a fluorescent result within 20 minutes

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



**510(k) Summary  
For  
K173670  
Celerity 20 Steam Incubator**

**Sponsor Facility**

STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060  
Phone: (440) 354-2600  
Fax No: (440) 357-9198

**Manufacturing Facility**

STERIS Corporation  
9325 Pinecone Drive  
Mentor, OH 44060  
Phone: (440) 392-7800  
Fax No: (440) 392-7896

Contact: Anthony Piotrkowski

Telephone: (440) 392-7437  
Fax No: (440) 357-9198  
e-mail: tony\_piotrkowski@steris.com

Submission Date: January 17, 2018

Premarket Notification Number: K173670

**STERIS TRADITIONAL 510(k) PREMARKET NOTIFICATION  
K173670 Celerity 20 Steam Incubator**

---

**1. Device Name**

Trade Name: Celerity 20 Steam Incubator

Common/usual Name: Incubator/Reader (accessory to Biological Indicator)

Device Classification: Class II

Classification Name: Indicator, Biological Sterilization Process  
(21 CFR 880.2800, FRC)

**2. Predicate Device**

K171587 - VERIFY Incubator for Assert VH2O2 Self-Contained Biological Indicator, (since renamed Celerity 20 HP Incubator)

**3. Description of Device**

The Celerity 20 Steam Incubator (Incubator) is an incubator/reader designed for use specifically with the Celerity 20 Steam Biological Indicator (BI), subject of a concurrent 510(k) currently under review. The incubator/reader provides a constant temperature range to allow for activation and outgrowth of *Geobacillus stearothermophilus* leading to production of a fluorescent moiety. The presence of an increasing fluorescence signal due to increasing concentrations of this fluorescent moiety in the BI is detected by the incubator/reader and indicates the presence of viable test microorganisms.

**4. Intended Use/ Indications for Use**

Use the Celerity 20 Steam Incubator to incubate and automatically read Celerity 20 Steam Biological Indicators at 57 °C for a fluorescent result within 20 minutes.

**5. Summary of Technical Characteristics**

A comparison of technical characteristics are summarized in **Table 5-1**.

**STERIS TRADITIONAL 510(k) PREMARKET NOTIFICATION  
K173670 Celerity 20 Steam Incubator**

**Table 5-1 Summary of Incubator Physical Description and Technological Properties**

Feature	Celerity Steam Incubator (proposed)	VERIFY Incubator (K71587) Predicate	Comparison
Intended Use	Use the Celerity 20 Steam Incubator to incubate and automatically read Celerity 20 Steam Biological Indicators at 57 °C for a fluorescent result within 20 minutes.	Use the VERIFY Incubator for Assert VH2O2 Self Contained Biological Indicators (Incubator) to incubate and automatically read VERIFY Assert VH2O2 Self-Contained Biological Indicators at 57 °C for a fluorescent result within 20 minutes.	Identical except for the BI they are meant to read.  Testing with the Celerity 20 Steam BI, included in the submission, demonstrates proper performance.
Basis of Readout	Photodiode detects fluorescence produced by enzymatic activity that results from growing biological indicator organisms	Photodiode detects fluorescence produced by enzymatic activity that results from growing biological indicator organisms	Same
Incubation Temperature Range	55 - 60 °C	55 - 60 °C	Same
Voltage Range	90 to 264 VAC with 12 VDC conversion.	90 to 264 VAC with 12 VDC conversion.	Same

**6. Summary of Nonclinical Tests**

Performance testing to demonstrate substantial equivalence to the predicate has been completed and is summarized in **Table 5-2** below.

**Table 5-2. Summary of Non-clinical Testing**

Test	Acceptance Criteria	Conclusion
Maintenance of Incubation Temperature	Maintain 55-60 °C for a minimum of 20 minutes (incubation time of Celerity 20 Steam BI)	PASS
Qualification testing with Celerity 20 Steam BI	20-minute fluorescent read meets >97% alignment with 7-day growth results per FDA guidance on reduced incubation time.	PASS

**STERIS TRADITIONAL 510(k) PREMARKET NOTIFICATION  
K173670 Celerity 20 Steam Incubator**

<b>Test</b>	<b>Acceptance Criteria</b>	<b>Conclusion</b>
Qualification testing with Celerity 20 Steam BI	<ul style="list-style-type: none"><li>• Pass testing with BI exposed to full cycle exposure and negative growth result in incubator</li><li>• Fail testing with BI exposed in abbreviated cycle exposure and positive result in incubator</li></ul>	PASS  PASS
Alarm, LED and Print function Test	Demonstrate proper function of alarms, LED and print outputs	PASS

**7. Conclusion**

The Celerity 20 Steam Incubator has met the established performance criteria. The conclusions drawn from the nonclinical tests performed demonstrate the subject device is as safe, as effective, and performs as well or better than the legally marketed predicate device VERIFY Incubator for Assert VH2O2 Self-Contained Biological Indicator cleared in K171587 (21 CFR 880.2800, Product code FRC).