



AmCAD BioMed Corporation
% Chiu S. Lin
President
Lin & Associates, LLC
5614 Johnson Avenue
BETHESDA MD 20817

August 31st, 2018

Re: K180006
Trade/Device Name: AmCAD-UT Detection 2.2
Regulation Number: 21 CFR 892.2050
Regulation Name: Picture Archiving And Communications System
Regulatory Class: Class II
Product Code: LLZ
Dated: August 1, 2018
Received: August 3, 2018

Dear Chiu Lin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



for

Robert A. Ochs, Ph.D.

Director

Division of Radiological Health

Office of In Vitro Diagnostics
and Radiological Health

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K180006

Device Name

AmCAD-UT® Detection 2.2

Indications for Use (Describe)

AmCAD-UT® Detection is a Windows-based computer-aided detection (CADe) device intended to assist the medical professionals in analyzing thyroid ultrasound images, acquired from FDA-cleared ultrasound systems, with user-selected regions of interest (ROI). After the initial review of the ultrasound images by the physicians, the device further provides detailed information with quantification and visualization of sonographic characteristics of thyroid nodules. The device is intended for use on ultrasound images of discrete thyroid nodules larger than 1cm, for which a biopsy recommendation is required.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary of Safety and Effectiveness

This 510(k) summary of safety and effectiveness information is submitted as part of the Premarket Notification in compliance with requirements of CFR Part 807, Subpart E and Section 807.92

5.1 Identification of Submitter:

Submitter: AmCad BioMed Corporation
Address: FL.5-2, NO.167, Fu Hsing N. RD., Taipei 105, Taiwan, R.O.C.
Phone: 886-2-2713-6227
Fax: 886-2-2514-0245
Contact: Jack Yang
Title: Vice President
Phone: 886-2-2713-6227 ext.358
Fax: 886-2-2514-0245
Email: jack.yang@amcad.com.tw
Manufacturer: AmCad BioMed Corporation

US Agent and Contact: Chiu S. Lin, Ph.D.
Lin & Associates, LLC
Address: 5614 Johnson Avenue
Bethesda, MD 20817
Phone: (O) 301-591-3895
E-mail: cslin@lin-associates.com

Date prepared: December 27, 2017

5.2 Identification of Product

Device Trade Name: AmCAD-UT[®] Detection 2.2
Common and Usual Name: Computer-Assisted Detection (CAdE) Device
Device Classification Name: Picture Archiving and Communications System
Regulation Number: 21 CFR 892.2050
Classification Product Code: LLZ
Classification Panel: Radiology

Classification: Class II
Manufacturer: AmCad BioMed Corporation

5.3 Predicate Device

This subject software medical device is substantially equivalent to the device listed below:

Model: AmCAD-UT[®] Detection, Version 2.0
Manufacturer: AmCad BioMed Corporation
510(k) Number: K122536

5.4 Device Description

AmCAD-UT[®] Detection 2.2 is a Windows-based computer-assisted detection (CADe) software application device designed to assist medical professionals in analyzing thyroid ultrasound images of user selected regions of interest (ROI).

After the initial review of thyroid ultrasound images by the physician, he/she can use AmCAD-UT[®] Detection to analyze thyroid images for further interpretation. The physician selects an ROI to define the initial boundary of the ROI. Once the ROI is confirmed, the physician may process the image for detection and quantification of sonographic characteristics (i.e., hyperechoic foci, echogenicity, texture, margin, orientation and anechoic areas) by AmCAD-UT[®] Detection. The device provides more detailed information with quantification and visualization of the sonographic characteristics of thyroid nodule that may assist physician in his/her complete interpretation.

The software application automatically generates reports given the user preference inputs (e.g., the nodule size, nodule location and shape, and the presence or absence of the sonographic characteristics) annotated during the image analysis process. A report form has been designed by AmCad to be consistent with the conventional clinical thyroid report form. An output of the report may be viewed and sent to paper printers or saved on the standalone PC or review station as PDF file.

5.5 Indications for Use

AmCAD-UT® Detection is a Windows-based computer-aided detection (CADe) device intended to assist the medical professionals in analyzing thyroid ultrasound images, acquired from FDA-cleared ultrasound systems, with user-selected regions of interest (ROI). After the initial review of the ultrasound images by the physicians, the device further provides detailed information with quantification and visualization of sonographic characteristics of thyroid nodules. The device is intended for use on ultrasound images of discrete thyroid nodules larger than 1cm, for which a biopsy recommendation is required.

5.6 Comparison with Predicate Devices

AmCAD-UT® Detection 2.2 is a computer-assisted detection (CADe) device which provides viewing and post-acquisition image processing and analysis of thyroid ultrasound images with user-selected regions of interest (ROI) and automatically generates reports from user inputs annotated during the image analysis process. This software medical device is substantially equivalent to the predicate device listed below:

Model: AmCAD-UT® Detection, Version 2.0
 Manufacturer: AmCad BioMed Corporation
 510(k) Number: K122536

The comparison as described in the following table:

	AmCAD-UT® Detection 2.2	AmCAD-UT® Detection 2.0
Manufacturer	AmCad BioMed Corp.	AmCad BioMed Corp.
510(k) Number	K180006	K122536
Device Common Name	Computer-Assisted Detection (CADe)	Computer-Assisted Detection (CADe)
Regulation Number	21 CFR 892.2050 - Class II	21 CFR 892.2050 - Class II

	AmCAD-UT® Detection 2.2	AmCAD-UT® Detection 2.0
Regulation Name	Picture archiving and communications system	Picture archiving and communications system
Product Code	LLZ	LLZ
Intended Use	AmCAD-UT® Detection 2.2 is intended to assist the medical professionals in analyzing thyroid ultrasound images of user-selected regions of interest (ROI). After the initial review of the ultrasound images by the physicians, the device further provides detailed information with quantification and visualization of sonographic characteristics of thyroid nodules.	AmCAD-UT® Detection 2.0 is intended to assist the medical professionals in analyzing thyroid ultrasound images of user-selected regions of interest (ROI). After the initial review of the ultrasound images by the physicians, the device further provides detailed information with quantification and visualization of sonographic characteristics of thyroid nodules.
Indications for Use	AmCAD-UT® Detection 2.2 is a Windows-based computer-aided detection (CADe) device intended to assist the medical professionals in analyzing thyroid ultrasound images, acquired from FDA-cleared ultrasound systems, with user-selected regions of interest (ROI). After the initial review of the ultrasound images by the physicians, the device further provides detailed information with quantification and visualization of sonographic characteristics of thyroid nodules. The device is intended for use on ultrasound	AmCAD-UT® Detection 2.0 is a Windows-based computer-aided detection (CADe) device intended to assist the medical professionals in analyzing thyroid ultrasound images of user-selected regions of interest (ROI). After the initial review of the ultrasound images by the physicians, the device further provides detailed information with quantification and visualization of sonographic characteristics of thyroid nodules. The device is intended for use on Philips HDI5000 images of discrete thyroid nodules larger than 1cm, for which a biopsy has

	AmCAD-UT® Detection 2.2	AmCAD-UT® Detection 2.0
	images of discrete thyroid nodules larger than 1cm, for which a biopsy recommendation is required.	been recommended. The device performance has been validated on images collected from Philips HDI5000 with a 5-12MHz multi-frequency probe.
Functional Capability of Image Processing	AmCAD-UT® Detection 2.2 analyzes the user-selected regions of interest (ROI) of thyroid ultrasound image for the detection and quantification of sonographic characteristics (hyperechoic foci, echogenicity, texture, margin, orientation and anechoic areas). The device further provides detailed information with visualization of sonographic characteristics of thyroid nodules.	AmCAD-UT® Detection 2.0 analyzes the user-selected regions of interest (ROI) of thyroid ultrasound image for the detection and quantification of sonographic characteristics (hyperechoic foci, echo-pattern, echo-texture and anechoic areas). The device further provides detailed information with visualization of sonographic characteristics of thyroid nodules.
Reading Paradigm	AmCAD-UT® Detection 2.2 is to provide quantification and visualization of sonographic characteristics after physicians' initial review of the images.	For use as "second reader" meaning that the function of AmCAD-UT® Detection 2.0 is to provide quantification and visualization of sonographic characteristics after physicians' initial review of the images.
Output Generated by the CAD Device	The image can be annotated with the detected sonographic characteristics and be recorded by the device. The software also automatically generates reports given the user preference inputs in the analysis process.	The image can be annotated with the detected sonographic characteristics and be recorded by the device. The software also automatically generates reports given the user preference inputs in the analysis process.
Type of Film to be Processed by the CAD Device	Digital ultrasound image	Digital ultrasound image

	AmCAD-UT® Detection 2.2	AmCAD-UT® Detection 2.0
Software Design	Based on Statistical Pattern Recognition and Quantification method	Based on Statistical Pattern Recognition and Quantification method
Ground Truth Establishment	The ground truth to be established for performance studies of the device includes the ROI, the presence of each sonographic characteristic, and the surgical pathology examination result.	The ground truth to be established for performance studies of the device includes the ROI, the presence of each sonographic characteristic, and the surgical pathology examination result.
Platform	Window-based	Window-based
Operating System	Standard PC or review station	Standard PC or review station
Clinical Application	Thyroid cancers	Thyroid cancers
Image Type	Ultrasound Image	Ultrasound Image
Image Format	DICOM3.0, Bitmap, JPEG	DICOM3.0, Bitmap, JPEG
ROI Quantification	Yes	Yes
Automatically Generating Report	Yes	Yes
Report Storage	Paper printers, Local disk	Paper printers, Local disk
Performance Testing Data to Support SE Determination	Results from standalone performance testing and clinical performance testing (MRMC study)	Standalone performance assessment and clinical MRMC study.

AmCAD-UT® Detection 2.2 is substantially equivalent to AmCAD-UT® Detection 2.0 that provides display and post-acquisition image analysis of ultrasound images assisting the physician in analyzing the ultrasound images of thyroid nodules. The standalone and clinical reader performance assessment results of AmCAD-UT® Detection 2.2 are substantially equivalent to the standalone performance assessment and clinical MRMC reader study of AmCAD-UT® Detection 2.0. Minor technological characteristics differences do not raise any new questions of safety and effectiveness.

Thus, AmCAD-UT® Detection 2.2 is substantially equivalent to the predicate device as a Computer-Assisted Detection (CADe) device intended to assist the physicians in clinical practice.

5.7 Performance Standards

No applicable FDA performance standards have been issued under the authority of Section 514.

5.8 Software

Software development for the AmCAD-UT® Detection 2.2 follows documented processes for software design, verification and validation testing. A risk assessment has been completed to identify potential design hazards that could cause an error or injury based on the use of the quantification results. Appropriate steps have been taken to control all identified risks for this type of image viewing and quantification device.

5.9 Summary of Performance Data to Support Substantial Equivalence

AmCad BioMed Corporation has conducted standalone and clinical reader performance studies to validate and assess the performance of the AmCAD-UT® Detection 2.2 for its intended use. The standalone studies evaluated the performance of the quantified sonographic characteristics (hyperechoic foci, echogenicity, texture, margin, orientation and anechoic areas) on images acquired from ultrasound systems of different brands and showed that the device is effective in detecting the sonographic characteristics of thyroid nodules acquired from various ultrasound systems.

The intended use of the AmCAD-UT® Detection 2.2 was validated in a clinical Multiple-Reader-Multiple-Case (MRMC) study. The results of the MRMC study demonstrated that the physician reading thyroid nodule sonography images with the assistance of AmCAD-UT® Detection 2.2 can enhance their ability in analyzing the sonographic characteristics and has led to a significant increase in effectiveness of making clinical judgment.

5.10 Conclusions

AmCAD-UT[®] Detection 2.2, being a computer assisted detection (CADe) software device, has the same intended use as the predicate device. In addition to generalizing the use of the device for ultrasound images acquired from various FDA-cleared ultrasound systems, two quantified and visualized sonographic characteristics, i.e. margin distinctness and tumor orientation, are added. On the user interface of the proposed device, gauge meters expressing the quantified values of sonographic characteristics are also added for convenience of the medical professionals' reading. The performance data demonstrates that the proposed device performs as safely and effectively as the predicate devices. Therefore, AmCAD-UT[®] Detection 2.2 is substantially equivalent to the predicate devices as the CADe device intended to assist the medical professionals in analyzing ultrasound images of thyroid nodules in clinical practice.