



September 28, 2018

Acessa Health, Inc.  
Kim Rodriguez  
President & CEO  
7004 Bee Cave Road, Bldg. 3, Suite 200  
Austin, TX 78746

Re: K181124  
Trade/Device Name: Acessa ProVu System  
Regulation Number: 21 CFR§ 884.4160  
Regulation Name: Unipolar endoscopic coagulator-cutter and accessories  
Regulatory Class: II  
Product Code: HFG, OEW, IYO, ITX  
Dated: August 31, 2018  
Received: August 31, 2018

Dear Kim Rodriguez:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

  
Sharon M. Andrews -S

for  
Benjamin R. Fisher, Ph.D.  
Director  
Division of Reproductive, Gastro-Renal,  
and Urological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K181124

Device Name

Acessa ProVu System

Indications for Use (Describe)

The Acessa ProVu System is indicated for use in percutaneous, laparoscopic coagulation and ablation of soft tissue, including treatment of symptomatic uterine fibroids under laparoscopic ultrasound guidance.

The Acessa ProVu System includes optional electromagnetic guidance for enhancing the ultrasonic image of the Acessa ProVu Handpiece and for predicting its future path on a computer monitor screen which also shows the ultrasound B-scan image.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(K) SUMMARY

In accordance with 21 CFR 807.87(h) and (21 CFR 807.92) the 510(k) Summary for the Acesa ProVu System is provided below.

	Primary	Secondary
Submitter's Name:	Acesa Health Inc.	Acesa Health Inc.
Submitter's Address:	7004 Bee Cave Rd, Bldg. 3, Suite 200, Austin, TX 78746	7004 Bee Cave Rd, Bldg. 3, Suite 200, Austin, TX 78746
Submitter's Telephone:	(877) 412-3828	(877) 412-3828
Submitter's Fax:	(925) 605-0327	(925) 605-0327
Contact Name:	Kim Rodriguez	Brian J. Bergeron
Date Summary was Prepared:	September 27, 2018	
Trade or Proprietary Name:	Acesa ProVu System	
Common or Usual Name:	Radiofrequency Ablation System	
Primary Product Code:	HFG (coagulator, laparoscopic, unipolar (and accessories))	
Regulation Name and Number:	Class II per 21 CFR 884.4160: Unipolar Endoscopic Coagulator-Cutter and Accessories	
Secondary Product Codes:	OEW, IYO, ITX	
Classification Panel:	Obstetrics and Gynecology	
Legally Marketed (Unmodified) Devices to Which Substantial Equivalence is Claimed:	<p>Primary Predicate: K132744 Acesa Guidance System</p> <p>Secondary Predicates: K121858 The Acesa System K163121 Telemed SmartUS Ultrasound System</p> <p>The predicate devices have not been subject to a design related recall.</p>	

**DESCRIPTION OF THE DEVICE SUBJECT TO PREMARKET NOTIFICATION:**

The Acesa ProVu System provides radiofrequency (RF) ablation, ultrasound visualization and guidance within a single console and includes additional accessories. The Acesa ProVu System consists of a console with push buttons to control menu functionality. The console contains the following hardware and electronic components:

1. RF Ablation System
2. Ultrasound Visualization System
3. Guidance System

The following accessories connect to the console:

- Dual Foot pedal (one for RF ablation and one for coagulation)
- Video Cable
- Power Cord
- Acesa Pads and Acesa Pad Cable
- Acesa ProVu Handpiece and Acesa ProVu Handpiece Cable
- Acesa Transducer, either:
  - Acesa ProVu Transducer used with Acesa ProVu Transducer Sleeve (embedded with sensor), or
  - Acesa ProVu Transducer with (embedded) sensor
- Acesa Table Top Field Generator or Acesa ProVu Planar Field Generator
- Monitor (hospital-owned accessory monitor, not provided by Acesa)

The Acesa ProVu System is designed to deliver up to 200 W of RF power at 460 kHz in three operational modes: Temperature Control, Power Control and Coagulation Mode. The graphical user interface (GUI) is displayed on a hospital-owned monitor. The system enables selection of operational parameters such as the mode of operation, the ablation time, the target temperature, and the power delivery level. With the Acesa ProVu Handpiece placed in the tissue to be ablated and its electrodes deployed, RF power can be turned on. The system parameters are continuously monitored and displayed on the monitor. If the measured parameters are outside the acceptable limits, the RF energy delivery automatically stops and a message appears on the graphical user interface. RF energy during an ablation or coagulation can also be stopped at any time by the user by pressing the foot pedal.

The Acesa ProVu System uses electromagnetic tracking technology to track the positions of the Acesa ProVu Transducer used with Acesa ProVu Transducer Sleeve (or the Acesa ProVu Transducer with embedded sensor) and the Acesa ProVu Handpiece shaft and draws virtual representations of them in their spatial relationship, so that a physician can predict the Acesa ProVu Handpiece shaft's future path in relation to the features in the ultrasound slice. The use of the guidance functionality is considered an optional accessory to procedures where ultrasound is currently used for visualization, such as the Acesa System procedure.

### INTENDED USE/INDICATIONS FOR USE

The Acesa ProVu System is indicated for use in percutaneous, laparoscopic coagulation and ablation of soft tissue, including treatment of symptomatic uterine fibroids under laparoscopic ultrasound guidance.

The Acesa ProVu System includes optional electromagnetic guidance for enhancing the ultrasonic image of the Acesa ProVu Handpiece and for predicting its future path on a computer monitor screen which also shows the ultrasound B-scan image.

The subject and predicate device have the same intended use.

### COMPARISON OF TECHNOLOGICAL CHARACTERISTICS

	<b>Acesa ProVu System</b> (K181124)	<b>Acesa Guidance System</b> (K132744)	<b>The Acesa System</b> (K121858)	<b>SmartUS EXT-1M/3M</b> (K163121)
<b>General Information</b>				
<b>Classification Regulation</b>	21 CFR §884.4160	21 CFR §884.4160	21 CFR §884.4160	21 CFR §892.1550
<b>Classification Name</b>	Unipolar endoscopic coagulator cutter and accessories	Unipolar endoscopic coagulator cutter and accessories	Unipolar endoscopic coagulator cutter and accessories	Ultrasonic pulsed doppler imaging system
<b>Product Code</b>	HFG - coagulator, laparoscopic, unipolar (and accessories)	HFG - coagulator, laparoscopic, unipolar (and accessories)	HFG - coagulator, laparoscopic, unipolar (and accessories)	IYN
<b>Secondary Product Codes</b>	OEW, IYO, ITX	OEW, IYO	N/A	ITX, IYO
<b>Product Class</b>	II	II	II	II
<b>Intended Users</b>	Health Care Professional	Health Care Professionals	Health Care Professionals	Health Care Professionals
<b>Classification Panel</b>	Obstetrics and Gynecology	Obstetrics and Gynecology	Obstetrics and Gynecology	Radiology
<b>System Accessories</b>	Field Generator Ultrasound Transducer RF Handpiece Dual foot pedal Single use pads Handpiece Cable Pad Cable Power Cord	Field Generator Ultrasound Transducer RF Handpiece Dual foot pedal Single use pads Handpiece Cable Pad Cable Power Cord	Field Generator Ultrasound Transducer RF Handpiece Dual foot pedal Single use pads Handpiece Cable Pad Cable Power Cord	
<b>RF Functionality</b>				
<b>Operating Principle</b>	Monopolar	Monopolar	Monopolar	
<b>Max Power Delivery</b>	200 W into 50-80 Ω	200 W into 50-80 Ω	200 W into 50-80 Ω	
<b>Frequency</b>	460 kHz	460 kHz	460 kHz	
<b>Wave Form</b>	Sinusoidal	Sinusoidal	Sinusoidal	
<b>Temperature Display Range</b>	15 - 125°C	15 - 125°C	15 - 125°C	
<b>Impedance Range</b>	0-511 Ohms	0-511 Ohms	0-511 Ohms	
<b>Duration Range</b>	0.1 to 12 minutes	0.1 to 12 minutes	0.1 to 12 minutes	
<b>Temperature Display</b>	Handpiece and Pad thermocouples	Handpiece and Pad thermocouples	Handpiece and Pad thermocouples	

	<b>Acessa ProVu System</b> (K181124)	<b>Acessa Guidance System</b> (K132744)	<b>The Acessa System</b> (K121858)	<b>SmartUS EXT-1M/3M</b> (K163121)
<b>Set Power Display</b>	Yes, coagulation only	Yes, coagulation only	Yes, coagulation only	
<b>Duration of Energy Delivery Displayed</b>	Yes	Yes	Yes	
<b>Control Modes</b>	Temperature, Manual, Coagulation	Temperature, Manual, Coagulation	Temperature, Manual, Coagulation	
<b>Power Level Controls</b>	Yes, coagulation only	Yes, coagulation only	Yes, coagulation only	
<b>Ablation Duration Controls</b>	Yes	Yes	Yes	
<b>Maximum Ablation Duration</b>	12 minutes at target time	12 minutes at target time	12 minutes at target time	
<b>Power Supply (PS)</b>	100-240 V, 50/60 Hz switchable PS	100-240 V, 50/60 Hz switchable PS	100-240 V, 50/60 Hz switchable PS	
<b>Peak Voltage</b>	230 V	230 V	230 V	
<b>Manual Mode Power Limit</b>	15 W	10-15 W	10-15 W	
<b>Number of Serial Communication Ports</b>	Single USB port	Single USB Port	Single USB Port	
<b>Pump used with Generator</b>	None	None	None	
<b>Display</b>	Hospital-provided Monitor	LCD/Touch Screen	LCD/Touch Screen	
<b>User Control Interface</b>	System settings controlled at the Console.	System settings controlled with the handpiece buttons and at the RF Generator LCD/Touch Screen.	System settings controlled with the handpiece buttons and at the RF Generator LCD/Touch Screen.	
<b>Guidance Functionality</b>				
<b>Tracking Technology</b>	Electromagnetic	Electromagnetic		
<b>Accuracy</b>	±10 mm	±10 mm		
<b>Views</b>	Provide a geometrically accurate “3D” representation of the Transducer and Handpiece. Display a 2D inset view of the ultrasound video, with the added feature of the target zone overlaid on it and the laparoscopic ultrasound transducer	Provide a geometrically accurate “3D” representation of the Transducer and Handpiece. Display a 2D inset view of the ultrasound video and the laparoscopic ultrasound transducer		
<b>Tracking System</b>	NDI V3	NDI V2		
<b>Ultrasound Functionality</b>				
<b>Ultrasound System</b>	Ultrasound beam formers and drivers	Ultrasound beam formers and drivers (Off the shelf Aloka 6 Ultrasound system K093488)	Ultrasound beam formers and drivers (Off the shelf Aloka 6 Ultrasound system K093488)	Ultrasound beam formers and drivers
<b>Ultrasound Transducer</b>	Acessa ProVu Transducer (Model Number: 7500 or 7700)	Aloka UST 5526L-75 Transducer (K040719)	Aloka UST 5526L-75 Transducer (K040719)	Aloka UST 5526L-75 Transducer (K040719)



The subject device is a combination system of the previously cleared devices that includes an RF ablation system (K121858), guidance control (K132744) and an OEM ultrasound component (K163121). As evidenced by the table above, there are differences in technological characteristics between the subject and predicate device. These differences do not raise different questions of safety and effectiveness.

#### SUMMARY OF PERFORMANCE DATA

To further support a determination of substantial equivalence, non-clinical bench testing was conducted to support the subject device. The specific types of non-clinical testing conducted are listed below.

- Biocompatibility Testing
- Use Life Testing of Transducers
- Electrical Safety and Electromagnetic Compatibility testing per:
  - ANSI/AAMI/ES 60601-1:2005 + A2:2010 + A1:2012
  - IEC 60601-2-2: 2017
  - IEC 60601-1-2: 2014
  - IEC 60601-2-37: 2015
- Software Validation and Verification Testing per FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices." The software for this device was considered as a "major" level of concern, since a failure or latent flaw in the software could directly result in serious injury or death to the patient or operator.
- Acesa ProVu Console and Guidance System Level Testing
- Benchtop Ablation Testing to validate that the ablations created by the subject device are substantially equivalent to the predicate.
- Acesa ProVu Handpiece Performance Testing to qualify and quantify the mechanical and functional properties of the Handpiece.
- Acesa Pads Performance Testing to qualify and quantify the mechanical and functional properties of the Pads.
- Acesa ProVu Transducer Performance Testing, including Acoustical and Thermal Measurement, in conformance with IEC 60601-2-37
- Acesa ProVu Transducer Sleeve Testing to qualify and quantify the mechanical and functional properties of the Transducer Sleeve.

#### CONCLUSION

The performance data demonstrate that the subject device is substantially equivalent to the predicate devices.