



August 7, 2018

Dong Guan Ou Kang Electronics CO., LTD  
% Field Fu  
Consultant  
Shenzhen Joyantech Consulting Co., Ltd.  
Room 1122, International Mayors Communication Centre  
Shenzhen, 511470 Cn

Re: K181234

Trade/Device Name: Self-adhesive Electrode  
Regulation Number: 21 CFR 882.1320  
Regulation Name: Cutaneous Electrode  
Regulatory Class: Class II  
Product Code: GXY  
Dated: May 9, 2018  
Received: May 9, 2018

Dear Field Fu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

**Michael J. Hoffmann -S**

for Carlos L. Peña, PhD, MS  
Director  
Division of Neurological  
and Physical Medicine Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

k181234

Device Name

Self-adhesive electrode

Indications for Use (Describe)

Self-adhesive electrode is intended for single use, conductive adhesive interface between the patient's skin and the electrical stimulation devices. Some common types of the stimulation devices include, but are not limit to, TENS (Transcutaneous Electrical Nerve Stimulation) and EMS (Electrical Muscular Stimulation). The electrode is for OTC (Over-The -Counter) or Prescription use, for single patient use only.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

Date of summary prepared: 08/06/2018

This summary of 510(K) safety and effectiveness information is submitted as required by requirements of SMDA and 21 CFR §807.92.

### 1. Administrative Information

<b>Date of summary prepared</b> <b>Manufacturer information</b>	Company title: Dong Guan Ou Kang Electronics CO.,LTD Company address: 3 <sup>rd</sup> floor, Xiqunli building, Shipai road, Yanwo village, Shipai town, Dongguan Phone: +86-0769-82686009 Fax: +86-0769-82686009 Contact person: Nianjun Liu E-mail: nianjun_liu@sohu.com
<b>Submission Correspondent</b> 	Name: Shenzhen Joyantech Consulting Co., Ltd Address: Room 1122, International Mayors Communication Centre, NO. 55 Shizhou middle road , Nanshan District, Shenzhen Contact person: Field Fu; Summer Wu; Elly Xv E-mail: <a href="mailto:summer@cefda.com">summer@cefda.com</a> ; <a href="mailto:elly@cefda.com">elly@cefda.com</a>

### 2. Device information

<b>Type of 510(k) submission:</b> <b>Trade name:</b> <b>Model:</b> <b>Classification Name:</b> <b>Review panel:</b> <b>Product Code:</b> <b>Device Class</b> <b>Regulation number:</b>	Traditional Self-adhesive electrode Lead wire type electrode and snap type electrode Electrode, Cutaneous Neurology GXY II 882.1320
---	--

### 3. Predicate Device Information

<b>Sponsor:</b> <b>Device:</b> <b>510(k) number:</b>	Shenzhen Mailuokang Technology Co., Ltd. Mailuokang Self-adhesive Electrode K152815
--	---

**4. Device description**

**1) Intended use/indications for use**

Self-adhesive electrode is intended for single use, conductive adhesive interface between the patient’s skin and the electrical stimulation devices. Some common types of the stimulation devices include, but are not limit to, TENS (Transcutaneous Electrical Nerve Stimulation) and EMS (Electrical Muscular Stimulation). The electrode is for OTC (Over-The -Counter) or Prescription use, for single patient use only.

**2) Device specification**

The specification of self-adhesive electrode is listed below.

Model	Shape	Product Size(mm)	Connector Size (hole diameter) (mm)
Non-woven lead wire type electrode	Rectangular and circular	Rectangular: 40×40;40×60;50×50;50×60;50×70;50×90;50×100;60×80;60×90;70×110;90×130 Circular: Φ30; Φ40; Φ50; Φ80	2.0, 2.5
Non-woven snap type electrode	Rectangular , circular, oval and hand shape	Rectangular: 40×40; 40×80; 50×50; 50×90; 50×100; 60×90 Circular: Φ30; Φ40; Φ50; Φ80 Oval: 43×30 Hand shape: 75×46	3.3, 3.6

**3) Device design**

Non-woven lead wire type electrode and non-woven snap type electrode have different appearance and construction. See Figure 1 and Figure 2 for the structure picture of non-woven lead wire type electrode and non-woven snap type electrode.




Figure 1 The structure of non-woven lead wire type electrode




Figure 2 The structure of non-woven snap type electrode

**4) Materials used**

Both non-woven lead wire type electrode and snap type electrode are composed of a non-woven fabric, a double side adhesive tape, conducting film, hydrogel and plastic film. The construction difference between non-woven lead wire electrode and non-woven snap electrode is lead wire and snap.

**5) Physical and performance characteristics of the device**

Electrical impedance: less than 3KΩ.

Current uniformity: the current deviation index of the electrode is less than 5%.

Adhesive performance: 30 minutes maximum duration use.

**6) Principle of Operation**

Self-adhesive electrode functions as a passive device by carrying an electrical signal from a stimulation device through the device cable and electrode lead wire or snap to the user skin. Electrical signal from an electrical stimulator is conducted to the self-adhesive electrode through a lead wire or snap; which is dispersed across the conductive film, then transmitted through the conductive adhesive hydrogel to the surface of the patient’s skin.

**7) Material of construction**

Component	Material
Non-woven fabric	100% polypropylene in white
Conducting film	Carbon paste
Hydrogel	Cross linked acrylic resin, polyhydric alcohol, electrolytic salt, additives and de-ionized water
Lead wire	PVC
Snap	Nickel-clad copper

**5. Technological characteristics of the subject device compared to the predicate device**

Items	Predicate Device (K152815)	Subject Device	Comparison
<b>General Characteristics</b>			

Items	Predicate Device (K152815)	Subject Device	Comparison
<b>Classification Name</b>	Cutaneous Electrode	Cutaneous Electrode	Same.
<b>Product Code</b>	GXY	GXY	Same.
<b>Regulation Number</b>	882.1320	882.1320	Same.
<b>Device Class</b>	Class II	Class II	Same.
<b>Intended Use/Indications for Use</b>	Mailuokang Self-adhesive Electrode is intended for use as reusable, conductive adhesive interface between the patient's skin and the electrical stimulation devices. Example electrical stimulation devices for current applications of the electrodes includes, but are not limited to, TENS (Transcutaneous Electrical Nerve Stimulation) and EMS (Electrical Muscle Stimulation). Mailuokang Self-adhesive Electrode is for prescription use and over-the-counter use, for single patient use only.	Self-adhesive electrode is intended for use as single use, conductive adhesive interface between the patient's skin and the electrical stimulation devices. Some common types of the stimulation devices include, but are not limit to, TENS (Transcutaneous Electrical Nerve Stimulation) and EMS (Electrical Muscular Stimulation). The electrode is for OTC (Over-The-Counter) or Prescription use, for single patient use only.	Similar. (Note 1)
<b>Intended User/Patient Population</b>	No specified population.	No specified population.	Same.
<b>Electrical Connection</b>	Lead wire or snap.	Lead wire or snap.	Same.
<b>OTC or prescription use</b>	OTC use and prescription use.	OTC use and prescription use.	Same.

Items	Predicate Device (K152815)	Subject Device	Comparison
<b>Reusable</b>	Yes	No.	Similar. (Note 1)
<b>Single Patient Use</b>	Yes	Yes	Same.
<b>Multiple Applications</b>	Yes	No.	Similar. (Note 1)
<b>Sterility Status</b>	Non-sterile	Non-sterile	Same.
<b>Technical Characteristics</b>			
<b>Electrode dimensions(cm)</b>	For lead wire style electrode: Rectangular: 3x7;4x4;5x5;5x6;5x7; 6x9;8x12;8x13 For snap style electrode: Rectangular: 9x4 Circle: Φ4; Φ5; Φ8.5 Oval: 4.5x3	For non-woven lead-wire type electrode: Rectangular: 4x4;4x6;5x5;5x6;5x7; 5x9;5x10;6x8;6x9; 7x11;9x13 Circular: Φ3; Φ4; Φ5; Φ8 For non-woven snap type electrode: Rectangular: 4x4; 4x8; 5x5; 5x9; 5x10; 6x9 Circular: Φ3; Φ4; Φ5; Φ8 Oval: 4.3x3 Hand shape: 7.5x4.6	Similar. (Note 2)
<b>Conductive surface area(cm<sup>2</sup>)</b>	12.56~104	7~117	Similar. (Note 2)
<b>Design Features</b>	Six basic components for Lead Wire Style Electrode: <ul style="list-style-type: none"> <li>• non-woven fabrics;</li> <li>• double sides adhesive tape;</li> <li>• conducting film;</li> <li>• hydrogel;</li> <li>• plastic film;</li> </ul>	Six basic components for Lead Wire Style Electrode: <ul style="list-style-type: none"> <li>• non-woven fabrics;</li> <li>• double sides adhesive tape;</li> <li>• conducting film;</li> <li>• hydrogel;</li> <li>• plastic film;</li> </ul>	Same.



Items	Predicate Device (K152815)	Subject Device	Comparison																		
	<ul style="list-style-type: none"> <li>carbon fiber wire.</li> </ul> <p>Six basic components for Snap Style Electrode:</p> <ul style="list-style-type: none"> <li>non-woven fabrics;</li> <li>double sided adhesive tape;</li> <li>conducting film;</li> <li>hydrogel;</li> <li>plastic film;</li> <li>snap.</li> </ul>	<ul style="list-style-type: none"> <li>carbon fiber wire.</li> </ul> <p>Six basic components for Snap Style Electrode:</p> <ul style="list-style-type: none"> <li>non-woven fabrics;</li> <li>double sided adhesive tape;</li> <li>conducting film;</li> <li>hydrogel;</li> <li>plastic film;</li> <li>snap.</li> </ul>																			
<b>Hydrogel Composition</b>	Glyrol, Polyacrylic acid, Water and Salt	Crosslinked acrylic resin, polyhydric alcohol, electrolytic salt, additives and de-ionized water.	Similar. (Note 43)																		
<b>Function and concentration of each ingredient in final formulation</b>	Not publicly available.	<table border="1"> <thead> <tr> <th data-bbox="785 1037 879 1238">Ingredient</th> <th data-bbox="879 1037 973 1238">Concentration (WT%)</th> <th data-bbox="973 1037 1075 1238">Function</th> </tr> </thead> <tbody> <tr> <td data-bbox="785 1238 879 1440">Crosslinked acrylic resin</td> <td data-bbox="879 1238 973 1440">15-25</td> <td data-bbox="973 1238 1075 1440">Cross-linking polymerization.</td> </tr> <tr> <td data-bbox="785 1440 879 1619">polyhydric alcohol</td> <td data-bbox="879 1440 973 1619">40-60</td> <td data-bbox="973 1440 1075 1619">Buffer and humectant.</td> </tr> <tr> <td data-bbox="785 1619 879 1821">electrolytic salt</td> <td data-bbox="879 1619 973 1821">2-10</td> <td data-bbox="973 1619 1075 1821">Ionic conducting media.</td> </tr> <tr> <td data-bbox="785 1821 879 1910">additives</td> <td data-bbox="879 1821 973 1910">&lt;0.5</td> <td data-bbox="973 1821 1075 1910">Stabilization.</td> </tr> <tr> <td data-bbox="785 1910 879 2020">de-ionized water</td> <td data-bbox="879 1910 973 2020">20-40</td> <td data-bbox="973 1910 1075 2020">Solution.</td> </tr> </tbody> </table>	Ingredient	Concentration (WT%)	Function	Crosslinked acrylic resin	15-25	Cross-linking polymerization.	polyhydric alcohol	40-60	Buffer and humectant.	electrolytic salt	2-10	Ionic conducting media.	additives	<0.5	Stabilization.	de-ionized water	20-40	Solution.	Note 3
Ingredient	Concentration (WT%)	Function																			
Crosslinked acrylic resin	15-25	Cross-linking polymerization.																			
polyhydric alcohol	40-60	Buffer and humectant.																			
electrolytic salt	2-10	Ionic conducting media.																			
additives	<0.5	Stabilization.																			
de-ionized water	20-40	Solution.																			

Items	Predicate Device (K152815)	Subject Device	Comparison
pH of hydrogel	Not publicly available.	6.5~7.5	Note 4
Conductivity of hydrogel (S/m)	Not publicly available.	50	Note 5
Resistivity ( $\Omega$ m)	Not publicly available.	0.02	Note 5
Volume of hydrogel ( $\text{cm}^3$ )	Not publicly available.	The volume of hydrogel varies by the shape of electrodes, ranging from 0.7 to 11.7.	Note 6
Maximum duration of use	30min.	30min.	Same.
Biocompatibility	Mailuokang's hydrogel has passed three biocompatibility tests: skin irritation, sensitization and cytotoxicity.	Self adhesive electrode has passed 3 biocompatibility tests: skin irritation, sensitization and cytotoxicity.	Same.
Performance Test	Measuring the electrode contact impedance and adhesion.	Size, electrical impedance and adhesion.	Same.
Conductive Surface Shapes	Various shapes (rectangular, oval, circle)	Various shapes (rectangular, oval, circular, hand shape)	Similar. (Note 2)

**Discussion of the similarities and differences:****Note 1:**

The difference between the subject device and the predicate device is that, the predicate device is intended for multiple uses, while the subject device is intended for single use. The difference does not adversely affect the safety and effectiveness of the device.

**Note 2:**

The electrode dimension and conductive surface of the predicate device is similar to the subject device. The difference does not adversely affect the safety and effectiveness of the device.

**Note 3:**

Both Mailuokang self-adhesive electrode and self-adhesive electrode have passed the biocompatibility test, electrode impedance test and adhesion test.

Therefore, the difference in hydrogel composition will not raise any safety or efficacy issues.

**Note 4:**

The conductive gel is called hydrogel, it is pre-applied to the electrode. The hydrogel pH of the predicate device is unknown, as the subject device has pass the biocompatibility tests performed per ISO 10993 series standards, the difference does not adversely affect the safety and effectiveness of the device.

**Note 5:**

The conductivity and impedance of the hydrogel in the predicate device is unknown. As both the electrical impedance of the predicate device and the subject device are less than  $3k\Omega$ , the difference does not adversely affect the safety and effectiveness of the device.

**Note 6:**

The difference in volume of hydrogel does not adversely affect the safety and effectiveness of the device.

**6. Brief discussion of the nonclinical tests**

Self-adhesive electrode has passed performance tests for size, electrical impedance, current uniformity test and adhesion. Adhesion test is performed according to IEC 60601-2-2:2009 Clause 201.15.101.7.

Self-adhesive electrode has passed 3 biocompatibility tests: cytotoxicity, skin irritation and skin sensitization. The biocompatibility tests are performed according to ISO 10993-1, 5 and 10.

The shelf life of self-adhesive electrode is verified according to ASTM F 1980-07.

**7. Brief discussion of clinical tests**

Not applicable.

**8. Other information (such as required by FDA guidance)**

No other information.

**9. Conclusions**

The subject device is substantially equivalent to Mailuokang Self-adhesive Electrode whose 510(k) number is K152815.