



October 12, 2018

OssDsign AB  
% David Weissburg  
Principal  
Weissburg Associates  
808 Williamson St., Suite 402  
Madison, Wisconsin 53703

Re: K181539  
Trade/Device Name: OSSDSIGN Cranioplug  
Regulation Number: 21 CFR 882.5250  
Regulation Name: Burr Hole Cover  
Regulatory Class: Class II  
Product Code: GXR  
Dated: September 10, 2018  
Received: September 11, 2018

Dear Mr. Weissburg:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

  
Matthew C. Krueger -S

for Carlos L. Peña, PhD, MS  
Director  
Division of Neurological  
and Physical Medicine Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K181539

Device Name

OSSDSIGN Cranioplug

Indications for Use (Describe)

OSSDSIGN Cranioplug is an implant intended to cover and plug holes drilled into the skull during surgery and to reattach cranial bone removed during surgery. These osseous defects are surgically created and are not intrinsic to the stability of the bony structure. The ceramic component of Cranioplug resorbs and is replaced with bone during the healing process. Cranioplug is indicated for use in adults and adolescents age 12 and older.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) SUMMARY

### K181539

1. 510(k) Owner Name and Address:  
OssDsign AB  
Virdings Allé 2  
SE 754 50, Uppsala  
Sweden  
Telephone: +46 (0) 18-55 39 93  
Email: info@ossdesign.com  
Contact: Ulrik Birgersson
2. Contact Person:  
David Weissburg  
Weissburg Associates  
808 Williamson St., Suite 402  
Madison, Wisconsin, 53703 USA
3. Date prepared: October 10, 2018
4. Trade Name: OSSDSIGN® Cranioplug
5. Regulation Description: Burr Hole Cover
6. Classification Name: Cover, Burr Hole (21 CFR 882.5250, Product Code GXR)
7. Class: 2
8. Predicate: K140309, Cranioplug, OssDsign
9. Device Description: OSSDSIGN Cranioplug003 is an osteoconductive calcium phosphate ceramic plug reinforced with a titanium mesh plate which together provide the mechanical performance, safety and efficacy properties. The fully cured calcium phosphate ceramic fills the void in the burr hole. The osteoconductive ceramic component of Cranioplug resorbs and is replaced with bone during the healing process. Cranioplug is sized to match standard 11mm and 14mm burr hole perforators.
10. Intended Use / Indications For Use: OSSDSIGN Cranioplug is an implant intended to cover and plug holes drilled into the skull during surgery and to reattach cranial bone removed during surgery. These osseous defects are surgically created and are not intrinsic to the stability of the bony structure. The ceramic component of Cranioplug resorbs and is replaced with bone during the healing process. Cranioplug is indicated for use in adults and adolescents age 12 and older.

11. Comparison of Technological Characteristics with the Predicate Device:

	OSSDSIGN Cranioplug003 (subject device)	OSSDSIGN Cranioplug001 (Predicate, K140309)
Intended Use	OSSDSIGN Cranioplug is an implant intended to cover and plug holes drilled into the skull during surgery and to reattach cranial bone removed during surgery. These osseous defects are surgically created and are not intrinsic to the stability of the bony structure. The ceramic component of Cranioplug resorbs and is replaced with bone during the healing process. Cranioplug is indicated for use in adults and adolescents age 12 and older.	OSSDSIGN Cranioplug is intended to cover and plug holes drilled into the skull during surgery and to reattach cranial bone removed during surgery. It is cleared for use for non-weight bearing applications in adults and adolescents age 12 and older.
Materials	Ti grade 2, proprietary calcium phosphate formulation (same as predicate)	Ti grade 3, proprietary calcium phosphate formulation
Re-sterilizable	No	No
Titanium thickness	0.4mm	0.5mm
Sizes (diameter)	11mm and 14mm	14mm
Ceramic hardening	Cured in manufacturer's facility	Cured in manufacturer's facility
Sterility on delivery	Sterile	Sterile
Shelf life	24 months	18 months

12. Testing vs. predicate: Cranioplug003 and its predicate were tested for bench performance characteristics and biologic performance. Direct comparison of performance, safety and effectiveness of Cranioplug003 and its predicate demonstrates that Cranioplug003 is substantially equivalent to its predicate in all characteristics. The table below provides a summary of some of the tests completed.

Test	Test Method Summary	Results
Energy absorption, flap reattachment	Device and predicate installed in anatomical model. Tested in Universal Testing Machine to 2mm displacement to determine deformation energy.	Subject device and predicate sustained deformation energy of 0.06 J.
Flap Fixation Dynamic Load	Flap fixated with 3 Cranioplugs installed in anatomic model. 8.5 kg head weight tested to simulate 50 years of sleep with hourly repositioning.	No flap deformation. All Cranioplug models tested.
Cytotoxicity	ISO elution method, ISO 10993-5, extracted in IX MEM at 37°C for 24 hours	No evidence of causing cell lysis or toxicity.
<i>in vivo</i> implantation	52-week sheep implantation study, ISO 10993-6	<i>in vivo</i> studies show biocompatibility, adequate resorption rate and osteoconduction.

13. Conclusions: Nonclinical tests demonstrate that Cranioplug003 is as safe and effective as its legally marketed predicate device.