



August 23, 2019

ACARE Technology Co., Ltd.
% Ming-Yie Jan
RusCert Technology Co., Ltd.
8F.-2, No.21, Sec. 6, Zhongxiao E. Rd., Nangang Dist.
Taipei City, Taiwan

Re: K182950
Trade/Device Name: Acare Suction Unit
Regulation Number: 21 CFR 878.4780
Regulation Name: Powered Suction Pump
Regulatory Class: Class II
Product Code: JCX
Dated: July 26, 2019
Received: July 29, 2019

Dear Ming-Yie Jan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Long Chen, Ph.D.
Assistant Director
DHT4A: Division of General Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

SECTION 4 INDICATIONS FOR USE STATEMENT

DEPARTMENT OF HEALTH AND HUMAN SERVICES
Food and Drug Administration

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2020

Indications for Use

See PRA Statement below.

510(k) Number *(if known)*

K182950

Device Name

Acare Suction Unit, model ASU-200

Indications for Use *(Describe)*

The device is to be used to remove fluids from the airway or respiratory system. The device creates a negative pressure (vacuum) that draws fluids through disposable tubing that is connected to a collection canister. The fluids are trapped in the collection canister for proper disposal. It is for use on the order of a physician only.

Type of Use *(Select one or both, as applicable)*

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

ACARE TECHNOLOGY CO., LTD.**510 (K) Summary****5.1 Submitter's Information**

Name: ACARE TECHNOLOGY CO., LTD.
Address: 5f-4, No.24, Wuquan 2nd Rd., Xinzhuang Dist.,
New Taipei City 242, Taiwan

Establishment Registration Number: 3006395447
Contact: Kevin Tsai / QA manager

Tel: + 886-2-2298-8170 #204
Fax: +886-2-22988560
E-mail: Kevin@acarecorp.com

Or

Linda Cheng / General Manager

Date Prepared: September 30, 2018

2. Device Name

Trade Name: Acare Suction Unit, model ASU-200
Common/usual Name: Acare Suction Unit, model ASU-200
Device Classification Names: Apparatus, Suction, Ward use, Portable, Ac-Powered
Panel: General & Plastic Surgery,
Product Code: JCX
Classification: Class II
Regulation Description: Powered Suction Pump
Regulation Number: §878.4780

5.3 Predicate Devices:

- (a) VacPlus Suction Unit, 510(k) Number: K103535
- (b) EMG Suction Unit, 510(k) Number: K112421

ACARE TECHNOLOGY CO., LTD.

5.4 Indication for Use

The device is to be used to remove fluids from the airway or respiratory system. The device creates a negative pressure (vacuum) that draws fluids through disposable tubing that is connected to a collection canister. The fluids are trapped in the collection canister for proper disposal. It is for use on the order of a physician only.

5.5 Device Description

The Acare Suction Unit, model ASU-200 is to be used to remove fluids from the airway or respiratory system. The device creates a negative pressure (vacuum) that draws fluids through disposable tubing that is connected to a collection canister. The fluids are trapped in the collection canister for proper disposal. The device is a rugged, compact and portable electrically-operated medical suction device, designed for multi-source powering from AC (External 100V– 240 V) or DC (External 12V DC) mains supply or internal rechargeable battery, and is ideally suited for field, transport and hospital/clinical use applications.

The Acare canister is an accessory to the Acare Pump and the removed body fluids are collected in it. The Acare Pump is an AC powered oil-less positive displacement pump capable of delivering up to 600 mmHg of vacuum to draw fluid and small particles. The Acare Pump has a plastic casing fitted with a power (ON/OFF) button, power indicator, one suction adjustment knob, a suction gauge and a canister holder. The Acare canister is a preassembled 1000 mL canister including canister tubing fitted with a microbial filter at the distal end. The canister is fitted with a removable lid for fluid containment and an overflow protection device. The Acare canister should be placed in the in-built canister holder and connect the Acare canister with the Acare Pump via the connectors of the canister tubing prior to use.

ACARE TECHNOLOGY CO., LTD.**5.6 Comparison of Technological Characteristics**

A summary of the technological characteristics of the device subject of this premarket notification in comparison to those of the predicate devices is included in **Table 5-1**.

Table 5-1. Summary of the Proposed and Predicate Devices Technological Characteristics

Feature	Proposed device	Predicate device (a)	Predicate device (b)	Comparison
Device	Acare Suction Unit, model ASU-200 (K182950)	VacPlus Suction Unit (K103535)	EMG Suction Unit (K112421)	
Classification Name	Apparatus, Suction, Ward use, Portable, Ac- Powered	Apparatus, Suction, Ward use, Portable, Ac- Powered	Apparatus, Suction, Ward use, Portable, Ac- Powered	Same
Classification / Product Code	Class II / JCX	Class II / JCX	Class II / JCX	Same
Indication for Use	The device is to be used to remove fluids from the airway or respiratory system. The device creates a negative pressure (vacuum) that draws fluids through disposable tubing that is connected to a collection canister. The fluids are trapped in the collection canister for proper disposal. It is for use on the order of a physician only.	The device is to be used to remove fluids from the treated Tracheotomy patient airway.	The device is to be used to remove body fluids from a patient's airway or respiratory system. It is for use on the order of a physician only.	Similarity
Product Type	AC-Powered suction pump	AC-Powered suction pump	AC-Powered suction pump	Same
Electrical requirement	AC 100-240V, 50/60Hz DC 18V	AC 110-240V, 50/60Hz DC 12V	AC 100-240V, 47-63Hz, DC14V,	Similarity
Power Consumption	35 VA	40W	38W	Similarity
Battery	Li-ion Battery: 14.8Vdc, 2.5Ah	Lead Acid battery, 12Vdc, 2.8Ah	Li-ion Battery: 12Vdc, 2.5Ah	Similarity
Indication of Vacuum Level	Analogue Vacuum Gauge/ 760mmHg/100kPa	Analogue Vacuum Gauge/ 760mmHg/100kPa	Analogue Vacuum Gauge 76cmHg/30inHg	Same
Maximum Vacuum (Pressure)	600 mmHg	600 mmHg	600 mmHg	Same
Maximum Flow Rate	24 LPM	24 LPM	20~25LPM	Same
Sound Level (Noise level)	Pass ISO 10079-1	Pass ISO 10079-1	Pass ISO 10079-1	Same

ACARE TECHNOLOGY CO., LTD.**Table 5-2.** (Continued) Summary of the Proposed and Predicate Devices Technological Characteristics

Feature	Proposed device	Predicate device (a)	Predicate device (b)	Comparison
Device	Acare Suction Unit, model ASU-200 (K182950)	VacPlus Suction Unit (K103535)	EMG Suction Unit (K112421)	
Type of Pump	Oil-less Piston pump	Oil-less Piston pump	Oil-less Piston pump	Same
Pressure Control	adjusted by tuning the knob of vacuum adjusting set	adjusted by tuning the knob of vacuum adjusting set	adjusted by tuning the knob of vacuum adjusting set	Same
Disposable Canister	Yes	Yes	Yes	Same
Connection tube	Yes	Yes	Yes	Same
Patient Tube	Yes	Yes	Yes	Same
Overflow protection	Yes	Yes	Yes	Same
Bacterial filter	Yes	Yes	Yes	Same
Operating Environment	Temp: 50°F~104°F (10°C ~40°C) Humidity:10%-90% RH	Temp: 50°F~104°F (10°C ~40°C) Humidity:10%-90% RH	Temp: 32°F~104°F (0°~40 °C) Humidity:< 95%RH	Similarity
Operating Atmospheric	700~1060hPa	700~1013hPa	700~1060hPa	Similarity
Storage Environment	Temp: 59°F~122°F (15°C ~50°C) Humidity:10%~ 90% RH	Temp: -15°C to 50°C (+5°F to +122°F) Humidity:10%~ 90% RH	Temp: 14°F~158°F (-10°C~70 °C) Humidity: < 95% RH	Similarity
Dimensions	19.0cm×11.0cm×17.5cm	35.2cm×20.6cm×19.2cm	38.0cm×25.0cm×26.0cm	Similarity
Weight	1.7 kg	3.3 kg	2.6 kg	Similarity

ACARE TECHNOLOGY CO., LTD.**Table 5-2.** (Continued) Summary of the Proposed and Predicate Devices Technological Characteristics

Feature	Proposed device	Predicate device (a)	Predicate device (b)	Comparison
Performance Testing	Acare Suction Unit, model ASU-200 (K182950)	VacPlus Suction Unit (K103535)	EMG Suction Unit (K112421)	
Electromagnetic compatibility (EMC) <i>IEC 60601-1-2</i>	Passed	Passed	Passed	Same
Basic safety and Essential performance <i>IEC 60601-1 (ANSI/AAMI ES 60601-1)</i>	Passed	Passed	Passed	Same
Electrically powered suction equipment - Safety requirements <i>ISO 10079-1</i>	Passed	Passed	Passed	Same
Degrees of protection provided by enclosures (IP Code) <i>IEC 60529</i>	IP 21	IP 21	IP 21	Same

ACARE TECHNOLOGY CO., LTD.

5.7 Summary of Non-Clinical Testing

The Acare Suction Unit, model ASU-200 has the identical intended use and indication for use as the predicate devices. The subject device was compared to the predicate device by testing the electrical specifications, vacuum range, flow rate, power options, and operation environment.

The results of the Acare Suction Unit, model ASU-200 validation studies demonstrate that the Suction Unit perform as intended. The results are summarized as follows:

- The electromagnetic compatibility (EMC) testing performed as described in IEC 60601-1-2 Edition 3: 2007-03, “Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests”. The results confirm that the testing results demonstrate of the Acare Suction Unit, model ASU-200 met the electromagnetic compatibility (EMC) encompasses both emissions (interference with electronic products) and immunity (interference with device performance created by emissions from other electronic products).
- The protection against electrical hazards from equipment, Protection against excessive temperatures and other hazards, hazards situations and fault conditions, construction of equipment and electromagnetic compatibility of equipment and systems, etc. testing performed as described in ANSI/AAMI ES 60601-1:2005/(R)2012 and C1:2009/(R)2012, “Medical electrical equipment – Part 1: General requirements for basic safety and essential performance”. The testing results demonstrate of the Acare Suction Unit, model ASU-200 met safety and essential performance testing acceptance criteria for each test mentioned above.
- The Needle Flame Test, Measurement of the patient leakage current, Measurement of the vibration and noise of equipment other than the low vacuum equipment, Liquid overflow/ spillage test, Air Leakage test, Interruption of the power supply, Abnormal operation and fault conditions, Environmental tests for equipment for field and/or transport use, Degree of collapse test, Airflow and Vacuum tests, and Resistance to implosion, etc. testing performed as described in ISO 10079-1:2009, “Medical suction equipment Part 1: Electrically powered suction equipment — Safety requirements”. The testing results demonstrate of the Acare Suction Unit, model ASU-200 met Safety requirements testing acceptance criteria for each test mentioned above.
- The classification of degrees of protection provided by enclosures for electrical equipment testing performed as described in IEC 60529-2013, “Degrees of protection provided by enclosures (IP Code)”. The testing results demonstrate of the Acare Suction Unit, model ASU-200 met IP 21 testing acceptance criteria for each test mentioned above.

8. CONCLUSION

The conclusions drawn from the-nonclinical tests demonstrate that the proposed device is as safe, as effective, and performs as well as the predicate device K103535 and K112421.