



January 2, 2020

NuvoAir AB
% Yolanda Smith
Consultant
Smith Associates
1468 Harwell Ave
Crofton, Maryland 21114

Re: K183089
Trade/Device Name: Air Next
Regulation Number: 21 CFR 868.1840
Regulation Name: Diagnostic Spirometer
Regulatory Class: Class II
Product Code: BZG
Dated: December 20, 2019
Received: December 30, 2019

Dear Yolanda Smith:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Michael Ryan
Division Director
DHT1C: Division of ENT, Sleep Disordered
Breathing, Respiratory and
Anesthesia Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices Office of Product
Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K183089

Device Name

Air Next

Indications for Use (Describe)

Air Next is intended to be used by:

Healthcare professionals trained to perform spirometry tests on patients of age ≥ 5 years old, > 10 kg and ≥ 110 cm.

Air Next is intended to perform basic lung function and spirometry testing. The actual diagnosis shall be done by a healthcare professional in hospital and clinical settings.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

SPONSOR

Company Name: NuvoAir AB
Company Address: 17 Riddargatan
Stockholm 114 57
Sweden

Telephone: +46 725538233
Contact Person: Vanessa Caldeira

Summary Preparation Date: December 31, 2019

DEVICE NAME

Trade Name: Air Next
Common/Usual Name: Spirometer
Classification Name: Diagnostic Spirometer
Regulation Number: 21 CFR 868.1840
Product Code: BZG
Device Class: Class II
Reviewing Panel: Anesthesiology

PREDICATE DEVICE

Legally Marketed Equivalent Device

Predicate	K Number	Brand Name	Company
Primary	K072979	Spirobank G	MIR Medical International Research
Reference	K061712	Spirobank II	MIR Medical International Research

DEVICE DESCRIPTION

Air Next is intended to perform basic lung function and spirometry testing. It measures parameters such as the forced expiratory volume in 1 sec (FEV1) and the forced vital capacity (FVC) in a forced expiratory maneuver. These measures can be used for detection, assessment and monitoring of diseases affecting the lung function, such as bronchial Asthma, COPD and Cystic Fibrosis.

Air Next is a hand-held spirometer, weighing only 75g and it is powered by 2 AAA alkaline 1.5V batteries. It consists of 3 main components: the Air Next device, the NuvoAir disposable turbine (delivered in one package) and the Air Next mobile application downloadable from Apple's and Google's Play Stores. The package includes one Air Next device, one NuvoAir disposable turbine through which the user can start using the Air Next, two AAA 1.5V alkaline batteries, one user manual and one cotton bag to carry the device. Dimensions of the Air Next are 98 x 62 x 26 mm, which can be seen in the engineering drawing in the attachments.

DEVICE INDICATIONS FOR USE

Air Next is intended to be used by:

Healthcare professionals trained to perform spirometry tests on patients of age ≥ 5 years old, > 10 kg and ≥ 110 cm.

Air Next is intended to perform basic lung function and spirometry testing. The actual diagnosis shall be done by a healthcare professional in hospital and clinical settings.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a dentist or physician.

COMPARISON OF TECHNICAL CHARACTERISTICS

Predicate Product Comparison Table

Feature	Subject Device	Primary Predicate K072979	Similarities and Differences
Manufacturer	NuvoAir	M.I.R. Medical International Research	
Product Code	BZG	BZG	Same
Regulation No,	21 CFR 868.1840	21 CFR 868.1840	Same
Indications for Use	<p>Air Next is intended to be used by: Healthcare professionals trained to perform spirometry tests on patients of age ≥ 5 years old, > 10 kg and ≥ 110 cm.</p> <p>Air Next is intended to perform basic lung function and spirometry testing. The actual diagnosis shall be done by a healthcare professional in hospital and clinical settings.</p>	<p>The Spirobank G spirometer is intended to be used by a physician or by a patient under the instruction of a physician or paramedic. The device is intended to test lung function and can make spirometry testing in people of all ages, excluding infants and neonates. It can be used in any setting.</p>	<p>Similar indications for the use. The only difference between subject device and predicate device is that the Indications for Use in the subject device includes also weight and height while Predicate device only uses age as a limitation. Subject device is not for use in home setting.</p>
Use Environment	Hospital and clinical settings.	All settings that include clinical and home use	Different Use is limited to clinical settings only.
Intended User	Trained health care professionals perform spirometry testing on Adults and pediatric patients over 5 years old	People of all ages, excluding infants and neonates.	Different
Prescription only	Yes	Yes	Same

Feature	Subject Device	Primary Predicate K072979	Similarities and Differences
Principle of Operation	<p>The Air Next is designed to work with NuvoAir disposable turbine. When performing a spirometry test, the user exhales into the turbine. The airflow generated is forcing a propeller to rotate inside the turbine. The Air Next registers the speed of the spinning propeller by counting the rotations with a digital infrared interruption sensor. The algorithm in the firmware inside the Air Next device then converts the rotations into airflow measured in liters per second.</p> <p>The device is also tested against B.T.P.S. (body temperature and pressure with saturated water vapor) conditions as prescribed by ATS guidelines and the results are well within range.</p>	<p>Spirobank G is a powerful and compact measurement device, intended for use by a respiratory specialist or by a suitably trained generalist.</p> <p>Two different types of turbine sensors can be used with the device, one is reusable and one is single-patient disposable. A disposable mouthpiece is required in order to connect a subject to the spirometer.</p> <p>The flow and volume measurement sensor is a digital turbine, based on the infrared interruption principle.</p> <p>Cleaning of the disposable turbine is not required, as it is supplied clean in a sealed plastic bag. It must be disposed of after use.</p> <p>Spirometry test interpretation is based on the Forced Vital Capacity (FVC) test and is based on the ATS standard.</p>	Same
Technical Feature			
Flow Sensor	Bidirectional turbine with infrared interruption	Bidirectional turbine with infrared	Same
Physical Characteristics			
Energy Type	Two AAA 1.5V Alkaline batteries	9V Alkaline battery	Different
Display	Touchscreen on smartphone or tablet	Touchscreen / LCD and membrane	Different
Operating Environment	T: min +10°C/max +40° C RH: min 10%/max 95 % ALT: max 2000 m	Temperature: MIN + 10 °C, MAX + 40 °C; Humidity: MIN 10% RH; MAX 95%RH	Same
Physical Size	Light hand-held device	Light hand-held device	Same

Feature	Subject Device	Primary Predicate K072979	Similarities and Differences
Weight	75g	160g	Different
Technical Specifications			
Volume accuracy	±3% of reading or ±0.050 L, whichever is greater	±3% of reading or ±0.050 L, whichever is greater	Same
Volume range	Up to 10 L	Up to 10 L	Same
Flow range	0-15L/s	0-16 L/s	Similar
Flow accuracy	±5% or 200 mL/s	±5% or 200 mL/s	Same
Flow resistance	<0.5 cmH2O/L/s	<0.5 cmH2O/L/s	Same
Connectivity	Bluetooth	Bluetooth or USB	Same - Bluetooth
Measured parameter			
Forced vital capacity	FVC	FVC	Same
Volume expired in the first	FEV0.75, FEV1, FEV3, FEV6	FEV0.75, FEV1, FEV3, FEV6	Same
Ratio between volume expired in a certain time period and FVC	FEV/FVC (FER) for 0.75 /1 /3 / 6	FEV/FVC (FER) for 0.75 /1 /3 / 6	Same
Peak expiratory flow	PEF	PEF	Same
Forced expiratory flow between the first 25% and 75% of the FVC	FEF25-75 (MEF)	FEF25-75 (MEF)	Same
Volume inspired in the first second of the test	FIV1	FIV1	Same
Forced inspiratory volume	FIVC	FIVC	Same
Peak inspiratory flow	PIF	PIF	Same

Feature	Subject Device	Primary Predicate K072979	Similarities and Differences
Forced inspiratory flow in the first 25% and 75% of the FIV	FIF25-75 (MIF25-75)	FIF25-75 (MIF25-75)	Same
Ratio between FEV1 and FEV6	FEV1/FEV6	FEV1/FEV6	Same
Forced expiratory flow at 50% of FVC divided by FVC	FEF50/FVC	FEF50/FVC	Same
Forced inspiratory flow in the first 25% and 75% of the FIV	FIF25-75 (MIF25-75)	FIF25-75 (MIF25-75)	Same
Ratio between FEV1 and FEV6	FEV1/FEV6	FEV1/FEV6	Same
Forced expiratory flow at 50% of FVC divided by FVC	FEF50/FVC	FEF50/FVC	Same
Force inspiratory flow at first second device forced inspiratory volume	FIV1/FIVC (FIR)	FIV1/FIVC (FIR)	Same
Forced expiratory time	FET	FET	Same
Maximum voluntary ventilation	N/A	MVV (ind)	Different
Device Use Design			
Used with PFT filter and / or mouthpiece	Mouthpiece Disposable	Mouthpiece Disposable or Reusable	Same for Disposable Turbine – See reference predicate chart below
Coaching	Display of flow-time curves, timer and effort during test	Display of volume-time and flow-volume curves	Different
Feedback on test quality	Yes - according to ATS guidelines	Yes – Traffic lights for interpretation	Similar
Graphic display	Yes- on App	Yes on device screen	Different
Safety Standards			

Feature	Subject Device	Primary Predicate K072979	Similarities and Differences
Applied Standards for safety	Electrical Safety Standard IEC 60601-1 Electro Magnetic Compatibility IEC 60601-1-2	Electrical Safety Standard IEC 60601-1 Electro Magnetic Compatibility IEC 60601-1-2	Same
Applied Standards for safety	Electrical Safety Standard IEC 60601-1 Electro Magnetic Compatibility IEC 60601-1-2	Electrical Safety Standard IEC 60601-1 Electro Magnetic Compatibility IEC 60601-1-2	Same
Biocompatibility	According to ISO 10993-1	According to ISO 10993- 1	Same
Meets ATS accuracy requirements	Yes	Yes	Same

Reference Device

The disposable turbine used in Air Next device is manufactured and purchased from MIR Medical International Research and is cleared under K061712. The disposable turbine is supplied in its final finished form and is identical to the reference predicate device.

PERFORMANCE TESTING

The following performance data were provided in support of the substantial equivalence determination.

Safety testing

Testing Standard	Test Description	Result
IEC 60601-1:2005 (Third Edition) + CORR. 1:2006 + CORR. 2:2007 + A1:2012	Medical electrical equipment – Part 1: general requirements for basic safety and essential performance	Pass
IEC 60601-1-6:2010, AMD1:2013	Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability	Pass
IEC 62366:2007	Medical devices – Application of usability engineering to medical devices	Pass
IEC 60601-1-2:2014 (ED. 4)	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests	Pass

Biocompatibility

Handle and Exterior Surface

The Handle and Exterior is designed to be a surface device for intact skin with limited contact A – limited (≤24 h). The following testing was performed:

Based on the device categorization outlined above, the nature of the material contact and the guidance given in ISO 10993-1, the tests deemed appropriate for the device materials were selected as shown in the table below.

Biological Effect	Test	Compliance Standard
Cytotoxicity	Mem Elution	ISO10993-5
Irritation or Intracutaneous Reactivity	Intracutaneous Injection Test	ISO10993-10
Irritation and Sensitization	Kligman Maximization test	ISO10993-10

Disposable Turbine

NuvoAir Turbines are disposable turbines intended for single use for a period of less than 24hr. NuvoAir has provided a declaration statement from the legal manufacturer, MIR, stating that the NuvoAir disposable turbine for Air Next Spirometer (ref. code 910007) is identical in its final finished form to the disposable turbine legally marketed and identified with code 910004 and that this turbine is in turn legally US-marketed Spirobank II (K061712).

Additional Testing Standards

- ATS/ERS 2005 guideline
- ISO 14971:2007 - Application of risk management to medical devices.
- IEC 62366:2007 + A1:2014. *Application of usability engineering to medical devices.*
- ISO 62304:2006/A1:2015: Medical device software. Software life-cycle processes
- EN ISO 26782: 2009, Anesthetic and respiratory equipment – Spirometers intended for the measurement of time forced expired volumes in humans
- EN ISO 23747:2015, Anesthesia and respiratory equipment – Peak expiratory flowmeters for the assessment of pulmonary function in spontaneously breathing humans
- FCC Part 15, Subpart B, C

Clinical testing was not performed with this device.

Substantial Equivalence Conclusion:

Based upon the foregoing performance testing and comparison to the legally marketed predicate devices for indications for use, technology, and performance we believe we have demonstrated that the Air Next is substantially equivalent and as safe and as effective as the predicate device. Furthermore, the differences versus the predicate device do not impact safety and effectiveness.