



August 16, 2019

Trudell Medical International  
c/o Paul Dryden  
Consultant  
725 Third Street  
London, N5V 5G4 CA

Re: K183108

Trade/Device Name: Combined VersaPAP device and Aerobika OPEP  
Regulation Number: 21 CFR 868.5690  
Regulation Name: Incentive Spirometer  
Regulatory Class: Class II  
Product Code: BWF  
Dated: July 18, 2019  
Received: July 19, 2019

Dear Paul Dryden:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

James Lee, Ph.D.  
Assistant Director  
DHT1C: Division of ENT, Sleep Disordered  
Breathing, Respiratory and  
Anesthesia Devices  
OHT1: Office of Ophthalmic, Anesthesia,  
Respiratory, ENT and Dental Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

**K183108**

Device Name

**Combined VersaPAP™ device and Aerobika® OPEP**

Indications for Use (Describe)

The Combined VersaPAP™ and Aerobika® OPEP device is intended for use as a Positive Airway Pressure (PAP) device and a Positive Expiratory Pressure (PEP) device. The combined device has the ability to provide supplemental oxygen when used with compressed oxygen. The combined device is for patients (ages 5 years and above) who are capable of following directions for Positive Airway Pressure Therapy and capable of generating exhalation flow of 10 lpm for 3 – 4 seconds. The combined device is a single patient, multiple use device intended to be used in a hospital environment under the supervision of a healthcare professional.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)  Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRStaff@fda.hhs.gov](mailto:PRStaff@fda.hhs.gov)

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

**510(k) Summary****Date Prepared** 16-Aug-19

Trudell Medical International  
725 Third Street  
London, ON N5V 5G4 CANADA  
Tel – 519-455-7060

**Official Contact:** Marianne Tanton – Director, Quality and Regulatory Affairs**Proprietary or Trade Name:** Combined VersaPAP™ device and Aerobika® OPEP**Common/Usual Name:** Incentive Spirometer**Classification Code/Name:** BWF – Incentive Spirometer  
21 CFR 868.5690  
Class II**Predicate Devices:** K150173 – Trudell – Aerobika® OPEP  
K173918 – Trudell – VersaPAP™ device**Reference Device:** K151689 – Hill-Rom – MetaNeb®**Device Description:**

The Aerobika® OPEP device and VersaPAP™ device can be used independently or combined. When combined for ease of use. The combined device creates oscillations during the exhalation phase (OPEP) to aid in secretion clearance and during the inhalation phase (PAP) helps to maintain positive airway pressure for the treatment and prevention of atelectasis.

**Indications for Use:**

The Combined VersaPAP™ and Aerobika® OPEP device is intended for use as a Positive Airway Pressure (PAP) device and a Positive Expiratory Pressure (PEP) device. The combined device has the ability to provide supplemental oxygen when used with compressed oxygen. The combined device is for patients (ages 5 years and above) who are capable of following directions for Positive Airway Pressure Therapy and capable of generating exhalation flow of 10 lpm for 3 – 4 seconds. The combined device is a single patient, multiple use device intended to be used in a hospital environment under the supervision of a healthcare professional.

Table 1 – Comparison to Predicates

	<b>Proposed Combined Aerobika® OPEP and VersaPAP™ device</b>	<b>Predicate Aerobika® OPEP K150173</b>	<b>Predicate VersaPAP™ device K173918</b>
<b>CFR Classification</b>	868.5690 BWF		
<b>Classification name</b>	Incentive Spirometer		
<b>Intended Use</b>	<p>Intended for use as a Positive Expiratory Pressure (PEP) device (exercise patient's lungs and to improve secretion clearance)</p> <p>Intended for use as a Positive Airway Pressure (PAP) device (Treatment and prevention of atelectasis)</p>	Intended for use as a Positive Expiratory Pressure (PEP) device.	Treatment and prevention of atelectasis
<b>Indications for Use</b>	<p>The Combined VersaPAP™ and Aerobika® OPEP device is intended for use as a Positive Airway Pressure (PAP) device and a Positive Expiratory Pressure (PEP) device. The combined device has the ability to provide supplemental oxygen when used with compressed oxygen. The combined device is for patients (ages 5 years and above) who are capable of following directions for Positive Airway Pressure Therapy and capable of generating exhalation flow of 10 lpm for 3 – 4 seconds. The combined device is a single patient, multiple use device intended to be used in a hospital environment under the supervision of a healthcare professional</p>	<p>The Aerobika® Oscillating Positive Expiratory Pressure device is intended for use as a Positive Expiratory Pressure (PEP) device. The Aerobika® Oscillating PEP device may also be used simultaneously with nebulized aerosol drug delivery. The device is intended to be used by patients capable of generating exhalation flow of 10 lpm for 3–4 seconds.</p>	<p>The VersaPAP™ device is indicated for the treatment and prevention of atelectasis. The VersaPAP™ device also has the ability to provide supplemental oxygen when used with compressed oxygen. The VersaPAP™ device is for patients (ages 5 years and above) who are capable of following directions for Positive Airway Pressure Therapy in a hospital environment. The VersaPAP™ device is a single patient, multiple use device intended to be used under the supervision of a healthcare professional.</p>
<b>Patient Population</b>	<p>Patients capable of generating exhalation flow of 10 lpm for 3–4 seconds.</p> <p>Patients (ages 5 years and above) who are capable of following directions.</p>	<p>Patients capable of generating exhalation flow of 10 lpm for 3–4 seconds.</p>	<p>Patients (ages 5 years and above) who are capable of following directions.</p>

	<b>Proposed Combined Aerobika® OPEP and VersaPAP™ device</b>	<b>Predicate Aerobika® OPEP K150173</b>	<b>Predicate VersaPAP™ device K173918</b>
<b>Environments</b>	Hospital settings	Hospital and clinical Home settings	Hospital settings
<b>Operating principles</b> <b>OPEP</b>	Oscillations generated by patient's exhaled breath. Adjustable exhalation resistance  Visual feedback via manometer	Oscillations generated by patient's exhaled breath. Adjustable exhalation resistance Visual feedback only when used with the manometer	
<b>PAP</b>	Venturi principle created by compressed gas source Supplemental oxygen Measures airway pressure with manometer		Venturi principle created by compressed gas source Supplemental oxygen Measures airway pressure with manometer
<b>Single patient, multi-use</b>	Components may be cleaned and reused by the same patient		

Table 2 - Comparison of the Subject vs. the Reference Hill-Rom MetaNeb (K151689)

	<b>Proposed Combined OPEP + PAP</b>	<b>Reference Hill-Rom - MetaNeb® K151689</b>
<b>Classification</b>	BWF 868.5690 Incentive Spirometer	NHJ 868.5905 Device, positive pressure breathing, intermittent (IPPB)
<b>Intended use</b>	Intended for use as a Positive Expiratory Pressure (PEP) device (exercise patient's lungs and to improve secretion clearance)  Intended for use as a Positive Airway Pressure (PAP) device (Treatment and prevention of atelectasis)	Intended for use as a Positive Expiratory Pressure (PEP) device.  Treatment and prevention of atelectasis
<b>Indications for Use</b>	The Combined VersaPAP™ and Aerobika® OPEP device is intended for use as a Positive Airway Pressure (PAP) device and a Positive Expiratory Pressure (PEP) device. The combined device has the ability to provide supplemental oxygen when used with compressed oxygen. The combined device is for patients (ages 5 years and above) who are capable of following directions for Positive Airway Pressure Therapy and capable of generating exhalation flow of 10 lpm for 3 – 4 seconds. The combined device is a single patient, multiple use device intended to be used in a hospital environment under the supervision of a healthcare professional	The MetaNeb® System is indicated for mobilization of secretions, lung expansion therapy, the treatment and prevention of pulmonary atelectasis, and also has the ability to provide supplemental oxygen when used with compressed oxygen.
<b>Patient Population</b>	Patients capable of generating exhalation flow of 10 lpm for 3–4 seconds.  Patients (ages 5 years and above) who are capable of following directions.	Patient 5 years and older
<b>Environments of Use</b>	Hospital	Hospital sub-acute facilities Nursing care Homecare

	<b>Proposed Combined OPEP + PAP</b>	<b>Reference Hill-Rom - MetaNeb® K151689</b>
<b>Mode of Operation</b>	OPEP Oscillations generated by patient's exhaled breath. Adjustable exhalation resistance Visual feedback via manometer  PAP Venturi vacuum created by compressed gas source Supplemental oxygen Measures airway pressure with manometer	Can deliver continuous expiratory pressure (CPEP) combined with medicated aerosol. Oscillations generated by device  Venturi from compressed gas source continuous  Supplemental oxygen Measures airway pressure with manometer
<b>Therapy</b>	OPEP PAP	CPEP CHFO Aerosol
<b>Patient Interface</b>	Mouthpiece	Disposable circuit includes connection for in-line nebulizer with mouthpiece
<b>Gas source</b>	Compressed air or oxygen	Pneumatic and air or oxygen
<b>Control settings</b>	OPEP resistance adjustable  PAP flow rate: 5-15 L/min.	Mode selection  Frequency selection for CHFO mode  Pressure adjustment for CPEP mode

**Predicate, Reference and Proposed Device Comparison and Substantial Equivalence Discussion:**

The proposed and predicates are similar in the following categories:

- Indications for Use
- Patient Population
- Technological Characteristics
- Materials in patient contact
- Shelf-life
- Cleaning methods

The Combined Aerobika® OPEP and VersaPAP™ device is substantially equivalent to the predicate and reference devices because:

**Indications –**

The proposed indications for use for secretions clearance and the treatment and prevention of pulmonary atelectasis, and also has the ability to provide supplemental oxygen when used with compressed oxygen are identical to the predicate and reference devices.

**Discussion –** There are no differences between the subject and predicate devices.

**Patient Population –**

The patient population is identical to the predicates.

**Discussion –** There are no differences between the subject and predicate devices.



**Environment of Use –**

The environments of use are similar to the predicates. Noting that the OPEP also includes the home setting when used individually, but the home setting is not intended for the subject device.

**Discussion** – There are no differences between the subject and predicate devices except the combined device will not be used in the home setting.

**Technology –**

The technology for generating OPEP and PAP is identical to the predicates. The technology of having a device which offers the combination of OPEP and PEP is similar to the reference device.

**Discussion** – There are no differences between the subject and predicate devices.

**Performance –**

The pressure was compared to the reference device. While the pressure of the subject device is less than that of the reference device, it is higher than that of the predicates when they are tested as standalone devices.

**Discussion** – The difference in performance is within the range of the separate predicate devices and the reference device and does not raise new safety or effectiveness concerns.

The frequency of the combined device is higher than that of the reference MetaNeb device, however, the measured frequency is comparable to the Aerobika® OPEP predicate device and is within acceptable clinical frequency efficacy range.

**Discussion** - The frequency of the subject device is comparable to the Aerobika® OPEP predicate device. As such, the frequency generated by the subject device does not add risk to the patient compared to other FDA cleared devices.

**Non-clinical Comparative Performance****Biocompatibility –**

The materials that are patient contacting and gas pathway are same as K150173 – Trudell – Aerobika® OPEP and K173918 – Trudell – VersaPAP™ device

**Discussion** – Additionally testing was performed, ISO 18562, which supported patient exposure in new configurations.

**Bench Testing -**

We performed comparative testing at adult and child breathing setting to the predicates and reference. The results demonstrate that the subject device is substantially equivalent to the predicate and reference devices.

**Discussion of Differences**

The only differences are:

- Environment is only hospital setting
  - Where the Aerobika® OPEP device when used individually can also be used in the home setting
- Certain performance categories differ when the devices are combined vs. when they are used separately.
  - These differences relate to pressure and frequency (oscillations) when they are combined. When compared to the predicates and reference device, these differences in pressure and frequency are within the performance specifications of the cleared predicates and reference device.

Any differences do not raise new or different concerns of safety and effectiveness.

**Substantial Equivalence Conclusion**

Based upon the comparative performance testing we have demonstrated that the proposed device compared to the predicate and reference devices can be found to substantially equivalent.