



January 4, 2019

BostonSight  
% Bret Andre  
Principal Consultant  
EyeReg Consulting, Inc.  
6119 Canter Ln  
West Linn, OR 97068

Re: K183175  
Trade/Device Name: BostonSight Scleral  
Regulation Number: 21 CFR 886.5916  
Regulation Name: Rigid Gas Permeable Contact Lens  
Regulatory Class: Class II  
Product Code: HQD  
Dated: November 13, 2018  
Received: November 16, 2018

Dear Bret Andre:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely yours,

**J Angelo Green**

for Malvina B. Eydelman, M.D.  
Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K183175

Device Name

BostonSight Scleral

### Indications for Use (Describe)

BostonSight Scleral daily wear contact lenses are indicated for the correction of refractive ametropia (myopia, hyperopia, astigmatism and presbyopia) in aphakic and non aphakic persons. Also, the lenses may be prescribed in eyes that require a rigid contact lens for the management of irregular corneal conditions such as keratoconus, pellucid marginal degeneration, keratoglobus, post-LASIK ectasia, or following penetrating keratoplasty or refractive surgery.

Furthermore, eyes suffering from certain ocular surface disorders may benefit from the physical protection, aqueous hydrated environment and the saline bath provided by scleral lens designs.

BostonSight Scleral daily wear contact lenses are indicated for therapeutic use for the management of irregular and distorted corneal surfaces where the subject:

1. cannot be adequately corrected with spectacle lenses
2. requires a rigid gas permeable contact lens surface to improve vision
3. is unable to wear a corneal rigid gas permeable lens due to corneal distortion or surface irregularities

Common causes of corneal distortion include, but are not limited to, corneal infections, trauma, tractions as a result of scar formation secondary to refractive surgery (e.g. LASIK or radial keratotomy) or corneal transplantation. Causes may also include corneal degeneration (e.g. keratoconus, keratoglobus, pellucid marginal degeneration, Salzmann's nodular degeneration) and corneal dystrophy (e.g., lattice dystrophy, granular corneal dystrophy, Reis-Bucklers dystrophy, Cogan's dystrophy).

The BostonSight Scleral daily wear contact lenses are also indicated for therapeutic use in eyes with ocular surface disease including, but not limited to, ocular Graft-versus-Host disease, Sjögren's syndrome, dry eye syndrome and Filamentary Keratitis, limbal stem cell deficiency (e.g. Stevens-Johnson syndrome, chemical radiation and thermal burns), disorders of the skin (e.g. atopy, ectodermal dysplasia), neurotrophic keratitis (e.g. Herpes simplex, Herpes zoster, Familial Dysautonomia), and corneal exposure (e.g. anatomic, paralytic) that might benefit from the presence of an expanded tear reservoir and protection against an adverse environment. When prescribed for therapeutic use for a distorted cornea or ocular surface disease, the BostonSight Scleral daily wear contact lenses may concurrently provide correction of refractive error.

The lenses may be disinfected using a chemical disinfection (not heat) system only.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**



## 510 (k) SUMMARY OF SAFETY AND EFFECTIVENESS

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

**The assigned 510(k) number is: 183175**

### **I. SUBMITTER**

Date Prepared: November 13<sup>th</sup>, 2018

Name: **BostonSight**  
Address: 464 Hillside Avenue, Suite 205  
Needham, MA 02494

Contact Person: Gene P. Guselli  
President and CEO  
Phone number: 781-726-7337

Consultant: Bret Andre  
EyeReg Consulting, Inc.  
6119 Canter Ln.  
West Linn, OR 97068  
Phone number: (503) 372-5226

### **II. DEVICE**

Trade Name: **BostonSight Scleral**

Common  
Name: Daily wear rigid gas permeable contact lens

Classification  
Name: Rigid gas permeable contact lens. (21 CFR 886.5916)

Regulatory  
Class: Class II

Product Code: HQD



## 510(K) Premarket Notification

~Reason for Submission~

Labeling Changes - Expanded Indications for Use

**III. PREDICATE DEVICE**

**BostonSight Scleral** daily wear contact lenses are substantially equivalent to the following predicate devices:

- “Boston XO® (hexafocon A) and Boston XO2® (hexafocon B) Rigid Gas Permeable Contact Lens”  
By Bausch + Lomb Inc.  
510(k) number; **K171404**
- “OPTIMUM GP (Roflufocon D, Roflurocon E) Daily Wear Contact Lens”  
By Contamac Ltd.  
510(k) number; **K180616**
- “BostonSight PD Prosthetic Device”  
By BostonSight  
510(k) number; **K161461**

**IV. DEVICE DESCRIPTION**

The **BostonSight Scleral** daily wear contact lenses are manufactured with a large diameter RGP lens design that vaults over the cornea and rests on the conjunctiva overlying the sclera. The **BostonSight Scleral** daily wear contact lenses are lathe cut and fabricated from one of the following fluoro-silicone acrylate rigid gas permeable (RGP) lens materials:

- roflufocon D supplied by Contamac Ltd.
- roflufocon E supplied by Contamac Ltd.
- oprifocon A supplied by Bausch and Lomb, Inc.
- hexafocon B supplied by Bausch and Lomb, Inc.

The **BostonSight Scleral** daily wear contact lenses may be shipped “dry” or “wet” in a polypropylene contact lens case. The primary container for shipping the **BostonSight Scleral** daily wear contact lenses is the Bausch & Lomb Frequent Replacement Contact Lens Case, with clearance under 510(k) K896685.

When shipped “wet”, The **BostonSight Scleral** daily wear contact lenses manufactured from material supplied by Bausch & Lomb, inc. may be packaged and shipped in Boston Advance Comfort Formula Conditioning Solution (K974466) or Boston SIMPLUS Multi-Action solution (K024289). The **BostonSight Scleral** daily wear contact lenses manufactured from material supplied by Contamac, Ltd. may be packaged and shipped “wet” in in the OPTIMUM by Lobob Cleaning and Disinfecting Storage solution, with clearance under 510(k) K014162.



## 510(K) Premarket Notification

The physical properties of the **BostonSight Scleral** daily wear contact lenses are as follows:

	<b>roflufocon D</b>	<b>roflufocon E</b>	<b>oprifocon A</b>	<b>hexafocon B</b>
<b>Refractive Index</b>	1.4333	1.4332	1.4230	1.4240
<b>Light Transmission (clear)</b>	>97%	>97%	>95%	>95%
<b>Light Transmission (tinted)</b>	>90%	>90%	>90%	>83%
<b>Specific Gravity</b>	1.166	1.155	1.24	1.19
<b>Oxygen Permeability (Dk) ISO/FATT Method</b>	100 x 10 <sup>-11</sup> (cm <sup>2</sup> /sec) (ml O <sub>2</sub> /ml x mm Hg @ 35°C)	125 x 10 <sup>-11</sup> (cm <sup>2</sup> /sec) (ml O <sub>2</sub> /ml x mm Hg @ 35°C)	85 x 10 <sup>-11</sup> (cm <sup>2</sup> /sec) (ml O <sub>2</sub> /ml x mm Hg @ 35°C)	141 x 10 <sup>-11</sup> (cm <sup>2</sup> /sec) (ml O <sub>2</sub> /ml x mm Hg @ 35°C)
<b>Visitint lenses contain one or more of the following color additives conforming to: 21 CFR Part 73 &amp; 74, Subpart D</b>	D & C Green No. 6, FD & C Red No. 17, CI Solvent Yellow 18	D & C Green No. 6, FD & C Red No. 17, CI Solvent Yellow 18	D&C Green No.6 and D&C Yellow No.10	D&C Green No. 6; C.I. Solvent Yellow No. 18; D&C Violet No. 2; D&C Red No. 17; C.I. Solvent Yellow No.18
<b>UV Light Blocking (UVB – 280nm – 315nm; UVA 316nm – 380nm)</b>	>98% UVB >95% UVA	>98% UVB >95% UVA	>95% UVB >97% UVA	>95% UVB >97% UVA
<b>Dynamic Receding Contact Angle</b>	3°	6°	56°	40°

The **BostonSight Scleral** daily wear contact lenses are available in the following lens parameters:

<b>Parameter</b>	<b>Range</b>	<b>Tolerance</b>
Base Curve	5.00mm to 9.00mm	± 0.05 mm
Center Thickness	0.05mm to 0.60mm	± 0.02 mm
Diameter	8.00mm to 26.00mm	± 0.10mm
Spherical Power	-25.00 D to +35.00 D (in 0.25D steps)	± 0.12 (0 to <= 5D) ± 0.18 (5 to <= 10.0D) ± 0.25 (10 to <= 15D) ± 0.37 (15 to <= 20D) ± 0.50 (over 20D)

## V. INDICATIONS FOR USE

**BostonSight Scleral** daily wear contact lenses are indicated for the correction of refractive ametropia (myopia, hyperopia, astigmatism and presbyopia) in aphakic and non aphakic persons. Also, the lenses may be prescribed in eyes that require a rigid contact lens for the management of irregular corneal conditions such as keratoconus, pellucid marginal degeneration, keratoglobus, post-LASIK ectasia, or following penetrating keratoplasty or refractive surgery.




---

 510(K) Premarket Notification
 

---

Furthermore, eyes suffering from certain ocular surface disorders may benefit from the physical protection, aqueous hydrated environment and the saline bath provided by scleral lens designs.

**BostonSight Scleral** daily wear contact lenses are indicated for therapeutic use for the management of irregular and distorted corneal surfaces where the subject:

1. cannot be adequately corrected with spectacle lenses
2. requires a rigid gas permeable contact lens surface to improve vision
3. is unable to wear a corneal rigid gas permeable lens due to corneal distortion or surface irregularities

Common causes of corneal distortion include, but are not limited to, corneal infections, trauma, tractions as a result of scar formation secondary to refractive surgery (e.g. LASIK or radial keratotomy) or corneal transplantation. Causes may also include corneal degeneration (e.g. keratoconus, keratoglobus, pellucid marginal degeneration, Salzmann's nodular degeneration) and corneal dystrophy (e.g., lattice dystrophy, granular corneal dystrophy, Reis-Bucklers dystrophy, Cogan's dystrophy).

The **BostonSight Scleral** daily wear contact lenses are also indicated for therapeutic use in eyes with ocular surface disease including, but not limited to, ocular Graft-versus-Host disease, Sjögren's syndrome, dry eye syndrome and Filamentary Keratitis, limbal stem cell deficiency (e.g. Stevens-Johnson syndrome, chemical radiation and thermal burns), disorders of the skin (e.g. atopy, ectodermal dysplasia), neurotrophic keratitis (e.g. Herpes simplex, Herpes zoster, Familial Dysautonomia), and corneal exposure (e.g. anatomic, paralytic) that might benefit from the presence of an expanded tear reservoir and protection against an adverse environment. When prescribed for therapeutic use for a distorted cornea or ocular surface disease, the **BostonSight Scleral** daily wear contact lenses may concurrently provide correction of refractive error.

The lenses may be disinfected using a chemical disinfection (not heat) system only.

## VI. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH PREDICATE DEVICE

The **BostonSight Scleral** daily wear contact lenses are substantially equivalent to the predicate devices (cleared under K171404, K180616, and K161461) in terms of the following:

- Intended use – daily wear contact lenses
- Indications for use
- Actions
- Classification – Lenses, Rigid Gas Permeable, Daily Wear (21 CFR 886.5916)
- FDA material group – group # 3 fluoro silicone acrylate
- USAN materials (roflucocon D, roflucocon E, oprifocon A and hexafocon B)
- Production method – lathe cut
- Final packaging and shipping



## 510(K) Premarket Notification

The following matrix illustrates the production method, lens function and material characteristics of the **BostonSight Scleral** daily wear contact lenses, as well as the predicate device.

### Substantial Equivalence Matrix

	BostonSight Scleral	Boston XO® and Boston XO2® RGP Contact Lenses	OPTIMUM GP Contact Lens	BostonSight PD Prosthetic Device
	Subject Device	Predicate Device (K171404)	Predicate Device (K180616)	Predicate Device (K161461)
<b>Classification</b>	Same as predicate	Class II Lenses, Rigid Gas Permeable, Daily Wear 21 CFR 886.5916	Class II Lenses, Rigid Gas Permeable, Daily Wear 21 CFR 886.5916	Class II Lenses, Rigid Gas Permeable, Daily Wear 21 CFR 886.5916
<b>Product Code</b>	Same as predicate	HQD	HQD	HQD
<b>FDA Group #</b>	Same as predicate	Group # 3 Fluoro Silicone Acrylate	Group # 3 Fluoro Silicone Acrylate	Group # 3 Fluoro Silicone Acrylate
<b>USAN</b>	roflufocon D, roflufocon E, oprifocon A, hexafocon B	hexafocon A, hexafocon B	roflufocon D, roflufocon E,	roflufocon D, roflufocon E, oprifocon A, hexafocon B
<b>Production Method</b>	Same as predicate	Lathe-Cut	Lathe-Cut	Lathe-Cut
<b>Actions</b>	Same as predicate	The contact lenses act as a refractive medium that focus light rays from near and distant objects on the retina. When placed on the eye for therapeutic use, the lens replaces or supports impaired ocular surface function.	The contact lenses act as a refractive medium that focus light rays from near and distant objects on the retina. When placed on the eye for therapeutic use, the lens replaces or supports impaired ocular surface function.	The contact lenses act as a refractive medium that focus light rays from near and distant objects on the retina. When placed on the eye for therapeutic use, the lens replaces or supports impaired ocular surface function.
<b>Intended Use</b>	Same as predicate	Daily Wear	Daily Wear	Daily Wear
<b>Indication for Use</b>	Same as predicate	Therapeutic (Scleral)	Therapeutic (Scleral)	Therapeutic (Scleral)
<b>Water Content (%)</b>	<1%	<1%	<1%	<1%
<b>UV Absorber Available</b>	Yes	Yes	Yes	Yes





## VII. PERFORMANCE DATA

### *~ Non-Clinical Studies ~*

The purpose of this application is to modify the labeling of previously FDA cleared RGP contact lenses/materials to include therapeutic indications for use. Non-clinical testing to demonstrate the safety and effectiveness of contact lenses manufactured from roflufocon D, roflufocon E, oprifocon A and hexafocon B materials has been addressed in previous applications.

### *~ Clinical Studies ~*

Clinical performance data to demonstrate the safety and effectiveness of contact lenses manufactured from roflufocon D, roflufocon E, oprifocon A and hexafocon B has been previously addressed.

## VIII. CONCLUSIONS

### Substantial Equivalence

Information presented in this Premarket Notification establishes that **BostonSight Scleral** daily wear contact lenses are as safe and effective as the predicate device when used in accordance with the labeled directions for use and for the proposed indications.

### Risks and Benefits

The risks of the subject device are the same as those normally attributed to the wearing of rigid gas permeable (RGP) daily wear contact lenses. The benefits to the patient are the same as those for other RGP contact lenses.