



March 1, 2019

Schwarzer CardioTek GmbH
% Melissa Walker
President and Chief Technology Officer
Graematter, Inc.
1324 Clarkson Clayton Ctr., #332
Ballwin, Missouri 63011

Re: K183266
Trade/Device Name: EP-TRACER System
Regulation Number: 21 CFR 870.1425
Regulation Name: Programmable Diagnostic Computer
Regulatory Class: Class II
Product Code: DQK
Dated: February 4, 2019
Received: February 4, 2019

Dear Melissa Walker:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K183266

Device Name

EP-TRACER System

Indications for Use (Describe)

The EP-TRACER System is an electrophysiology measurement system used to acquire, filter, digitize, amplify, display, and record signals obtained during electrophysiological studies and related procedures.

The system allows the user to monitor, display and record the signals. The system incorporates a stimulator intended to be used for diagnostic cardiac stimulation during electrophysiological testing of the heart.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Section 5: 510(k) Summary per 21CFR §807.92

EP TRACER Summary

Submitter's information

Contact: Tatjana Nuber
Schwarzer CardioTek GmbH
Im Zukunftspark 3
74076 Heilbronn / Germany

Contact: Melissa Walker
Graematter, Inc.
1324 Clarkson Clayton Ctr #332
Ballwin, MO 63011
Phone: 636-405-7498
Date: 28 Feb 2019

Device/ classification name

Device Name:
The EP-TRACER is a Class 2 device (product code DQK).

Classification/Common name:

- Programmable diagnostic computer, 21 CFR §870.1425.

The marketed device(s) to which substantial equivalence is claimed:

- EP-TRACER System, K161245.
-

Device description

The EP-TRACER system is a computerized electrophysiology measurement system designed for both regular and experimental EP studies.

The EP-TRACER is comprised of these major components:

1. EP-TRACER hardware – Amplifier/stimulator
 2. EP-TRACER Software – Software pre-installed
-

Indications for use

“The EP-TRACER System is an electrophysiology measurement system used to acquire, filter, digitize, amplify, display, and record signals obtained during electrophysiological studies and related procedures.

“The system allows the user to monitor, display and record the signals. The system incorporates a stimulator intended to be used for diagnostic cardiac stimulation during electrophysiological testing of the heart”.

Continued on next page

510(k) Summary per 21CFR §807.92, Continued

Technological characteristics The table below lists the technological characteristics for both the new and predicate devices. Modifications to the EP TRACER for V2.2 are limited to user interface modifications.

Device Characteristic	Predicate Device EP-TRACER V2.0	New Device EP-TRACER V2.2
Amplifier Dimensions (with integrated stimulator) WxDxH	38 channels: 28x27x7 in cm 70/102 channels: 28x27x12 in cm	Same
Environmental Specifications		
Temperature Operating	+10°C to +30°C	Same
Temperature Transport/Storage	-29°C to +66°C	Same
Humidity Operating	20 - 80 % rH (non-condensing)	Same
Humidity Transport/Storage	< 95 % rH (non-condensing)	Same
Power Specifications		
Power Requirements	100 - 240 V AC, 50 - 60 Hz	Same
Power Input	38 channels: +5 V, 0.3 A & +12 V, 0.9 A 70 channels: +5 V, 0.3 A & +12 V, 1.5 A 102 channels: +5 V, 0.3 A & +12 V, 2 A	Same
Design		
Sampling and Hold	Each channel sampled prior the acquisition	Same
Sampling Rate	1 kHz	Same
CMRR	> 100 dB	Same
Input Impedance	Typical 20 MΩ	Same
Leakage Current		
Patient Source	< 10 μA	Same
Patient Sink	< 10 μA	Same
Patient Sink (measured at patient leads under single fault conditions)	< 50 μA	Same
Chassis Leakage	< 100 μA	Same
ECG Input		
Outputs	12 lead ECG produced	Same
High Pass Filter	0.05 Hz, 0.2 Hz	Same
Low Pass Filter	150 Hz	Same
RF Filtering	All inputs	Same
Gain	Between 0 and 255 mm/mV - continuous	Same
Saturation Recovery	< 1 sec (manual reset)	Same
Notch Filter	Power line (50/60 Hz)	Same

Device Characteristic	Predicate Device EP-TRACER V2.0	New Device EP-TRACER V2.2
Dynamic Range	±5 mV	Same
Baseline Correction	±300 mV	Same
Input/ Output		
Inputs 32/38 channels	20 intracardiac channels, 6 auxiliary channels, 12 ECG channels	Same
Inputs 64/70 channels	52 intracardiac channels, 6 auxiliary channels, 12 ECG channels	Same
Inputs 102 channels	84 intracardiac channels, 6 auxiliary (pressure) channels, 12 ECG channels	Same
Outputs	No outputs	Same
Switching	Each channel can be either bipolar or unipolar with manual switching	Same
High Pass Filter	0.05 Hz, 0.2 Hz, 40 Hz, 80 Hz	Same
Low Pass Filter	350 Hz	Same
RF Filtering	All inputs	Same
Gain	Between 0 and 255 mm/mV – continuous	Same
Saturation Recovery	< 1 s (manual reset)	Same
Notch Filter	Power line (50/60 Hz)	Same
Dynamic Range	±5 mV	Same
Baseline Correction	±300 mV	Same
Stimulator (integrated into device)		
Isolated Stimulus Channels	2	Same
Pulse Amplitude		
Range	0 - 25.5 mA into 1000 Ω load	Same
Increment	0.1 ms	Same
Accuracy	±0.15 ms	Same
Pulse Duration		
Range	0.1 - 9.9 ms	Same
Increment	0.1 ms	Same
Accuracy	±0.15 ms	Same
Inter-Stimulus Interval (ISI)		
Range	10 - 9999 ms	Same
Range (Burst)	10 - 9999 ms	Same
Increment	10 ms	Same
Sequential Delay (AV)		
Range	11 - 250 ms	Same
Increment	1 ms	Same
Programmed Protocols	Preprogrammed protocols: BASIC 1 (induction) BASIC 2 (termination) BASIC 3 (backup) ACUTE Multi-Sx	Same

Device Characteristic	Predicate Device EP-TRACER V2.0	New Device EP-TRACER V2.2
Programmed Protocols	Pace Automatic mode Wenkebach mode User defined protocols	Same
Programmable Protocol Key	10	Same
Number of Extra-Stimuli	5 (S2 - S6)	Same
Sensing (ECG Synchronization)		
Automatic or Manual Trigger Setting – Sensitivity	Internal from any surface or intra-cardiac channel	Same
Automatic or Manual Trigger Setting – Trigger lockup (refractory time)	5 - 5000 ms	Same
Automatic or Manual Trigger Setting – ECG Delay	5 - 5000 ms	Same
Additional Outputs	No	Same
Power Source	Integrated with amplifier 38 channels: +5 V, 0.3 A & +12 V, 0.9 A 70 channels: +5 V, 0.3 A & +12 V, 1.5 A 102 channels: +5 V, 0.3 A & +12 V, 2 A	Same
Pacing Channels		
Isolated Channels	(i) atrial and (ii) ventricular and (iii) emergency fixed pace output to atrium and ventricle	Same
Circuit Isolation	Compliant with IEC 60601-1, Type CF, 5 kV, common & differential mode	Same
Computer Controlled Stimulus Pulses		
Current	0 - 25.5 mA into 1000 Ω load	Same
Current Steps	0.1 mA	Same
Accuracy	± 0.1 mA	Same
Pulse duration	Pulse width 0.1 - 9.9 ms, steps of 0.1 ms	Same
Accuracy	± 0.1 mA	Same
Load Impedance	1000 Ω	Same
Max. Output Voltage	25 V	Same
Inter-Stimulus Intervals		
S1 Range	10 - 9999 ms	Same
Stability	Quartz computer clock, ± 30 parts per million at +25°C	Same
Extra-Stimuli	5 (S2 - S6)	Same
Coupling Interval	30 - 9990 ms	Same
Accuracy	± 10 ms	Same
Protocol Automation		
Auto decrement/increment	Yes	Same

Device Characteristic	Predicate Device EP-TRACER V2.0	New Device EP-TRACER V2.2
Backup Manually Controlled Stimulation	Use external backup stimulator	Same
Emergency Backup Pacing	Use external backup stimulator	Same
Compliance with Standards		
Standards	EN 60601-1 EN 60601-1-2	Same
Device Directive	European Union Medical Device Directive (CE Marked)	Same
Environmental/ Electrical Specifications		
Operating Temperature	+10°C to +30°C	Same
Storage Temperature	-29°C to +66°C	Same
Operating Humidity	20 - 80 % rH (non-condensing)	Same
Storage Humidity	< 95 % rH (non-condensing)	Same
Max. Current Draw	15 A/115 V, 7 A/230 V	Same
Chassis Leakage Current	< 100 µA	Same
Advanced Features		
Display Ablation Parameters	Connection to RF ablation generator(s)	Same
Certification		
MDD Device Class	Class IIb	Same
IEC 60601-1	Certified	Same
EMC Compliance	Certified	Same
CE Marking	Certified CE 0197	Same
US Regulations	K161245	Current submission
System Dimensions		
W x D x H	Customer Option Mobile Cart: 100x80x170 cm Mobile Desk: 125x80x150 cm Control Desk: customer defined Laptop: customer defined	Same

Continued on next page

510(k) Summary per 21CFR §807.92, Continued

Software – user interface convenience features added The EP TRACER V2.2 included software modifications affecting the user interface including improvements in the display and organization of the system’s outputs. There were no modifications to the hardware of the new device.

The table below shows the software differences between the EP-TRACER V2.2 vs. the predicate device EP-TRACER V2.0.

Device Characteristic	Predicate Device EP-TRACER V2.0	New Device EP-TRACER V2.2
User Interface Convenience Features		
Pressure mode display	Pressure data is displayed in a full-screen window which is separate from the main trace display window.	In addition to the full-screen pressure mode available in V2.0, pressure channels can be displayed in a “minimized” pressure mode display within the main trace display window.
Calipers	Horizontal calipers can be displayed	Both horizontal and vertical calipers can be displayed
One-click report generation	Enables transfer of procedural data in a limited format	Enables transfer of procedural data to a customizable, preconfigured template.
Clock	Real-time clock	Added a procedure clock, displaying procedure time.
New Icons		
Split-screens	Executed with keystrokes	Icon added for existing feature
Pressure mode	Only one pressure mode display	Toggles between full-screen and minimized pressure mode display
Stimulator power settings	Executed with keystrokes	Icon added for existing feature
Add comment	Executed with keystrokes	Icon added for existing feature
Add mark	Executed with keystrokes	Icon added for existing feature
Take screenshot	Executed with keystrokes	Icon added for existing feature
Toggle stimulator protocol/parameters windows	Executed with keystrokes	Icon added for existing feature
Software compatibility		
Operating system	Windows 7 (32 bit)	Windows 7 (64 bit) and Windows 7 (32 bit)

Continued on next page

510(k) Summary per 21CFR §807.92, Continued

Performance testing

Completion of all verification and validation activities demonstrated that the subject devices meet their predetermined design and performance specifications. Verification activities performed confirmed that the modifications to the predicate device did not adversely affect the safety and effectiveness of the modified device.

Testing for the EP-Tracer V 2.2 shows conformance to the following harmonized standards:

- EN ISO 13485:2012
- EN ISO 14971:2012
- EN 60601-1:2006+A1: 2013
- EN 60601-1-2:2007+ AC: 2010
- EN 60601-1-6:2010
- EN 60601-2-27:2006
- EN 60601-2-34:2011
- EN 10993-1:2009
- EN 980:2008
- EN 1041:2008
- EN 62366: 2008
- EN 62304: 2006

Conclusion

Both the subject and predicate devices have the same intended use, indications for use, and operate using the same fundamental scientific technology.

Based on the technical characteristics and the results of the performance testing, the EP TRACER V2.2 is substantially equivalent to V2.0.
