



January 3, 2018

STERIS Corporation  
Jennifer Nalepka  
Senior Regulatory Affairs Specialist  
5960 Heisley Road  
Mentor, Ohio 44060

Re: K183297

Trade/Device Name: Vis-U-All Low Temperature Sterilization Pouches/Tubing  
Regulation Number: 21 CFR 880.6850  
Regulation Name: Sterilization Wrap  
Regulatory Class: Class II  
Product Code: FRG  
Dated: November 26, 2018  
Received: November 27, 2018

Dear Jennifer Nalepka:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

  
**Elizabeth F. Claverie -S**

For Tina Kiang, Ph.D.  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K183297

Device Name

Vis-U-All Low Temperature Sterilization Pouches/Tubing

Indications for Use (Describe)

The Vis-U-All Low Temperature Sterilization Pouches/Tubing are sterilization containment pouches for use by health care providers to enclose:

- medical devices in a single or double pouch configuration
- trays\* containing medical devices in a single or double pouch configuration
- small items requiring surface sterilization in a single pouch configuration within a tray\*

to be sterilized in the Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO® Low Temperature Sterilization Systems. The pouches maintain the sterility of the enclosed devices until used.

When used to enclose medical devices, the pouches are intended to contain the devices in such a manner as to leave a minimum of one inch between the devices and seal on all sides. When used to enclose a tray, the tray must fit loosely within the pouch.

Intended Sterilization Cycles and Intended Pouch Loads when Medical Devices are:

- Directly pouched
- Placed inside of a tray\* and the tray pouched\*

V-PRO 60 & s2 Lumen Cycle

- Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors
- Non-lumened devices including non-lumened rigid and semi-rigid endoscopes
- Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations:
  - single or dual lumen devices
    - $\geq 0.77$  mm internal diameter (ID) and  $\leq 410$  mm in length
  - triple lumen devices
    - $\geq 1.2$  mm ID and  $\leq 275$  mm in length
    - $\geq 1.8$  mm ID and  $\leq 310$  mm in length
    - or
    - $\geq 2.8$  mm ID and  $\leq 317$  mm in length

V-PRO 60 & s2 Non Lumen Cycle

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.

V-PRO 60 & s2 Flexible Cycle

Load 1: One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:

- single or dual lumen device with lumens that are  $\geq 1$  mm ID and  $\leq 990$  mm in length

Load 2: Non-lumened devices including non-lumened rigid, semi-rigid, and flexible endoscopes and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes, with the following dimensions:

- $\geq 2$  mm ID and  $\leq 400$  mm in length
- $\geq 0.76$  mm ID and  $\leq 233$  mm in length
- $\geq 1.0$  mm ID and  $\leq 254$  mm in length

---

#### V-PRO s2 Fast Cycle

- Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes, and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors.
- Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations:
  - Single or dual channeled devices with stainless steel lumens that are  $\geq 0.77$  mm ID and  $\leq 410$  mm in length
  - Triple channeled devices with stainless steel lumens that are either:
    - $\geq 1.2$  mm ID and  $\leq 275$  mm in length
    - $\geq 1.8$  mm ID and  $\leq 310$  mm in length
    - or
    - $\geq 2.8$  mm ID and  $\leq 317$  mm in length

#### V-PRO 1, 1 Plus, maX & maX 2 Lumen Cycle

- Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors
- Non-lumened devices including non-lumened rigid and semi-rigid endoscopes
- Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations:
  - single or dual lumen devices:
    - $\geq 0.77$  mm internal diameter (ID) and  $\leq 410$  mm in length
  - triple lumen devices:
    - $\geq 1.2$  mm ID and  $\leq 275$  mm in length
    - $\geq 1.8$  mm ID and  $\leq 310$  mm in length
    - or
    - $\geq 2.8$  mm ID and  $\leq 317$  mm in length

#### V-PRO 1, 1 Plus, maX & maX 2 Non Lumen Cycle

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors

#### V-PRO 1, 1 Plus, maX & maX 2 Flexible Cycle

Load 1: Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load.

The flexible endoscopes may contain either:

- a single lumen that is  $\geq 1$  mm ID and  $\leq 1050$  mm in length
- or two lumens with:
  - one lumen that is  $\geq 1$  mm ID and  $< 990$  mm in length
  - and the other lumen that is  $\geq 1$  mm ID and  $\leq 850$  mm in length

Load 2: Non-lumened instruments including instruments with diffusion-restricted areas such as the hinged portion of forceps or scissors.

#### V-PRO maX 2 Fast Non Lumen Cycle

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.

\*Trays must be legally marketed for use in the V-PRO Low Temperature Sterilization Systems and contain a vent surface area to tray volume ratio  $\geq 0.135$  in-1 with the maximum number of instrument organizers installed.

---

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

---

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

---

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
*PRASStaff@fda.hhs.gov*

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*



**510(k) Summary**  
**for**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

**Sponsor Facility**

STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060  
Phone: (440) 354-2600  
Fax No: (440) 357-9198

Contact: Jennifer Nalepka  
Senior Regulatory Affairs Specialist

Telephone: (440) 392-7458  
Fax No: (440) 357-9198  
e-mail: [jennifer\\_nalepka@steris.com](mailto:jennifer_nalepka@steris.com)

Submission Date: December 30, 2018

**1. Device Name**

Trade Name: Vis-U-All Low Temperature Sterilization Pouches/Tubing

Device Classification: Class II

Common/Usual Name: Sterilization pouch

Classification Name: Sterilization wrap

Classification Number: 21 CFR 880.6850

Product Code: FRG

**2. Predicate Device**

Vis-U-All Low Temperature Sterilization Pouches/Tubing (K172749)

**3. Description of Device**

The proposed Vis-U-All Low Temperature Sterilization Pouches/Tubing is identical to the predicate and is a Tyvek/plastic film sterilization containment pouch designed for devices to be sterilized by the health care provider in V-PRO Low Temperature Sterilization Systems. As is the predicate device, the proposed device is available as a self seal pouch, a heat seal pouch, or heat seal tubing. Available sizes and configurations are shown in **Table 5-1**.

**Table 5-1.** Sizes and Configurations of Vis-U-All Low Temperature Sterilization Pouches/Tubing

Type	Size (inches unless specified)
Heat Seal Pouch	3 x 7
	4 x 9
	4 x 12
	4 x 22
	6 x 10
	8 x 12
	10 x 15
	12 x 18
Self Seal Pouch	3 x 7
	4 x 9
	4 x 12
	4 x 22
	6 x 10
	8 x 12

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

---

	10 x 15
	12 x 18
	8 x 21*
	8 x 27*
	9 x 27*
	11 x 22*
	12 x 27*
Tubing	3" x 100'
	4" x 100'
	6" x 100'
	9" x 100'
	14" x 100'

\* Additional pouch sizes

The purpose of this submission is to add claims for:

- Use with devices having stainless steel lumens of the following configurations:
  - $\geq 2.0$  mm internal diameter (ID) and  $\leq 400$  mm in length
  - $\geq 1.0$  mm ID and  $\leq 254$  in length
  - $\geq 0.76$  mm ID and  $\leq 233$  mm in length
 in the V-PRO 60 and V-PRO s2 Sterilizer Flexible Cycle
- Use with diffusion restricted spaces in addition to stainless steel and titanium when used in the V-PRO 60 and V-PRO s2 Sterilizer Non Lumen and Flexible Cycles.
- Use in the V-PRO s2 Sterilizer Fast Lumen Cycle
- Addition of the V-PRO s2 Sterilizer to the indications for use in the Lumen Cycle of the V-PRO 60 Sterilizer. The V-PRO s2 Sterilizer Lumen Cycle is performed identically to the previously cleared cycle listed in **Table 5-2**. The updated claims being sought for the V-PRO 60/s2 Sterilizer Non Lumen and Flexible Cycles as well as the V-PRO s2 Sterilizer Fast Cycle are the subject of a concurrent submission (K182528).
- Notify FDA of an extension of the self seal pouch line to include five additional pouch models whose sizes were optimized to contain the smaller PRO-LITE Sterilization Trays that has recently been implemented.



**Table 5-2: V-PRO Cycles Provided V-PRO 60 and s2 Sterilizers**

Cycle	Most Recent Premarket Notification Clearance	V-PRO 60	V-PRO s2
Lumen Cycle	K172319	X	X
Non Lumen Cycle	K172319 and Concurrent submission (K182827)	X	X
Flexible Cycle		X	X
Fast Cycle	Concurrent Submission (K182827)		X

\* Shaded cell indicates that the V-PRO Sterilizer does not include the identified cycle.

**NOTE: V-PRO Sterilizer cycles are pre-set and thus the name of the cycle describes the cycle conditions.**

**4. Intended Use/ Indications for Use**

The Vis-U-All Low Temperature Sterilization Pouches/Tubing are sterilization containment pouches for use by health care providers to enclose:

- medical devices in a single or double pouch configuration
- trays\* containing medical devices in a single or double pouch configuration
- small items requiring surface sterilization in a single pouch configuration within a tray\*

to be sterilized in the Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO® Low Temperature Sterilization Systems. The pouches maintain the sterility of the enclosed devices until used.

When used to enclose medical devices, the pouches are intended to contain the devices in such a manner as to leave a minimum of one inch between the devices and seal on all sides. When used to enclose a tray, the tray must fit loosely within the pouch.

Intended Sterilization Cycles	Intended Pouch Load when Medical Devices are: <ul style="list-style-type: none"> <li>• Directly Pouched or</li> <li>• Placed Inside of a Tray* and the Tray* Pouched</li> </ul>
V-PRO 60 & s2 Lumen Cycle	<ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm internal diameter (ID) and <math>\leq 410</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> <li>or</li> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul>

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

Intended Sterilization Cycles	Intended Pouch Load when Medical Devices are: <ul style="list-style-type: none"> <li>• Directly Pouched or</li> <li>• Placed Inside of a Tray* and the Tray* Pouched</li> </ul>
V-PRO 60 & s2 Non Lumen Cycle	Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.
V-PRO 60 & s2 Flexible Cycle	<p><u>Load 1:</u> One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:</p> <ul style="list-style-type: none"> <li>○ single or dual lumen device with lumens that are <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> </ul> <p><u>Load 2:</u> Non-lumened devices including non-lumened rigid, semi-rigid, and flexible endoscopes and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes, with the following dimensions:</p> <ul style="list-style-type: none"> <li>○ <math>\geq 2</math> mm ID and <math>\leq 400</math> mm in length</li> <li>○ <math>\geq 0.76</math> mm ID and <math>\leq 233</math> mm in length</li> <li>○ <math>\geq 1.0</math> mm ID and <math>\leq 254</math> mm in length</li> </ul>
V-PRO s2 Fast Cycle	<ul style="list-style-type: none"> <li>• Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes, and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors.</li> <li>• Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations: <ul style="list-style-type: none"> <li>• Single or dual channeled devices with stainless steel lumens that are <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> <li>• Triple channeled devices with stainless steel lumens that are either: <ul style="list-style-type: none"> <li><math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li><math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> or <ul style="list-style-type: none"> <li><math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul>
V-PRO 1, 1 Plus, maX & maX 2 Lumen Cycle	<ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 527</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> or <ul style="list-style-type: none"> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul>
V-PRO 1 Plus, maX & maX 2 Non Lumen Cycle	Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.
V-PRO maX & maX 2 Flexible Cycle	<p><u>Load 1:</u> Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load. The flexible endoscopes may contain either:</p>

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

Intended Sterilization Cycles	Intended Pouch Load when Medical Devices are: <ul style="list-style-type: none"> <li>• Directly Pouched or</li> <li>• Placed Inside of a Tray* and the Tray* Pouched</li> </ul>
	<ul style="list-style-type: none"> <li>• a single lumen that is <math>\geq 1</math> mm ID and <math>\leq 1050</math> mm in length</li> <li>• or two lumens with: <ul style="list-style-type: none"> <li>• one lumen that is <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> <li>• and the other lumen that is <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>Load 2:</u> Non-lumened instruments including instruments with diffusion-restricted areas such as the hinged portion of forceps or scissors.</p>
V-PRO maX 2 Fast Non Lumen Cycle	Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.

\* Trays must be legally marketed for use in the V-PRO Low Temperature Sterilization Systems and contain a vent surface area to tray volume ratio  $\geq 0.135$  in<sup>-1</sup> with the maximum number of instrument organizers installed.

**5. Description of Safety and Substantial Equivalence**

The proposed and predicate devices are single use sterilization pouches for use in V-PRO Sterilizers. **Table 5-3** summarizes the difference between the proposed device and predicate device cleared under K172749.

**Table 5-3.** Technical Comparison to the Predicate.

Characteristic	Proposed	Predicate K172749	Comparison
Intended Use	<p>The Vis-U-All Low Temperature Sterilization Pouches/Tubing are sterilization containment pouches for use by health care providers to enclose:</p> <ul style="list-style-type: none"> <li>• medical devices in a single or double pouch configuration</li> <li>• trays* containing medical devices in a single or double pouch configuration</li> <li>• small items requiring surface sterilization in a single pouch configuration within a tray* to be sterilized in the:</li> <li>• Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO® Low Temperature Sterilization Systems</li> </ul>	<p>The Vis-U-All Low Temperature Sterilization Pouches/Tubing are sterilization containment pouches for use by health care providers to enclose:</p> <ul style="list-style-type: none"> <li>• medical devices in a single or double pouch configuration</li> <li>• trays* containing medical devices in a single or double pouch configuration</li> <li>• small items requiring surface sterilization in a single pouch configuration within a tray* to be sterilized in the Lumen, Non Lumen, Flexible and Fast Non Lumen Cycles of the V-PRO Low Temperature Sterilization Systems.</li> </ul>	<p>Similar</p> <p>Data is presented in this submission to support safe and effective use of the pouches in the additional cycles and load configurations.</p> <p>Combined V-PRO 1, 1 Plus and maX Sterilizer Lumen, Non Lumen and Flexible Cycle claims with the V-PRO maX 2 Lumen, Non Lumen and</p>

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

	<p>The pouches maintain the sterility of the enclosed devices until used.</p> <p>When used to enclose medical devices, the pouches are intended to contain the devices in such a manner as to leave a minimum of one inch between the devices and seal on all sides. When used to enclose a tray, the tray must fit loosely within the pouch.</p> <p><u>V-PRO 60 &amp; s2 Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm internal diameter (ID) and <math>\leq 410</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> <li>or</li> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul> <p><u>V-PRO 60 &amp; s2 Non Lumen Cycle</u></p> <p>Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-</p>	<p>The pouches maintain the sterility of the enclosed devices until used.</p> <p>When used to enclose medical devices, the pouches are intended to contain the devices in such a manner as to leave a minimum of one inch between the devices and seal on all sides. When used to enclose a tray, the tray must fit loosely within the pouch.</p> <p><u>V-PRO 60 Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> <li>or</li> <li>• <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul> <p><u>V-PRO 60 Non Lumen Cycle</u></p> <p>Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted</p>	<p>Flexible Cycle claims.</p>
--	--	--	-------------------------------

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

	<p>restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO 60 &amp; s2 Flexible Cycle</u>  <u>Load 1:</u> One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:</p> <ul style="list-style-type: none"> <li>• single or dual lumen device with lumens that are <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> </ul> <p><u>Load 2:</u> Non-lumened devices including non-lumened rigid, semi-rigid, and flexible endoscopes and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes, with the following dimensions:</p> <ul style="list-style-type: none"> <li>• <math>\geq 2</math> mm ID and <math>\leq 400</math> mm in length</li> <li>• <math>\geq 0.76</math> mm ID and <math>\leq 233</math> mm in length</li> <li>• <math>\geq 1.0</math> mm ID and <math>\leq 254</math> mm in length</li> </ul> <p><u>V-PRO s2 Fast Cycle</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes, and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors.</li> <li>• Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations:</li> </ul>	<p>spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO 60 Flexible Cycle</u>  One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:</p> <ul style="list-style-type: none"> <li>• single or dual lumen device with lumens that are <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> </ul>	
--	---	---	--

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

	<ul style="list-style-type: none"> <li>• Single or dual channeled devices with stainless steel lumens that are <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> <li>• Triple channeled devices with stainless steel lumens that are either:  <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length  <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length  or</li> <li>• <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> <p><u>V-PRO 1, 1 Plus, maX &amp; maX 2 Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 527</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length or</li> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul> <p><u>V-PRO 1, 1 Plus, maX &amp; maX2 Non Lumen Cycle</u>  Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened</p>	<p><u>V-PRO 1, 1 Plus &amp; maX Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 527</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length or</li> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul> <p><u>V-PRO 1 Plus &amp; maX Non Lumen Cycle</u>  Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened</p>	
--	---	--	--

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

	<p>devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO maX and maX 2 Flexible Cycle</u>  <u>Load 1:</u> Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load.  The flexible endoscopes may contain either:</p> <ul style="list-style-type: none"> <li>• a single lumen that is <math>\geq 1</math> mm ID and <math>\leq 1050</math> mm in length</li> <li>• or two lumens with: <ul style="list-style-type: none"> <li>• one lumen that is <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> <li>• and the other lumen that is <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>Load 2:</u> Non-lumened instruments including instruments with diffusion-restricted areas such as the hinged portion of forceps or scissors.</p> <p><u>V-PRO maX 2 Fast Non Lumen Cycle</u>  Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p>*Trays must be legally marketed for use in the V-PRO Low Temperature Sterilization Systems and</p>	<p>devices with stainless steel diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO maX Flexible Cycle</u>  <u>Load 1:</u> Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load.  The flexible endoscopes may contain either:</p> <ul style="list-style-type: none"> <li>• a single lumen that is <math>\geq 1</math> mm ID and <math>\leq 1050</math> mm in length</li> <li>• or two lumens with: <ul style="list-style-type: none"> <li>• one lumen that is <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> <li>• and the other lumen that is <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>Load 2:</u> Non-lumened instruments including instruments with diffusion-restricted areas such as the hinged portion of forceps or scissors.</p> <p><u>V-PRO maX 2 Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations:</li> </ul>	
--	--	--	--

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

	<p>contain a vent surface area to tray volume ratio <math>\geq 0.135 \text{ in}^{-1}</math> with the maximum number of instrument organizers installed.</p>	<ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77 \text{ mm ID}</math> and <math>\leq 527 \text{ mm}</math> in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2 \text{ mm ID}</math> and <math>\leq 275 \text{ mm}</math> in length</li> <li>▪ <math>\geq 1.8 \text{ mm ID}</math> and <math>\leq 310 \text{ mm}</math> in length</li> <li>or</li> <li>▪ <math>\geq 2.8 \text{ mm ID}</math> and <math>\leq 317 \text{ mm}</math> in length</li> </ul> </li> </ul> <p><u>V-PRO maX 2 Non Lumen Cycle</u>  Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO maX 2 Flexible Cycle</u>  <u>Load 1:</u> Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load.  The flexible endoscopes may contain either:</p> <ul style="list-style-type: none"> <li>• a single lumen that is <math>\geq 1 \text{ mm ID}</math> and <math>\leq 1050 \text{ mm}</math> in length</li> <li>• or two lumens with: <ul style="list-style-type: none"> <li>• one lumen that is <math>\geq 1 \text{ mm ID}</math> and <math>\leq 990 \text{ mm}</math> in length</li> <li>• and the other lumen that is <math>\geq 1 \text{ mm ID}</math> and <math>\leq 850 \text{ mm}</math> in length</li> </ul> </li> </ul> <p><u>Load 2:</u> Non-lumened instruments including instruments with diffusion-</p>	
--	---	--	--



**STERIS Traditional 510(k) PREMARKET NOTIFICATION  
Vis-U-All Low Temperature Sterilization Pouches/Tubing**

		<p>restricted areas such as the hinged portion of forceps or scissors.</p> <p><u>V-PRO maX 2 Fast Non Lumen Cycle</u> Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p>* Trays must be legally marketed for use in the V-PRO Low Temperature Sterilization Systems and contain a vent surface area to tray volume ratio <math>\geq 0.135 \text{ in}^{-1}</math> with the maximum number of instrument organizers installed.</p>	
Device Features	<ul style="list-style-type: none"> <li>▪ Chevron end of pouches for ease of opening</li> <li>▪ Chemical process indicator for EO</li> </ul>	<ul style="list-style-type: none"> <li>▪ Chevron end of pouches for ease of opening</li> <li>▪ Chemical process indicator for EO</li> </ul>	Same
Maintenance of Sterility	1 year	1 year	Same
Materials of Construction	Tyvek and plastic	Tyvek and plastic	Same
Types	Self Seal, Heat Seal, Tubing	Self Seal, Heat Seal, Tubing	Same
Chemical Indicator	Ethylene Oxide Process Chemical Indicator Printed on both sides of Tyvek	Ethylene Oxide Process Chemical Indicator Printed on both sides of Tyvek	Same

**Table 5-4** summarizes the testing of the Vis-U-All Low Temperature Sterilization Pouches/Tubing to demonstrate that the proposed pouch is qualified for use in the V-PRO s2 Fast Cycle and for extended sterilization claims in the V-PRO 60 and s2 Sterilizers Non Lumen and Flexible Cycles and is as safe, as effective, and performs the same as the predicate device.

**STERIS Traditional 510(k) PREMARKET NOTIFICATION**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

**Table 5-4. Performance Test Summary**

Test	Acceptance Criteria	Conclusion	
<u>Effective Sterilant Penetration into Pouches</u> (including pouched trays and, if applicable, pouches placed within a tray):	V-PRO s2 Sterilizer Fast Cycle	Worst case stainless steel lumen test articles shall be reproducibly sterilized under worst case ½ cycle conditions for the Fast Cycle.	PASS
	V-PRO s2 Sterilizer Fast Cycle	Diffusion-restricted spaces shall be reproducibly sterilized under worst case ½ cycle conditions for Fast Cycle.	PASS
	V-PRO s2 Sterilizer Fast Cycle	Contact points between medical devices and tray accessories shall be reproducibly sterilized under worst case ½ cycle conditions for the Fast Cycle.	PASS
	V-PRO s2 Sterilizer Fast Cycle	Worst case surface sterilization test articles shall be reproducibly sterilized under worst case ½ cycle conditions for the Fast Cycle.	PASS
	V-PRO 60 & s2 Sterilizer Non Lumen Cycle, diffusion-restricted spaces claims	Diffusion-restricted spaces shall be reproducibly sterilized under worst case ½ cycle conditions for the Non Lumen Cycle.	PASS
	V-PRO 60 & s2 Sterilizer Flexible Cycle, diffusion-restricted spaces claims	Diffusion-restricted spaces shall be reproducibly sterilized under worst case ½ cycle conditions for the Flexible Cycle	PASS
V-PRO 60 & s2 Sterilizer Flexible Cycle extended claims	Stainless steel test articles shall be reproducibly sterilized under worst case ½ cycle conditions for the Flexible Cycle.	PASS	
<u>Pouch Integrity: Physical and Microbial Barrier Properties</u>	Tensile Strength	Pouch material tensile strength will show no statistical difference between processed and unprocessed samples.	PASS
	Whole Package Integrity (Burst)	Pouch burst strength will show no statistical difference between processed and unprocessed pouches.	PASS
	Seal Strength	Pouch seal strength will show no statistical difference between processed and unprocessed pouches.	PASS
	Microbial Retention	Tyvek microbial retention will show no statistical difference between processed and unprocessed pouches.	PASS
Maintenance of Package Integrity	Packaged instruments shall remain sterile through event related and real time studies.	PASS	
Aeration: Hydrogen Peroxide Residuals	Hydrogen peroxide residuals on the pouch will be reduced to acceptable levels for dermal contact in Fast Non Lumen Cycle with 3 minute aeration.	PASS	

**8. Conclusion**

Based on the intended use, technological characteristics and non-clinical performance data, the subject device is as safe, as effective and performs at least as well as the legally marketed predicate device (K172749), Class II (21 CFR 880.6850), product code FRG.