



Shenzhen Leading Perfection Technology Co., Ltd
% Rain Yip
Registered Engineer
Feiying Drug & Medical Consulting Technical Service Group
Rm. 3005, Area B, Bldg.1, Southward Ruifeng Business Center,
Guimiao Road
Shenzhen, Guangdong 518000 CN

Re: K183674

Trade/Device Name: Electronic Pulse Stimulator, Model S3
Regulation Number: 21 CFR 882.5890
Regulation Name: Transcutaneous Electrical Nerve Stimulator For Pain Relief
Regulatory Class: Class II
Product Code: NUH. NGX
Dated: March 30, 2019
Received: September 9, 2019

Dear Rain Yip:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Vivek Pinto, PhD
Director (Acting)
DHT5B: Division of Neuromodulation
and Physical Medicine Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K183674

Device Name

Electronic Pulse Stimulator (Model: S3)

Indications for Use (Describe)

Mode 1 (PMS):

To be used for stimulating healthy muscles in order to improve and facilitate muscle performance.

Mode 2 (TENS):

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities (arm), lower extremities (leg), abdomen and bottom due to strain from exercise or normal household work activities

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

"510(K) Summary" as required by 21CFR 807.92.

Date: 2019-01-25

I. Submitter

Shenzhen Leading Perfection Technology Co., Ltd
4th Floor, J Building, San Yi Dui Industrial Park, Zhoushi Road, Jiuwei, Hangcheng Street,
Bao'an District, Shenzhen, China

Tel.: +86 755 2999 1449

Tax: +86 755 2371 3490

Simon Hou (General Manager)

Tel: +86 134 8075 8801

Email: info@leadingind.com

II. Device

Type of 510(k): Traditional

Common Name: Powered muscle stimulator

Trade Name: Electronic Pulse Stimulator

Model: S3

Classification Name: Transcutaneous electrical nerve stimulator for pain relief

Review Panel: Neurology

Regulatory Class: II

Product Code: NUH, NGX

Regulation Number: 21 CFR 882.5890

III. Predicate Device

<u>Applicant</u>	<u>Predicate Device</u>	<u>510(k) Number</u>	<u>Approval Date</u>
Philips Consumer Lifestyle	<u>(Primary):</u> PulseRelief	K151035	July 21,2015
HIVOX BIOTEK INC.	HIVOX OTC Electrical Stimulator /SEM44	K171803	November 29, 2017
Shenzhen OSTO Technology Company Limited	Health Expert Electronic Stimulator, Model AST-300C and AST-300D	K133929	November 12, 2014

IV. Device Description

Electronic Pulse Stimulator is a product that adopts modern electronic science and technology to delivers electric pulses generated to the user's skin through the electrodes. It has two modes of different pulse frequencies, mode 1 uses Electrical Muscle Stimulation (EMS) technology for muscle training and mode 2 uses Transcutaneous Electrical Nerve Stimulation (TENS) technology to temporary relieve pain. The device includes operating elements of Power ON/OFF button, intensity increase"+" button, and intensity decrease"-" button and Mode selection button, and could be attached and detached to the electrode pad through the snap-on connector.

In additional, according to different simulation-needed bodies, users can choose the four types of electrode pad based on their own situation. Electronic Pulse Stimulator is suitable for the simulation-needed bodies including hips, arms and legs, abdomen, and feet.

V. Indications for Use

Mode 1 (PMS): To be used for stimulating healthy muscles in order to improve and facilitate muscle performance.

Mode 2 (TENS): To be used for temporary relief of pain associate with sore and aching muscles in the shoulder, waist, back, neck, upper extremities (arm), lower extremities (leg), abdomen and bottom due to strain from exercise or normal household work activities.

VI. Comparison of Technological Characteristics With the Predicate Devices

The Electronic Pulse Stimulator is substantially equivalent to the predicate devices based on intended use, design, specifications and performance.

The Electronic Pulse Stimulator raises no safety or efficacy concerns when compared to the following predicate devices.

Information for predicate device was obtained from publicly available sources, including the 510(k) Summary and device instruction manual. A technical comparison to the predicate is provided below:

<u>Comparison Elements</u>	<u>Subject Device</u>	<u>(Primary) Predicate Device K151035</u>	<u>Predicate Device K171803</u>	<u>Predicate Device K133929</u>	<u>Comparison</u>
Trade Name/Model	Electronic Pulse Stimulator/S3	PulseRelief	HIVOX OTC Electrical Stimulator/SEM44	Health Expert Electronic Stimulator, Model AST-300C and AST-300D	/
510(k) Number	To be assigned	K151035	K171803	K133929	/
Manufacturer	Shenzhen Leading Perfection Technology Co., Ltd	Philips Consumer Lifsty	HIVOX BIOTEK INC.	Shenzhen OSTO Technology Company Limited	/
Regulation Number	21CFR 882.5890	21CFR 882.5890 and 890.5850	21CFR 882.5890	21 CFR 882.5890	SE
Device Class	II	II	II	II	SE
Product Code	NUH, NGX	NUH, NGX	NUH, NGX	NUH, NGX	SE <u>NOTE 1</u>
Indication for Use/Intended Use	Mode 1(PMS): To be used for stimulating healthy muscles in order to improve and facilitate muscle performance. Mode 2(TENS): To be used for temporary relief of pain associate with sore and aching muscles in the shoulder, waist, back, neck, upper extremities (arm), lower	The OTC TENS/EMS stimulator PulseRelief is designed to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities(arm) and	SEM44 (EMS): The device is designed to be used for stimulate healthy muscles in order to improve and facilitate muscle performance. SEM44 (TENS): The device is designed to be used for temporary relief	PMS (Mode 1~8) It is intended to stimulate healthy muscles in order to improve and facilitate muscle performance. TENS (Mode 9~25) To be used for temporary relief of pain associated with sore and aching muscles in	SE <u>NOTE 1</u>

	extremities (leg), abdomen and bottom due to strain from exercise or normal household work activities				
--	--	--	--	--	--

<u>Comparison Elements</u>	<u>Subject Device</u>	<u>(Primary) Predicate Device K151035</u>	<u>Predicate Device K171803</u>	<u>Predicate Device K133929</u>	<u>Comparison</u>
		lower extremities (leg) due to strain from exercise or normal household work activities. It should be applied to normal, healthy, dry and clean skin of adult patients, and is to be used for stimulate healthy muscles in order to improve and facilitate muscle performance.	of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities (arm), lower extremities (leg), abdomen and bottom due to strain from exercise or normal household work activities.	the shoulder, waist, back, back of the neck, arm, leg, and foot due to strain from exercise or normal household work activities by applying current to stimulate nerve.	
Shape of electrodes	Rectangle	Square	Unknown	Rectangle	<u>SE</u>
Anatomical application	Arm, leg, abdomen, hip, feet	Shoulder, waist, back, upper extremities (arm), and lower extremities (leg)	Shoulder, waist, back, neck, upper extremities (arm), lower extremities (leg), abdomen and bottom	Shoulder, waist, back, back of the neck, arm, leg, and foot	<u>SE</u>
Applicable People	Adults	Adults	Adults	Adults	SE
Location for Use	OTC	OTC	OTC	OTC	SE
Basic Unit Characteristics					

<u>Comparison Elements</u>		<u>Subject Device</u>	<u>(Primary) Predicate Device K151035</u>	<u>Predicate Device K171803</u>	<u>Predicate Device K133929</u>	<u>Comparison</u>
Power Supply		USB rechargeable battery, 3.7V	Li-ion 3.7V 500mAH	4.5V (batteries, 3×1.5V AAA)	Adaptor Input: 100-240Vac, 50-60Hz, 0.1A Output: 5Vdc, 1A Unit Input: 5Vdc, 1A	SE <u>NOTE 2</u>
Method of Line Current Isolation		N/A	N/A	N/A	Type BF Applied Part	SE
Patient Leakage Current	Normal Condition	< 0.1	<10μA	N/A	AC: 54.5μA, DC: 0.5μA	SE
	Signal Fault Condition	N/A	<50μA	N/A	AC:120.0μA, DC: 0.6μA	SE
Number of Output Modes		2 (1 TENS and 1 EMS)	EMS: 5 TENS:15	EMS: 35 TENS:15	EMS: 8 TENS:17	SE <u>NOTE 2</u>
Number of Output Channels		1	1	2	2	SE <u>NOTE 2</u>
Synchronous or Alternating?		N/A	N/A	Synchronous	Synchronous	SE
Method of Channel Isolation		N/A	N/A	By electrical circuit and software	Voltage Transform Isolation “BODY▼” and “ BODY▼” buttons for body channel, “ SOLE▲” and “SOLE▼” buttons for feet channel	SE

<u>Comparison Elements</u>		<u>Subject Device</u>	<u>(Primary) Predicate Device K151035</u>	<u>Predicate Device K171803</u>	<u>Predicate Device K133929</u>	<u>Comparison</u>
Software /Firmware /Microprocessor Control?		Yes	Yes	Yes	Yes	SE
Automatic Overload Trip?		Yes	Yes	Yes	No	SE
Automatic No-Load Trip?		Yes	Yes	Yes	No	SE
Automatic Shut OFF?		Yes	Yes	Yes	Yes	SE
Patient Override Control?		Yes	Yes	Yes	Yes	SE
Indicat or Display	On/Off Status?	Yes	Yes	Yes	Yes	SE
	Low Battery?	Yes	Yes	Yes	No	SE
	Voltage/ Current Level?	No	Yes	Yes	Yes	SE
Timer Range		25 minutes	1~59minutes and continuous	5~100minutes	25 minutes	SE <u>NOTE 2</u>
Compliance with Voluntary Standards?		IEC60601-1-2 IEC60601-1 IEC60601-11	IEC60601-1-2 IEC60601-1 IEC60601-2-10	IEC60601-1-2 IEC60601-1 IEC60601-2-10	IEC60601-1-2 IEC60601-1 IEC60601-2-10	SE

<u>Comparison Elements</u>		<u>Subject Device</u>	<u>(Primary) Predicate Device K151035</u>	<u>Predicate Device K171803</u>	<u>Predicate Device K133929</u>	<u>Comparison</u>
		IEC60601-2-10 ISO10993-5/10	ISO10993-5/10	ISO10993-5/10	ISO10993-5/10	
Compliance with 21CFR898		Yes	Yes	Yes	Yes	SE
Weight		14.5g(with battery)	62g (excl. electrodes)	89g (including belt clip, without batteries) 123g (including belt clip, and batteries)	2Kg (Without accessories)	SE <u>NOTE 3</u>
Dimensions		Φ4.2cm×1.38cm	54×54×14mm (unit)	132 x 63 x 29.5mm(including belt clip)	428mm x 428.8mm x 185mm	SE <u>NOTE 3</u>
Housing Materials and Construction (Main unit)		ABS	ABS, PC	ABS	ABS	SE
Output Specifications						
Waveform		Symmetrical rectangular	Biphasic, rectangular	Biphasic, square	Pulsed, symmetric, Biphasic, Rectangular,	SE <u>NOTE 4</u>
Maxim um Output Voltage (Volts,	@500Ω	40±20%	31±20%	100±10 % (50±10 %(Vp))	44±10%	SE <u>NOTE 4</u>
	@2kΩ	64±20%	69±20%	180±10 % (90±10 %(Vp))	80±10%	SE <u>NOTE 4</u>

Comparison Elements		Subject Device	(Primary) Predicate Device K151035	Predicate Device K171803	Predicate Device K133929	Comparison
Vp)	@10kΩ	76±20%	70±20%	250±10 % (125±10 %(Vp))	112±10%	SE <u>NOTE 4</u>
Maximum Output Current (mA)	@500Ω	80±20%	62±20%	200±10 % (100±10 %(Ip))	88±10%	SE <u>NOTE 4</u>
	@2kΩ	32±20%	34±20%	90±10 % (45±10 %(Ip))	40±10%	SE <u>NOTE 4</u>
	@10kΩ	7.6±20%	7±20%	25±10 % (12.5±10 %(Ip))	11.2±10%	SE <u>NOTE 4</u>
Pulse Width		200μs	TENS: 60~350μs EMS: 150~350μs	50-450μs	120μs	SE <u>NOTE 4</u>
Frequency (Hz)		Mode 1: 4-35Hz Mode 2: 4-41Hz	TENS: 1~100Hz EMS: 45~65Hz	1-150Hz	77.3Hz	SE <u>NOTE 4</u>
Maximum Phase Charge (μC) @500Ω		22.31	TENS: 1.6~6.8 EMS: 4.7~10.9	0.045	12.78μC @ 500Ω	SE <u>NOTE 4</u>
Maximum Current Density (mA/cm²) @500Ω		Hips pad: 0.079 Arms and legs pad: 0.104 Abdomen pad: 0.111 Feet pad: 0.016	TENS: 0.002~0.045 EMS: 0.019~0.037	0.667	0.235mA/cm² @ 500Ω	SE <u>NOTE 4</u>
Maximum Power Density (W/cm²) @500Ω		Hips pad: 0.0001 Arms and legs pad: 0.0002 Abdomen pad: 0.0002	(average) TENS: 0.24~1.69 EMS: 0.62~1.15	0.0046	1.38mW/cm² @ 500Ω	SE <u>NOTE 4</u>

<u>Comparison Elements</u>		<u>Subject Device</u>	<u>(Primary) Predicate Device K151035</u>	<u>Predicate Device K171803</u>	<u>Predicate Device K133929</u>	<u>Comparison</u>
		Feet pad: 0.00003				
Burst Mode	Pulses per burst	N/A	TENS: 5, 7 EMS: N/A	3	--	SE
	Bursts per second	N/A	TENS: 1, 2, 3 EMS: N/A	2/60Hz	--	SE
	Burst duration (seconds)	N/A	TENS: 62.5~87.5 EMS: N/A	36ms	--	SE
	Duty Cycle [Line(b) x Line (c)]	N/A	TENS: 6.3%~19% EMS: N/A	36ms/390ms	--	SE
ON Time (seconds)		1s	TENS: N/A EMS: --intensity ramp-up: 2~4sec --constant intensity: 5~6sec --intensity ramp-down: 1~2sec	2s	0.6s	SE
OFF Time (seconds)		4s	TENS: N/A EMS: N/A	2s	0.6s	SE

Comparison in details:

NOTE 1: The "Intended Use" of the subject device is the substantially equivalent to the predicate devices selected, all have adopted TENS (Transcutaneous Electrical Nerve Stimulation) and EMS (Electrical Muscle Stimulation) technologies to the treatment of temporary relief of pain associate with sore and aching muscles and/or the muscle performance improvement of healthy muscles respectively. Meanwhile, the

"Product Code" of the subject device is the same as the predicate devices. This define does not affect the intended use or normal use of the subject device.

NOTE 2: Although the "Power Supply", "Number of Output Modes", "Number of Output Channels", and "Timer Range" of the subject device are a little different from the predicate devices, they are all complaint with the requirements of IEC 60601-1-2 and IEC 60601-1 standard. So these differences will not raise any safety or effectiveness issue.

NOTE 3: Although the "Weight" and "Dimensions" of the subject device are different from the predicate devices, these differences are insignificant in the terms of safety or effectiveness.

NOTE 4: Although the "Waveform", "Maximum Output Voltage", "Maximum Output Current", "Pulse Duration", "Frequency", "Maximum Phase Charge", "Maximum Current Density", and "Maximum Power Density" of the subject device are a little different from the predicate devices, they all comply with the requirements of IEC 60601-1 and IEC 60601-2-10 standard, as well as the Guidance for Powered Muscle Stimulator for Muscle Conditioning. So these differences will not raise any safety or effectiveness issue.

CONCLUSION: The subject device Electronic Pulse Stimulator is substantially equivalent to the predicate devices.

VII. Performance Data

The following performance data were provided in support of the substantial equivalence determination.

1) Biocompatibility Testing

The direct components of applied parts contacting with the users in the Electronic Pulse Stimulator is as the following item:

Component Name	Material of Component	Body Contact Category	Contact Duration
Electrode pads used for hips, arms and legs, and abdomen (hips pad, arms and legs pad, and abdomen pad)	Conductive hydrogel	Surface-contacting skin	Max. 25 minutes, less than 24 hours
Electrode pads used for feet (feet pad)	PU fabric (contained sliver paste and conductive oil)	Surface-contacting skin	Max. 25 minutes, less than 24 hours

The electrode pads are directly purchased from qualified supplier which has been tested and passed the In Vitro Cytotoxicity, Skin Sensitization, and Skin Irritation test in accordance with the ISO10993-1 standards. So we have reason to believe that the electrode pads are safe for the users. For details, please refer to "Biocompatibility Discussion" in this application.

2) Electrical and EMC Safety

Electrical safety and EMC safety testing was performed to, and passed, the following standards:

- IEC 60601-1 Medical electrical equipment –Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-11 Medical electrical equipment –Part 1-11: General requirements for basic safety and essential performance –Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- IEC 60601-2-10 Medical electrical equipment -- Part 2-10: Particular requirements for the basic safety and essential performance of nerve and muscle stimulators
- IEC 60601-1-2 Medical electrical equipment –Part 1-2: General requirements for basic safety and essential performance –Collateral standard: electromagnetic compatibility – Requirements and tests

In addition to the compliance of voluntary standards:

- The software verification has been carries out according to the FDA Guidance for the Content of Premarket Submission for Software Contained in Medical Devices.
- The waveform test report has also been conducted to verify the output specifications of the subject device according to the FDA Guidance for Powered Muscle Stimulator 510(k)s
- The adhesion test report has been conducted to verify the maximum duration of use for the electrodes used by the subject device according to the requirements of AAMI EC 12_2000(R) 2010-Section 5.4.

-
- The dispersion and shelf life test report has been conducted to verify the current dispersion and shelf-life of the electrodes used by the subject device in the expiration date according to the requirements of the FDA guidance –Shelf Life of Medical Device and ASTM F 1980-07 standard.

Summary

Based on the above performance as documented in this application, Electronic Pulse Stimulator was found to have a safety and effectiveness profile that is similar to the predicate devices.

VIII. Conclusions

The subject device Electronic Pulse Stimulator is to be concluded substantial equivalent to its predicate devices.