



November 5, 2019

SunBio, Inc.
% Stuart Goldman
Sr. Consultant, RA/QA
Emergo by UL
2500 Bee Cave Road; Bldg. 1, Suite 300
Austin, Texas 78746

Re: K190144
Trade/Device Name: MucoPEG
Regulatory Class: Unclassified
Product Code: LFD
Dated: August 9, 2019
Received: August 13, 2019

Dear Stuart Goldman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Srinivas Nandkumar, Ph.D.
Acting Director
DHT1B: Division of Dental Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K190144

Device Name

MucoPEG™

Indications for Use (Describe)

MucoPEG™ is indicated to relieve the symptoms and discomfort of dry mouth, refresh, moisturize/hydrate, clean, soothe oral irritation, and lubricate oral dryness.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K190144

510(k) Summary

MucoPEG™

1. Submission Sponsor

SunBio, Inc.
95 Sanbon-ro
Gunpo-si
Gyeonggi-do
South Korea 15849

Contact: Ms. Sun S. Kim
Title: Executive Director, Clinical Research
Phone: (925) 876-0439 (USA); 82-31-423-5467 (S. Korea)

2. Submission Correspondent

Emergo Global Consulting, LLC
2500 Bee Cave Road
Building 1, Suite 300
Austin, TX 78746

Contact: Stuart R. Goldman
Title: Senior Consultant RA/QA
Phone: (512) 327-9997

3. Date Prepared

March 13, 2019

4. Device Identification

Trade Name: MucoPEG™
Common Name: Artificial Saliva
Classification Name: Pre-Amendment
Regulation Number: Pre-Amendment
Product Code: LFD
Class: Unclassified
Classification Panel: Dental

5. Legally Marketed Predicate Device

Predicate Device – Hydris™ Oral Rinse (K163029; Product Code LFD)

6. Indications for Use Statement

MucoPEG™ is indicated to relieve the symptoms and discomfort of dry mouth, refresh, moisturize/hydrate, clean, soothe oral irritation, and lubricate oral dryness.

7. Device Overview

MucoPEG™ is an artificial saliva to relieve xerostomia symptoms and discomfort. The main ingredient of MucoPEG™ is a polyethylene glycol (PEG) derivative. The additives to MucoPEG™ are mint (for flavor) and sodium bicarbonate (for pH maintenance when dissolved in water). MucoPEG™ is formulated as a powder and packaged in a single-use packet.

8. Substantial Equivalence Discussion

The substantial equivalence comparison of MucoPEG™ to Hydris™ with respect to its intended use and technological characteristics (ingredients, conditions of use, physical properties and biocompatibility) is shown in **Table 5-1** and provides detailed information regarding the basis for the determination of substantial equivalence between the subject and predicate devices. Two additional reference devices, Neutrasal® (K093642) and Coseal® (P030039), are being added to this 510(k) for MucoPEG™ to incorporate technological features of the subject device that are not contained in the predicate device Hydris™.

Table 5-1 – Substantial Equivalence Comparison

Attributes	Predicate Device	Subject Device	Similarities / Differences
Regulatory Information			
Device Name	Hydris™ Oral Rinse	MucoPEG™	-
Manufacturer	Colgate-Palmolive	SunBio	-
510(k) #	K163029	K190144	-
Product Code	LFD	LFD	Same
Regulation	Pre-Amendment	Pre-Amendment	Same
Class	Unclassified	Unclassified	Same
Review Panel	Dental	Dental	Same
Intended Use	Temporary relief of xerostomia (i.e., dry mouth).	Temporary relief of xerostomia (i.e., dry mouth).	Same
Indications for Use	To relieve the symptoms and discomfort of dry mouth, refresh, moisturize/hydrate, clean, soothe oral irritation, and lubricate oral dryness.	To relieve the symptoms and discomfort of dry mouth, refresh, moisturize/hydrate, clean, soothe oral irritation, and lubricate oral dryness.	Same
Ingredients			
Solvent	H ₂ O	H ₂ O	Same
Humectants/Moisturizers	Glycerin, Propylene Glycol (PG)	Polyethylene Glycol (PEG) Tetra Succinimidyl Glutarate	Similar. The PEG used in the subject device and the PG used in the predicate device provide the same mechanical mechanism of action and are alkylene oxide by-products sourced from the petroleum refining process, and then purified to USP grade.

Buffers	Disodium Phosphate, Sodium Phosphate	Sodium Bicarbonate	Similar. These different buffering agents are generally recognized as safe (GRAS) and provide similar pH levels for the subject and predicate devices.
Flavors	Mint	Mint	Same
Thickeners	Cellulose Gum, Xanthan Gum, Carbomer	None	Thickeners are not used in the subject device.
Surfactants	Poloxamer 407	None	Surfactants are not used in the subject device.
Preservatives	Cetylpyridinium Chloride, Sodium Benzoate	None	Preservatives are not used in the subject device.
Colorants	FD&C Blue 1	None	Colorants are not used in the subject device.
Sweeteners	Sorbitol, Sodium Saccharin, Sucralose	None	Sweeteners are not used in the subject device.
Conditions of Use			
Area of Use	Oral cavity	Oral cavity	Same
Dosage Form	Oral rinse	Oral rinse	Same
Dosage (per use)	20 mL (4 teaspoons)	1g dissolved in 20 mL water.	Same
Applications/Day	Up to 2 times daily.	As often as needed.	Similar
Method of Use	Ready to use liquid.	Powder to be mixed with H ₂ O prior to use.	Similar, both are in liquid form when used.
Environment of Use	Home	Home	Same
Prescription/OTC	OTC	R _x	The subject device is sold by prescription to maintain its prescribed storage condition and product integrity while it remains at the pharmacy prior to being sold to the end user.
Packaging Unit	8/16/33.8 fl-oz bottles	30 (1g) packets / box Dissolution Bottle	Similar. The differences between the way the subject and predicate devices are packaged is a function of their physical state.
Physical Properties (Performance Testing)			
Appearance	Viscous liquid, when applied	Viscous liquid, when applied	Same

pH	5.91	7.29	Similar. The subject and predicate devices have pH values near the pH of saliva (5.3 - 7.8).
Viscosity (Pa-s)	0.0126	0.0019	Similar. The values are not significantly different when compared to the viscosity of the saliva (0.0078 Pa-s).
Surface Tension (dyn/cm)	30.24	49.28	Similar. The values are not significantly different when compared to the surface tension of saliva (58 dyn/cm).
Moisturization/Hydration of Cell	12.0±3.38%	71.9±47.26%	Similar. The difference is not statistically significantly different in protecting epithelial cells from dryness (p-value = 0.093).
Sterility	Non-Sterile	Non-Sterile	Similar.
Shelf-Life	24 months	18 months	Similar. The validated shelf-life of the subject device is 18 months when stored as instructed.
Storage Condition	Not specified	-20±5 °C (-4 ± 9°F)	The subject device is stored under the prescribed storage condition to maintain the integrity (stability) of the product.
Biocompatibility Testing			
Biocompatibility	Conforms with ISO 10993-1 <ul style="list-style-type: none"> • cytotoxicity • sensitization • irritation 	Conforms with ISO 10993-1 <ul style="list-style-type: none"> • cytotoxicity • sensitization • irritation 	Same

9. Non-Clinical Performance Data

As part of demonstrating substantial equivalence of MucoPEG™ to the predicate device, SunBio submitted the subject device for testing in accordance with the applicable parts of the following voluntary standards:

- ISO 10993-5 (cytotoxicity)
- ISO 10993-10 (irritation)
- ISO 10993-10 (sensitization)

SunBio also performed various physical properties tests on MucoPEG™ in a side-by-side manner against the predicate device for certain performance characteristics that are relevant to artificial saliva products. These tests included:

- Appearance
- pH
- Viscosity
- Surface Tension
- Moisturization/Hydration of cultured epithelial layer cells

10. Statement of Substantial Equivalence

MucoPEG™ artificial saliva has the same intended use and indications for use, and similar technological features as Hydris™ Oral Rinse. Any minor differences in the ingredients used to make the subject device when compared to the predicate device have been successfully evaluated by SunBio through side-by-side comparative testing of the physical properties of the two devices, as well as by biocompatibility testing on the subject device, such that the information submitted to the FDA demonstrates that the subject device, when compared to the predicate device, does not raise any new questions of safety and effectiveness. MucoPEG™ artificial saliva, as designed and manufactured by SunBio, has been determined to be substantially equivalent to Hydris™ Oral Rinse.