



Durr Dental SE
% Mr. Daniel Kamm
Principal Engineer
Kamm & Associates
8870 Ravello Ct
NAPLES FLORIDA 34114

April 1, 2019

Re: K190629

Trade/Device Name: DBSWIN and VISTAEASY Imaging Software
Regulation Number: 21 CFR 892.2050
Regulation Name: Picture archiving and communications system
Regulatory Class: Class II
Product Code: LLZ
Dated: March 4, 2019
Received: March 12, 2019

Dear Mr. Daniel Kamm:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Thalia Mills, Ph.D.
Director
Division of Radiological Health
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K190629

Device Name
DBSWIN and VISTAEASY Imaging
Software

Indications for Use (Describe)

DBSWIN and VistaEasy imaging software are intended for use by qualified dental professionals for windows based diagnostics. The software is a diagnostic aide for licensed radiologists, dentists and clinicians, who perform the actual diagnosis based on their training, qualification, and clinical experience. DBSWIN and VistaEasy are clinical software applications that receive images and data from various imaging sources (i.e., radiography devices and digital video capture devices) that are manufactured and distributed by DÜRR Dental and Air Techniques. It is intended to acquire, display, edit (i.e., resize, adjust contrast, etc.) and distribute images using standard PC hardware. In addition, DBSWIN enables the acquisition of still images from 3rd party TWAIN compliant imaging devices (e.g., generic image devices such as scanners) and the storage and printing of clinical exam data, while VistaEasy distributes the acquired images to 3rd party TWAIN compliant PACS systems for storage and printing.

DBSWIN and VistaEasy software are not intended for mammography use.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

**510(k) Summary, DÜRR DENTAL SE,
DBSWIN and VISTAEASY Imaging Software
K19 0629**

This 510(k) is being submitted in accordance with the requirements of 21 CFR §807.92.

1. Date Summary Prepared:

March 5, 2019

2. Submitter's Identification:

DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Deutschland/Germany
Phone: + 49 (0) 7142 70 5-0
Fax: + 49 (0) 7142 705-500
E-Mail info@duerr.de Internet:
www.duerrdental.com

Responsible Individual Contact/Prepared by:
Oliver Lange, Head Director of Quality Management
DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Deutschland/Germany
Phone: 49-7142-705-190
Email: lange.o@duerr.de
oliver.lange@duerrdental.com

U.S. Agent Contact:
Vincent Vega
Air Techniques, Inc.
1295 Walt Whitman Road
Melville, NY 11747
Phone: 516 2145512
Fax: 516 4336171
Email: Vincent.Vega@airtechniques.com

3. Device:

Trade /Proprietary Name:	DBSWIN and VISTAEASY Imaging Software
Device:	Medical Imaging Software
Regulation Description:	Picture archiving and communications system.
Regulation Medical Specialty:	Radiology
Review Panel	Radiology
Product Code	LLZ
Regulation Number	892.2050

4. Predicate Device:

Legally Marketed Predicate Device Information:

510(k) Number:	K161444
Manufacturer:	"DÜRR DENTAL SE".
Trade /Proprietary Name	DBSWIN and VistaEasy Imaging Software
Device:	Medical Imaging Software
Regulation Description:	Picture archiving and communications system.
Regulation Medical Specialty:	Radiology
Review Panel	Radiology
Product Code	LLZ
Regulation Number	892.2050
Device Class	2

5. Device Description:

DBSWIN and VistaEasy imaging software is an image management system that allows dentists to acquire, display, edit, view, store, print, and distribute medical images. DBSWIN and VistaEasy software runs on user provided PC-compatible computers and utilize previously cleared digital image capture devices for image acquisition. VistaEasy is included as part of DBSWIN. It provides additional interfaces for Third Party Software. VistaEasy can also be used by itself, as a reduced feature version of DBSWIN.

6. Indications for Use

DBSWIN and VistaEasy imaging software are intended for use by qualified dental professionals for windows based diagnostics. The software is a diagnostic aide for licensed radiologists, dentists and clinicians, who perform the actual diagnosis based on their training, qualification, and clinical experience. DBSWIN and VistaEasy are clinical software applications that receive images and data from various imaging sources (i.e., radiography devices and digital video capture devices) that are manufactured and distributed by DÜRR Dental and Air Techniques. It is intended to acquire, display, edit (i.e., resize, adjust contrast, etc.) and distribute images using standard PC hardware. In addition, DBSWIN enables the acquisition of still images from 3rd party TWAIN compliant imaging devices (e.g., generic image devices such as scanners) and the storage and printing of clinical exam data, while VistaEasy distributes the acquired images to 3rd party TWAIN compliant PACS systems for storage and printing.

DBSWIN and VistaEasy software are not intended for mammography use.

7. Summary of the technological characteristics of the device compared to the predicate device:

DBSWIN and VistaEasy imaging software by DÜRR DENTAL SE are two software components that are substantially equivalent to DBSWIN and VistaEasy Imaging Software K161444 software applications that have identical indications for use, functionality, performance, and features as shown in the following comparison Table.

Comparison of the Technological Characteristics

Descriptive Information	DBSWIN and VistaEasy Imaging Software K161444	DBSWIN and VistaEasy Imaging Software (Modified)
Indications for Use	DBSWIN and VistaEasy imaging software are intended for use by qualified dental professionals for windows based diagnostics. The software is a diagnostic aide for licensed radiologists, dentists and clinicians, who perform the actual diagnosis based on their training, qualification, and clinical experience. DBSWIN and VistaEasy are clinical software applications that receive images and data from various imaging sources (i.e., radiography devices and digital video capture devices) that are manufactured and distributed by DÜRR Dental and Air Techniques. It is intended to acquire, display, edit (i.e., resize, adjust contrast, etc.) and distribute images using standard PC hardware. In addition, DBSWIN enables the acquisition of still images from 3rd party TWAIN compliant imaging devices (e.g., generic image devices such as scanners) and the storage and printing of clinical exam data, while VistaEasy distributes the acquired images to 3rd party TWAIN compliant PACS systems for storage and printing. DBSWIN and VistaEasy software are not intended for mammography use.	SAME, unchanged.
Patient Management	YES	YES
Image Management	YES	YES
Acquisition Sources		
X-ray (i.e., Phosphor Plate, Digital Panoramic)	YES	YES
Laser Fluorescence Caries Detection Aid	YES	YES
Video	YES	YES

Descriptive Information	DBSWIN and VistaEasy Imaging Software K161444	DBSWIN and VistaEasy Imaging Software (Modified)
Photos	YES	YES
Documents	YES	YES
Import*	YES	YES
Display Images	YES	YES
Save/Store Images*	YES	YES
Produce Reports*	YES	YES
Print/Export Images*	YES	YES
Enhance Images		
Brightness	YES	YES
Contrast	YES	YES
Colorize*	YES	YES
Crop	YES	YES
Rotate	YES	YES
Zoom In/Out	YES	YES
Invert*	YES	YES
Sharpen*	YES	YES
Measure*	YES	YES
Over/Under Exposure	YES	YES
Annotate*	YES	YES
Run on standard PC compatible computers	YES	YES

Descriptive Information	DBSWIN and VistaEasy Imaging Software K161444	DBSWIN and VistaEasy Imaging Software (Modified)
Supported Devices	VistaCam iX Proof CamX Spectra VistaCam iX Cam CamX Elara CamX Polaris VistaCam iX HD CamX Triton HD CamX Luna HD VistaCam iX Macro VistaPano S ProVecta S-Pan VistaPano S CEPH ProVecta S-PAN CEPH ScanX Swift ScanX Duo ScanX Classic ScanX Ortho ScanX Intraoral	Same, plus these newly integrated devices: ScanX Swift View VistaScan Nano ScanX Classic View CamX Triton HD Proxi CamX Triton HD Spectra
Computer operating systems	Microsoft Windows 7, 32 bit (from Home Premium), Microsoft Windows 7, 64 bit (from Home Premium), Microsoft Windows 8, 64 bit Microsoft Windows 10, 64-bit Microsoft Windows Server 2012	Microsoft Windows 7, 32-bit (from Home to Premium) Microsoft Windows 7, 64-bit (from Home to Premium) Microsoft Windows 8.1, 64-bit Microsoft Windows 10, 64-bit Microsoft Windows Server 2012 Microsoft Windows Server 2016
CPU	≥ Intel Pentium IV compatible, 1.4 GHz, ≥ Intel Core i3	≥ Intel Pentium IV compatible, 1.4 GHz
RAM	≥ 1GB (2GB recommended), ≥ 4 GB	≥ 1GB (2GB recommended)
Drive	DVD-ROM	No change
Hard Disk	Workstation (without database) ≥50 GB	Workstation (without database) ≥50 GB The memory requirements of the database depend on the number of images taken at the office. (Camera image: Approx. 1 MB, X-ray image: Approx. 2 MB – 10 MB)
Data Backup	Daily data back up	No change
Interface	Ethernet ≥ 100 Mbit	No change
Diagnostic Monitor	SVGA ≥ 17", ≥ 1024 x 768 pixel, 24/32 bit color depth	No change
Resolution /Graphics	≥ 1024 x 768	≥ 1024 x 768 Depth of color 32-bit, 16.7 million colors
Run on standard PC compatible computers	Yes	No change

* Not available on VistaEasy

8. Discussion of Similarities and Differences:

The list of compatible operating systems has changed. We no longer support Windows XP and Windows 8. We now support these additional operating systems: Microsoft Windows 8.1, 64-bit, Microsoft Windows 10, 64-bit Microsoft Windows Server 2016.

We have added Cybersecurity precautions and information.

Compatibility with these additional image sources has been added: ScanX Swift View; VistaScan Nano; ScanX Classic View; CamX Triton HD Proxi; CamX Triton HD Spectra.

Firmware Files for DÜRR Dental SE and Air Techniques scanner devices has been updated. UDI information display has been added.

9. Non-Clinical Data and Performance Testing

DBSWIN/VistaEasy was developed in compliance with the harmonized standard of IEC 62304 for medical device software life cycle requirements.

DBSWIN product has been in sales and distribution in the European dental market for over 19 years serving and performing the same intended use, functionality, and hardware compatibility interfaces with 3rd party software.

Bench testing, effectiveness, and functionality were successfully conducted and verified between DBSWIN and VistaEasy, and image capture devices.

DBSWIN is DICOM compliant.

Risk Analysis based design development and design reviews were conducted.

Full functional software cross check testing was performed.

Labeling has been updated to assure compliance with the FDA Guidance: *Content of Premarket Submissions for Management of Cybersecurity in Medical Devices Guidance for Industry and Food and Drug Administration Staff Document Issued on: October 2, 2014*

10. Clinical Data: Not required for a finding of substantial equivalence.

11. Conclusion:

In accordance with the Federal Food, Drug and Cosmetic Act, 21 CFR Part 807 and based on the similarity to the predicate device in terms of technology, performance and indications for use, DÜRR DENTAL SE concludes that the DBSWIN and VistaEasy Imaging Software is substantially equivalent to the predicate device as described herein.

The minor device modifications to DBSWIN/VistaEasy do not alter the fundamental scientific technology of the predicate device and summary level information is adequate to assess the modifications. The verification testing demonstrates that the device continues to meet its performance specifications and the results of the testing did not raise new issues of safety or effectiveness. Therefore, the modified DBSWIN/VistaEasy can be found substantially equivalent to the predicate device as cleared in K161444.