



June 25, 2019

Roshana Ahmed  
President  
Quaras, LLC  
2101 Camino Rey  
Fullerton, California 92833

Re: K190784

Trade/Device Name: BUHLMANN fCAL turbo  
Regulation Number: 21 CFR 866.5180  
Regulation Name: Fecal calprotectin immunological test system  
Regulatory Class: Class II  
Product Code: NXO  
Dated: March 25, 2019  
Received: March 27, 2019

Dear Roshana Ahmed:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Doug Jeffery, Ph.D.  
Acting Deputy Division Director  
Division of Immunology and Hematology Devices  
OHT7: Office of In Vitro Diagnostics and Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K190784

Device Name

BÜHLMANN fCAL turbo

Indications for Use (Describe)

The BÜHLMANN fCAL turbo is an in vitro diagnostic assay intended for the quantitative measurement of fecal calprotectin, a neutrophilic protein that is a marker of intestinal mucosal inflammation, in human stool. The BÜHLMANN fCAL turbo aids in the diagnosis of inflammatory bowel disease (IBD), specifically Crohn's disease (CD) and ulcerative colitis (UC) and aids in the differentiation of IBD from irritable bowel syndrome (IBS) in conjunction with other laboratory and clinical findings.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

### I. Submitter

BÜHLMANN Laboratories AG  
Baselstrasse 55  
Schönenbuch CH-4124  
Switzerland  
Phone: +41 61 487 12 50

Contact Person: Laura Zurbrügg  
Date Prepared: June 20, 2019

### II. Device

Device Proprietary Name:	BUHLMANN fCAL <sup>®</sup> turbo
Common or Usual Name:	Fecal calprotectin immunological test system
Classification Name:	Calprotectin, Fecal
Regulation Number:	21 CFR 866.5180
Product Code:	NXO
Device Classification	II

### III. Predicate Device

Substantial equivalence is claimed to the following device:

- BÜHLMANN fCAL<sup>®</sup> ELISA, K181012, BÜHLMANN Laboratories AG

### IV. Device Description

The BÜHLMANN fCAL<sup>®</sup> turbo, a particle-enhanced turbidimetric immunoassay (PETIA), is performed using patient stool extracts collected without preservatives. Calprotectin within the sample extract mediates immunoparticle agglutination; sample turbidity is proportional to calprotectin concentration. The detected light absorbance allows quantification of calprotectin concentration via interpolation of an established calibration curve. The assay is validated for use on clinical chemistry analyzers such as the Roche cobas<sup>®</sup> c501/c502 platforms.

The BÜHLMANN fCAL<sup>®</sup> turbo Reagent Kit is to be used in conjunction with the BÜHLMANN fCAL<sup>®</sup> turbo Calibrator Kit and BÜHLMANN fCAL<sup>®</sup> turbo Control Kit, which are available separately.

## V. Indications for Use

The BÜHLMANN fCAL® turbo is an in vitro diagnostic assay intended for the quantitative measurement of fecal calprotectin, a neutrophilic protein that is a marker of intestinal mucosal inflammation, in human stool. The BÜHLMANN fCAL® turbo aids in the diagnosis of inflammatory bowel disease (IBD), specifically Crohn's disease (CD) and ulcerative colitis (UC) and aids in the differentiation of IBD from irritable bowel syndrome (IBS) in conjunction with other laboratory and clinical findings.

## VI. Comparison of Technological Characteristics

The BÜHLMANN fCAL® turbo and the predicate device share the following characteristics:

- measurement of human fecal calprotectin in human stool;
- use of a quantitative platform;
- use of manual weighing extraction method for stool samples; and
- clinical decision thresholds.

The BÜHLMANN fCAL® turbo is technologically different from the predicate device as follows:

- assay method;
- use of an alternate, automated, detection method;
- use of polyclonal antibodies; and
- broader measuring range.

The tables below compare key technological features between the subject and predicate device.

### Technological comparison

#### Comparison of Assay

Similarities		
	BÜHLMANN fCAL® turbo	BÜHLMANN fCAL® ELISA (K181012)
Analyte	Human fecal calprotectin (MRP8/14)	Human fecal calprotectin (MRP8/14)
Assay format	Quantitative	Quantitative
Specimen type	Human stool	Human stool
Extraction Method	Manual Weighing (1:50 dilution in Extraction Buffer)	Manual Weighing (1:50 dilution in Extraction Buffer)
Clinical Decision Thresholds	Normal: < 80 µg/g Gray-zone/borderline: 80 - 160 µg/g Elevated: > 160 µg/g	Normal: < 80 µg/g Gray-zone/borderline: 80 - 160 µg/g Elevated: > 160 µg/g

<b>Differences</b>		
	<b>BÜHLMANN fCAL® turbo</b>	<b>BÜHLMANN fCAL® ELISA (K181012)</b>
Method	PETIA	ELISA
Automation	Automated	Not automated
Solid phase	Polystyrene nanoparticles (beads)	96-well polystyrene microtiter plate
Detection method	Automated clinical chemistry analyzer read at 546 nm	Microtiter plate reader read at 450 nm
Analyte-specific antibody components	Polyclonal antibodies against human calprotectin coated on polystyrene beads	Capture antibodies: monoclonal antibodies against human calprotectin coated on microtiter plates
		Detection antibodies: monoclonal antibodies against human calprotectin conjugated to HRP
Measuring range	Direct measuring range: 30 - 2000 µg/g Measuring range with automatic dilution: 30 – 10,000 µg/g	30 – 1800 µg/g

### Comparison of Calibrators

	<b>BÜHLMANN fCAL® turbo</b>	<b>BÜHLMANN fCAL® ELISA</b>
Indications for Use	The BÜHLMANN fCAL® turbo Calibrator Kit is intended for use with the BÜHLMANN fCAL® turbo Reagent Kit for the determination of fecal calprotectin levels in extracted stool samples. Comprised of six (6) calibrators, each calibrator establishes a point of reference for the working curve that is used to calculate test results from patient samples.	N/A
Method	BÜHLMANN fCAL® turbo	BÜHLMANN fCAL® ELISA
Analyte	Native human calprotectin Source: human granulocyte extract	Native human calprotectin Source: human serum
Calibrators	6 levels: Target values: 0, 50, 200, 500, 1000, 2000 µg/g	5 levels: 4, 12, 40, 120, and 240 ng/mL
Conversion factor	N/A	7.5
Value assignment:	Calibrator values assigned using a value transfer protocol for each calibrator lot. Values are indicated in the QC datasheet.	Nominal values: 30 µg/g, 90 µg/g, 300 µg/g, 900 µg/g, 1800 µg/g

	<b>BÜHLMANN fCAL® turbo</b>	<b>BÜHLMANN fCAL® ELISA</b>
Configuration	Available as a separate BÜHLMANN fCAL® turbo Calibrator Kit.	Included within the BÜHLMANN fCAL® ELISA kit.

**Comparison of Controls**

	<b>BÜHLMANN fCAL® turbo</b>	<b>BÜHLMANN fCAL® ELISA</b>
Indications for Use	The BÜHLMANN fCAL® turbo Control Kit, comprised of a high and low control, is intended for use with the BÜHLMANN fCAL® turbo Reagent Kit, for quality control, in the determination of fecal calprotectin levels in extracted stool samples.	N/A
Method	BÜHLMANN fCAL® turbo	BÜHLMANN fCAL® ELISA
Analyte	Native human calprotectin Source: human granulocyte extract	Native human calprotectin Source: human serum
Levels	2 (low and high)	2 (low and high)
Physio-chemical characteristics	Ready to use	Ready to use
Configuration	Available as a separate BÜHLMANN fCAL® turbo Control Kit	Included within the BÜHLMANN fCAL® ELISA kit

Discussion

As seen above, differences between the subject and predicate device include the assay and detection methods, antibodies, and reportable measuring range. These technological differences do not create new questions of safety and effectiveness and the differences are addressed by the performance studies identified below.

## VII. Performance Data

### A. Clinical Thresholds

Calprotectin Concentration	Interpretation	Follow-Up
< 80 µg/g	Normal	None
80 µg/g ≤ x ≤ 160 µg/g	Gray-zone/Borderline	Retest within 4 – 6 weeks
> 160 µg/g	Elevated	Retest as appropriate

### B. Precision

#### Single-Site Repeatability Study Results:

ID	Mean [µg/g]	n	Within-run (Repeatability)		Between-run		Between-day		Within-laboratory	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
S01	42.9	80	3.6	8.3%	1.2	2.7%	1.1	2.5%	3.9	9.1%
S02	98.4	80	2.5	2.6%	1.8	1.8%	2.2	2.2%	3.7	3.8%
S03	166.5	80	4.3	2.6%	0.8	0.5%	1.9	1.2%	4.8	2.9%
S04	267.6	80	3.9	1.4%	2.5	0.9%	1.8	0.7%	5.0	1.9%
S05	642.0	80	20.1	3.1%	14.9	2.3%	0.0	0.0%	25.1	3.9%
S06	1414.2	80	19.6	1.4%	11.1	0.8%	3.5	0.2%	22.8	1.6%
S07	3251.4	80	40.8	1.3%	21.4	0.7%	19.7	0.6%	50.1	1.5%
S08	5405.6	80	40.2	0.7%	56.6	1.0%	34.5	0.6%	77.5	1.4%

#### Multi-Site Reproducibility Study Results:

ID	Mean [µg/g]	n	Within-run (Repeatability)		Between-day		Between-site		Total Precision	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
S01	47.2	75	3.6	7.6	2.4	5.0	0.0	0.0	4.3	9.1
S02	91.1	75	3.5	3.8	3.5	3.8	2.8	3.1	5.7	6.2
S03	185.4	75	5.1	2.7	2.7	1.4	5.5	3.0	7.9	4.3
S04	276.9	75	6.4	2.3	4.5	1.6	9.7	3.5	12.5	4.5
S05	674.5	75	12.9	1.9	1.2	0.2	22.8	3.4	26.3	3.9
S06	1519.6	75	25.3	1.7	17.8	1.2	62.3	4.1	69.6	4.6
S07	3343.8	75	54.6	1.6	35.6	1.1	100.0	3.0	119.4	3.6
S08	5475.6	75	72.1	1.3	35.8	0.7	154.2	2.8	173.9	3.2



*Lot-to-Lot Precision Study Results:*

ID	Mean [µg/g]	n	Within-Run (Repeatability)		Between-Day		Between-Lot		Total Precision	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
S1	45.2	75	3.22	7.1%	1.36	3.0%	3.70	8.2%	5.09	11.3%
S2	86.4	75	3.69	4.3%	1.19	1.4%	5.66	6.6%	6.86	7.9%
S3	175.8	75	5.04	2.9%	0.29	0.2%	9.90	5.6%	11.11	6.3%
S4	263.9	75	7.55	2.9%	0.00	0.0%	9.98	3.8%	12.52	4.7%
S5	647.4	75	15.47	2.4%	0.00	0.0%	15.28	2.4%	21.74	3.4%
S6	1460.7	75	33.66	2.3%	11.64	0.8%	41.01	2.8%	54.32	3.7%
S7	3234.5	75	71.23	2.2%	8.90	0.3%	130.29	4.0%	148.76	4.6%
S8	5303.1	75	97.42	1.8%	11.18	0.2%	163.87	3.1%	190.97	3.6%

*Extraction Reproducibility Study Results:*

Sample	n	Mean (µg/g)	Within-Run (Repeatability)		Between-extraction		Between-day		Between-operator		Total Precision	
			SD (µg/g)	%CV	SD (µg/g)	%CV	SD (µg/g)	%CV	SD (µg/g)	%CV	SD (µg/g)	%CV
S1	80	47.7	2.8	5.9	1.1	2.4	0.7	1.5	1.4	2.9	3.4	7.2
S2	80	72.3	3.8	5.2	3.9	5.4	4.2	5.8	0.0	0.0	6.8	9.5
S3	80	96.1	3.8	3.9	2.2	2.3	1.4	1.4	0.0	0.0	4.6	4.8
S4	80	170.6	4.0	2.4	2.5	1.5	8.7	5.1	0.0	0.0	9.9	5.8
S5	80	277.0	3.7	1.4	27.9	10.1	10.0	3.6	11.0	4.0	31.8	11.5
S6	80	421.1	9.8	2.3	5.9	1.4	15.3	3.6	0.0	0.0	19.1	4.5
S7	80	573.9	5.4	0.9	39.5	6.9	0.0	0.0	0.0	0.0	39.9	6.9
S8	80	1387.4	39.1	2.8	75.1	5.4	159.9	11.5	0.0	0.0	180.9	13.0
S9	80	3264.9	87.2	2.7	236.2	7.2	256.9	7.9	0.0	0.0	359.7	11.0
S10	80	3330.4	89.8	2.7	92.4	2.8	75.7	2.3	0.0	0.0	149.4	4.5

C. Linearity

Study procedures were performed using two dilution series. For each dilution series, a stool specimen extract with concentration above the anticipated upper limit of the analytical measuring range was combined with a stool specimen extract with concentration below the anticipated lower limit of the analytical measuring range, in various mixing ratios covering the range; each dilution was tested in 4 replicates. Results of the linear regression analyses are presented in the table below.

Best	Measuring Range [µg/g]	Linear regression parameters		
		Intercept (95% C.I.)	Slope (95% C.I.)	R <sup>2</sup>
1	37.6 – 12,216.0	5.7 (1.6, 16.9)	1.057 (1.044, 1.075)	0.9983
2	33.5 – 13,339.5	3.8 (-0.4, 13.3)	1.031 (1.014, 1.042)	0.9984

The data supports the following claims for analytical measuring range:

Direct analytical measuring range: 30 – 2000 µg/g

Measuring range with automatic dilution: 30 – 10,000 µg/g

D. High Dose Hook Effect

No high dose hook effect at theoretical concentrations up to 45,715 µg/g.

E. Accuracy/Recovery

Sample No	7226	7228	7238	7236	7244	7234	7246
Baseline result [µg/g]	44.10	65.45	116.43	138.48	230.88	510.78	1076.33
Expected post-spike result [µg/g]	101.04	122.39	173.37	195.42	458.65	738.55	1304.10
Observed post-spike result [µg/g]	94.55	114.53	170.23	186.93	453.10	753.18	1309.28
Total recovery [%]	93.6%	93.6%	98.2%	95.7%	98.8%	102.0%	100.4%

F. Analytical Sensitivity

Results of the analytical sensitivity studies support a claimed direct measuring range of 30 – 2000 µg/g and a measuring range of 30 – 10,000 µg/g with automatic dilution.

LoB = 16.7 µg/g

LoD = 23.7 µg/g

LoQ = 30 µg/g

G. Interfering Substances

Study procedures were performed using stool specimen extracts with the following approximate calprotectin concentrations: 30 µg/g, 100 µg/g, 300 µg/g, and 550 µg/g. The following analytes, pharmaceuticals, and nutritional supplements did not interfere with the BÜHLMANN fCAL® turbo:

Trade name	Active component	Solvent	Concentration mg/50 mg stool
gyno-Tardyferon	Iron (II) sulfate	HCl/NaOH	0.11
Prednisone	Prednisone	DMSO	0.31
Imurek	Azathioprine	DMSO	0.19
Salofalk	Mesalamine; 5-ASA	DMSO	5.21
Agopton	Lansoprazole	Dimethylformamide	0.18
Asacol	Mesalamine; 5-ASA	DMSO	2.50
Vancocin	Vancomycin	H <sub>2</sub> Odd	2.00
Sulfamethoxazole	Sulfamethoxazole	DMSO	1.60
Trimethoprim	Trimethoprim lactate	DMSO/Exbuffer	0.35
Ciproxine	Ciprofloxacin	solvent from manufacturer/H <sub>2</sub> Odd	1.25
Vitamin E	DL-α Tocopherol Acetate	H <sub>2</sub> O + Tween	0.30
Bion 3	Multivitamin	HCl/NaOH	1.06
Hemoglobin	Hemoglobin	H <sub>2</sub> Odd	1.25

The following enteropathological microorganisms did not interfere with the BÜHLMANN fCAL® turbo when added to stool extracts at the given cell counts:

Microorganism	Concentration (cfu/mL)
<i>Escherichia coli</i>	3.3 x 10 <sup>7</sup>
<i>Salmonella enterica subsp. enterica</i>	9.0 x 10 <sup>7</sup>
<i>Klebsiella pneumoniae subsp. pneumonia</i>	5.3 x 10 <sup>7</sup>
<i>Citrobacter freundii</i>	12.9 x 10 <sup>7</sup>
<i>Shigella flexneri</i>	5.0 x 10 <sup>7</sup>
<i>Yersinia enterocolitica subsp. enterocolitica</i>	9.8 x 10 <sup>7</sup>

*H. Method Comparison*

A total of 248 clinical study samples were tested using the BÜHLMANN fCAL<sup>®</sup> turbo and the predicate device (BÜHLMANN fCAL<sup>®</sup> ELISA assay); valid results within the linear measuring range for both assays were obtained for 220 of these samples. Results were analyzed by Passing-Bablok regression analysis.

<b>Slope (95% CI)</b>	<b>Intercept (µg/g) (95% CI)</b>	<b>Bias at 80 µg/g (95% CI)</b>	<b>Bias at 160 µg/g (95% CI)</b>	<b>Correlation <i>r</i></b>
1.025 (0.990, 1.058)	-4.5 (-8.7, 0.3)	-3.1% (-7.2%, 0.5%)	-0.3% (-2.4%, 2.7%)	0.972

Frequency counts of BÜHLMANN fCAL<sup>®</sup> turbo test results and corresponding BÜHLMANN fCAL<sup>®</sup> ELISA assay results within each of the diagnostic ranges of these tests are provided below.

		<b># in BÜHLMANN fCAL<sup>®</sup> ELISA assay range (µg/g)</b>			
		<b>&lt; 80</b>	<b>80 - 160</b>	<b>&gt; 160</b>	<b>Total</b>
<b># in fCAL turbo range (µg/g)</b>	<b>&lt; 80</b>	84	10	0	94
	<b>80 - 160</b>	8	41	6	55
	<b>&gt; 160</b>	0	7	92	99
<b>Total</b>		92	58	98	248

Estimates of positive percent agreement (PPA) and negative percent agreement (NPA) between the BÜHLMANN fCAL<sup>®</sup> turbo results and corresponding BÜHLMANN fCAL<sup>®</sup> ELISA assay results, using both sets of assay cutoffs, with respect to IBD subjects, IBS subjects, other GI subjects, normal subjects, and all subjects combined are shown in the table below.

Subgroup	Metric	Estimate	95% C.I.
IBD	PPA (lower cutoff)	68/70 = 97.1%	[90.1%, 99.7%]
	NPA (lower cutoff)	6/7 = 85.7%	[42.1%, 99.6%]
	PPA (upper cutoff)	52/56 = 92.9%	[82.7%, 98.0%]
	NPA (upper cutoff)	19/21 = 90.5%	[69.6%, 98.8%]
IBS	PPA (lower cutoff)	28/32 = 87.5%	[71.0%, 96.5%]
	NPA (lower cutoff)	31/33 = 93.9%	[79.8%, 99.3%]
	PPA (upper cutoff)	12/13 = 92.3%	[64.0%, 99.8%]
	NPA (upper cutoff)	49/52 = 94.2%	[84.1%, 98.8%]
Other GI	PPA (lower cutoff)	20/21 = 95.2%	[76.2%, 99.9%]
	NPA (lower cutoff)	16/16 = 100%	[79.4%, 100%]
	PPA (upper cutoff)	13/14 = 92.9%	[66.1%, 99.8%]
	NPA (upper cutoff)	23/23 = 100%	[85.2%, 100%]
Normal	PPA (lower cutoff)	30/33 = 90.9%	[75.7%, 98.1%]
	NPA (lower cutoff)	31/36 = 86.1%	[70.5%, 95.3%]
	PPA (upper cutoff)	15/15 = 100%	[78.2%, 100%]
	NPA (upper cutoff)	52/54 = 96.3%	[87.3%, 99.5%]
All subjects	<b>PPA (lower cutoff)</b>	<b>146/156 = 93.6%</b>	<b>[88.5%, 96.9%]</b>
	<b>NPA (lower cutoff)</b>	<b>84/92 = 91.3%</b>	<b>[83.6%, 96.2%]</b>
	<b>PPA (upper cutoff)</b>	<b>92/98 = 93.9%</b>	<b>[87.1%, 97.7%]</b>
	<b>NPA (upper cutoff)</b>	<b>143/150 = 95.3%</b>	<b>[90.6%, 98.1%]</b>

I. Clinical Sensitivity/Specificity

IBD vs. IBS:

Borderline Values Considered Positive		Clinical Diagnosis		Total
		IBD	IBS	
BÜHLMANN fCAL® turbo	Positive	123	31	154
	Negative	12	99	111
	<b>Total</b>	135	130	265
Sensitivity = 91.1%; 95% C.I. (85.0%, 95.3%)				
Specificity = 76.2%; 95% C.I. (67.9%, 83.2%)				
PPV = 79.9%; 95% C.I. (72.7%, 85.9%)				
NPV = 89.2%; 95% C.I. (81.9%, 94.3%)				

Borderline Values Considered Negative		Clinical Diagnosis		Total
		IBD	IBS	
BÜHLMANN fCAL® turbo	Positive	108	16	124
	Negative	27	114	141
	<b>Total</b>	135	130	265
Sensitivity = 80.0%; 95% C.I. (72.3%, 86.4%)				
Specificity = 87.7%; 95% C.I. (80.8%, 92.8%)				
PPV = 87.1%; 95% C.I. (79.9%, 92.4%)				
NPV = 80.9%; 95% C.I. (73.4%, 87.0%)				

*IBD vs. non-IBD:*

Borderline Values Considered Positive		Clinical Diagnosis		Total
		IBD	Non-IBD	
BÜHLMANN fCAL® turbo	Positive	123	52	175
	Negative	12	150	162
	<b>Total</b>	135	202	337
Sensitivity = 91.1%; 95% C.I. (85.0%, 95.3%)				
Specificity = 74.3%; 95% C.I. (67.7%, 80.1%)				
PPV = 70.3%; 95% C.I. (62.9%, 76.9%)				
NPV = 92.6%; 95% C.I. (87.4%, 96.1%)				

Borderline Values Considered Negative		Clinical Diagnosis		Total
		IBD	IBS	
BÜHLMANN fCAL® turbo	Positive	108	30	138
	Negative	27	172	199
	<b>Total</b>	135	202	337
Sensitivity = 80.0%; 95% C.I. (72.3%, 86.4%)				
Specificity = 85.1%; 95% C.I. (79.5%, 89.8%)				
PPV = 78.3%; 95% C.I. (70.4%, 84.8%)				
NPV = 86.4%; 95% C.I. (80.9%, 90.9%)				

*J. Expected Values/Reference Range:*

Stool samples were obtained from 141 apparently healthy normal adults (> 21 years of age) with no symptoms or signs of gastrointestinal disease. The test results were categorized by the assay cut-offs below.

	Calprotectin level by BÜHLMANN fCAL turbo			Total
	< 80 µg/g	80 - 160 µg/g	> 160 µg/g	
Number of subjects (%)	106 (75.2%)	18 (12.8%)	17 (12.1%)	141 (100%)

**VIII. Conclusion**

The information provided above supports that the BÜHLMANN fCAL<sup>®</sup> turbo is as safe and effective as the predicate device. Although differences exist between the subject and predicate device, verification and validation studies support that these differences do not raise any new questions of safety and effectiveness. Therefore, it is concluded that the BÜHLMANN fCAL<sup>®</sup> turbo is substantially equivalent to the predicate device.

-----**This space intentionally left blank**-----