



May 10, 2019

Exact Imaging, Inc.
% Mr. Mark Job
Responsible Third Party Official
Regulatory Technology Services, LLC
1000 Westgate Drive
Suite 510k
SAINT PAUL MN 55114

Re: K190995

Trade/Device Name: ExactVu™ High Resolution Micro-Ultrasound System
Regulation Number: 21 CFR 892.1550
Regulation Name: Ultrasonic pulsed doppler imaging system
Regulatory Class: Class II
Product Code: IYN, IYO, IYX, OIJ
Dated: April 12, 2019
Received: April 16, 2019

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Thalia T. Mills, Ph.D.
Director
Division of Radiological Health
OHT7: Office of In Vitro Diagnostics
and Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K190995

Device Name

ExactVu™ High Resolution Micro-Ultrasound System

Indications for Use (Describe)

The ExactVu micro-ultrasound system is intended for use by qualified medical professionals for diagnostic ultrasound imaging or fluid flow analysis of the human body. The indications for use (clinical applications) are:

- Small Organ
- Transrectal
- Abdominal

The system may be used with patients of all ages, but is not designed for pediatric or fetal use.

The system is contraindicated for direct cardiac application and for ophthalmic use or any application that causes the acoustic beam to pass through the eye.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Diagnostic Ultrasound Indications for Use Form – ExactVu™ High Resolution Micro-Ultrasound System

System	ExactVu™ High Resolution Micro-Ultrasound System						
Transducer	N/A						
Intended Use	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
	Mode of Operation						
Clinical Application	B (2D Mode)	M	PWD	CWD	Color Doppler	Combined (specify)	Other (specify)
Ophthalmic							
Fetal							
Abdominal	N				N (3)		N (2)
Intra-operative (Abdominal organs and vascular)							
Intra-operative (Neuro.)							
Laparoscopic							
Pediatric							
Small Organ (prostate)	P						P, 1
Neonatal Cephalic							
Adult Cephalic							
Transrectal	P						P, 1
Trans-vaginal							
Trans-urethral							
Trans-esoph. (non-Card.)							
Musculo-skel. (Convent.)							
Musculo-skel. (Superfic.)							
Intra-luminal							
Other (spec.)							
Cardiac Adult							
Cardiac Pediatric							
Trans-esophageal (card.)							
Other (spec.)							
Peripheral vessel							
Other (spec.)							
Dermatology							

N= new indication; P= previously cleared ; E= added under this appendix

Items marked "P" were previously cleared by 510(k) number K180636.

Additional Comments:

1. Includes imaging to assist in the placement of needles for prostate biopsy procedures.
2. Includes imaging to assist in the placement of needles for kidney biopsy procedures.
3. ExactVu supports simultaneous color flow imaging with B-Mode.

Diagnostic Ultrasound Indications for Use Form – Ultrasound Indications for Use Form – EV29L™ High Resolution Transrectal Side-fire Transducer

System	ExactVu™ High Resolution Micro-Ultrasound System						
Transducer	EV29L						
Intended Use	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
	Mode of Operation						
Clinical Application	B (2D Mode)	M	PWD	CWD	Color Doppler	Combined (specify)	Other (specify)
Ophthalmic							
Fetal							
Abdominal							
Intra-operative (Abdominal organs and vascular)							
Intra-operative (Neuro.)							
Laparoscopic							
Pediatric							
Small Organ (prostate)	P						P, 1
Neonatal Cephalic							
Adult Cephalic							
Transrectal	P						P, 1
Trans-vaginal							
Trans-urethral							
Trans-esoph. (non-Card.)							
Musculo-skel. (Convent.)							
Musculo-skel. (Superfic.)							
Intra-luminal							
Other (spec.)							
Cardiac Adult							
Cardiac Pediatric							
Trans-esophageal (card.)							
Other (spec.)							
Peripheral vessel							
Other (spec.)							
Dermatology							

N= new indication; P= previously cleared ; E= added under this appendix

Items marked "P" were previously cleared by 510(k) number K180636.

Additional Comments:

1. Includes imaging to assist in the placement of needles for prostate biopsy procedures.

Diagnostic Ultrasound Indications for Use Form – EV9C™ Transrectal End-fire Transducer

System	ExactVu™ High Resolution Micro-Ultrasound System						
Transducer	EV9C						
Intended Use	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
	Mode of Operation						
Clinical Application	B (2D Mode)	M	PWD	CWD	Color Doppler	Combined (specify)	Other (specify)
Ophthalmic							
Fetal							
Abdominal							
Intra-operative (Abdominal organs and vascular)							
Intra-operative (Neuro.)							
Laparoscopic							
Pediatric							
Small Organ (prostate)	P						P, 1
Neonatal Cephalic							
Adult Cephalic							
Transrectal	P						P, 1
Trans-vaginal							
Trans-urethral							
Trans-esoph. (non-Card.)							
Musculo-skel. (Convent.)							
Musculo-skel. (Superfic.)							
Intra-luminal							
Other (spec.)							
Cardiac Adult							
Cardiac Pediatric							
Trans-esophageal (card.)							
Other (spec.)							
Peripheral vessel							
Other (spec.)							
Dermatology							

N= new indication; P= previously cleared ; E= added under this appendix

Items marked "P" were previously cleared by 510(k) number K180636.

Additional Comments:

1. Includes imaging to assist in the placement of needles for prostate biopsy procedures.

Diagnostic Ultrasound Indications for Use Form – EV5C™ Abdominal Transducer

System	ExactVu™ High Resolution Micro-Ultrasound System						
Transducer	EV5C						
Intended Use	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
	Mode of Operation						
Clinical Application	B (2D Mode)	M	PWD	CWD	Color Doppler	Combined (specify)	Other (specify)
Ophthalmic							
Fetal							
Abdominal	N				N (2)		N (1)
Intra-operative (Abdominal organs and vascular)							
Intra-operative (Neuro.)							
Laparoscopic							
Pediatric							
Small Organ (prostate)							
Neonatal Cephalic							
Adult Cephalic							
Transrectal							
Trans-vaginal							
Trans-urethral							
Trans-esoph. (non-Card.)							
Musculo-skel. (Convent.)							
Musculo-skel. (Superfic.)							
Intra-luminal							
Other (spec.)							
Cardiac Adult							
Cardiac Pediatric							
Trans-esophageal (card.)							
Other (spec.)							
Peripheral vessel							
Other (spec.)							
Dermatology							

N= new indication; P= previously cleared ; E= added under this appendix

Additional Comments:

1. Includes imaging to assist in the placement of needles for kidney biopsy procedures.
2. ExactVu supports simultaneous color flow imaging with B-Mode.

Diagnostic Ultrasound Indications for Use Form – Needle guide for use with EV29L™ High Resolution Transrectal Side-fire Transducer

System	ExactVu™ High Resolution Micro-Ultrasound System						
Transducer	EV29L						
Intended Use	Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:						
	Mode of Operation						
Clinical Application	B (2D Mode)	M	PWD	CWD	Color Doppler	Combined (specify)	Other (specify)
Ophthalmic							
Fetal							
Abdominal							
Intra-operative (Abdominal organs and vascular)							
Intra-operative (Neuro.)							
Laparoscopic							
Pediatric							
Small Organ (prostate)	P						P, 1
Neonatal Cephalic							
Adult Cephalic							
Transrectal	P						P, 1
Trans-vaginal							
Trans-urethral							
Trans-esoph. (non-Card.)							
Musculo-skel. (Convent.)							
Musculo-skel. (Superfic.)							
Intra-luminal							
Other (spec.)							
Cardiac Adult							
Cardiac Pediatric							
Trans-esophageal (card.)							
Other (spec.)							
Peripheral vessel							
Other (spec.)							
Dermatology							

N= new indication; P= previously cleared ; E= added under this appendix

Items marked "P" were previously cleared by 510(k) number K180636.

Additional Comments:

1. Provides a mechanical means for performing needle / instrument guided procedures with the use of the diagnostic ultrasound transducer, EV29L.

Date: May 3, 2019

510(k) Summary

K190995

This summary of safety and effectiveness is provided as part of this Premarket Notification in compliance with 21 CFR, Part 807, Subpart E, Section 807.92.

1) Submitter's name, address, telephone number, contact person:

Exact Imaging, Inc.
7676 Woodbine Avenue, Unit 15
Markham, ON L3R 2N2
Canada

Corresponding Official: Randy AuCoin
President & CEO

Address: Exact Imaging, Inc.
7676 Woodbine Avenue, Unit 15
Markham, ON L3R 2N2, Canada

E-mail: raucoin@exactimaging.com

Telephone: (905) 415-0030

Facsimile: (905) 415-0031

Establishment Registration Number:

Exact Imaging, Inc.
Registration No. 3012402886
Owner/Operator No. 10053728

Date prepared: May 3, 2019

2) Name of the device, including the trade or proprietary name if applicable, the common or usual name, and the classification name, if known:

Common/ Usual Name

Diagnostic Ultrasound System with Accessories

Proprietary Name

ExactVu™ High Resolution Micro-Ultrasound System

Classification Names

Name	CFR Number	Product Code
Ultrasonic pulsed echo imaging system	892.1560	IYO
Diagnostic Ultrasound Transducers	892.1570	ITX
Biopsy Needle Guide	892.1560	OIJ
Ultrasonic pulsed doppler imaging system	892.1550	IYN

Table 1: ExactVu Classifications

Classification Panel

Radiology

3) Identification of the predicate or legally marketed device:

- ExactVu High Resolution Micro-Ultrasound System (K180636) (Primary Predicate)
- Ultrasound Scanner Flex Focus 1202 by BK Medical (K132335) (Reference Device) (used as a Reference Device to support differences in technological characteristics)

4) Device Description:

The ExactVu high resolution micro-ultrasound system is intended for use by qualified medical professionals for diagnostic ultrasound imaging, data processing and guidance of puncture and biopsy. The ExactVu System comprises transducers responsible for ultrasound signal generation and recording, needle guides and a main unit that controls the transducers, processes the acoustic data, and processes and displays images.

5) Intended Use /Indications for use:

The ExactVu micro-ultrasound system is intended for use by qualified medical professionals for diagnostic ultrasound imaging or fluid flow analysis of the human body. The indications for use (clinical applications) are:

- Small Organ
- Transrectal
- Abdominal

The system may be used with patients of all ages, but is not designed for pediatric or fetal use.

The system is contraindicated for direct cardiac application and for ophthalmic use or any application that causes the acoustic beam to pass through the eye.

6) Technological Characteristics:

The ExactVu High Resolution Micro-Ultrasound System ("ExactVu") that is the subject of this submission is substantially equivalent to its predicate device. The same fundamental scientific technology and general ultrasound principles that are used for ExactVu are also used on the predicate device, and both are Track 3 devices.

Primary differences between the subject and predicate device are support for an abdominal transducer, the addition of color flow imaging modes and the addition of a DICOM service class.

Comparisons of relevant significant features are presented in two tables below:
Table 2 for a comparison between the Subject Device and the Predicate Device and Table 3 for a comparison between the Subject Device and the Reference Device.

Feature	Subject Device: ExactVu 2.5	Predicate Device: ExactVu 2.0 (K180636)
Manufacturer	Exact Imaging Inc.	Exact Imaging Inc.
Intended use / Indications for Use	The ExactVu micro-ultrasound system is intended for use by qualified medical professionals for diagnostic ultrasound imaging or fluid flow analysis of the human body. The indications for use (clinical applications) are: • Small Organ • Transrectal • Abdominal The system may be used with patients of all ages, but is not designed for pediatric or fetal use. The system is contraindicated for direct cardiac application and for ophthalmic use or any application that causes the acoustic beam to pass through the eye.	The ExactVu High Resolution Micro-Ultrasound System is intended for use by qualified medical professionals for diagnostic ultrasound imaging or fluid flow analysis of the human body. The indications for use (clinical applications) are: • Small Organ (prostate) • Transrectal The system may be used with patients of all ages, but is not designed for pediatric or fetal use.
Modes of Operation	B-mode (2-D Grayscale Imaging, Transverse, CFM and combinations - Color Doppler and Power Doppler Mode)	B-mode (2-D Grayscale Imaging, Transverse, CFM and combinations
MI Indication	Yes	Yes
TI Indication	Yes	Yes

Feature	Subject Device: ExactVu 2.5	Predicate Device: ExactVu 2.0 (K180636)
Electrical Safety	TUV, IEC 60601-1	TUV, IEC 60601-1
510(k) Track	Track 3	Track 3
Center Frequency Range	EV9C: 6.5 MHz EV29L: 22.5 MHz EV5C: 3.5 MHz	EV9C: 6.5 MHz EV29L: 22.5 MHz
Transducer Types	Endocavity Linear Endocavity Curved Array EV5C: Abdominal Curved Array	Endocavity Linear Endocavity Curved Array
Measurement	Manual measurements are available on the ultrasound system, including distance, area approximations and volume calculations.	Manual measurements are available on the ultrasound system, including distance, area approximations and volume calculations.
#Transmit Channels	128 channels	128 channels
#Receive Channels	128 channels	128 channels
DICOM	DICOM Storage SCU Version 3.0 DICOM Storage Commitment SCU Version 3.0 DICOM Query Retrieve SCU Version 3.0 ExactVu reads MRI study data in DICOMDIR format, where MRI markup uses DICOM GSPS (Grayscale Softcopy Presentation State)	DICOM Storage SCU Version 3.0 DICOM Storage Commitment SCU Version 3.0 ExactVu reads MRI study data in DICOMDIR format, where MRI markup uses DICOM GSPS (Grayscale Softcopy Presentation State)

Table 2: Technological Characteristics Comparison Between Subject Device and Predicate Device

Feature	Subject Device: ExactVu 2.5	Reference Device: Ultrasound Scanner Flex Focus 1202 (K132335)
Manufacturer	Exact Imaging Inc.	BK Medical
Intended Use / Indication for Use	The ExactVu micro-ultrasound system is intended for use by qualified medical professionals for diagnostic ultrasound imaging or fluid flow analysis of the human body. The indications for use (clinical applications) are:	Ultrasound scanner and transducers for B, Tissue and Contrast Harmonic Imaging, M, PWD, CWD, Color Doppler, Vector Flow Imaging and combined mode imaging. Signal analysis and display.

Feature	Subject Device: ExactVu 2.5	Reference Device: Ultrasound Scanner Flex Focus 1202 (K132335)
	• Small Organ • Transrectal • Abdominal The system may be used with patients of all ages, but is not designed for pediatric or fetal use. The system is contraindicated for direct cardiac application and for ophthalmic use or any application that causes the acoustic beam to pass through the eye.	Guidance of biopsy needles, geometrical measurements and calculation of parameters. An optional 3-DI unit can reconstruct a series of 2-DI images into a single 3-DI volume and display this on the screen. An optional Vector Flow Imaging (VFI) module: Color Flow Mapping (CFM) imaging mode with the ability to visualize both the axial and the transverse velocity. An optional RF wireless function with the ability to wireless transmit for printing and archive connectivity purpose. Clinical applications: Abdominal, Cardiac, Fetal, Intraoperative, Neurosurgery, Obstetrics, Pediatrics, Transrectal, Small organs, Transvaginal, Musculoskeletal
Abdominal Transducer	Abdominal Curved Array (transducer model EV-5C) Center Frequency: 3.5 MHz Maximum depth: 180 mm Global Maximum Outputs/Worst Case Settings • 2D mode ◦ Ispta.3: 11.93 mW/cm2 ◦ TI Type: TIS ◦ TI Value: <1.0 ◦ MI: 1.04 ◦ Ipa.3@MI_Max: 117.8 W/cm2 • CFI mode ◦ Ispta.3: 42.84 mW/cm2 ◦ TI Type: TIS ◦ TI Value: 1.00 ◦ MI: 0.91	Abdominal Curved Array (transducer model 8823) Center Frequency: 1.8-6 MHz Maximum depth: 171-260mm depending on settings Global Maximum Outputs/Worst Case Settings • Ispta.3: 720 mW/cm2 • TI Type: TIC • TI Value: 6.0 • MI: 1.9 • Ipa.3@MI_Max: 367 W/cm2 Patient contact material: not available

Feature	Subject Device: ExactVu 2.5	Reference Device: Ultrasound Scanner Flex Focus 1202 (K132335)
	○ Ipa.3@MI_Max: 79.29 W/cm2 Patient contact material: Epotek 301, PEEK LSG, Polyurethane, Rexolite, Silicon	
Modes of Operation	Color flow imaging maps fluid flow velocities over the 2D image using color. ExactVu provides two color flow imaging modes ("CFI Modes"): • Color Doppler Mode • Power Doppler Mode	• Color mode (CFM, color flow mapping, color Doppler) ultrasound displays color-coded real-time information about direction and velocity of flow in the tissues. • Power mode displays color-coded information about the amount of flow but not the direction.
DICOM	DICOM Storage SCU Version 3.0 DICOM Storage Commitment SCU Version 3.0 DICOM Query Retrieve SCU Version 3.0 ExactVu reads MRI study data in DICOMDIR format, where MRI markup uses DICOM GSPS (Grayscale Softcopy Presentation State)	DICOM, Print, Modality Worklist, Storage Commitment, Modality Performed Procedure Steps (MPPS)

Table 3: Technological Characteristics Comparison Between Subject Device and Reference Device

7) Determination of Substantial Equivalence:

Summary of Non-Clinical Tests:

The ExactVu System has been evaluated for electrical, thermal, mechanical and EMC safety. Additionally, cleaning/disinfection, biocompatibility, and acoustic output have been evaluated, and the device has been found to conform to applicable mandatory medical device safety standards. Assurance of quality was established by employing the following elements of product development: Design Phase Reviews, Risk Assessment, Requirements Development, System and Software Verification, Hardware Verification, Risk Mitigation Verification, Clinical Evaluation.

All patient contact materials are biocompatible and are materials that are already used in other legally marketed devices or meet 10993-1.

The ExactVu System is designed to comply with the following voluntary standards:

Reference No	Recognition No	Title
AAMI/ANSI/ISO 10993-1	2-220	ISO 10993-1:2009, Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process
IEC 60601-1	19-4	AAMI / ANSI ES60601-1:2005/(R) 2012 and A1:2012, c1:2009/(R) 2012 and a2:2010/(R) 2012. Medical electrical equipment - part 1: general requirements for basic safety and essential performance (IEC 60601-1:2005, mod)
IEC 60601-1-2	19-1	IEC 60601-1-2:2014, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests (Edition 3)
IEC 60601-2-18	9-91	IEC 60601-2-18: Edition 3.0 2009-08, Medical Electrical Equipment - Part 2-18: Particular Requirements For The Basic Safety And Essential Performance Of Endoscopic Equipment.
IEC 60601-2-37	12-209	IEC 60601-2-37:2015, Particular Requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment
NEMA UD 2-2004	12-105	Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment

Table 4: Standards with which ExactVu Complies

Summary of Clinical Tests:

The ExactVu System, transducers and needle guides, subject of this submission, did not require clinical studies to support the determination of substantial equivalence. An analysis of clinical literature, adverse event and recall data is provided in the clinical evaluation justifies the safety and efficacy of the functionality of the ExactVu device without pre-market clinical investigation.

8) Conclusion:

Intended uses and other key features are consistent with traditional clinical practice and FDA guidance. The ExactVu device and predicate device conform to applicable

electrical medical device safety standards with compliance verified through independent evaluation and meet FDA requirements for Track 3 devices, have biosafety equivalence and are manufactured using a comparable quality system that has been updated to meet ISO 13485:2016 and the requirements for MDSAP. Exact Imaging, Inc. believes that the ExactVu system is substantially equivalent with regard to safety and effectiveness to the predicate device.