



September 5, 2019

BroadMaster Biotech, Corp.
% Dr. Ke-Min Jen
Contact Person
Chinese-European Industrial Research Society
No. 58, Fu-Chiun St
Hsin-Chu City, 30067
Taiwan

Re: K191004

Trade/Device Name: Advocate Non-Contact Infrared Thermometer, Model EF001S
Regulation Number: 21 CFR 880.2910
Regulation Name: Clinical Electronic Thermometer
Regulatory Class: Class II
Product Code: FLL
Dated: July 19, 2019
Received: August 6, 2019

Dear Dr. Ke-Min Jen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmnm.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Alan Stevens

Assistant Director

DHT3C: Division of Drug Delivery and
General Hospital Devices,
and Human Factors

OHT3: Office of Gastrorenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K191004

Device Name
Advocate Non-Contact Infrared Thermometer, EF001S

Indications for Use (Describe)

Advocate Non-Contact Infrared Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in clinic as reference. It is intended for measuring human body temperature of all ranges of people by detecting infrared heat from the forehead.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



BroadMaster Biotech, Corp.

1F., 2F., No. 91, Xiyuan Rd., Zhongli Dist., Taoyuan City
32057, Taiwan (R.O.C.)
Tel: +886-3-451-7600
Fax: +886-3-451-9500
Website: www.broadmaster-biotech.com

510(k) Summary (Per 21 CFR 807.92)

510(k) number: K191004

Trade or proprietary name	Advocate Non-Contact Infrared Thermometer, model EF001S
Common Name	Digital Thermometer
Classification Name	Clinical Electronic Thermometer 21 CFR 880.2910
Class	II
Panel	80 General Hospital
Product Code	FLL
Owner/Operator	BroadMaster Biotech, Corp. 1F., 2F., No. 91, Xiyuan Rd., Zhongli Dist., Taoyuan City 32057, Taiwan (R.O.C.) Tel: +886-3-451-7600 Fax: +886-3-451-9500 Website: www.broadmaster-biotech.com
Date prepared	July 19, 2019
Application Correspondent	Dr. Jen, Ke-Min ROC Chinese-European Industrial Research Society No. 58, Fu-Chiun St., Hsin-Chu City, 30067, Taiwan Tel: +886-3-5208829 Fax: +886-3-5209783 Email: ceirs.jen@msa.hinet.net
Predicate Device	Manufacturer: BroadMaster Biotech, Corp. Product name: Advocate Non-Contact Infrared Thermometer Model No: EF001A 510(k) number: K180355


● **Indications for Use:**

Advocate Non-Contact Infrared Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in clinic as reference. It is intended for measuring human body temperature of all ranges of people by detecting infrared heat from the forehead.

- **Descriptions:**

The Advocate Non-Contact Infrared Thermometer, model EF001S, measures temperatures of people by detecting the infrared energy radiated directly from the forehead without physical contact. The device is composed of a Probe of metals with infrared sensor inside to detect the infrared energy, an LCD Display, a SCAN button to start measuring temperatures, a Power switch button to switch on or off the device, a Body/Ambient button to switch between two measuring modes, and an Enclosure of ABS, , as shown in the following diagram.

The device has the following features: one-second measuring time, measuring Body or Ambient temperature, 12-memory recalls, °F/°C unit switchable, over range message (Hi/Lo), low battery indication, auto display for the last reading when power is on, auto shut-off when the device is idle for 60 seconds and voice function. When the device starts, “Please measure” in English or “Mira la temperatura” in Spanish will be heard. When completes, the result will be heard in addition to the data display. The Advocate Non-Contact Infrared Thermometer, model EF001S, including the voice function, is not intended for use by the visually impaired individuals.



Components of EF001S

- **Principle Operation**

The Advocate® Non-Contact Infrared Thermometer measures temperatures of people by detecting the infrared energy radiated directly from the forehead without physical contact. As soon as the distance between the probe and the forehead is within 2 -3.94 inches(5-10 cm), the IR radiation sensor is activated, and the measurement will be taken instantly by detection of the infrared heat.



BroadMaster Biotech, Corp.

1F., 2F., No. 91, Xiyuan Rd., Zhongli Dist., Taoyuan City
32057, Taiwan (R.O.C.)
Tel:+886-3-451-7600
Fax:+886-3-451-9500
Website: www.broadmaster-biotech.com

● **Substantial Equivalence Comparison Table**

Comparison Items	Predicate Device	Subject Device	Remarks
Applicant	<i>BroadMaster Biotech, Corp.</i>	<i>BroadMaster Biotech, Corp.</i>	Same
Proprietary (Trade) Name	<i>Advocate Non-Contact Infrared Thermometer</i>	<i>Advocate Non-Contact Infrared Thermometer</i>	Same
Model	EF001A	EF001S	Series model
510(k) No.	K180355	K191004	--
Classification Name and Regulation Number	thermometer, electronic, clinical 21 CFR 880.2910	<i>thermometer, electronic, clinical 21 CFR 880.2910</i>	Same
Product Code	FLL	FLL	Same
Device Class	II	II	Same
Intended uses	Advocate Non-Contact Infrared Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in clinic as reference. It is intended for measuring human body temperature of all ranges of people by detecting infrared heat from the forehead.	<i>Advocate Non-Contact Infrared Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in clinic as reference. It is intended for measuring human body temperature of all ranges of people by detecting infrared heat from the forehead.</i>	Same
Intended users	Lay user and professional	<i>Lay user and professional</i>	Same
Measurement method	Infrared radiation detection	<i>Infrared radiation detection</i>	Same
Measurement mode	Forehead measurement mode	<i>Forehead measurement mode</i>	Same
Measuring range	Body measurement mode: 89.6°F to 109.4°F(32°C - 43 °C)	<i>Body measurement mode: 89.6°F-109.4°F (32°C - 43 °C)</i>	Same
Display resolution	0.1°F/ 0.1°C	<i>0.1°F / 0.1 °C</i>	Same
C/F unit switchable	Yes	Yes	Same
Measuring accuracy	<i>Body measurement mode: ±0.5 °F / 0.3 °C (93.2 °F -94.8 °F, 34.0 °C-34.8°C) ±0.4 °F/0.2°C (95.0 °F -107.6 °F, 35.0°C-42.0 °C) ±0.5°F/0.3 °C (107.8 °F-109.4 °F, 42.1°C-43.0 °C)</i>	<i>Body measurement mode: ±0.5 °F / 0.3 °C (93.2 °F -94.8 °F, 34.0 °C-34.8°C) ±0.4 °F/0.2°C (95.0 °F -107.6 °F, 35.0°C-42.0 °C) ±0.5°F/0.3 °C (107.8 °F-109.4 °F, 42.1°C-43.0 °C)</i>	Same
Display	LCD display	<i>LCD display</i>	Same
Measurement distance	2-3.94 inch (5-10 cm)	<i>2-3.94 inch (5-10 cm)</i>	Same
Memory set	12 sets	<i>12 sets</i>	Same
Power source	Two 1.5V AAA alkaline batteries	<i>Two 1.5V AAA alkaline batteries</i>	Same
Low battery indication	Yes	Yes	Same



BroadMaster Biotech, Corp.

1F., 2F., No. 91, Xiyuan Rd., Zhongli Dist., Taoyuan City
 32057, Taiwan (R.O.C.)
 Tel:+886-3-451-7600
 Fax:+886-3-451-9500
 Website: www.broadmaster-biotech.com

Degree of protection	IP20	IP20	Same
Operating condition	50.0 °F- 104.0 °F (10.0 °C-40.0 °C) ≤ 80% RH	50.0 °F- 104.0 °F (10.0 °C-40.0 °C) ≤ 80% RH	Same
Storage condition	-13.0 °F-131 °F (-25.0 °C-55.0 °C) ≤ 95% RH	-13.0 °F-131 °F (-25.0 °C-55.0 °C) ≤ 95% RH	Same
Cleaning method	The thermometer enclosure and probe can be cleaned and disinfected by 70% alcohol.	The thermometer enclosure and probe can be cleaned and disinfected by 70% alcohol.	Same
Human -contacting materials	Enclosure of red & black ABS, LCD Lens of PMMA and Probe of Metals	Enclosure of red & black ABS, LCD Lens of PMMA and Probe of Metals	Same
Biocompatibility	EN ISO 10993-5:2009 & ISO 10993-10:2010	Complying EN ISO 10993-1:2009	Same
Software	Software validation report	Revised software validation report	Different
Voice function	No	Yes, extra 5 voice hardware: U3, U3A, C10 (10 uF), C13(0.1 uF) and Speaker	Different
EMC	IEC 60601-1-2:2014 FCC 47 CFR Part 18, Subpart B	IEC 60601-1-2:2014 FCC 47 CFR Part 18, Subpart B	Same
Electrical Safety	IEC 60601-1: 2005/A1:2012	ANSI AAMI ES60601-1:2005 IEC 60601-1: 2005/A1:2012	Same
Performance	ASTM E1965-98(2016) ISO 80601-2-56:2017	*Declaration of conformity to ASTM E1965-98(2016) *Declaration of conformity to ISO 80601-2-56:2017	Same
PCB	FR4 PCB	FR4 PCB	Same
Materials	Patient contacting materials is ABS (device housing / handle and button).	Patient contacting materials is ABS (device housing / handle and button).	Same
MCU	HYCON HY11P13 8-Bit RISC-like Mixed Signal Microcontroller Embedded 4x20 LCD Driver Low Noise Amplifier 18-Bit ADC	HYCON HY11P13 8-Bit RISC-like Mixed Signal Microcontroller Embedded 4x20 LCD Driver Low Noise Amplifier 18-Bit ADC	Same
Sensor	GE THERMOPILE IR SENSOR ZTP-148SR	GE THERMOPILE IR SENSOR ZTP-148SR	Same
LCD	TN LCD	TN LCD	Same
Speaker	Buzzer	Speaker	Different

● **Comparison discussion**

Because the subject device is almost identical to the predicate device, except for the voice function, we will discuss the differences raised by the existence of voice function for the subject device.

The software is different due to the addition of the voice function. However, software validation was performed and this difference of software did not raise any new or different questions of safety and effectiveness for the subject device.

The extra 5 hardware components used by EF001S due to addition of the voice function are voice ICs:U3, U3A, C10 (10 uF), C13 (0.1 uF) and Speaker. Thus, we subjected the subject



BroadMaster Biotech, Corp.

1F., 2F., No. 91, Xiyuan Rd., Zhongli Dist., Taoyuan City
 32057, Taiwan (R.O.C.)
 Tel:+886-3-451-7600
 Fax:+886-3-451-9500
 Website: www.broadmaster-biotech.com

device EF001S to meeting the provisions of Safety and EMC standards, i.e., ANSI AAMIES60601-1:2005, IEC 60601-1-2:2014, and FCC 47 CFR Part 18, Subpart B. This voice function is not intended for use by vision-impaired patients, and it is an auxiliary function for the lay users and medical professionals. So the voice function brings some convenience to the lay users and professionals. Regarding the extra 5 hardware components used by the subject device, there is no different or new safety and effectiveness questions.

● **Non-Clinical Testing**

Testing name	Referenced standard	Summary result	Verdict
Electric safety testing	ANSI AAMIES60601-1:2005 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance FDA recognition number: 19-4	The subject device complies with the applicable requirements set forth in the referenced electric safety standard, ANSI AAMIES60601-1:2005.	Pass
EMC testing	IEC60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances- Requirements and tests, FDA recognition number: 19-8	The subject device complies with the applicable requirements set forth in the referenced EMC standard, IEC 60601-1-2:2014.	Pass
	FCC 47 CFR Part 18: Industrial, Scientific, And Medical Equipment, Subpart B:Applications and Authorizations	The subject device complies with the applicable requirements set forth in the referenced EMC standard, FCC 47 CFR Part 18.	Pass
Performance testing	ISO 80601-2-56: 2017. Medical electrical equipment – Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement FDA recognition number: 6-403	The subject device complies with the applicable requirements set forth in the referenced performance standard, ISO 80601-2-56:2017.	Pass
Biocompatibility testing	EN ISO 10993-1:2018 Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process FDA recognition number: 2-258	The subject device EF001S has the same human-contacting materials as predicate device EF001A. There is no need to proceed the biocompatibility evaluation according to EN ISO 10993-1:2018.	Pass



BroadMaster Biotech, Corp.

1F., 2F., No. 91, Xiyuan Rd., Zhongli Dist., Taoyuan City
 32057, Taiwan (R.O.C.)
 Tel:+886-3-451-7600
 Fax:+886-3-451-9500
 Website: www.broadmaster-biotech.com

Revised Software Validation Report	Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices, issued on May 11,2005	The software contained in the subject device complies with the applicable requirements set forth in the referenced guidance document, "Guidance for the Content of Premarket Submissions for Software Contained, issue done May 11, 2005.	Pass
------------------------------------	--	---	------

● **Clinical Testing**

Name of clinical testing	Referenced standard	Summary of testing	Patient population (age groups, number of subjects)	Verdict
EF001 Clinical Accuracy	*ASTM E1965-98(2016) Standard Specification for Infrared Thermometers for Intermittent Determination of Patient Temperature FDA recognition number: 6-125	The methods and criteria of EF001 Clinical Test had been clinically assessed to meet the requirements of clinical accuracy per the referenced standards.	40 subjects in each age group, infants (0-1 year), children (1-5 years) and adults (>5 years) (Total 120 subjects)	Pass

● **Conclusion**

Non-clinical performance and clinical tests were conducted on the subject device and all tests met specified criteria. Based on the information provided in this submission the subject device, Advocate Non-Contact Infrared Thermometer, EF001S, is substantially equivalent to the predicate device, EF001A.