



February 21, 2020

Coloplast A/S  
Lykke Forchhammer  
Director of RA  
Holtedam 1  
Humblebaek, 3050 Dk

Re: K191536

Trade/Device Name: Biatain Silicone Ag  
Regulatory Class: Unclassified  
Product Code: FRO  
Dated: January 16, 2020  
Received: January 21, 2020

Dear Lykke Forchhammer:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see

<https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Kimberly M. Ferlin, Ph.D.  
Assistant Director (Acting)  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K191536

Device Name

Biatain Silicone Ag

Indications for Use (Describe)

Biatain® Silicone Ag is intended to provide a moist wound environment and exudate management of acute and chronic wounds.

Biatain® Silicone Ag is indicated for the management of exuding leg and foot ulcers, pressure ulcers, diabetic foot ulcers, superficial and partial thickness burns, donor sites, and, traumatic and post-operative wounds.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 5. TRADITIONAL 510(K) SUMMARY

**Submitted by:** Coloplast A/S  
Holtedam 1  
3050 Humlebaek  
Denmark

**Contact Person:** Lykke Forchhammer  
Director of RA Innovation and Market Expansion  
Coloplast A/S  
Holtedam 1  
3050 Humlebaek  
Denmark  
Telephone: 0045-49111271  
Email: [dklfo@coloplast.com](mailto:dklfo@coloplast.com)

**Date of Summary:** 21 February 2020

**Trade or Proprietary Name:** Biatain<sup>®</sup> Silicone Ag

**Common or Usual Name:** Dressing, Wound, Drug

**Classification and CFR:** Unclassified

**Product Code:** FRO

**Review Panel:** General and Plastic Surgery

**Predicate Device:** K100218, Biatain<sup>®</sup> Ag (Coloplast A/S)

**Reference Device:** K120828, Biatain<sup>®</sup> Silicone Ag (Coloplast A/S)

**Device Description:**

The device is available as sterile, single use only. The device is intended to provide a moist environment and exudate management of acute and chronic wounds. The device maintains a moist wound providing the optimal environment for wound healing.

The device consists of:

- Polyurethane (PU) top film – a 25 µm semi-permeable barrier, printed with grey dots.
- Absorbing pad - composed of a 3.0 mm PU foam containing an anti-bacterial silver complex (silver sodium hydrogen zirconium phosphate).
- Silicone adhesive bi-layer - composed of a middle PU film with a silicone adhesive gel. The PU top film and the middle PU film are laminated together at the border.
- Protective film – composed of three or five parts; the center part and the remaining side parts.

The device:

- contains a silver compound (silver sodium hydrogen zirconium phosphate). In vitro testing demonstrated that Biatain® Silicone Ag has an antimicrobial effect in the dressing against the following strains of 3 gram positive bacteria, 3 gram negative bacteria, 1 yeast, and 1 mold: *Staphylococcus aureus*, *Enterococcus faecalis*, *Streptococcus pyogenes*, *Proteus mirabilis*, *Klebsiella pneumoniae*, *Escherichia coli*, *Aspergillus brasiliensis* and *Meyerozyma guilliermondii*. The dressing sustains antimicrobial activity for up to 7 days.

The device may be used in hospitals, healthcare facilities, and home care.

Sterilized using ethylene oxide (EO).

Prescription Use Only.

**Indications for Use:**

Biatain<sup>®</sup> Silicone Ag is intended to provide a moist wound environment and exudate management of acute and chronic wounds.

Biatain<sup>®</sup> Silicone Ag is indicated for the management of exuding leg and foot ulcers, pressure ulcers, diabetic foot ulcers, superficial and partial thickness burns, donor sites, and, traumatic and post-operative wounds.

**Technological  
Characteristics:**

The subject device is substantially equivalent to the predicate device based upon the information below.

.

	<b>Subject Device Biatain® Silicone Ag</b>	<b>Predicate Device Biatain® Ag</b>	<b>Reference Device: Biatain® Silicone Ag</b>
<b>510(k) Submitter</b>	Coloplast A/S	Same	Same
<b>510(k) Number</b>	K191536	K100218	K120828
<b>Regulatory Class</b>	Unclassified	Same	Same
<b>Classification Name</b>	Dressing, Wound, Drug	Same	Same
<b>Product Code</b>	FRO	Same	Same
<b>Intended Use</b>	Biatain® Silicone Ag is intended to provide a moist wound environment and exudate management.	Same	Same

	<b>Subject Device</b> <b>Biatain® Silicone Ag</b>	<b>Predicate Device</b> <b>Biatain® Ag</b>	<b>Reference Device:</b> <b>Biatain® Silicone Ag</b>
<b>Indications for Use</b>	<p>Biatain® Silicone Ag is intended to provide a moist wound environment and exudate management of acute and chronic wounds.</p> <p>Biatain® Silicone Ag is indicated for the management of exuding leg and foot ulcers, pressure ulcers, diabetic foot ulcers, superficial and partial thickness burns, donor sites, and, traumatic and post-operative wounds.</p>	<p>Biatain Ag foam dressings are indicated for wounds with moderate to high amounts of exudate, including leg ulcers and Category II-IV pressure ulcers with delayed healing due to bacteria, or where there is a risk of infection. Biatain Ag foam dressings may be used for second-degree burns, donor sites, postoperative wounds and skin abrasions with delayed healing due to bacteria, or where there is a risk of infection. Biatain Ag non-adhesive foam dressings are additionally indicated for diabetic foot ulcers. Biatain Ag foam dressings may reduce odor caused by micro-organisms in the wound. Biatain Ag foam dressings may be used to support moist wound healing on patients who are in treatment for a local or systemic infection under the discretion of a healthcare professional. Depending on the prognosis of the wound, Biatain Ag foam dressings may be used throughout the healing process to provide protection for the indicated types of wounds. Biatain Ag foam dressings are suitable for use in combination with compression therapy. Biatain Ag foam filler/cavity is indicated for deep exuding wounds, including leg ulcers, pressure ulcers, diabetic foot ulcers and acute wounds.</p>	<p>Biatain® Silicone Ag Foam Dressings are indicated for the management of moderately to highly exuding leg ulcers and pressure sores. The dressing can also be used for 2nd degree burns, donor sites, postoperative wounds, and skin abrasions.</p>
<b>Prescription Device</b>	Yes	Same	Same
<b>Single Use</b>	Yes	Same	Same
<b>Sterile Device</b>	Yes	Same	Same
<b>Method of Sterilization</b>	EO SAL 10 <sup>-6</sup>	Irradiated SAL 10 <sup>-6</sup>	Same (EO) SAL 10 <sup>-6</sup>

	<b>Subject Device Biatain® Silicone Ag</b>	<b>Predicate Device Biatain® Ag</b>	<b>Reference Device: Biatain® Silicone Ag</b>
<b>Recommended Wear Time</b>	Up to 7 days. Repeated use of the device should not exceed 30 days	Same	Same
<b>Shelf Life</b>	2-years	3-years	3-years
<b>Silver (Ag)</b>	Silver Compound (Ionic) [silver sodium hydrogen zirconium phosphate]	Same	Same
<b>Silver Content</b>	0.95 mg/cm <sup>2</sup>	Same	Same
<b>Silver Release (within the dressing)</b>	24 h: 207 – 367 µg/cm <sup>2</sup> 48 h: 350 – 589 µg/cm <sup>2</sup> 72 h: 439 – 726 µg/cm <sup>2</sup> 96 h: 500 – 813 µg/cm <sup>2</sup> 168 h: 591 – 937 µg/cm <sup>2</sup>	Same	Same

	<b>Subject Device</b> <b>Biatain® Silicone Ag</b>	<b>Predicate Device</b> <b>Biatain® Ag</b>	<b>Reference Device:</b> <b>Biatain® Silicone Ag</b>
<b>Antimicrobial Effectiveness within the Dressing</b>	<p>Fulfills both <math>\geq 4</math> log reduction (AATCC 100) as described below:</p> <p>Fulfills <math>\geq 4</math> log reduction (AATCC 100) for antimicrobial effect against 3 gram positive bacteria, 3 gram negative bacteria, 1 yeast, and 1 mold:  <i>Staphylococcus aureus</i>, <i>Enterococcus faecalis</i>, <i>Streptococcus pyogenes</i>, <i>Proteus mirabilis</i>, <i>Klebsiella pneumoniae</i>, <i>Escherichia coli</i>, <i>Aspergillus brasiliensis</i> and <i>Meyerozyma guilliermondii</i>.</p>	<p>Effective antimicrobial activity within the dressing.</p>	<p>Effective antimicrobial activity within the dressing.</p>
<b>Primary Packaging</b>	Paper/Polyester Film	<p>Biatain® Ag Non-adhesive: White Polyester Film/Transparent Polyester Film</p> <p>Biatain® Ag Adhesive: Foil</p>	Tyvek/Transparent Foil

	<b>Subject Device</b> <b>Biatain® Silicone Ag</b>	<b>Predicate Device</b> <b>Biatain® Ag</b>	<b>Reference Device:</b> <b>Biatain® Silicone Ag</b>
<b>Top Film</b>	Vapor permeable polyurethane (PU) top film	Same	Same
<b>Pigment: Top Film Colorant</b>	CI Pigment Black 7 Titanium dioxide Aluminum hydroxide	Black iron oxide Titanium dioxide C.I Pigment Black 7	Black iron oxide Titanium dioxide
<b>Pigment: Top Film Dots</b>	C.I. Pigment White 6 C.I. Pigment Black 7 Titanium chelates Aluminum oxide	Not Applicable	Not Applicable
<b>Pigment: Top Film Logo</b>	CI Pigment Black 7	Not Applicable	Not Applicable
<b>Absorbaent Polyurethane Foam with Silver Compound</b>	Polyurethane (PU) foam with silver complex (silver sodium hydrogen zirconium phosphate)	Same	Same and contains hotmelt Pressure Sensitive Adhesive (PSA)
<b>Polyurathane Foam Thickness</b>	3.0 mm PU foam with silver complex	3.0 mm and 4.4 mm PU foam with silver complex	Same (3.0 mm PU foam with silver complex)
<b>Adhesive Bilayer</b>	Silicone adhesive bi-layer across the device with perforations across the absorbent pad	Biatain® Ag Non-adhesive Dressings have no adhesive layer. Biatain® Ag Adhesive: Hydrocolloid adhesive border	Perforated Silicone adhesive bi-layer on the border
<b>Protective Film (Release Liner)</b>	Polypropylene (PP) 3 and 5-piece design	Polyethylene (PE) 2 and 3-piece design	Polyethylene (PE) 3-piece design
<b>Device Sizes and Shapes.</b>	Square: 7.5 x 7.5 cm 10 x 10 cm 12.5 x 12.5 cm 15 x 15 cm 17.5 x 17.5 cm	Square: 10 x 10 cm (non-adhesive) 12.5 x 12.5 cm (adhesive) 15 x 15 cm (non-adhesive) 18 x 18 cm (adhesive) 20 x 20 cm (non-adhesive)	Square: 7.5 x 7.5 cm 10 x 10 cm 12.5 x 12.5 cm

	<b>Subject Device Biatain® Silicone Ag</b>	<b>Predicate Device Biatain® Ag</b>	<b>Reference Device: Biatain® Silicone Ag</b>
	Rectangle: 10 x 20 cm 10 x 30 cm	Rectangle: 5 x 7 cm 10 x 20 cm	Not Applicable
	Heel: 18 x 18 cm	Heel: 19 x 20 cm	Not Applicable
	Sacral: 15 x 19 cm 25 x 25 cm	Sacral: 23 x 23 cm	Not Applicable

<p><b>Performance Data:</b></p>	<p>Performance testing for the Biatain® Silicone Ag device was conducted according to applicable sections of voluntary standards in order to document the following properties of the subject device. The proposed changes do not impact the performance specifications:</p> <ul style="list-style-type: none"> <li>• Real Time and Accelerated Aged shelf life testing [per ASTM F1980-07 (2011)]</li> <li>• Packaging transportation and integrity testing per ASTM D4169-16</li> <li>• Biocompatibility according to ISO 10993-1 (2009)</li> <li>• Free swelling absorptive capacity, wear time, border permeability, and Fluid Handling Capacity (absorption and moisture vapour transmission rate) per EN 13726-1:2002/AC2003</li> <li>• Waterproofness per EN 13726-3:2002</li> <li>• Dynamic Friction of the PU Film per EN ISO 8295 1 2004</li> <li>• Conformability to body per EN 13726-4:2003</li> <li>• Ease of release - peel adhesion per ASTM D3330/D3330M-04 Method A</li> <li>• Antimicrobial effectiveness per modified ATCC 100:2012</li> </ul> <p>All tests passed the pre-determined acceptance criteria.</p>
<p><b>Clinical Data</b></p>	<p>No clinical data was required to support substantial equivalence.</p>
<p><b>Substantial Equivalence Conclusion:</b></p>	<p>Based on the same intended use with no new patient population or wound type, similar technological characteristics and materials, and performance testing, Coloplast believes the proposed Biatain® Silicone Ag is substantially equivalent to the predicate device Biatain® Ag (K100218), and the reference device, Biatain® Silicone Ag (K120828).</p>