



August 16, 2019

Konica Minolta, Inc.
% Mr. Russell D. Munves
Official Correspondent
Storch Amini PC
140 East 45th Street
25th Floor
NEW YORK NY 10017

Re: K191645
Trade/Device Name: SKR 4000
Regulation Number: 21 CFR 892.1680
Regulation Name: Stationary x-ray system
Regulatory Class: Class II
Product Code: MQB
Dated: June 17, 2019
Received: June 19, 2019

Dear Mr. Munves:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For

Thalia T. Mills, Ph.D.
Director
Division of Radiological Health
OHT7: Office of In Vitro Diagnostics
and Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K191645

Device Name

SKR 4000

Indications for Use (Describe)

The SKR 4000 is indicated for use in generating radiographic images of human anatomy. It is intended to replace a radiographic film/screen system in general purpose diagnostic procedures.

The SKR 4000 is not indicated for use in mammography, fluoroscopy, and angiography applications.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



K191645

510(k) Summary

Submitter's Name: KONICA MINOLTA, INC.

Address: 1 Sakura-machi,
Hino-shi, 191-8511 Japan

Contact: Tsutomu Fukui

Telephone: +81 42 589 8429

Date: July 26, 2019

Trade Name: SKR 4000

Model No: P-41

Common Name: Digital Radiography

Regulation Name / Number: Stationary x-ray system / 21 CFR 892.1680

Regulatory Class: Class II

Product Code(s): 90-MQB

Predicate Device(s): K160810 - ViZionDR + Wireless (Viztek LLC)
Regulation Name: Stationary x-ray system
Regulation Number: 21CFR 892.1680
Product Codes: 90-MQB
K182688 - SKR 3000 (KONICA MINOLTA, INC.)
Regulation Name: Stationary x-ray system
Regulation Name: 21CFR 892.1680
Product Codes: 90-MQB, 90-LLZ

Device Description:

The SKR 4000 consists of the DR Detector (P-41), Battery, Battery Charger, AC Adapter, AC Power Cable and DC Power Cable. It is intended for use replacing a radiographic film/screen system in general-purpose diagnostic procedures of human anatomy.



The SKR 4000 performs radiography imaging of the human body using an X-ray FPD that outputs a digital signal, which is then input into an image processing device.

The SKR 4000 can communicate with the image processing device through the Wireless LAN (IEEE802.11a/b/g/n/ac and FCC compliant). The WPA2-PSK (AES) encryption is adopted for a security of wireless connection.

The SKR 4000 is designed to comply with the following standards; AAMI/ANSI ES 60601-1 “Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance”, IEC 60601-1-2 “Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests”, and ISO 10993-1 “Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process”.

Indications for Use:

The SKR 4000 is indicated for use in generating radiographic images of human anatomy. It is intended to replace a radiographic film/screen system in general-purpose diagnostic procedures.

The SKR 4000 is not indicated for use in mammography, fluoroscopy and angiography applications.

Comparison Table:

	Subject Device	Predicate Device 1 (PD1)	Predicate Device 2 (PD2)
	SKR 4000	ViZionDR + Wireless	SKR 3000
510(K) Number	This Submission	K160810	K182688
Indications for Use	The SKR 4000 is indicated for use in generating radiographic images of human anatomy. It is intended to replace radiographic film/screen system in general-purpose diagnostic procedures. The SKR 4000 is not indicated for use in	Intended for digital image capture use in general radiographic examinations, wherever conventional screen-film systems may be used, excluding fluoroscopy, angiography and mammography. ViZion DR + Wireless allows imaging of skull, chest,	The SKR 3000 is indicated for use in generating radiographic images of human anatomy. It is intended to replace radiographic film/screen system in general-purpose diagnostic procedures. The SKR 30000 is not indicated for use in



	Subject Device	Predicate Device 1 (PD1)	Predicate Device 2 (PD2)
	SKR 4000	ViZionDR + Wireless	SKR 3000
510(K) Number	This Submission	K160810	K182688
	mammography, fluoroscopy, and angiography applications.	shoulders, spine, abdomen, pelvis, and extremities.	mammography, fluoroscopy and angiography applications.
Specification			
Detection method	Same as PDs	Indirect conversion method	Indirect conversion method
Scintillator	Same as PDs	CsI (Cesium Iodide)	CsI (Cesium Iodide)
Image area size	345.6×420.0mm (2,304 x 2,800 pixels)	343.2mm x 420mm (2,288 x 2,784 pixels)	348.95×425.25mm (3,488 x 4,256pixels)
Pixel size	Same as PD1	150 µm	100 µm / 200 µm / 400 µm
A/D conversion	Same as PD2	14 bit	16 bit
MTF (1.0 cycle/mm)	More than 0.60	0.83 at 0.5 (1/mm) (Better)	More than 0.53
DQE (1.0 cycle/mm)	More than 30%@0.3mR	0.55 at 0.5 (1/mm) (Better)	More than 49% @ 1mR More than 39% @ 0.1mR
Electrical			
Power source	Same as PD1	AC Line and/or Removable Battery	AC Line and/or Built-in Battery
Battery type	Same as PD1	Lithium Ion Battery	Lithium Ion Capacitor
Mechanical			
External dimensions	384(W)×460(D)×15(H)mm	384(W)×460(D)×15.2(H)mm	384(W)×460(D)×15(H)mm
Weight	3.6kg (with battery)	3.82 kg (with battery)	2.6 kg
IPX	IPX1	IPX0	IPX6
Communication	Wireless	Wireless/Wired	Wireless/Wired

Summary of Technological Characteristics

Compared to Predicate Devices:

The SKR 4000 employs the same fundamental scientific technologies as the PDs (K160810 and K182688). The indications for use is identical to the PD1 (K182688), and the other summary of comparisons of technological characteristics for both systems is provided below;

Operational principles and designing:

The SKR 4000 has same scientific technologies and operational principal as the predicate devices (K160810 and K182688). All the verification activities required by the specification and the risk analysis for the SKR 4000 were performed and the results demonstrated that the predetermined acceptance criteria were met.

Performance test:

The performance tests according to the “Guidance for the Submission of 510(k)s for Solid State X-ray Imaging Devices” and the other verification and validation including the items required by the risk analysis for the SKR 4000 were performed and the results demonstrated that the predetermined acceptance criteria were met.

The concurrence study is not necessary because the proposed modifications are categorized as the type of unnecessary to provide additional clinical data in accordance with the SSXI guidance. Besides, the results of risk management did not require clinical studies to demonstrate the substantial equivalency of the proposed device.

Safety:

The system is in conformance with the standards described above, which are same standards to those of predicate device.

The risk analysis for the SKR 4000 has been conducted in accordance with ISO 14971 “Medical devices -- Application of risk management to medical devices”. The risks associated with all the identified hazards were reduced to acceptable level by the risk control measures as shown in risk assessment record.

Biocompatibility:

The all patient contact materials for human body surface are evaluated under ISO 10993 “Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process” and determined as acceptable for this usage. The proposed P-41 and the predicate devices achieve the same acceptance level for biocompatibility.

Relevant FDA Guidance

The SKR 4000 was developed in accordance with the following FDA guidance:

- Guidance for the Submission of 510(k)s for Solid State X-ray Imaging Devices (September 2016)
- Radio Frequency Wireless Technology in Medical Devices (August 2013)
- Appropriate Use of Voluntary Consensus Standards in Premarket Submissions for Medical Devices (September 2018)
- Refused to Accept Policy for 510(k)s (February 2019)
- Format for Traditional and Abbreviated 510(k)s (August 2005)
- Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices (May 2005)



Conclusion:

The SKR 4000 has the same intended use and indications for use, technological characteristics, and principal operations. The technological differences raised no new issues of safety or effectiveness as compared to its predicate devices. Performance tests demonstrate that the SKR 4000 performs according to specifications and functions as intended. Therefore, the SKR 4000 is substantially equivalent to its predicate devices (K160810 and K182688).