



February 20, 2020

Cadwell Industries, Inc.
Mr. James Blevins
Product Manager - Sleep
909 N. Kellogg Street
Kennewick, Washington 99336

Re: K192624

Trade/Device Name: ApneaTrak
Regulation Number: 21 CFR 882.1400
Regulation Name: Electroencephalograph
Regulatory Class: Class II
Product Code: OLV, MNR
Dated: January 7, 2020
Received: January 21, 2020

Dear Mr. James Blevins:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see

<https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jay Gupta
Assistant Director
DHT5A: Division of Neurosurgical,
Neurointerventional
and Neurodiagnostic Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K192624

Device Name

ApneaTrak

Indications for Use (Describe)

The Cadwell ApneaTrak device is intended for home sleep testing, including the acquisition of physiological and environmental data. The recorded signals are then transmitted to a PC so that they can be viewed. ApneaTrak is intended for use on patients older than 2 years of age.

ApneaTrak is intended for use in hospitals, sleep centers and other sleep testing environments, including the patient's home. ApneaTrak is intended to be used when prescribed by a qualified healthcare provider for use on patients suspected of sleep disorders, specifically Sleep Disordered Breathing (SDB) and requires review by qualified medical personnel. ApneaTrak is NOT intended to perform automatic diagnosis.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Submitter: Cadwell Industries, Inc.
909 N. Kellogg Street
Kennewick, Washington 99336
509-735-6481

Contact Person: James Blevins
Email: jamesbl@cadwell.com

Date Prepared: September 20, 2019

Trade Name: ApneaTrak

Regulation Name: Standard Polysomnograph with Electroencephalograph

Regulation Number: 21 CFR 882.1400

Regulatory Classification: Class II

Product Codes: OLV, MNR

Classification Panel: Neurology

Predicate Device: Primary Predicate:
Zmachine Synergy Sleep Monitoring System from Consolidated Research of Richmond, Inc.
Product Code: OLV, OMC, MNR
510(k) Number: K172986
Reference Predicate:
Nox T3 Sleep Recorder from Nox Medical
Product Code: MNR
510(k) Number: K082113

Device Description: Cadwell's ApneaTrak is a sleep diagnostic system consisting of: (1) acquisition hardware that can acquire, record, store, and transfer up to 3 channels of ExG (including EEG, EMG, ECG, and EOG signals) data, 2 respiratory effort channels, 1 thermistor channel, 1 pressure channel, 1 snore channel, and 1 oximetry channel; (2) a host electronic device (typically a PC) capable of running the software as well as charging

and interfacing with the acquisition device; and (3) software that allows for device configuration and data download.

ApneaTrak is connected, by a clinical user, to a host device via USB cable for initialization. After initialization and having been given instruction on correct clinical use of the device, ApneaTrak is then used by the patient at home. The device acquires and stores physiological and/or environmental data to onboard memory. After use, the device is returned to the clinical user, who connects the device to the host PC. The software downloads and stores data from the device in European Data Format (EDF).

Indications for Use:

The Cadwell ApneaTrak device is intended for home sleep testing, including the acquisition of physiological and environmental data. The recorded signals are then transmitted to a PC so that they can be viewed. ApneaTrak is intended for use on patients older than 2 years of age.

ApneaTrak is intended for use in hospitals, sleep centers and other sleep testing environments, including the patient's home. ApneaTrak is intended to be used when prescribed by a qualified healthcare provider for use on patients suspected of sleep disorders, specifically Sleep Disordered Breathing (SDB) and requires review by qualified medical personnel. ApneaTrak is NOT intended to perform automatic diagnosis.



510(k) Summary

Technology The ApneaTrak employs the same technological characteristics as the predicate devices.

Comparison:

	Zmachine Synergy (K172986, Primary Predicate)	Nox T3 Sleep Recorder (K082113, Reference Predicate)	Cadwell ApneaTrak (Proposed Device)	Discussion of Differences
General				
Product Codes	OLV, OMC, MNR	MNR	OLV, MNR	Equivalent
Classification Regulation	21 CFR 882.1400	21 CFR 868.2375	21 CFR 882.1400	Equivalent
Population	Adults	Greater than 2 years of age	Greater than 2 years of age	Differences in patient population do not raise concerns of safety or effectiveness.
Prescription Use	Yes	Yes	Yes	Equivalent
Intended Environments	Home or clinical environment	Home or clinical environment	Home or clinical environment	Equivalent
Indication for Use	The Zmachine Synergy is an EEG and respiratory signal recorder. The device is intended for use by adult patients in the home or clinical environment, under the direction of a qualified	The Nox T3 device is intended for ambulatory recording of physiological signals during sleep. The recorded signals are then downloaded to a PC where the signals can be viewed and analyzed by use of the Nox T3 application	The Cadwell ApneaTrak device is intended for home sleep testing, including the acquisition of physiological and environmental data. The recorded signals are then transmitted to a PC so that they can be viewed. ApneaTrak is	The indications for use between the proposed device and the predicate have no substantive differences.



510(k) Summary

	healthcare practitioner, to aid in the diagnosis of sleep disorders.	(Noxturnal). The Nox T3 system is indicated for use in patients greater than 2 years of age. The Nox T3 system is NOT intended for any patient monitoring or automatic diagnosis. The intended environments are hospitals, institutions, sleep centers, sleep clinics, or other test environments, including patient's home.	intended for use on patients older than 2 years of age. ApneaTrak is intended for use in hospitals, sleep centers and other sleep testing environments, including the patient's home. ApneaTrak is intended to be used when prescribed by a qualified healthcare provider for use on patients suspected of sleep disorders, specifically Sleep Disordered Breathing (SDB) and requires review by qualified medical personnel. ApneaTrak is NOT intended to perform automatic diagnosis.	
Physical				
Case Material	ABS Plastic	ABS Plastic	ABS Plastic	Equivalent
Dimension	61mm (2.40") W 120mm (4.72) H 24mm (0.94") D	80mm (2.5") W 111mm (4.9") H 18.5mm (0.8") D	73mm (2.87") W 115mm (4.52") H 25mm (0.98") D	Differences in physical dimension do not raise concerns of safety or effectiveness.
Weight	86g	280g	143.5g	Differences in physical dimension do not raise concerns of safety or effectiveness.
Functional				
Acquisition Units	One unit	One unit	One unit	Equivalent

Number of Channel Inputs	6	7	8	The proposed device has 2 additional EXG channel inputs. These additional EXG channels are composed of the same hardware and are controlled in the same manner as the other EXG channels. As such, whether the overall number of EXG channels is 2 or 3, the risk profile of the device does not change. Therefore, this difference does not raise concerns of safety or effectiveness.
Recording Time	Up to 30 hours	Up to 24 hours	Up to 24 hours	Differences in recording time do not raise concerns of safety or effectiveness.
Data Storage	On-board storage	On-board storage	On-board storage	Equivalent
Data Interface (PC)	USB	USB	USB	Equivalent
Connections to Patient	RIP belt for respiratory effort Probes or Flexi wrap for oximetry Plastic tubing and cannula for pressure sensing	RIP belt for respiratory effort Probes or Flexi wrap for oximetry Plastic tubing and cannula for pressure sensing	RIP or PVDF belts for respiratory effort Probes for oximetry Plastic cannula for pressure sensing Plastic snore microphone	All accessories are previously cleared, widely available medical devices. Differences do not raise questions of safety and effectiveness.



510(k) Summary

Display Type	LEDs on device for signal check and device status	LCD on device for signal check and device status	LEDs on device for signal check and device status	Equivalent
Signals and Sensors				
Signals Recorded	Respiratory Effort Body position Activity Oxygen Saturation Pulse Airflow Snore EEG	Respiratory Effort (Abdomen and Thorax) Body position Activity Oxygen Saturation Pulse Nasal/mask pressure Airflow Snore EEG, EOG, EMG, ECG Respiratory sound	Respiratory Effort (Abdomen and Thorax) Body position Activity Oxygen Saturation Pulse Nasal/mask pressure Airflow (Pressure and Thermal) Snore EEG, EOG, EMG, ECG	The addition of an independent snore sensor (both sensors derive snore from the cannula as well) and additional ExG inputs (3 vs 1) do not raise questions of safety or effectiveness as the channels are just multiplied to meet clinical user's needs. Nasal/mask pressure is the same as pressure-based airflow and would be equivalent in the predicate despite not being listed.
Sensor Technology	Solid state pressure sensor Solid state position/activity sensor Respiratory effort sensors (RIP Technology) Oximetry Zmachine EEG Technology	Solid state pressure sensor Solid state position/activity sensor Respiratory effort sensors (RIP technology) Oximetry Microphone Gold cup electrodes Ag/AgCL electrodes	Solid state pressure sensor Solid state position/activity sensor Respiratory effort sensors (RIP or PVDF technology) Oximetry Snore Microphone Gold cup electrodes Ag/AgCL electrodes	Equivalent ApneaTrak uses common ExG technology using gold or Ag/AgCL cup electrodes. The predicate device uses a snap style connector. Electrode connection styles do not raise questions of safety or effectiveness as they are



510(k) Summary

				both used commonly in the industry.
Power and Isolation				
Power Source	Internally powered using li-ion rechargeable battery	Internally powered using AA disposable battery	Internally powered using li-ion rechargeable battery	Equivalent
Patient Isolation	Device has no galvanic connections to mains during operation as it is a battery-operated device Not possible to connect auxiliary devices to the device	Device has no galvanic connections to mains during operation as it is a battery-operated device Not possible to connect auxiliary devices to the device	Device has no galvanic connections to mains during operation as it is a battery-operated device Not possible to connect auxiliary devices to the device	Equivalent
Transmitter				
RF Data transfer	None	Bluetooth wireless technology	Bluetooth Transmitter Equipped, but disabled	Bluetooth transmitter has been included as part of the ApneaTrak Hardware but is disabled by firmware/software at this time and is not available to the user. Bluetooth functionality is subject to future 510(k) clearance. Equivalent currently.
Recorded Data				
Access to recorded data	Recorded data is stored in the device. When the device is connected to a PC via USB cable the device provides access to its internal memory.	Recorded data is stored in the device. When the device is connected to a PC via USB cable the device provides access to its internal memory.	Recorded data is stored in the device. When the device is connected to a PC via USB cable the device provides access to its internal memory.	Equivalent



510(k) Summary

Recorded data format	Each channel of recorded data is stored in an individual file of the GCS2 data format.	All channels of recorded data are stored in proprietary data format.	All channels of recorded data are stored in EDF data format.	GCS2 appears to be the predicate device's proprietary data format. EDF is non-proprietary. Data format does not raise new concerns of safety or effectiveness.
Device Initialization	Yes	Yes	Yes	Equivalent
Data Download	Yes	Yes	Yes	Equivalent

Electrical Safety:

The ApneaTrak was tested for safety and essential performance in accordance with the following safety standards:

- IEC 60601-1:2005 + CORR. 1:2006 + CORR. 2:2007 + A1:2012
- AAMI ES60601-1:2005 +C1+A2 [R2012]
- IEC 60601-2-40:2016
- IEC 60601-2-26:2012
- ISO 80601-2-61:2017
- IEC 60601-1-11:2015
- IEC 60601-1-6:2013
- IEC 62304:2006 + A1:2015

Test results indicate that the ApneaTrak complies with the applicable standards.

Electromagnetic Disturbances:

The ApneaTrak was tested for performance in accordance with the following standards:

- IEC 60601-1-2:2014

Test results indicate that the ApneaTrak complies with the applicable standards.

Performance Testing:

The ApneaTrak was tested in accordance with internal software requirements, system requirements, and usability requirements as a result of the risk analysis External tests have been completed for electrical safety and EMC as indicated above. In addition, bench testing results are summarized below to demonstrate the ApneaTrak complies with its predetermined specifications and accurately captures each biologically relevant signal being collected by the predicate devices.

Test	Test Method Summary	Results
ExG	ExG functionality of the subject device is validated by complying with essential performance requirements from the following standards: IEC 60601-2-26 and IEC 60601-2-40.	All tests results demonstrate compliance with the standards.
Pulse Oximetry	The pulse oximetry functionality of the subject device is validated by complying with ISO 80601-2-61 Particular requirements for pulse oximeters	All tests results demonstrate compliance with the standard.

Respiratory Effort	A known oscillating input signal was injected into the respiratory channel of the subject device. The input and output data were plotted and quantitatively compared.	Passing result based on high measure of equivalence between input and output signals.
Airflow - Pressure	A known oscillating input signal was input to the airflow pressure channel of the subject device. The input and output data were plotted and quantitatively compared.	Passing result based on high measure of equivalence between input and output signals.
Airflow - Thermal	A known oscillating input signal was input to the airflow thermal channel of the subject device. The input and output data were plotted and quantitatively compared.	Passing result based on high measure of equivalence between input and output signals.
Snore	A known oscillating input signal was input to the snore channel of the subject device. The input and output data were plotted and quantitatively compared.	Passing result based on high measure of equivalence between input and output signals.

All tests confirm the ApneaTrak meets the requirements for both the external tests and performance bench tests. Cadwell Industries, Inc has determined the ApneaTrak is Substantially Equivalent to the predicate devices.

Conclusion:

Verification and validation activities were conducted to establish the performance and safety characteristics of the ApneaTrak. The results of these activities demonstrate that the ApneaTrak is as safe, as effective, and performs as well as or better than the predicate devices.

Therefore, the ApneaTrak is considered substantially equivalent to the predicate devices.