



October 22, 2019

RAY CO., Ltd
c/o Changhwan Lee
RA Manager
332-7, Samsung 1-ro
Hwaseong-si, Gyeonggi-do 18380
REPUBLIC OF KOREA

Re: K192737
Trade/Device Name: RCT800
Regulation Number: 21 CFR 892.1750
Regulation Name: Computed Tomography X-Ray System
Regulatory Class: Class II
Product Code: OAS
Dated: September 25, 2019
Received: September 27, 2019

Dear Changhwan Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For

Thalia T. Mills, Ph.D.
Director
Division of Radiological Health
OHT7: Office of In Vitro Diagnostics
and Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K192737

Device Name

RCT800

Indications for Use (Describe)

RCT800 is CBCT and panoramic x-ray imaging system with cephalometric.

Which is intended to radiographic examination of the dento-maxillofacial, sinus, TMJ, Airway for diagnostic support for adult and pediatric patients. And a model scan is included as an option.

Cephalometric image is also includes wrist to obtain carpus images for growth and maturity assessment for orthodontic treatment.

The device is to be operated and used by dentists or other legally qualified health care professionals.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

K192737

1. 510(k) Summary

The 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR Part 807.92.

2. Date: September 25, 2019

3. Administrative Information

APPLICANT RAY Co.,Ltd

ADDRESS #332-7, Samsung 1-ro, Hwaseong-si, Gyeonggi-do, 18380, Korea

Manufacturer RAY Co.,Ltd
332-7, Samsung 1-ro, Hwaseong-si, Gyeonggi-do, 18380, Korea
TEL : +82-31-605-1000
FAX : +82-2-6280-5534

Contact Person e-mail : ch0406.lee@raymedical.co.kr

4. Device Information**Device Name**

Trade/Proprietary Name: RCT800
Common Name: Dental panoramic/tomography and cephalometric x-ray system

Classification

Classification Name: Computed tomography x-ray system
Regulation Number : 21 CFR 892.1750
Class : II
Product code : OAS
Panel : Radiology

5. Predicate device

Parameter	Predicate Device-1
Device Name	RCT800
Manufacturer	RAY Co., Ltd

510(K) Number	K182805
Classification name	Computed tomography x-ray system
Regulation number	892.1750
Primary product code	OAS

7. Device Description

System purpose RCT800 is 3D computed tomography for scanning hard tissues like bone and teeth. By rotating the c-arm which is embedded with high voltage generator all-in-one x-ray tube and a detector on each end, CBCT images of dental maxillofacial is attained by recombining data from the same level that are scanned from different angle.

Panoramic image scanning function for attaining image of whole teeth, cephalometric scanning option for attaining cephalic image, and Model Scan option for attaining dental model CBCT image are included.

8. Indication for use

RCT800 is CBCT and panoramic x-ray imaging system with cephalometric.

Which is intended to radiographic examination of the dento-maxillofacial, sinus, TMJ, Airway for diagnostic support for adult and pediatric patients. And a model scan is included as an option.

Cephalometric image is also includes wrist to obtain carpus images for growth and maturity assessment for orthodontic treatment.

The device is to be operated and used by dentists or other legally qualified health care professionals.

9. Patient population

The patient population can be the possible person who can be taken X-ray diagnostic radiation exposure.

There is no restriction for ethnic group, Gender, weight, health, or condition.

We recommend patients for x-ray diagnostic radiation exposure to be over 5 years old.

10. Comparison with predicate device

The following table provides the summary of the technological characteristics of RCT800 compared to the predicate device.

Parameter	Proposed Device	Predicate Device	
Manufacturer	RAY Co., Ltd.	RAY Co., Ltd.	
Device name	RCT800	RCT800	
510(K) Number	(Traditional 510K)	K182805 (Special 510K)	
Common Name	Dental panoramic/tomography and cephalometric x-ray system	Dental panoramic/tomography and cephalometric x-ray system	
Indications for use	RCT800 is CBCT and panoramic x-ray imaging system with cephalometric. Which is intended to radiographic examination of the dento-maxillofacial, sinus, TMJ, Airway for diagnostic support for adult and pediatric patients. And a model scan is included as an option. Cephalometric image is also includes wrist to obtain carpus images for growth and maturity assessment for orthodontic treatment. The device is to be operated and used by dentists or other legally qualified health care professionals.	RCT800 is CBCT and panoramic x-ray imaging system with cephalometric. Which is intended to radiographic examination of the dento-maxillofacial, sinus, TMJ, Airway for diagnostic support for adult and pediatric patients. And a model scan is included as an option. Cephalometric image is also includes wrist to obtain carpus images for growth and maturity assessment for orthodontic treatment. The device is to be operated and used by dentists or other legally qualified health care professionals.	
Mode of Operation	Same as predicate device #1	Continuous operation with intermittent, stated permissible loading	
3D technology	Same as predicate device #1	CBCT Cone beam Computed Tomography	
Performance Specification	Same as predicate device #1	1) CBCT Computed tomography - Patient - Dental Model Scan(Optional) 2) Panoramic 3) Cephalometric(optional) - One shot type - Scan type	
Functional Option	Same as predicate device #1	Base CT+PANO Option(CEPH) CT + PANO + SCAN CEPH CT + PANO + One shot(One shot, Standard Type) CT + PANO + One shot(One shot, Large Type).	
Detector Type	CT	Same as predicate device #1	FXDD-0606CA
		Same as predicate device #1	FXDD-1012CHA
	PANO	Same as predicate device #1	FXDD-0606CA
		Same as predicate device #1	FXDD-1012CHA

	Ceph (Scan)	Same as predicate device #1	XID-C24DC
	Ceph (One shot)	Same as predicate device #1	1717SCC
		Same as predicate device #1	PaxScan 2530C
Exposure switch Type		Same as predicate device #1	"Deadman" Button type
Main Components		Same as predicate device #1	Ceph Apparatus
		Same as predicate device #1	Vertical Carriage
		Same as predicate device #1	Rotator
		Same as predicate device #1	X-RAY Generator
		Same as predicate device #1	X-ray tube
		Same as predicate device #1	High Frequency Generator
		Same as predicate device #1	Column
		Same as predicate device #1	Touch monitor (panel)
		Same as predicate device #1	Detector - CT FXDD-0606CA FXDD-1012CHA - PANO FXDD-0606CA FXDD-1012CHA - Ceph XID-C24DC(Scan) 1717SCC(One shot, Large Size) PaxScan 2530C(One shot, Standard Size)
		Same as predicate device #1	Chinrest
		Same as predicate device #1	Head rest
		Same as predicate device #1	Automatic Collimator
		Same as predicate device #1	Exposure switch
		Same as predicate device #1	Emergency stop switch
	Same as predicate device #1	Console PC set	
Automatic Collimator		Same as predicate device #1	CT exams Panoramic exams Cephalometric exams
Display Type		Same as predicate device #1	TFT LCD type(Normally black) *1280x800 pixel
Class		Same as predicate device #1	Class I with type B applied parts according to IEC 60601-1
Focal size		Same as predicate device #1	Patient 0.5 Model scan 0.04 (Optional)
Field of View(CT)		Same as predicate device #1	FXDD-0606CA : Max.160x100 mm FXDD-1012CHA : Max. 200x200 mm

X-ray Voltage(Patient)	60~100kVp	60~90kVp	
X-ray Current(Patient)	Same as predicate device #1	4~17mA	
X-ray Voltage(Model Scan, Optional)	Same as predicate device #1	50~80kVp	
X-ray Current(Model Scan, Optional)	Same as predicate device #1	0.4~0.7mA	
Total Filtration	Same as predicate device #1	Min. 2.8 mm Al equivalent	
Detector Pixel size	CT	Same as predicate device #1	FXDD-0606CA: 119 μ m
		Same as predicate device #1	FXDD-1012CHA: 124 μ m
	PANO	Same as predicate device #1	FXDD-0606CA: 119 μ m
		Same as predicate device #1	FXDD-1012CHA: 124 μ m
	Ceph (Scan)	Same as predicate device #1	XID-C24DC: 100 μ m
	Ceph(One shot)	Same as predicate device #1	1717SCC: 127 μ m
Same as predicate device #1		PaxScan 2530C: 139 μ m	
Magnification	CT	Same as predicate device #1	FXDD-0606CA: 1.44(Patient) 1.91(Model Scan)
		Same as predicate device #1	FXDD-1012CHA: 1.44(Patient) 1.91(Model Scan)
	PANO	Same as predicate device #1	FXDD-0606CA: 1.3
		Same as predicate device #1	FXDD-1012CHA: 1.3
	Ceph (Scan)	Same as predicate device #1	XID-C24DC: 1.11
	Ceph(One shot)	Same as predicate device #1	1717SCC: 1.13
Same as predicate device #1		PaxScan 2530C: 1.12	
Scan time	CT : below 20sec(Patient) CT : below 180sec(Model Scan)		CT : below 14sec(Patient) CT : below 180sec(Model Scan)
	Same as predicate device #1		Pano : below 14sec
	Ceph[Scan type] : below 20sec		Ceph[Scan type] : below 19sec
	Same as predicate device #1		Ceph[One shot type]: below 2sec
Format compatible	Same as predicate device #1	DICOM 3.0 Format compatible	
Image Viewing Software	Same as predicate device #1	RayScan (Cleared under K182805)	
Image acquisition	Same as predicate device #1	Giga-Ethernet Network	
Total Height	Same as predicate device #1	Max 2,296mm	

Weight	Same as predicate device #1	<p>1) Computed Tomography(CT) + Panoramic(PANO)=189kg(416.6lb) \pm 10%</p> <p>2) Computed Tomography(CT) + Panoramic(PANO) + Ceph (Scan type)= 219kg (482.8lb) \pm 10%</p> <p>3) Computed Tomography(CT) + Panoramic(PANO) + Ceph (One shot type, installed in Standard size)= 217kg (478.4lb) \pm 10%</p> <p>4) Computed Tomography(CT) + Panoramic(PANO) + Ceph (One shot type, installed in Large size) 212kg (467.3lb) \pm 10%</p>
Type of installation	Same as predicate device #1	Wall or floor mount
Patient position	Same as predicate device #1	Standing / Wheelchair
Applicable Standards	Same as predicate device #1	<p>IEC 60601-1</p> <p>IEC 60601-1-3</p> <p>IEC 60601-2-63</p> <p>IEC 60601-1-2</p>

The product is principally just the same as in the previous 510(k) #K182805.

The complete of differences of the subject device to the predicate device is as follows

- The maximum X-ray voltage of the tube has been changed from 90kV to 100kV.
- The irradiation time of CT and Scan Ceph has been changed.
- Face Scan mode has been added in option. Face Scan is 3D Face Photo to capture 3D facial picture to visualize soft tissue. This is the same feature as CS Face Scan on the CS 9600 (K118136) and ProFace on the Planmeca ProMax 3D Max (K160506).

12. Safety and Effectiveness Information

RCT800 system described in this 510(k) is similar to the predicate device in terms of indications for use, materials, safety characteristics, and X-ray source.

The following information further substantiates the substantial equivalence between the subject device and predicate device.

The fundamental technological characteristics of the subject and predicate device are similar.

The imaging modes are similar; PANO, CEPH (Optional), CBCT, Model Scan All viewing software programs have been cleared with previous 510k submissions; RAYSCAN(K182805).

The sponsor tested the subject device in a laboratory and provided a non-clinical performance report. The same test protocol was used to test the performance of the subject and the predicate device for comparison. The sponsor certifies that adequate design and development controls (according to 21 CFR 820.30) were in place for manufacturing the subject device.

The differences are as follows.

- The maximum X-ray voltage of the tube has been changed from 90kV to 100kV.
- The irradiation time of CT and Scan Ceph has been changed.

Electrical, mechanical and environmental safety testing according to standard of IEC 60601-1: 2005/AMD1:2012(3.1 Edition), IEC 60601-1-3: 2008/AMD1:2013(Second Edition), IEC 60601-1-6:2010(Third Edition) and IEC 60601-2-63: 2012(first Edition) were performed. EMC testing was conducted in accordance with the standard IEC 60601-1-2: 2014(Edition 4.0).

The software of RCT800 saves the patient and image data and offers an inquiry function, in addition, supports the image generate function intended to obtain images using the RCT800 equipment and various sensors for diagnosis. And that has been validated according to FDA "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" and "Content of Premarket Submissions for Management of Cybersecurity in Medical Devices" to assure substantial equivalence. The software for this device was considered as a "moderate" level of concern, since a failure or latent flaw in the software would not directly result in serious injury or death to the patient or operator. As a result, we identified the level of concern associated with new device and provided documentation consistent with that level. Based on our risk analysis of software, the difference does not affect its safety and effectiveness.

Bench testing was conducted according to FDA Guidance "Format for Traditional and Abbreviated 510(k)s, section 18, Performance Testing – Bench"

Bench testing is used to assess whether or not the parameter measured required for describing functionalities related to imaging properties of the dental X-ray device and patient dosage satisfies the designated tolerance.

Performance (Imaging performance) testing was conducted according to standard of IEC 61223-3-4 and IEC 61223-3-5.

All test results were satisfactory.

Non-clinical considerations were conducted in accordance with FDA Guidance "Guidance for the submissions of 510(k)'s for Solid State X-ray Imaging Devices". Because the subject device used the same detector as the predicate device, there is no significant difference between the two devices as a result of non-clinical testing.

Clinical considerations were conducted according to FDA Guidance "Format for Traditional and Abbreviated 510(k)s, section 20".

Clinical images were provided these images were not necessary to establish substantial equivalence based on the modifications to the device but they provide further evidence in addition to the laboratory performance data to show that the complete system works as intended.

The features of RCT800 have been clinically tested and approved by two licensed practitioners/clinicians.

The clinical imaging samples are collected from the all detector on propose device at the 2 offices where the predicate device is installed on clinical consideration report for the clinical test images. These images were gathered from the all detector installed with RCT800 on any protocols with random patient age, gender, and size. A licensed practitioner reviewed the sample clinical images and found them to be of acceptable quality for the intended use.

12. Conclusions

In accordance with the Federal Food, Drug and Cosmetic Act, 21 CFR Part 807 and based on the information provided in this premarket notification. RAY Co., Ltd. concludes that the newly RCT800 is safe and effective and substantially equivalent to predicate device as described herein.