



December 6, 2019

Nu Eyne Co., Ltd  
% Dongha Lee  
Regulatory Affairs Consultant  
KMC, Inc.  
Room no. 904, 27, Digital-ro 27ga-gil, Guro-gu  
Seoul, 08375 Kr

Re: K192773

Trade/Device Name: Allive  
Regulation Number: 21 CFR 882.5891  
Regulation Name: Transcutaneous Electrical Nerve Stimulator To Treat Headache  
Regulatory Class: Class II  
Product Code: PCC,  
Dated: September 26, 2019  
Received: September 30, 2019

Dear Dongha Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Timothy Marjenin  
Assistant Director  
DHT5B: Division of Neuromodulation  
and Physical Medicine Devices  
OHT5: Office of Neurological  
and Physical Medicine Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K192773

Device Name

ALLIVE (MODEL: ALLIVE01, AD-01)

Indications for Use (Describe)

The Allive is indicated for the prophylactic treatment of episodic migraine in patients 18 years of age or older.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

---

## 510(k) SUMMARY

This summary of 510(k) –safety and effectiveness information is being submitted in accordance with requirements of 21 CFR Part 807.92.

Date: SEPTEMBER 26, 2019

### 1. INFORMATION

#### 1.1 Submitter Information

- Submitter Name: Nu Eyne Co., Ltd.
- Address  
: #403, Seoul Biohub, 117-3 Hoegi-ro, Dongdaemun-gu, Seoul, 02455, Korea
- Telephone Number: +82-2-6953-8120      ▪ Fax: +82-303-347-0017
- Email: dongseong.lee@nueyne.com

#### 1.2 Contact Person

- Name: DongHa Lee (Consultant / KMC, Inc.)
- Address: Room no. 904, 27, Digital-ro 27ga-gil, Guro-gu, Seoul, 08375, Korea
- Telephone Number: +82-70-8965-5554      ▪ Fax: +82-2-2672-0579
- E-mail: dhlee@kmcerti.com

### 2. DEVICE INFORMATION

2.1 Trade Name / Proprietary Name: ALLIVE (Model: ALLIVE01, AL-2HP)

2.2 Common Name: Transcutaneous electrical nerve stimulator to treat headache

2.3 Classification Name: Stimulator, Nerve, Electrical, Transcutaneous, For Migraine

2.4 Product Code: PCC

2.5 Classification Regulation: 21CFR 882.5891

2.6 Device Class: Class II (Special Control)

2.7 Classification Panel: Neurology

---

### 3. PREDICATE DEVICE

Predicate Device	
Manufacturer	CEFALY Technology
Device Name (Trade Name)	CEFALY®
510(k) Number	K160237

### 4. SUBJECT DEVICE DESCRIPTION



The Allive device is a transcutaneous electrical nerve stimulator (TENS) that is applied to the forehead using a self-adhesive electrode positioned over the upper branches of the trigeminal nerve bilaterally. It is intended to stimulate the upper branches of the trigeminal nerve in order to reduce the frequency of migraine attack.

### 5. INTENDED USE

The Allive is indicated for the prophylactic treatment of episodic migraine in patients 18 years of age or older.

## 6. SUBSTANTIAL EQUIVALENCE

Items	Subject Device	Predicate Device	Comparison Result
<b>Manufacturer</b>	Nu Eyne Co., Ltd.	CEFALY Technology	Different
<b>Device</b>	Stimulator, Nerve, Electrical, Transcutaneous, For Migraine	Stimulator, Nerve, Electrical, Transcutaneous, For Migraine	Same
<b>Trade/Device Name</b>	Allive	CEFALY®	Different
<b>K Number</b>	None	K160237	Different
<b>Regulation Number</b>	21CFR 882.5891	21CFR 882.5891	Same
<b>Regulation Description</b>	Transcutaneous electrical nerve stimulator to treat headache	Transcutaneous electrical nerve stimulator to treat headache	Same
<b>Regulatory Class</b>	Class II	Class II	Same
<b>Product Code</b>	PCC	PCC	Same
<b>Definition</b>	Used to apply an electrical current to a patient's cranium through electrodes placed on the skin.	Used to apply an electrical current to a patient's cranium through electrodes placed on the skin.	Same
<b>Review Panel</b>	Neurology	Neurology	Same
<b>Physical State</b>	Electrical stimulation unit with leads and cutaneous electrodes.	Electrical stimulation unit with leads and cutaneous electrodes.	Same
<b>Technical Method</b>	Applies an electrical current through electrodes on patient's skin.	Applies an electrical current through electrodes on patient's skin.	Same

<b>Target Area</b>	Afferent cranial nerves.	Afferent cranial nerves.	Same
<b>Intended use</b>	The Allive is indicated for the prophylactic treatment of episodic migraine in patients 18 years of age or older.	The Cefaly® device is indicated for the prophylactic treatment of episodic migraine in patients 18 years of age or older.	Same
<b>Picture</b>			Different
<b>Power Source</b>	1rechargeable Lipo 3.7 V battery	1rechargeable Lipo 3.7 V battery	Same
<b>Channels</b>	1	1	Same
<b>Computerized</b>	Yes	Yes	Same
<b>S/W provided</b>	1 fixed program	1 fixed program	Same
<b>Constant current</b>	Yes	Yes	Same
<b>Constant voltage</b>	No	No	Same
<b>Max output current</b>	16mA	16mA	Same
<b>Max output voltage (2kOhm)</b>	32 Volts	32 Volts	Same
<b>Patient Override Control Method</b>	On/Off button	On/Off button	Same

<b>Max Leakage Current</b>		None (battery operated)	None (battery operated)	Same
<b>Electrode</b>		Self-adhesive	Self-adhesive	Same
<b>Indicator display: Unit functioning</b>		Yes	Yes	Same
<b>Low battery indicator</b>		Yes	Yes	Same
<b>Standards</b>		IEC 60601-1 IEC 60601-1-2 IEC 60601-1-6 IEC 60601-1-11 IEC 60601-2-10 IEC 62366-1 IEC 62304 ISO 10993-1 ISO 10993-5 ISO 10993-10	IEC 60601-1 IEC 60601-1-2 IEC 60601-1-6 IEC 60601-1-11 IEC 60601-2-10 IEC 62366-1 IEC 62304 ISO 10993-1 ISO 10993-5 ISO 10993-10	Same
<b>Device</b>	<b>Timer Setting</b>	Yes	Yes	Same
	<b>Weight</b>	110.2g	12g	Different
	<b>Dimensions</b>	160mm x 66.5mm x 55mm	55mm x 35mm x 15mm	Different
	<b>Expected Service Life</b>	5 years	7 years	Different
	<b>IP Classification</b>	IP22	IP22	Same
	<b>Electrical Protection</b>	Type BF	Type BF	Same
<b>Battery</b>	<b>Battery Type</b>	Lithium ion Battery	Lithium ion Battery	Same
	<b>Expected Service Life</b>	7 years (300Cycles of complete discharge)	7 years (300Cylces of complete discharge)	Same
	<b>Maximum input voltage (USB connector)</b>	5.25 Vdc	5.25 Vdc	Same
<b>Electrode</b>	<b>Dimensions</b>	94mm (W) x 30mm (H)	94mm (W) x 30mm (H)	Same
	<b>Expected Service Life</b>	20 times	20 times	Same



<b>Special Requirement b.3) in accordance with 21CFR 882.5891</b>				
<b>Waveform</b>		Biphasic	Biphasic	Same
<b>Maximum output voltage (V)</b>	<b>500 ohms</b>	8	8	Same
	<b>2,000 ohms</b>	32	32	Same
	<b>10,000 ohms</b>	65	60	Same
<b>Maximum output current (mA)</b>	<b>500 ohms</b>	16	16	Same
	<b>2,000 ohms</b>	16	16	Same
	<b>10,000 ohms</b>	6.5	6	Same
<b>Pulse duration (µS)</b>		505	505	Same
<b>Frequency Range</b>		60Hz fixed	60Hz fixed	Same
<b>Net charge (µC) per pulse</b>		0	0	Same
<b>Duration of the primary (depolarizing) phase (µS)</b>		250	250	Same
<b>Standby duration between the two phase (µS)</b>		5	5	Same
<b>Maximum current density (mA/cm<sup>2</sup>, r.m.s.)</b>	<b>500 ohms</b>	2.37	2.37	Same
<b>Maximum average power density (W/cm<sup>2</sup>)</b>	<b>500 ohms</b>	0.000017	0.000017	Same
<b>Maximum average current (Average absolute value, mA)</b>	<b>500 ohms</b>	0.48	0.48	Same
<b>Other Technical Item</b>				
<b>Wave Shape</b>		Rectangular Full compensated Symmetrical	Rectangular Full compensated Symmetrical	Same

<b>Phase rise time</b>		2 $\mu$ S	2 $\mu$ S	Same
<b>Phase decay time</b>		2 $\mu$ S	2 $\mu$ S	Same
<b>Maximum Phase Charge (<math>\mu</math>S) at 500 ohms</b>		4	4	Same
<b>Modulation option</b>	<b>Amplitude</b>	0 to 16 mA Fixed	0 to 16 mA Fixed	Same
	<b>Frequency</b>	60Hz fixed	60Hz Fixed	Same
	<b>Duration</b>	250 $\mu$ S	250 $\mu$ S	Same
<b>Ramp Modulation</b>	<b>Ramp up</b>	14min	14min.	Same
	<b>Ramp down</b>	45s	45s	Same
<b>Steady time</b>		6min	6min	Same
<b>Session Duration</b>		20min 45s	20min 45s.	Same
<b>Amplitude modulation</b>		Amplitude is adjusted by the user	Amplitude is adjusted by the user	Same
<b>Material</b>				
<b>Device housing materials</b>		Plastic ABS	Plastic ABS	Same
<b>Electrode top layer</b>		Polyethylene terephthalate	Polyurethane	Different
<b>Electrode intermediate layer</b>		Conductive silver carbon ink	Carbon with silver coating	Different
<b>Electrode bottom layer</b>		Uv curable dielectric ink	x	Different
<b>Central pin</b>		N/A	N/A	Same
<b>Metallic parts for magnetic attraction</b>		Tinplate	Tinplate	Same

## 7. NON-CLINICAL DATA

### 7.1 Safety Test

#### 1) Biocompatibility

The biocompatibility tests were performed to protect patients from undue risks arise from biological hazards associated with materials of manufacture and final device. The tests were performed in accordance with the following standards and FDA Guidance - Use of International Standard ISO 10993-1, "Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process".

No.	Test Items	Standards
1	Cytotoxicity	ISO 10993-5:2009
2	Sensitization	ISO 10993-10:2010
3	Intracutaneous Reactivity Test	ISO 10993-10:2010

#### 2) Electrical Safety and EMC

The electrical safety tests were performed to protect patients from undue risks arise from any hazards associated with final device. The tests were performed in accordance with the following standards.

No.	Test Items	Standards
1	General requirement for basic safety and essential performance	<ul style="list-style-type: none"> <li>IEC 60601-1:2005+A1: 2012 (AAMI/ANSI ES 60601-1: 2005+A1: 2012)</li> </ul>
2	General requirement for safety - Electromagnetic disturbances	<ul style="list-style-type: none"> <li>IEC 60601-1-2:2014</li> </ul>
3	General requirement for safety - Medical electrical equipment used in the home healthcare environment	<ul style="list-style-type: none"> <li>IEC 60601-1-11:2015 and</li> <li>FDA Guidance ("Design Considerations for Devices Intended for Home Use")</li> </ul>
4	Particular requirement for safety – Nerve and muscle stimulators	<ul style="list-style-type: none"> <li>IEC 60601-2-10:2012/Amd1:2016</li> </ul>

### 7.2 Performance Test

The following tests were performed to assess effectiveness of performance of the device. The tests were performed in accordance with following standards.

No.	Test Items	Standards
1	Particular requirement for safety – Nerve and muscle stimulators	IEC 60601-2-10:2012/Amd1:2016
2	Technical Test	Manufacturer's SOP

### 7.3 Usability V&V

The following tests were performed to assess effectiveness of usability of the device. The test was performed in accordance with following standards.

No.	Test Items	Standards
1	General requirement for safety - Usability	<ul style="list-style-type: none"> <li>• IEC 60601-1-6:2013</li> <li>• IEC 62366-1:2015 and</li> <li>• FDA Guidance (“Applying Human Factors and Usability Engineering to Medical Devices”)</li> </ul>

### 7.4 Software

The following tests were performed to assess effectiveness of software of the device. The test was performed in accordance with following standards.

No.	Test Items	Standards
1	General requirement for safety - Programmable electrical medical systems (PEMS)	<ul style="list-style-type: none"> <li>• IEC 62304:2006/Amd1:2015</li> <li>• FDA Guidance (“Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices”)</li> </ul>

---

## 8. CLINICAL DATA

Although clinical performance data are required to demonstrate that the device is safe and effective as a treatment for headache in the indicated patient population on the special control according to clause b.6) of 21CFR882.5981, we consider that the subject device (Allive) is not applied with clause b.6) of 21CFR882.5981.

Although there are not the clinical performance data of the subject device, we prepare the clinical evaluation report by using the collected clinical data of the predicate device.

In this evaluation report, the subject device is safe and effective as a treatment for headache in the indicated patient population.

The clinical evaluation was performed in accordance with following standards.

No.	Items	Standards
1	Clinical Evaluation	<ul style="list-style-type: none"><li>MEDDEV 2.7.1. rev.4</li></ul>

---

## **9. CONCLUSION**

Under the comparing substantial equivalence between the subject device and the predicate device, there are the same points such as general information, some technical and material information. Although there are some differences, the safety and performance test reports are supported to the safety and effectiveness of the subject device.

In this regard, we conclude that the subject device is substantially equivalent to the predicate device.