

SYSMEX CORPORATION OF AMERICA

6699 Wildlife Way Long Grove IL 60047-9596 (800) 379-7639 (847) 726-3505 Facsimile

JUN 2 5 2002

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: KON1241.

1. Submitted by:	Sysmex Corporation of America 6699 Wildlife Way Long Grove, IL 60047-9596 Phone: (847) 726-3675; FAX: (847) 726-3559 Contact person: Nina Gamperling Date prepared: April 18, 2002
2. Name of Device:	Trade or proprietary name: Sysmex® XT-Series
	Common name: Automated Hematology Analyzer
	<u>Classification name</u> : Automated Differential Cell Counter 21 CFR 864.5220
3. Predicate Device:	Sysmex® XE-2100
4. Device Description:	The XT-Series is an automated hematology analyzer which consists of four principle units: (1) Main Unit which aspirates, dilutes, mixes, and analyzes whole blood samples; (2) Sampler Unit which supplies samples to the Main Unit automatically; (3) IPU (Information Processing Unit) which processes data from the Main Unit and provides the operator interface with the system; (4) Pneumatic Unit which supplies pressure and vacuum from the Main Unit. Additional information on the XT is presented in the following table.
5. Intended Use:	The Sysmex® XT-Series is intended for <i>in vitro</i> diagnostic use in the clinical laboratory as a multi-parameter hematology analyzer.
6. Substantial	The following table compares the XT-Series with the XE-
equivalence-similarities	2100.
and differences	
7. Clinical Performance	Studies were performed to evaluate the equivalency of the
Data:	XT-Series to the XE-2100. Results indicated equivalent performance.
8. Conclusions:	The performance data demonstrated substantial equivalence.
	Dage ()



SYSMEX CORPORATION OF AMERICA

6699 Wildlife Way Long Grove IL 60047-9596 (800) 379-7639 (847) 726-3505 Facsimile

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS (continued) Comparison Table to Predicate Device

Features	XT-Series	XE-2100	SF-3000
(Submission #)	A1-Series	(K992875)	(K950508)
FDA Clearance		3-Nov-99	3-Nov-95
Intended Use	Automated blood cell	Automated blood cell	Automated blood cell
intended Ose	analyzer	analyzer	analyzer
Comple Trme	Whole blood	Whole blood	Whole blood .
Sample Type			
Sample Volume	150μL- Cap piercer	200μL- Cap piercer	270μL- Cap piercer
	85μL -Manual	130μL -Manual	170μL -Manual
	40μL-Capillary dilution	40μL-Capillary dilution	40μL-Capillary dilution
Performance	Same	Proven performance in	Proven performance in
		FDA submission	FDA submission
Parameters	WBC, Neut%/#,	WBC, Neut%/#,	WBC, Neut%/#,
	Lymph%/#, Mono%/#,	Lymph%/#, Mono%/#,	Lymph%/#, Mono%/#,
and the second s	Eos%/#, Baso%/#,	Eos%/#, Baso%/#,	Eos%/#, Baso%/#,
	RBC, HGB, HCT,	NRBC%/#, RBC, HGB,	RBC, HGB, HCT,
	MCV, MCH, MCHC,	HCT, MCV, MCH,	MCV, MCH, MCHC,
	RDW-CV, RDW-SD,	MCHC, RDW-CV,	RDW-CV, RDW-SD,
	RET%/#, IRF, HFR*,	RDW-SD, RET%/#,	PLT, PDW, MPV, P-
	MFR*, LFR*, PLT,	IRF, HFR*, MFR*,	LCR.
	MPV, PDW*, P-LCR*,	LFR*, PLT, MPV,	
	PCT*	PDW*, P-LCR*, PCT*	
	*Not reportable in USA	*Not reportable in USA	
Reagents	Cellpack,	Cellpack, Cellsheath	Cellpack, Sulfolyzer,
	Stromatolyser-FB,	Stromatolyser-FB,	StromatolyserFD(I),
	Stromatolyser-4DL,	Stromatolyser-4DL,	StromatolyserFD(II),
	Stromatolyser-4DS,	Stromatolyser-4DS,	Stromatolyser-FB
	Sulfolyser,	Stromatolyser, NR,	
	Ret-Search II	Stromatolyser-IM,	
		Sulfolyser,	
		Ret-Search II	
Principles	RBC, PLT: DC	RBC, PLT: Sheath-flow	RBC, PLT: DC
	detection method,	DC detection method,	detection method,
	WBC: Flow Cytometry	WBC: Flow Cytometry	WBC: Flow using
	method using	method using	semiconductor laser
	semiconductor laser	semiconductor laser	
	detection method	detection method	
	HGB: SLS-Hgb method	HGB: SLS-Hgb method	HGB: SLS-Hgb method
Dimensions	630x520x720	711x706x912	600x580x450
(HxWxD) (mm)			
Weight (kg)	59	93	60
QC System	L-J: 20 Files with 300	L-J: 10 Files with 300	L-J: 12 Files with 180
	points per file	points per file	points per file
Bar Code	Yes	Yes	Yes
No. of Test / Hr	Approx 80	Approx 113-150	Approximately 80

- Page 10

DEPARTMENT OF HEALTH & HUMAN SERVICES



JUN 2 5 2002

Food and Drug Administration 2098 Gaither Road Rockville MD 20850

Ms. Nina M. Gamperling, MBA, MT (ASCP), RAC Manager, Regulatory Affairs Sysmex Corporation of America 6699 Wildlife Way Long Grove, Illinois 60047-9596

Re: k021241

Trade/Device Name: Sysmex® XT-Series Regulation Number: 21 CFR § 864.5220

Regulation Name: Automated differential cell counter

Regulatory Class: II Product Code: GKZ Dated: June 17, 2002 Received: June 18, 2002

Dear Ms. Gamperling:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsma/dsmamain.html".

Sincerely yours,

Steven I. Gutman, M.D., M.B.A.

Steven Butman

Director

Division of Clinical

Laboratory Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

INDICATIONS FOR USE STATEMENT

510(k) Number (if known): KD2134)	
Device Name: XT-Series, Automated Hematology Analyzer	
Indications For Use:	
The Sysmex® XT-Series is intended for <i>in vitro</i> diagnostic use in the clinical laboratory as a multi-parameter hematology analyzer.	
•	
(PLEASE DO NOT WRITE BELOW THIS LINE- CONTINUE ON ANOTHER PAGE IF NEEDED)	
Concurrence of CHRD, Office of Device Evaluation (ODE)	-
Prescription Use OR Over-The-Counter Use	
Division Sign-Off)	