

MAR 19 2003

## 510(K) Summary ChromoCheck™ Protein C

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K023990

**Submitters Name & Address:** Precision Biologic Incorporated  
900 Windmill Road, Suite 100  
Dartmouth, Nova Scotia B3B 1P7  
Canada

**Contact Name:** Stephen L. Duff – Director of New Business  
Development  
Phone: 902-468-6422 ext. 224  
Fax: 902-468-6421  
Email: sduff@precisionbiologic.com

**Preparation Date:** November 18, 2002

**Device Name & Classification:** ChromoCheck™ Protein C  
Common Name: Protein C chromogenic assay  
Classification Name: Test, Quantitative factor  
deficiency  
Regulatory Class II, 81 GGP

**Predicate Device:** Chromogenix AB/Instrumentation Laboratory  
Taljegardsgatan 3  
S-431 53 Molndal  
Sweden, SW

**Device Description:** ChromoCheck™ Protein is a chromogenic assay  
consisting of a synthetic substrate and Protein C  
Activator.

**Device Intended Use:** ChromoCheck™ Protein C is intended for use as a  
chromogenic assay for the quantitative  
determination of Protein C activity in citrated human  
plasma.

**Comparison to Predicate Device:**

Parameter	ChromoCheck™ Protein C	Coamatic Protein C
Intended Use	Test, Quantitative factor deficiency	Test, Quantitative factor deficiency
Analytes	Protein C activity	Protein C activity
Component Reagent Matrices	<b>Reagent 1:</b> Chromogenic substrate in distilled water matrix <b>Reagent 2:</b> Protein C activator in a distilled water matrix	<b>Reagent 1:</b> Chromogenic substrate in distilled water matrix <b>Reagent 2:</b> Protein C activator in a distilled water matrix
Format	Lyophilized	Lyophilized
Packaging	4 x Protein C Activator (0.65 IU) 4 x Substrate (4 mg) (Reconstituted volume – <b>2.5 mL</b> ) 4 x Protein C Activator (0.65 IU) 4 x Substrate (4 mg) (Reconstituted volume – <b>5.0 mL</b> )	2 x Protein C Activator (1.2 IU) 2 x Substrate S-2366 (6 mg) (Reconstituted volume – <b>7.2 mL</b> )

**Comments on Substantial Equivalence:**

It is the opinion of Precision BioLogic Inc. that **ChromoCheck™ PC** is substantially equivalent to **Coamatic Protein C** (K922201), manufactured by Chromogenix AB, and currently marketed in the United States by Instrumentation Laboratory. This opinion is based on the following similarities:

1. Both products are intended for use in the quantitative determination of Protein C activity in citrated human plasma
2. Both devices are based on synthetic chromogenic substrates
3. Both devices contain a synthetic chromogenic substrate and Protein C Activator and are reconstituted with distilled water
4. Both devices present results as a % activity of Protein C
5. Both devices are offered in a lyophilized format

**Correlation:**

Two lot numbers of ChromoCheck™ Protein C were compared to Coamatic Protein C in a correlation study using a mix of 60 individual normal and pathological patient samples. The following correlation was achieved:

Correlation parameter	ChromoCheck Lot 1	ChromoCheck Lot 2
Y-intercept	0.632	-1.308
Slope	0.979	1.025
R <sup>2</sup>	0.990	0.994

**Conclusion:** ChromoCheck™ Protein C is substantially equivalent to Coamatic Protein C.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Mr. Stephen L. Duff  
Director of New Business Development  
Precision BioLogic Inc.  
900 Windmill Road, Suite 100  
Dartmouth, Nova Scotia  
Canada B3B 1P7

MAR 19 2003

Re: k023990  
Trade/Device Name: ChromoCheck™ Protein C  
Regulation Number: 21 CFR § 864.7290  
Regulation Name: Quantitative Factor Deficiency Test  
Regulatory Class: II  
Product Code: GGP  
Dated: February 26, 2003  
Received: February 27, 2003

Dear Mr. Duff:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 –

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive style with a large initial 'S'.

Steven I. Gutman, M.D., M.B.A.  
Director  
Office of In Vitro Diagnostic Device Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

## Indications for Use Statement

510(k) Number: K023990

Device Name: ChromoCheck™ Protein C

### Indications for Use:

ChromoCheck™ Protein C is intended for use as an *in vitro* chromogenic assay for the quantitative determination of Protein C activity in citrated human plasma.

*Josephine Bantata*

Sign-Off  
Clinical Laboratory Devices  
Number K023990

Prescription Use