

August 12, 2020

Tyece Limited % Maria Griffin Senior Consultant mdi Consultants, Inc. 55 Northern Blvd, Ste 200 Great Neck, New York 11021

Re: K200838

Trade/Device Name: Tyece OTC TENS Device EM26, EM27, EM28, and EM29 Regulation Number: 21 CFR 882.5890 Regulation Name: Transcutaneous electrical nerve stimulator for pain relief Regulatory Class: Class II Product Code: NUH Dated: May 27, 2020 Received: June 1, 2020

Dear Maria Griffin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Vivek Pinto, PhD Director DHT5B: Division of Neuromodulation and Physical Medicine Devices OHT5: Office of Neurological and Physical Medicine Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K200838

Device Name

Tyece TENS Device Models: EM26, EM27, EM28, EM29

Indications for Use (Describe)

The Tyece OTC TENS Device EM26 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the shoulder and neck due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM27 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the lower extremities (leg) due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM28 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the upper extremities (arm) due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM29 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household and work activities.

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	⊠ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



Unit 803, Block A, Po Lung Centre 11 Wang Chiu Road, Kowloon Bay, Hong Kong Tel: (852) 2349 7456 . Fax (852) 2349 9166 www.tyece.com

510k Summary

The assigned 510(k) number is: **K200838**

 <u>Submitter's identifications:</u> Tyece Limited Unit 803, Block A, Po Lung Centre, 11 wang Chiu Road, Kowloon Bay, Hong Kong Contact: Parshid Falahati Telephone: +852 23497456 Email: <u>parshid@tyece.com</u>

Date of Summary Preparation: July 28, 2020

2. <u>Device</u>:

Trade Name: Tyece OTC TENS Device EM26, EM27, EM28 and EM29

Common Name: TENS Device Classification Name: Transcutaneous Electrical Nerve Stimulator for Pain Relief Product code: NUH Regulation class: II Regulation Number: 21 CFR 882.5890

3. Predicate Device

Trade/Device Name: LT3060 510(k) Number: K130802

4. <u>Reference Device</u>

Trade/Device Name: Savia OTC TENS Model EM-38 510(k) Number: K113321

5. <u>Device Description</u>:

The Tyece OTC TENS Device EM26 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the shoulder and neck due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM27 is intended for use by healthy adults for temporary relief of

pain associated with sore and aching muscles in the lower extremities (leg) due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM28 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the upper extremities (arm) due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM29 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM26, EM27, EM28 and EM29 powered by 4.5V (3 x 1.5V AAA /Alkaline batteries), is similar to the predicate device, LT3060 (K130802), with the following features:

- a. It is a portable single-channel, battery operated Transcutaneous Electrical Nerve Stimulator stimulation system. The predicate is a dual-channel TENS device.
- b. It contains 4 programs, similar to the predicate device which has 12 programs
- c. The output strength is adjustable from 0-110mA, similar to the predicate device at 0-96mA, via regulated voltage.
- d. The LCD display is provided for the indication of operation status including operation mode, output program mode, output intensity, time to cut-off, and battery low indication.

6. Intended Use:

The Tyece OTC TENS Device EM26 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the shoulder and neck due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM27 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the lower extremities (leg) due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM28 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the upper extremities (arm) due to strain from exercise or normal household and work activities.

The Tyece OTC TENS Device EM29 is intended for use by healthy adults for temporary relief of pain associated with sore and aching muscles in the upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household and work activities.

7. <u>Technological Comparison to Predicate Devices:</u>

Tyece OTC TENS Device EM26, EM27, EM28 and EM29 has similar technological characteristics to the predicate device in product design, material, energy source type, main program modes and the main output waveform etc. Through detailed calculation comparison

of stimulation output energy, we found that the output levels for the subject device and the predicate device are very similar and within acceptable ranges as specified in the FDA guidance. We believe that the differences between the two devices do not affect the determination of substantial equivalence.

8. <u>Discussion of Non-clinical Tests Performed for Determination of Substantial Equivalence are</u> <u>as follows</u>:

Non-clinical testing was performed in order to validate the design according to the company's specified design requirements, and to assure conformance with the following voluntary design standards. The device passes all testing performed below:

- IEC 60601-1: 2005 + AM1 (2012) Medical Devices Part 1: General Requirement for Safety Report
- IEC 60601-1-2:2014 Medical Devices Part 2: General Requirements for Safety Harmonized Standards: Electromagnetic Compatibility –Test and Requirement Report
- IEC 60601-1-11 Medical Electrical Equipment Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment Report
- IEC 60601-2-10 Medical electrical equipment Part 2-10: Particular requirements for the basic safety and essential performance of nerve and muscle stimulators Report
- Software Validation

The Tyece OTC TENS Device, Models EM26, EM27, EM28 and EM29, by applying the same stimulus parameters has the same intended use as the cleared predicate device, LT3060 (K130802).

A summary of the technological characteristics of Models EM26, EM27, EM28 and EM29 compared to the predicate device is given below:

Basic Unit Characteristics:

Item	Subject Device	Predicate Device	Discussion of differences
510(k) Number	K200838	K130802	Similar
Device Name, Model	Tyece OTC TENS Device EM26, EM27, EM28 and EM29	LT3060	Similar
Manufacturer	Savia Electronics (Shenzhen) Co. Ltd.	Shenzhen Dongdixin Technology Co., Ltd	Different but does not adversely impact safety and effectiveness of subject device

Intended Use	The Tyece OTC TENS Device EM26 is	Device:	Same
	intended for use by	The device is	
	healthy adults for	designed to be used	
	temporary relief of	for temporary relief	
	pain associated with	of pain associated	
	sore and aching	with sore and	
	muscles in the	aching muscles in	
	shoulder and neck	the shoulder, waist,	
	due to strain from	back, neck, upper	
	exercise or normal	extremities (arm),	
	household and work	and lower	
	activities.	extremities (leg) due	
	The Tyece OTC TENS	to strain from	
	Device EM27 is	exercise or normal	
	intended for use by	household work	
	healthy adults for	activities.	
	temporary relief of		
	pain associated with		
	sore and aching		
	muscles in the lower		
	extremities (leg) due		
	to strain from exercise		
	or normal household		
	and work activities.		
	The Tyece OTC TENS		
	, Device EM28 is		
	intended for use by		
	healthy adults for		
	temporary relief of		
	pain associated with		
	sore and aching		
	muscles in the upper		
	extremities (arm) due		
	to strain from exercise		
	or normal household		
	and work activities.		
	The Tyece OTC TENS		
	Device EM29 is		
	intended for use by		
	healthy adults for		
Prescription or OTC	OTC	отс	Same
Power Source	4.5V (3 x 1.5V AAA,	Battery	Different but does not
	type	powered, d.c.	adversely impact safety
	LR03 Alkaline	9.0V, one 6F22	and effectiveness of
	batteries)	battery	subject device

· · · · ·			Γ
-Method of Line Current Isolation	No line connection	No line connection	Same
	possible when	possible when	
	connected to body	connected to body	
-Patient Leakage Current	Not applicable, no	0.61uA	Different but does not
- Normal condition	line connection, no	0.68uA	adversely impact safety and
 Single fault condition 	AC charger		effectiveness of subject
	connection or		device
	operation.		
	Connection method		
	does		
	not allow AC		
	charger		
	connection to		
	patient		
Average DC current through			Same
electrodes when device is on	0μΑ	0μΑ	
but no pulses are being applied			
Number of Output Models	TENS	TENS	Same
Number of Output Channels	1	2	Different but does not
			adversely impact safety and
			effectiveness of subject
			device
-Synchronous or Alternating?	Synchronous	Alternating	Different but does not
			adversely impact safety and
			effectiveness of subject
-Method of Channel Isolation	Single Channel	Dual Channel	device Different but does not
	Single Channel	Dual Channel	adversely impact safety and
			effectiveness of subject
			device
Regulated Current or Regulated	Constant Voltage	Constant Voltage	Same
Software/Firmware/Micropr	YES	YES	Same
ocessor Control?			
Automatic Overload Trip?	YES	YES	Same
Automatic No-Load Trip?	YES	YES	Same
Automatic Shut Off?	YES	YES	Same
Patient Override Control?	YES	YES	Same
Indicator Display:			
- On/Off Status?	YES	YES	Same
- Low Battery?	YES	YES	
 Voltage/Current Level? Timer Range (minutes) 	YES (Voltage) 30 mins - Program A,	YES (Voltage) 1 - 60 mins	Different but does not
	25 mins - Programs B,		adversely impact safety and
	C, D		effectiveness of subject
			device
Compliance with Voluntary Standards?	60601-1 & 60601-1-2	60601-1 & 60601-1-	Same
		2	
Compliance with 21 CFR 898?	YES	YES	Same
	1	1	1

Weight (lbs., oz.)	EM26 - Approx. 0.6 lb (280g) with batteries EM27 - Approx. 0.34 lb (155g) with batteries EM28 - Approx. 0.39 lb (175g) with batteries EM29 - Approx. 0.46 lb (210g) with batteries	Approx. 0.28 lb with batteries	Different but does not adversely impact safety and effectiveness of subject device
Dimensions (in.) [W x H x D]	3.4 x 1.4 x 2.3 in	4.5 x 2.55 x 0.9 in	Different but does not adversely impact safety and effectiveness of subject device
Housing Materials and Construction	ABS	ABS	Same
Electrode Cable Length	1.2m	1.2m	Same

Output Specifications

Attribute		Subject Device:	Predicate Device:	Discussion of Difference
		Tyece OTC TENS	LT3060 (K130802)	
		Device, Models		
		EM26, EM27, EM28		
		and EM29		
Waveform		Symmetrical	Symmetrical	Same
		biphasic	biphasic	
Shape		Rectangular	Rectangular	Same
Channels		Single	Dual	Different but does not
				adversely impact safety
				and effectiveness of
				subject device
Maximum Output Voltage	@500Ω	54.4Vp-p	48Vp-p	Different but does not
(Volts)				adversely impact safety
(±20%)				and effectiveness of
	@2kΩ	112Vр-р	114Vp-p	subject device
	@10Ω	198Vp-p	115Vp-p	
Maximum Output Current	@500Ω	108.8mAp-p	96mAp-p	Different but does not
(mA)				adversely impact safety
(±20%)				and effectiveness of
	@2kΩ	56mAp-p	57mAp-p	subject device
	@10Ω	19.8mAp-p	11.5mAp-p	
Amplitude (mA)	<u> </u>	0-110	0-96	Different but does not
Frequency (Hz)		4-110	1-150	adversely impact safety
Pulse Width (µS)		60-220µs per phase	50-300µs per phase	and effectiveness of
				subject device

For multiphasic waveforms only	Symme phases	?	No multiphasic waveforms	Yes	Different but does not adversely impact safety and effectiveness of subject
	Phase [Duration		50-300uS	device.
Power ON indicator		LCD	LCD	Same	
Net Charge (µC pe	er pulse) (if zero	0.00246 μC@500Ω	0 μC@500Ω	Different but does not
state method of achieving zero net				adversely impact safety	
charge)					and effectiveness of
Maximum Phase (Charge (µ	C)	0.02398	0.0288	subject device
@500Ω					
Maximum Current	t Density	(mA/cm²)	0.52	1.15	
@500Ω					
Maximum Power	Density		0.033	0.373	
(W/ cm²) @500Ω	(using sm	allest			
electrode conduct	ive area)				
RMS Voltage (RMS	SV)	@500Ω	5.5V	5.2V	Different but does not
(±20%)				adversely impact safety	
		@2kΩ	16.3V	14.2V	and effectiveness of
		@10Ω	27.2V	20.2V	subject device
RMS Current (RM (±20%)	SA)	@500Ω	11mA	10mA	
		@2kΩ	7.35mA	8.85mA	
		@10Ω	2.02mA	1.88mA	
Burst Mode		•			Different but does not
- Pulses per burst		110	N/A	adversely impact safety	
- Burst per	second		2		and effectiveness of
- Burst duration (seconds)		0.5s		subject device	
- Duty cycle: Line (b) x Line (c)		1			
ON Time (seconds)		N/A	N/A	Same	
OFF Time (seconds)		N/A	N/A		
Additional Features (specify, if applicable)		Not applicable	Not applicable	Same	

Summary for the technology comparison

The Tyece OTC TENS Device EM26, EM27, EM28 and EM29 has similar technological characteristics to the predicate device in product design, material, energy source type, main program modes and the main output waveform etc. Through detailed calculation comparison of stimulation output energy, we found that the output levels for the subject device and the predicate device are very similar and within acceptable ranges as specified in the FDA guidance. We believe that the differences between the two devices do not affect the determination of substantial equivalence.

Conclusion

The Tyece OTC TENS Device EM26, EM27, EM28 and EM29 has the same intended use and similar technological characteristics as the cleared predicate device. Moreover, verification and validation tests demonstrate that the differences in the submitted models maintain the same safety and effectiveness as that of the cleared device. Therefore, we conclude that the subject device is substantially equivalent to the predicate device.