



June 18, 2020

Immucor GTI Diagnostics, Inc.  
Allison Stray  
Sr. Director QA/RA  
20925 Crossroads Circle  
Waukesha, Wisconsin 53185

Re: K201311

Trade/Device Name: PF4 IgG assay  
Regulation Number: 21 CFR 864.7695  
Regulation Name: Platelet Factor 4 Radioimmunoassay  
Regulatory Class: Class II  
Product Code: LCO  
Dated: May 15, 2020  
Received: May 18, 2020

Dear Allison Stray:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Takeesha Taylor-Bell  
Chief  
Division of Immunology  
and Hematology Devices  
OHT7: Office of In Vitro Diagnostics  
and Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

Device Name

PF4 IgG

Indications for Use (Describe)

PF4 IgG assay is a qualitative screening assay for the detection of heparin associated IgG antibodies in human serum or plasma. The presence of heparin-associated antibodies are commonly found in patients with Heparin Induced Thrombocytopenia (HIT).

IVD

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## Section 5.0 510(k) Summary

### 1. Submitter

Owner's Name: Immucor GTI Diagnostics, Inc.  
Address: 20925 Crossroads Circle, Waukesha, WI 53186  
Phone: (262) 754-1000  
Fax: (262) 754-9831  
Establishment Registration Number: 2183608  
Name of Contact Person: Allison Stray (astray@immucor.com)  
Date Summary Prepared: May 15, 2020

### 2. Name of Device

Device Name: PF4/PVS Complex IgG Antibody Detection Assay  
Proprietary Name: PF4 IgG  
Common Name: Platelet Factor 4 Assay  
Regulation Number: 21CFR 864.7695  
Classification: Class II  
Product Code: LCO  
Establishment Registration Number: 2183608

### 3. Name of Predicate Device for Claiming Equivalence

PF4 IgG (K071781)

### 4. Description of Device

PF4 IgG assay is an enzyme linked immunosorbent assay (ELISA). The PF4 IgG ELISA is intended to detect IgG antibodies in human serum or plasma that react with Platelet Factor 4 (PF4) when it is complexed to heparin or other polyanionic compounds. The PF4 IgG kit contains all of the reagents necessary to perform the assay.

Antibodies that react with PF4 when it is complexed to heparin are found in some patients undergoing heparin therapy. The presence of these antibodies has been shown to be associated with heparin induced thrombocytopenia Type II (Type II HIT). Heparin Induced Thrombocytopenia (HIT) is an adverse reaction that can occur in patients that are undergoing heparin therapy. HIT occurs in up to 3% of patients that are receiving heparin and can result from treatment with either low molecular weight heparin (LMWH) or unfractionated heparin (UFH). HIT is usually associated with thrombocytopenia and in some cases the more severe complications of arterial or venous thrombosis. The pathophysiology of HIT is linked to the formation of antibodies which are specific for epitopes formed when heparin binds to PF4. These antibodies can subsequently bind to platelets via the Fc $\gamma$ R11a receptor resulting in platelet activation.

It is now known that antibodies which bind to PF4/heparin complexes also bind to PF4 when it is complexed with other polyanionic compounds such as polyvinyl sulfonate (PVS). In the PF4 IgG assay, a complex of PF4/PVS, which has been immobilized in the microwells serve as a target for the binding of antibodies associated with Type II HIT.

#### Test Procedure

Patient sample is added to microwells coated with platelet factor 4 (PF4) complexed to polyvinyl sulfonate (PVS). If an antibody recognizing a site on PF4:PVS is present, binding will occur. Unbound antibodies are then washed away. An alkaline phosphatase labeled anti-human globulin reagent (Anti-IgG) is added to the wells and incubated. The unbound Anti-IgG is washed away and the substrate PNPP (p-nitrophenyl phosphate) is added. After a 30-minute incubation period, the reaction is stopped with Stopping Solution. The optical density of the color that develops is measured in a spectrophotometer.

Engineering drawings or schematics are not applicable to the PF4 IgG Assay. The device consists of reagents and drawings are not pertinent to describe the device.

#### **5. Intended Use**

PF4 IgG is a qualitative screening assay for the detection of heparin associated IgG antibodies in human serum or plasma.

#### **6. Substantial Equivalence**

This Special 510(k) is for a Positive Serum Control raw material change and other changes the prior company did not file. The Positive Serum Control raw material modification in this submission is from human serum containing Bovine Albumin and 0.1% NaN<sub>3</sub> to human serum containing Bovine Albumin and 0.1% NaN<sub>3</sub> plus a chimeric (mouse Fab/Human Fc) monoclonal antibody in the PF4 IgG positive control reagent.

The modification of the PF4 IgG and need for this submission is due to supply of the predicate PF4 IgG positive control material. The current supply of the PF4 IgG positive control material is depleting and will not be available in the future.

Two (2) additional changes, the prior company, did not file are being included. Those changes include: 1) Addition of plasma as a sample type, 2) Stopping Solution changed from 3M NaOH to a buffered EDTA solution. The prior company did document a 510(k) Letter to File for both the sample type change and Stopping Solution change. Their rationale for not filing with the Food and Drug Administration included the changes did not have any significant impact to the safety or effectiveness of the device.

There were no design or functionality changes. All elements of the Design History File (DHF) for the proposed device were used as possible and design outputs were created for the proposed device.

The Similarities between the Predicate and Candidate Device:

- The intended use has not changed
- The indications for use have not changed
- The technology (ELISA) and assay steps have not changed
- Reportable results have not changed
- The packaging configuration has not changed
- There are no design or functionality changes between the predicate and candidate device
- All reagents are identical with the exception of the Positive Serum Control and the Stopping Solution

The Differences between the Predicate and Candidate Device:

- The predicate device only allowed for Serum as the specimen collection type. The candidate device allows for without anticoagulant (Serum) and ACD or Sodium Citrate (Plasma) as the specimen collection type.
- The formulation of the Stopping Solution for PF4 IgG is different than the predicate device. The PF4 IgG Stopping Solution is a buffered EDTA solution which is ready for use, whereas the predicate device utilizes 3M NaOH.



- The candidate device contains Positive Serum Control as Human serum containing Bovine Albumin and 0.1% NaN<sub>3</sub>. The predicate device contains Positive Serum Control as Human serum containing Bovine Albumin and 0.1% NaN<sub>3</sub> plus chimeric (mouse Fab/Human Fc) monoclonal antibody.
- The candidate device has changes to the IFU as a result of the specimen collection type change and Stopping Solution type change.

The table below provides the comparison between the current version of PF4 IgG and the proposed device.

#	Features / Characteristics	Predicate Device	Candidate Device	Comments
1	Intended Use	PF4 IgG is a qualitative screening assay for the detection of heparin associated IgG antibodies in human serum.	PF4 IgG is a qualitative screening assay for the detection of heparin associated IgG antibodies in human serum or plasma.	The intended use of the qualitative assay to detect IgG antibodies has not changed. The revised verbiage adds plasma as a specimen type.
2	Indications for Use	PF4 IgG is a qualitative screening assay for the detection of heparin associated IgG antibodies in human serum. The presence of heparin-associated antibodies are commonly found in patients with Heparin Induced Thrombocytopenia (HIT).	PF4 IgG is a qualitative screening assay for the detection of heparin associated IgG antibodies in human serum or plasma. The presence of heparin-associated antibodies are commonly found in patients with Heparin Induced Thrombocytopenia (HIT).	N/A
3	Technology	ELISA with a colorimetric measurement system	Same	N/A
4	Specimen Collection	Serum	Without anticoagulant (Serum) and ACD or Sodium Citrate (Plasma)	Prior company made this change in 2010
5	Reportable Results	Qualitative assay; results are reported as positive or negative	Same	N/A
6	Packaging Configuration	13 and 45 Test kits	Same	N/A



#	Features / Characteristics	Predicate Device	Candidate Device	Comments
7	Reagents			
	Microwell Strips	Immobilized PF4/PVS Complex	Same	N/A
	Concentrated Wash Solution	10X Tris Buffer, NaCl, Tween 20, 1% NaN <sub>3</sub>	Same	N/A
	Specimen Diluent	Phosphate Buffered Saline, 0.05% NaN <sub>3</sub>	Same	N/A
	Substrate Buffer	Diethanolamine and magnesium chloride, 0.02% NaN <sub>3</sub>	Same	N/A
	Substrate	PNPP (crystalline powder)	Same	N/A
	Stopping Solution	3 M NaOH	Buffered EDTA Solution	Prior company made this change in 2010
	Positive Serum Control	Human serum containing Bovine Albumin and 0.1% NaN <sub>3</sub>	Human serum containing Bovine Albumin and 0.1% NaN <sub>3</sub> plus Chimeric (mouse Fab/Human Fc) monoclonal antibody	N/A
	Negative Serum Control	Human serum containing 0.1% NaN <sub>3</sub>	Same	N/A
	Conjugate	Goat anti-human IgG conjugated to alkaline phosphatase enzyme	Same	N/A
8	Instructions for Use (IFU)	PF4 IgG Direction Insert	No changes to the IFU for the Positive Control Serum raw material modification. IFU changes for the addition of plasma, new Stopping Solution and minor changes are included in Section 13.0, Proposed Labeling.	N/A

## 7. Company Name Change

Genetic Testing Institute Inc. was purchased by Gen-Probe Inc. in December 2010 and the company name was changed to Gen-Probe GTI Diagnostics, Inc. Immucor purchased Gen-Probe GTI Diagnostics Inc. in March 2013. The company name was subsequently changed to Immucor GTI Diagnostics, Inc.

## 8. Non-Clinical Studies Supporting Safety and Effectiveness of the Candidate Device

The details of each of the following studies are covered in Sections 14.0 (Sterilization and Shelf Life) and 18.0 (Performance Testing Bench). Only a brief summary of these studies is provided in this section. The manufacturing facility is in conformance with the design control procedure requirements as specified in 21 CFR 820.30 and the records are available for review.

### 8.1 Modification, Positive Serum Control Raw Material Change

#### Process Validation

Three validation batches of PF4 Positive Serum Control were manufactured for this process validation.

The three validation batches were tested per the PF4 IgG Instructions For Use (IFU-303290.IFUEN, Rev. E) and all acceptance criteria were met.

Based upon the data generated, the PF4 Positive Serum Control manufacturing process consistently delivers product meeting established acceptance criteria and is acceptable for use in the production of the PF4 IgG assay.

#### Precision Study

Testing was performed by 3 operators and five assays were conducted by each operator.

Testing was completed using a single lot of the PF4 IgG kit. Each of the three validation lots of the Positive Serum Control were tested concurrently on each of the 5 assays (days).

All results met the IFU acceptance criteria and were 100% positive, demonstrating 100% agreement within run and between run. All pre-defined Acceptance Criteria were met.

### 8.2 Modification, Addition of Plasma as a Sample Type

#### Comparison of Serum and Plasma (ACD or Sodium Citrate) Sample Types

Samples were obtained from Florida Hospital to perform an internal study. 174 frozen serum and plasma samples were used for the study and testing was performed according to the PF4 IgG Instructions For Use.

98.85% Agreement (172/174) was met between the qualitative results obtained with serum and plasma. A linear regression was performed comparing the OD values obtained from serum and ACD or Sodium Citrate (plasma) samples for PF4 IgG. The regression analysis of plasma versus serum for the PF4 IgG assay demonstrated very good correlation between OD values ( $R^2 = 99.3\%$ ). In addition, there was very good agreement between the OD values (slope = 0.9476 and intercept = -0.009914).

Equivalent results using serum or plasma (ACD or Sodium Citrate) were obtained in the PF4 IgG assay. Both sample types are acceptable for use in the PF4 IgG assay.

#### Comparison of ACD and Sodium Citrate (plasma) to Serum

The objective of the internal study was to demonstrate plasma samples collected in ACD or Sodium Citrate give equivalent results to serum.

Blood was collected from 24 individuals known to be PF4: heparin antibody negative in ACD (yellow top), Serum (red top) or Sodium Citrate (blue top) tubes. The serum or plasma were separated and equal volume aliquots of each were prepared.

A high level reacting anti:PF4 heparin antibody sample was used for the spiking study. A dilution (1:100) intended to give a medium positive result and a dilution (1:1000) intended to give a low positive result were chosen. The samples were also run without spiking to ensure each sample showed negative reactivity.

The results for all samples tested gave equivalent results to the serum samples. The use of ACD or Sodium Citrate (plasma) samples in the PF4 IgG assay gives equivalent results to those of serum samples.

### **8.3 Modification, Stopping Solution Change**

#### **Lot to Lot Comparison**

A lot to lot comparison study was performed to demonstrate the new Stopping Solution (ESS) performs equivalent to the predicate Stopping Solution (3 M NaOH) meeting established acceptance criteria.

Three lots of ESS (referred to as the test solution) were manufactured according to SOP-292 and tested per QAP-017A. The predicate Stopping Solution, 3 M NaOH, was used as a reference in the study.

All acceptance criteria were met demonstrating the new Stopping Solution (ESS) is acceptable to use in the PF4 IgG assay.

#### **Process Validation**

Three lots of the new Stopping Solution (ESS) were manufactured according to SOP-292 and tested with the PF4 IgG assay according to the IFU.

Well-to-well variation and expected qualitative results were evaluated to determine if the new Stopping Solution (ESS) was equivalent to the predicate 3 M NaOH Stopping Solution. The predicate Stopping Solution, 3 M NaOH, was used as the reference in the study.

For Well-to-well variation, all wells met acceptance criteria. For qualitative results, the finished product QC panel was tested with all four lots of Stopping Solution. All 21 positive samples and 24 negative samples met their expected qualitative result on all three lots of the new Stopping Solution and the one lot of 3M NaOH Stopping Solution. All assays met the QC requirements in the IFU.

The data demonstrate the new Stopping Solution (ESS) can be used in the PF4 IgG assay.

### **8.4 Shelf Life**

The PF4 IgG assay has a shelf life of 24 months. See Section 14.0, Sterilization and Shelf Life for stability studies performed to support this premarket submission. All acceptance criteria were met.

### **Conclusion**

Based on the non-clinical performance and stability studies, the data demonstrate the modifications to the PF4 IgG device does not present new issues of safety and effectiveness.