



Scientia Vascular LLC  
Amy McManus  
Regulatory Affairs Manager  
3487 West 2100 South Suite 100  
West Valley City, UT 84119

July 29, 2020

Re: K201760  
Trade/Device Name: Zoom 14 Guidewire  
Regulation Number: 21 CFR 870.1330  
Regulation Name: Catheter Guide Wire  
Regulatory Class: Class II  
Product Code: DQX, MOF  
Dated: June 23, 2020  
Received: June 29, 2020

Dear Amy McManus:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Naira Muradyan, Ph.D.  
Assistant Director (Acting)  
DHT5A: Division of Neurosurgical,  
Neurointerventional  
and Neurodiagnostic Devices  
OHT5: Office of Neurological  
and Physical Medicine Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K201760

Device Name  
Zoom 14 Guidewire

Indications for Use (Describe)

The Zoom 14 Guidewire is intended for general vascular use within the neuro and peripheral vasculatures to introduce and position catheters and other interventional devices. The guidewire is not intended for use in the coronary vasculature.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



**510(K) SUMMARY**

**(PER 21 CFR 807.92)**

**SCIENTIA VASCULAR LLC**

**Special 510(K): Device Modification  
Zoom 14 Guidewire**

510(k) Sponsor: Scientia Vascular LLC  
3487 West 2100 South Suite 100  
West Valley City, UT 84119  
Tel (888) 385.9016

Contact Person: Amy McManus,  
Regulatory Affairs Manager  
Tel: (888) 385.9016  
E-mail: amcmanus@scientiavascular.com

Date Prepared: June 10, 2020

Subject Device Information:

---

Trade Name:	Zoom 14 Guidewire
Common Name:	Guidewire
Classification Name	Catheter Guide Wire per 21 CFR 870.1330
Primary Product Code:	DQX
Secondary Product Code	MOF
Predicate Device:	Aristotle 14 Guidewire (K173235)

## **DEVICE DESCRIPTION**

The Zoom 14 Guidewire is a modification of Scientia Vascular's Aristotle 14 Guidewire. It is a 0.014" diameter steerable guidewire with a shapeable tip to aid in accessing vasculature. The guidewire is supplied sterile and is for single use only. It is provided in a range of stiffness profiles: support and extra support. The product is provided in lengths of 200cm or 300cm.

The distal portion of the guidewire tip includes a radiopaque platinum wire marker coil to facilitate fluoroscopic visualization. The guidewire has a hydrophilic polymer coating on the distal portion and a polytetrafluoroethylene (PTFE) coating on the proximal portion to reduce friction during manipulation in vessels.

The guidewire is provided with a shaping mandrel, an introducer (to aid with the insertion of the guidewire into a catheter hub and/or a hemostasis valve) and a torque device (to attach to the proximal portion to facilitate gripping and manipulation of the guidewire during use). The mandrel, introducer and torque accessory devices are included to facilitate use of the guidewire and are not intended to contact the patient's body.

The Zoom 14 Guidewire is substantially equivalent with respect to technological characteristics, design and materials to the previously cleared Aristotle 14 Guidewire.

## **INDICATION FOR USE**

The Zoom 14 Guidewire is intended for general vascular use within the neuro and peripheral vasculatures to introduce and position catheters and other interventional devices. The guidewire is not intended for use in the coronary vasculature.

This is the same indication for use as previously cleared for the Aristotle 14 Guidewire, K173235.

## **INTENDED USE**

The Zoom 14 Guidewire is intended for use by a physician to help introduce and position catheters or other interventional devices within the neuro and peripheral vasculature.

## **TECHNOLOGICAL CHARACTERISTICS**

The Zoom 14 Guidewire has the following similarities to the previously cleared Aristotle 14 Guidewire:

- Both devices have the same indicated use,
- Both devices have the same intended use
- Both devices use the same operating principle,
- Both devices incorporate the same basic guidewire design,
- Both devices incorporate the same materials, and

- Both devices are packaged and sterilized using the same materials and processes.

The technological characteristic comparison of the subject and predicate device is summarized in Table 1 below.

<b>Table 1: Comparison between Subject &amp; Predicate Device Technological Characteristics:</b>		
<b>Note: Areas of difference will be indicated in RED</b>		
Characteristic	Subject Device Zoom 14 Guidewire	Predicate Aristotle 14 Guidewire (K173235)
Anatomical Location	Neuro and peripheral vasculature	Neuro and peripheral vasculature
Dimensions	<i>Max O.D.: 0.014" (0.36mm)</i> <i>Length: 200cm and 300 cm</i>	<i>Max O.D.: 0.014" (0.36mm)</i> <i>Length: 200cm and 300 cm</i>
Core Wire	Stainless Steel	Stainless Steel
Distal Tip	Shapeable <i>Length: 35cm</i> <i>Material: Nitinol</i>	Shapeable <i>Length: 35cm</i> <i>Material: Nitinol</i>
Stiffness Profiles	Range from <b>extra support (stiff) to support (less stiff)</b>	Range from <b>support (stiff) to soft (flex)</b>
Coatings	<i>Distal End: Hydrophilic</i> <i>Proximal End: PTFE</i>	<i>Distal End: Hydrophilic</i> <i>Proximal End: PTFE</i>
Radiopaque Marker	Radiopaque marker at distal tip	Radiopaque marker at distal tip
Centering Coil	One (1) Centering coil	One (1) Centering coil
Guidewire Introducer (Accessory)	Provided with each guidewire	Provided with each guidewire
Torque Device (Accessory)	Provided with each guidewire	Provided with each guidewire
Sterilization Method	100% Ethylene Oxide (EO)	100% Ethylene Oxide (EO)

The subject device, Zoom 14 Guidewire, has one technological characteristic difference when compared to the predicate device. This differing technological characteristic does not result in new materials used or new questions of safety or effectiveness, nor does it result in new risks for the subject device. Testing of the subject device has been performed, demonstrating the safety and effectiveness of the device.

## **NON-CLINICAL PERFORMANCE TESTS**

### **Biocompatibility**

The materials used in the manufacture of the subject device Zoom 14 Guidewire are identical to those used in the manufacturing of the predicate device Aristotle 14 Guidewire, as seen in the table below.

The Zoom 14 Guidewire uses the same manufacturing assembly equipment and processes as the Aristotle 14 Guidewire. The Zoom 14 Guidewire and the Aristotle 14 Guidewire differ in the stiffness profile, with no overall dimensional or material differences. Since the Zoom 14 is identical to the Aristotle 14 in this respect, with no new materials or processes, the product has already been tested for biocompatibility and chemical analysis.

<b>Characteristic</b>	<b>Subject Device Zoom 14 Guidewire</b>	<b>Aristotle 14 Guidewire</b>
Overall Dimensions	NO CHANGE	O.D.: 0.014" (0.36mm) Length: 200cm to 300cm
Core Wire	NO CHANGE	Stainless Steel
Distal Tip	NO CHANGE	Length: 35cm Material: Nitinol
Coatings	NO CHANGE	Distal End: Hydrophilic Proximal End: PTFE
Radiopaque Marker	NO CHANGE	1 radiopaque marker at distal tip
Centering Coil	NO CHANGE	1 centering coil
Sterilization Method	NO CHANGE	100% Ethylene Oxide (EO)

### **Sterilization**

100% EO is used to sterilize the device to achieve a SAL of at least  $10^{-6}$ . The Zoom 14 Guidewire and the Aristotle 14 are similar with regards to all device and packaging characteristics that could affect the ability to sterilize the devices.

### **Functional Testing**

Functional testing on the subject device was performed after conducting a risk assessment in accordance with *EN ISO 14971:2012 Medical Devices – Risk Management*. This testing was performed in accordance with ISO 11070:2014 *Sterile single-use intravascular introducers, dilators and guidewires* and the FDA Guidance Document *Coronary, Peripheral and Neurovascular Guidewires – Performance Tests and Recommended Labeling (2019)*. Table 2 summarizes the functional tests performed and test results obtained to demonstrate substantial equivalence to the predicate device:

**Table 2: Summaries of Functional Tests Conducted to Support this Premarket Notification for the Modified (Subject) Device, Zoom 14 Guidewire**

Test	Test Method Summary	Results
Visual Inspection	Tests per ISO 11070: Visual inspection per engineering drawings.	The Zoom 14 Guidewires met testing acceptance criteria.
Dimensional Verification	Tests per ISO 11070: Dimensional inspection per engineering drawings.	The Zoom 14 Guidewires met testing acceptance criteria.
Tensile Testing	Tensile testing per ISO 11070.	The Zoom 14 Guidewires met testing acceptance criteria
Column Buckling	The force required to buckle the tip of the guidewire at 5mm, 10mm, 20mm.	The Zoom 14 Guidewires met testing acceptance criteria
Flexing	Tests per ISO 11070: Inspection for defects and damage or flaking of the coating after flexing.	The Zoom 14 Guidewires met testing acceptance criteria
Fracture	Tests per ISO 11070: Inspection for fracture, loosening, or failure after wrapping around mandrel.	The Zoom 14 Guidewires met testing acceptance criteria
Torqueability	Measurement of torque response (average input to output lag) in an anatomical model.	The Zoom 14 Guidewires met testing acceptance criteria
Torque Strength	Torque turns to failure in an anatomical model.	The Zoom 14 Guidewires met testing acceptance criteria
Model Evaluation	Anatomical model designed to simulate the tortuous anatomy of the neurovasculature used for simulated use testing	The Zoom 14 Guidewires met testing acceptance criteria

**CONCLUSION:**

Scientia Vascular, LLC has presented information in this premarket notification supporting the Zoom 14 Guidewire as substantially equivalent with respect to technological characteristics and indications for use to the predicate device.