

September 30, 2020

Gi Supply, Inc. Erika Parry Manager, Quality and Regulatory 5069 Ritter Road Suite 104 Mechanicsburg, Pennsylvania 17055

Re: K202376

Trade/Device Name: EverLift Submucosal Lifting Agent Regulation Number: 21 CFR 876.1500 Regulation Name: Endoscope And Accessories Regulatory Class: Class II Product Code: PLL Dated: September 1, 2020 Received: September 2, 2020

Dear Erika Parry:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Shanil P. Haugen, Ph.D.
Assistant Director
DHT3A: Division of Renal, Gastrointestinal, Obesity and Transplant Devices
OHT3: Office of GastroRenal, ObGyn, General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known) K202376

Device Name EverLift Submucosal Lifting Agent

Indications for Use (Describe)

EverLift Submucosal Lifting Agent is indicated for use in gastrointestinal endoscopic procedures for submucosal lift of polyps, adenomas, early-stage cancers, or other gastrointestinal lesions prior to excision with a snare or other appropriate endoscopic device.

Type of Use (Select one or both, as applicable)	
🛛 Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

Uver-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

> Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

I. SUBMITTER [Per 807.92(a)(1)]

Sponsor/Manufacturer GI Supply 5069 Ritter Road Suite 104 Mechanicsburg, PA 17055 USA Phone: (800)-451-5797

<u>Contact Person</u> Erika Parry Manager, Quality and Regulatory Phone: (717)-562-7580 Email: <u>e.parry@gi-supply.com</u>

Date Prepared

August 19, 2020 Amended September 29, 2020

II. <u>DEVICE</u> [Per 807.92(a)(2)]

Device Trade/Proprietary Name	EverLift™ Submucosal Lifting Agent	
Device Common or Usual Name	Submucosal Injection Agent	
Regulation Number	21 CFR 876.1500	
Regulation Name	Endoscope and Accessories	
Regulatory Class	Class II	
Submission Type	Special 510(k)	
Product Code	PLL	
Classification Panel	Gastroenterology/Urology	

III. PREDICATE DEVICE [Per 807.92(a)(3)]

The EverLift[™] Submucosal Lifting Agent (10mL) [subject device] is substantially equivalent in terms of its intended use to the claimed predicate device, the 5mL EverLift[™] Submucosal Lifting Agent (5mL) (K191923) with respect to device design, fundamental technology, physical characteristics, performance and intended use.

The device design and fundamental technology of the subject device are nearly identical to that of the predicate device. The formulation and composition of the primary packaging materials are identical. There are slight differences in the dimensions of the subject device and predicate device which do not impact safety or effectiveness. Both the subject and predicate devices are delivered sterile and are indicated for single-use only. The Sterilization Method and Shelf Life for the subject device are identical as that for the predicate device.

Predicate	The EverLift™ Submucosal Lifting Agent (10mL) [subject device] is substantially equivalent to	
Device	the following predicate device manufactured by GI Supply:	
	 EverLift[™] Submucosal Lifting Agent (5mL) (K191923) 	

IV. DEVICE DESCRIPTION [Per 807.92(a)(4)]

The GI Supply EverLift[™] Submucosal Lifting Agent (10mL) is a prefilled plastic syringe with attached plunger rod containing 10mL of lifting agent. The syringe has a luer lock connection capable of interfacing with a standard, commercially available endoscopic injection needle.

EverLift[™] Submucosal Lifting Agent is an injectable liquid composition for use as a submucosal injection agent during endoscopic mucosal resection (EMR), endoscopic mucosal dissection (ESD), and polypectomy procedures in the gastrointestinal tract. The device is intended for use in endoscopic resection procedures in the upper and the lower gastrointestinal tract, including the esophagus, the stomach, the small intestine, the colon, the sigmoid colon, and the rectum, as a submucosal injectable agent during the removal of polyps, adenomas, early-stage cancers, and other pathological lesions by EMR, ESD, or polypectomy.

EverLift[™] Submucosal Lifting Agent is injected into the submucosal layer by means of a standard, commercially available endoscopic injection needle, which is inserted into the working channel of the endoscope. The composition, when injected, creates a cushion in situ by lifting the gastrointestinal mucosa from the submucosal layer, allowing the endoscopist to perform an easy and safe resection procedure (EMR, ESD, or polypectomy).

V. INTENDED USE / INDICATIONS FOR USE [Per 807.92(a)(5)]

Intended Use / Indications for Use

EverLift[™] Submucosal Lifting Agent is indicated for use in gastrointestinal endoscopic procedures for submucosal lift of polyps, adenomas, early-stage cancers, or other gastrointestinal lesions prior to excision with a snare or other appropriate endoscopic device.

VI. <u>COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE</u> [Per 807.92(a)(6)]

The EverLift[™] Submucosal Lifting Agent (10mL) [subject device] is substantially equivalent to the GI Supply EverLift[™] Submucosal Lifting Agent (5mL) (K191923) [predicate device] based on the same indication for use as the predicate device and the similar or identical functional and performance characteristics of the subject device when compared to the predicate device. The differences between the subject device and predicate device do not raise different issues of safety and effectiveness.

The detailed substantial equivalence comparison of the similarities and differences between the GI Supply EverLift[™] Submucosal Lifting Agent (10ml) [subject device] and the GI Supply EverLift[™] Submucosal Lifting Agent (5mL) (K191923) [predicate device] is provided in the table below.

Regulatory			Similarities /
Information	[Subject Device]	[Predicate Device]	Differences
Manufacturer	GI Supply	GI Supply	Same
Device Trade			
or	EverLift [™] Submucosal Lifting Agent	EverLift™ Submucosal Lifting Agent	
Proprietary	(10mL)	(5mL)	
Name			

Regulatory Information	[Subject Device]	[Predicate Device]	Similarities / Differences
510(k)	K202376	K191923	
Number	N202370		
Device Class	Class II	Class II	Same
Device			
Classification	Endoscope and Accessories	Endoscope and Accessories	Same
Name			
Device	Submussed Injection Agent	Submussed Injection Agent	Sama
Namo	Submucosal Injection Agent	Submucosal Injection Agent	Same
Product Code	DII	DLI	Same
Regulation			Same
Number	21 CFR 876.1500	21 CFR 876.1500	Same
	Design Features and Ca	apabilities of the Device	
	EverLift™ Submucosal Lifting Agent	EverLift™ Submucosal Lifting Agent	
	is indicated for use in	is indicated for use in	
	gastrointestinal endoscopic	gastrointestinal endoscopic	
Indications	procedures for submucosal lift of	procedures for submucosal lift of	
for Use	polyps, adenomas, early-stage	polyps, adenomas, early-stage	Same
101 032	cancers, or other gastrointestinal	cancers, or other gastrointestinal	
	lesions prior to excision with a	lesions prior to excision with a	
	snare or other appropriate	snare or other appropriate	
	endoscopic device.	endoscopic device.	-
Intended Use	Same as Indications for Use	Same as Indications for Use	Same
Prescription			
or Over-the-	Prescription Use	Prescription Use	Same
Environment	Hospital / Clinic	Hospital / Clinic	Same
Sterile	Sterile	Sterile	Same
Single-Use	Single-Use Only	Single-Use Only	Same
	Design	Features	
	Each individual 10mL COC syringe is	Each individual 5mL COC syringe is	Different. The
	packed inside a thermoformed tray.	packed along with three (3) patient	differences in
	The tray is sealed with a labeled lid	labels inside a propionate tube with	packaging do
	stock with three (3) patient labels	a friction fit press on closure. Ten	not raise
	affixed to the lid stock. Ten (10) of	(10) tubes are packed into a shelf	different
	the sealed trays are packed into a	box on top of a copy of the IFU. The	questions of
Packaging	shelf box with a copy of the IFU	shelf box is shipped in a regular die	safety or
	placed between the trays and the	cut and scored 200-pound B-flute	effectiveness.
	shelf box wall. The shelf box is	kraft shipper box sealed with a	A Distribution
	shipped in a regular die cut and	shipper box label.	Study was
	scored 200-pound B-flute kraft		pertormed on
	shipper box sealed with a shipper		the new
	DOX IADEI.		packaging
			configuration.

Regulatory Information	[Subject Device]	[Predicate Device]	Similarities / Differences
Images of Device			
Sterilization Method	Sterilized by moist-heat (steam) sterilization	Sterilized by moist-heat (steam) sterilization	Same
Shelf Life	2-Year Shelf Life	2-Year Shelf Life	Same
Composition	Each steam sterilized syringe contains 10mL of lifting agent with the following ingredients: • Water • Hydroxyethyl Cellulose • Glycerin • Methylene Blue • Benzyl Alcohol • Sodium Phosphate • Potassium phosphate	Each steam sterilized syringe contains 5mL of lifting agent with the following ingredients: • Water • Hydroxyethyl Cellulose • Glycerin • Methylene Blue • Benzyl Alcohol • Sodium Phosphate • Potassium Phosphate	Same concentration; however, volume of lifting agent provided differs.

VII. <u>SUMMARY OF PERFORMANCE DATA AND PERFORMANCE TEST CONCLUSIONS</u> [Per 807.92(b)(1)(2)(3)]

The determination of substantial equivalence is based on an assessment of non-clinical performance data. To verify that the device design meets its functional and performance requirements, the EverLift[™] Submucosal Lifting Agent (10mL) underwent performance testing. Biocompatibility testing was performed only on the 5mL variant as it represented worst-case surface area to volume ratio.

The EverLift[™] Submucosal Lifting Agent was developed under GI Supply's risk management process in accordance with ISO 14971:2012 - Medical Devices-Application of Risk Management to Medical Devices. The identified risks were adequately mitigated and verified by means of non-clinical performance testing.

Summary of Performance Testing

A series of non-clinical performance testing was conducted on the subject device. Please refer to the table below for a summary of all non-clinical performance testing conducted on the subject device in support of substantial equivalence.

Test Description	Conforming Standard(s)	Conclusion in Support of Substantial Equivalence
Syringe Tip Cap	None	Syringe tip cap removal force testing was performed on
Removal Force		the subject device using the same method and
Testing		acceptance criteria as the predicate device. Results of
		the testing confirmed that the syringe tip cap removal
		force is equivalent between the subject and predicate
		devices, thereby demonstrating substantial equivalence.

Test Description	Conforming Standard(s)	Conclusion in Support of Substantial Equivalence
Product Color	None	Product color testing was performed on the subject
Testing		device using the same method and acceptance criteria as
		the predicate device. Results of the testing confirmed
		that product color is identical between the two devices,
		thereby demonstrating substantial equivalence.
Injection Flow	None	Injection flow rate testing was performed on the subject
Rate Testing		device using the same method and acceptance criteria as
		the predicate device. Results of the testing confirmed
		that the injection flow rate is equivalent between the
		subject and predicate devices, thereby demonstrating
		substantial equivalence.
Graduation	ISO 7886-1:2017, Sterile	Graduation mark tolerance testing was performed on
Marking	hypodermic needles for single	the subject device using the same method as the
Tolerance	use – Part 1: Syringes for manual	predicate device, and acceptance criteria as defined in
	use	ISO 7886-1. Results of the testing confirmed that the
		marking tolerance meets the acceptance criteria defined
		In ISO 7886-1, thereby demonstrating substantial
Containan	Neze	equivalence.
Container	None	Container closure integrity testing was performed on the
Integrity Testing		criteria as the predicate device. Pocults of the testing
integrity resting		confirmed that the container closure integrity is identical
		between the subject and predicate devices, thereby
		demonstrating substantial equivalence
Tray Lid Peel	None	The secondary packaging of the subject device required
Force Testing		a minor design change to accommodate the larger 10ml
		syringe. Tray lid peel force test results post-distribution
		simulation met acceptance criteria, confirming that the
		differences in secondary packaging do not raise different
		questions of safety or effectiveness.
Sterility	ISO 17665-1, Sterilization of	A sterilization validation was performed on the subject
Assurance Level	healthcare products – Moist	device using the same sterilization process and
	heat – Part 1: Requirements for	sterilization parameters as the predicate device.
	the development, validation and	Validation results confirmed that the Sterility Assurance
	routine control of a sterilization	Level (SAL) is identical for both devices (10 ⁻⁶), thereby
	process for medical devices	demonstrating substantial equivalence.
LAL Testing	USP <85>, Bacterial Endotoxins	LAL testing was performed on the subject device using
		the same method and acceptance criteria as the
	USP <161>, Transfusion and	predicate device. Results of the testing confirmed that
	Infusion Assemblies and Similar	the LAL level is identical between the subject and
	Medical Devices	predicate devices, thereby demonstrating substantial
		equivalence. Additionally, LAL testing will be performed
Distribution and	D7296 16 ASTM Standard	Dost distribution to tind product release.
Distribution and	D7380-10, ASTIVI Stalluard	Post-distribution testing was performed on the subject
Distribution	of Packages for Single Parcel	the predicate device. Results of the testing confirmed
Testing	Delivery Systems	that package integrity is identical between the subject
		and predicate devices, thereby demonstrating
		substantial equivalence.

Test Description	Conforming Standard(s)	Conclusion in Support of Substantial Equivalence
Accelerated and	F1980-16, ASTM Standard Guide	Post-accelerated aging testing was performed on the
Post-	for Accelerated Aging of Sterile	subject device using the same method and acceptance
Accelerated	Barrier Systems for Medical	criteria as the predicate device. Results of the testing
Aging Studies	Devices	was identical between the subject and predicate devices,
		thereby demonstrating substantial equivalence.

VIII. CONCLUSIONS

Based on the results of the non-clinical performance testing conducted on the subject device, it has been concluded that the proposed EverLift[™] Submucosal Lifting Agent is as safe and effective and performs as well as the legally marketed predicate device (K191923). The similar indications for use, intended use, technological characteristics, and performance characteristics for the proposed EverLift[™] Submucosal Lifting Agent have been assessed to be substantially equivalent to the predicate device, and any differences do not raise different issues of safety and effectiveness when compared to the predicate device. Therefore, the EverLift[™] Submucosal Lifting Agent is substantially equivalent to the predicate device.