SUMMARY OF SAFETY AND EFFECTIVENESS DATA (SSED)

I. <u>GENERAL INFORMATION</u>

Device Generic Name: Multipurpose Defibrillator Pads (Manual or Automated External Defibrillator Pad)

Device Trade Name: HeartSync Multifunction Disposable Single-Use AED Defibrillator Pads

Device Procode: MKJ

Applicant's Name and Address: Graphic Controls dba Nissha Medical Technologies/ Vermed/ Biomedical Innovations 400 Exchange Street, Buffalo, NY 14204

Date(s) of Panel Recommendation: None

Premarket Approval Application (PMA) Number: P200007

Date of FDA Notice of Approval: 6/23/2023

II. **INDICATIONS FOR USE**

These are multifunction pads, and can be used with automatic or manual defibrillators for monitoring, pacing, cardioversion, as well as defibrillation. These indications are consistent with current AHA Guidelines.

For Automatic External Defibrillators:

(Compatible Model AEDs: Physio Control: LifePak-15, LifePak-20/20e, LifePak -1000; Zoll Medical: R-Series, X-Series; Cardiac Science: PowerHeart AED G3 Plus, PowerHeart AED G3 Pro).

When used with an external defibrillator, these electrode pads are for treating patients in cardiopulmonary arrest who are:

- Unconscious,
- Not breathing spontaneously
- Without circulation (without a pulse).

The pads are single use and intended to be used in conjunction with an external defibrillator to monitor and deliver defibrillation energy to the patient. The pads are used

on patients over 8 years of age or greater than 55 pounds. The pads are intended for short term use (less than 8 hours).

DO NOT DELAY THERAPY IF YOU ARE NOT SURE OF EXACT AGE OR WEIGHT.

For Manual Defibrillators:

Manual Defibrillators can be used for monitoring, pacing, cardioversion, as well as defibrillation.

When used for defibrillation, these electrode pads are for treating patients in cardiopulmonary arrest who are:

- Unconscious,
- Not breathing spontaneously
- Without circulation (without a pulse).

The pads are single use and intended to be used in conjunction with an external defibrillator to monitor and deliver defibrillation energy to the patient. The pads are used on patients greater than 10 kg or 22 pounds. The pads are intended for short term use (less than 24 hours).

DO NOT DELAY THERAPY IF YOU ARE NOT SURE OF EXACT AGE OR WEIGHT.

III. <u>CONTRAINDICATIONS</u>

Contraindications for use should be followed as per the compatible AED/ defibrillation unit.

IV. WARNINGS AND PRECAUTIONS

The warnings and precautions can be found in the HeartSync Defibrillation Electrodes labeling.

V. <u>DEVICE DESCRIPTION</u>

The HeartSync electrodes are intended to be used in external pacing, cardioversion, defibrillation, and monitoring applications as a non-sterile, disposable device for single patient use only. The electrodes provide the conductive interface between the defibrillator and/or the external transcutaneous (noninvasive) cardiac pacemaker and the patient's skin. When a patient requires defibrillation, cardioversion or external pacing, these electrodes will be applied to the patient and connected to the instrument. The pads are packaged in airtight pouches either with the leads inside the package or leads outside the package. The package is designed to keep the pads clean and prevent them from drying out. See Figure 1 below.

These are single use, non-sterile, self-stick defibrillator electrodes packaged in pairs. The effective electrode area is 207.6 cm2. They are packaged with a connector to match the

specific defibrillator. The patient contact material is a conductive adhesive hydrogel. Each electrode consists of foam backing, conductor plate (either tin or carbon), layer of electrically conductive and adhesive hydrogel, protective release liner, insulated wire connected to the conductor plate, and a connector specific to the defibrillator manufacturer.



Figure 1



See Table 1 below for a list of HeartSync Defibrillation Models and the compatible devices.

Table 1

OEM	Defibrillators	HeartSync Models
Physio Control	LifePak-15	T100AC-Physio
	LifePak- 20/20e	T100LOAC-Physio
	LifePak-1000	C100AC-Physio
		C100LOAC-Physio
Zoll Medical	R-Series	T100AC-Zoll
	X-Series	T100LOAC-Zoll
		T100AC-Zoll-10
		C100LOAC-Zoll
		C100AC-Zoll
Cardiac Science	PowerHeart AED G3 Plus	T100-CS
	PowerHeart AED G3 Pro	

VI. <u>ALTERNATIVE PRACTICES AND PROCEDURES</u>

There are several other alternatives for the correction of external pacing, defibrillation, and monitoring. Alternative products (including those marketed by the AED device

original equipment manufacturer can be utilized. Traditional CPR can be considered in the absence of an AED. Each alternative has its own advantages and disadvantages.

VII. MARKETING HISTORY

The HeartSync Defibrillation Electrodes have been marketed in the United States as seen in Table 2 below.

Device Name	Applicant	510(K) Number	Decision Date
Heart Sync	Heart Sync, Inc	K131550	12/06/2013
Heart Sync	Heart Sync, Inc	K131494	09/16/2013
Heart Sync	Heart Sync, Inc	K20536	08/06/2012
pediatric Physio			
AED Pad			
Heart Sync	Heart Sync, Inc	K081442	09/24/2008
Pediatric, Model			
Ped-100			
Heart Sync,	Heart Sync, Inc	K080421	02/29/2008
Models C-100 and			
T-100			

Table 2 510(k) References

HeartSync Defibrillation Electrodes are marketed in Europe, Canada, Middle East, Latin America and the United States and have not been withdrawn from any market due to concerns regarding safety or effectiveness.

Single use defibrillator pads have been in distribution since 1994 and have been the subject of numerous 510(k) submissions. As a result of FDA Final Order [Docket no. FDA-2013-N-0234], single use defibrillator pads (accessories) for AED use are now required to have premarket approval under product code MKJ as a Class III device. Prior to Feb 3, 2015, defibrillator pads were marketed as Class II devices requiring a 510(k). HeartSync previously submitted and marketed single use defibrillator pads under the 510(k)s described above.

VIII. POTENTIAL ADVERSE EFFECTS OF THE DEVICE ON HEALTH

Below is a list of the potential adverse effects (e.g., complications) associated with the use of the device.

- There is a risk of failure to defibrillate a patient in cardiac arrest if the connections to the pad fail.
- There is a risk of electrical shock to the person administering the AED therapy.
- There is a risk of burn to the patient at the electrode site if the skin-pad resistance is high.
- There is a slight risk of skin irritation at the electrode site, mitigated by biocompatibility testing.

IX. <u>SUMMARY OF NONCLINICAL STUDIES</u>

A. Laboratory Studies

Tort		Standard/A acortan ac	Degulta
lest	Objective	Standard/Acceptance	Results
Biocompatibility	To determine acceptable use	ISO 10993-1, ISO	Pass
	on patient skin through	10993-5, and ISO	
	cytotoxicity studies,	10993-10	
	sensitization studies, and skin		
	irritation studies per ISO		
	standard		
	Patient Contacting Materials:		
	Hydrogel Foam Adhesives		
	and Foam		
	und i ouni.		
	Cytotoxicity: Hydrogel –		
	Direct Contact		
	Catatonicity From and From		
	Cytotoxicity: Foam and Foam		
	Adhesives – ISO Agarose		
	Overlay Method.		
	Sensitization: Closed Patch		
	Sensitization Study in Guinea		
	Pigs		
	1 155		
	Skin Irritation: ISO Skin		
	Irritation Study in Rabbits		
Packaging			
Shelf Life	To verify established shelf life	ASTM F1980	Shelf Life
	and function.		
		ASTM F1886-08R13	
	Samples were both real-time		
	aged and aged using accelerated		
	methods to 3.5 years then	Must meet criteria for	
	subjected to Bubble Emissions,	each post-test after	
	Pouch Seal Strength, and Visual	aging	
	inspection.		

Table 3

Test	Objective	Standard/Acceptance	Results
Pouch Integrity:	To determine pouch ability to	ASTM F2096-11	Pouch Integrity:
Bubble Emissions	maintain hydrogel moisture		Bubble Emissions
	and pressure within pouch.		
	Uses externally applied pressure	Must not have steady	
	to detect escaping air.	stream of bubbles.	
Pouch Integrity:	To determine pouch ability to	Pressure must be greater	Pouch Integrity:
Pressure Decay	maintain hydrogel moisture	than 7.0in H2O.	Pressure Decay
	and pressure within pouch.		
	Uses internally applied pressure		
	(13.8 in H2O) and is measured		
	after 60s.		
Pouch Seal	To determine pouch ability to	ASTM F088-09	Pass
Strength	maintain hydrogel moisture		
	and pressure within pouch.		
		Peel is 1.0lbf/min or	
		greater	
	Pulls two sides of pouch away		
	180° angle (ASTM Technique		
	C).		
Transit Test	To ensure pouches maintain	ASTM D4169-16	Pass
	seal and internal integrity after		
	shipping and handling.	ASTM D4332-14	
	Samples are pre-conditioned	ASTM F1886-08R13	
	then subjected to Distribution		
	cycle for Handling, Vehicle		
	Stacking, Loose-Load Vibration,	Must meet criteria for	
	Vibration Concentrated Impact	each sub-test after	
	and Second Handling tests.	aging.	
	Samples are then tested for		
	Bubble Emissions, Pouch Seal		
	Strength, and Visual Inspection.		D
Storage Test	To ensure pouches and	ASIM D4332-14	Pass
	hydrogel maintain properties	ASTM F1886-08R13	
	alter storage at temperature		
	extremes.		
		Must most suitaris for	
	_	each sub-test after	
	After storage at extreme low	aging.	
	temperature, tested for Bubble		
	and Visual Inspection.		

Test	Objective	Standard/Acceptance	Results
Mechanical			
Adhesion Peel	To verify the hydrogel adheres	ASTM D6252	Pass
Strength	to the patient.		
	Samples are pulled at a rate of 12in/min at a 90° angle to a sliding stainless-steel plate while the force is measured.	Average peak force shall be greater than 0.5ozf/min	
Wearability	To verify that electrodes	AAMI EC12-	Pass
	remain adhered to the patient for up to 24 hours.	2000+(R2020)	
	participants for over 24 hours.	Samples are still adhered after 24 hours	
Connector Compatibility	To ensure the connector properly connects to the generator and transmits energy. Perform measurements against drawings and perform tolerance analysis with Root Sum Squared compared with the mating connector.	Measurements do not preclude connector from functioning to retain receptacle and electrical contact.	Pass
Cable Anchorage	To ensure cable can withstand	ANSI/AAMI/IEC	Pass
	loading off-axis without loss of function.	60601-2- 4:2010/A1:2018	
	Electrode and connector tested at termination of cable. The connection material is oscillated with a 5N load applied at angles of 45° and 90° and a rate of 30 cycles/min. Samples are then subjected to Electrical Performance tests.	Cable has not worked loose and passes post- test acceptance criteria.	
First Article Inspections	To ensure materials meet design specifications for size and manufacturer testing processes.	Acceptance per each component's specifications.	Pass
	material.		

Test	Objective	Standard/Acceptance	Results
Electrode Area	To ensure an area large enough	ANSI/AAMI/IEC	Pass
	to disperse energy without	60601-2-	
	burning the patient.	4:2010/A1:2018	
	Measure dimensions of the		
	effective area of the electrode to		
	calculate the area.	Area of each individual	
		electrode shall be	
		greater than 50cm ² and the sum of the pair shall	
		be greater than 150 cm^2 .	
Electrical			-
Oscilloscope Wayafarm Tasting	To compare the waveform	AAMI TIR62	Pass
wavelorm resting	parameters with the waveform		
	of the OEM electrodes.		
		Energy must be	
		within 15%, or 3J of	
	Electrodes are shocked at	the OEM, whichever	
	impedance from 25Ω to 200Ω	1s greater.	
	(when applicable) at maximum	All other parameters	
	energy or highest energy	were nearly identical to	
	protocol and at lowest energy or	the OEM and met	
	energy protocol.	acceptance criteria	
Surface & Volume	To ensure sufficient resistivity	Resistivity must be	Pass
Resistivity	to disperse the energy across	within $\pm 5\%$ of the	
	the effective area of the	published value for the	
	electrode.	gel resistivity.	
	Resistive measurements taken on		
	a sample of gel material.		-
Electrical	To ensure adequate energy	ANSI/AAMI/IEC	Pass
Defibrillator	dispersion in a timely manner.	60601-2-	
Recovery		4:2010/A1:2018	
	Potential measured across		
	shock and performed before and	Potential shall not	
	after pacing for 60min.	4s Pre-Pace: 500mV	
		60s Pre-Pace: 400mV	
		4s Post-Pace: 1000mV	
		60s Post-Pace: 750mV	

Test	Objective	Standard/Acceptance	Results
Electrical	Offset voltage to ensure proper	ANSI/AAMI/IEC	Pass
Performance:	baselining of measurement	60601-2-	
Direct Current	equipment.	4:2010/A1:2018	
Offset			
	Measured across electrodes after		
	1 min stabilization through a		
	$10M\Omega$ resistor.	Offset Voltage shall not	
Flectrical	To ensure proper measurement	ANSI/AAMI/IEC	Dass
Performance	of impedance of the patient for	ANSI/AAMI/ILC	1 455
Alternating	later in the state of the state of the state	4.2010/41.2019	
Current Signal	determination of shock	4:2010/A1:2018	
Impedances (Small	delivery energy.		
10Hz & 30kHz	Using frequencies of 10Hz and		
and Large)	30kHz, impressed current of	AC Small 10Hz:	
	100µA measured peak-to-peak	3000Ω	
	across electrode pair (AC		
	Small), measures the impedance.	AC Small 30kHz:	
	Test is repeated with 50Ω load	10Ω	
	and measured at 360J of		
	Defibrillator shock (AC Large).	AC Large: 5Ω	
	All tests are repeated after 60min		
Energy Discharge	To ensure electrodes can	Total energy shall not	Pass
Verification	withstand 50 defibrillation	decrease by more than	1 455
	shocks	10% of the 3^{rd} shock on	
	SHOCKS.	any of the 50 shocks.	
	Total energy of the shock is		
	measured over 50 repeated		
	shocks at maximum energy of		
Desing	the defibrillator.	$T_{ha} 0.50/0.00/$	Deag
Confirmation	10 ensure electrodes can	Confidence/Population	Pass
Commination	perform per published pacing	of the lower limit on	
	charts.	total energy delivered	
		across all shocks shall	
		not be below -10% of	
	At both 20ms and 40ms pulse	the selected energy for	
	widths, samples were tested at	each of the	
	various parameters per the	configurations	
	HeartSync Pacing Guidance	(Conductor and Pulse	
	charts, then subjected to three	width).	
	(3) defibrillation shocks at		
	maximum energy.		

Test	Objective	Standard/Acceptance	Results
Dielectric Strength	To ensure cable can withstand	ANSI/AAMI/IEC	Pass
Testing	high-voltage shocks and	60601-2-	
	protect the user.	4:2010/A1:2018	
		т 1 (1111	
	Perform H1-Pot test on samples	Leakage current shall be	
	of all connector lead types.	less than 15 μA.	

B. <u>Animal Studies</u> There were no animal studies conducted.

C. Additional Studies

Table 4				
Test	Objective	Standard/Acceptance	Results	
Labeling		· · · · ·		
GHTF Label and	To ensure compliance to	GHTF/SG1/N70:2011	Pass	
Instructions for	GHTF/SG1/N70:2011			
Use		ISO 15223-1:2021		
	Compare GHTF requirements	All labels and IFUs must		
	against labels and IFU	meet standard		
FDA Device	To ensure compliance to 21	21 CFR 820.120	Pass	
Labeling 120	CFR 820.120.			
_				
		All labels and IFUs must		
	Compare CFR requirements	meet standard		
D . 001 I 1 1	against labels and IFU	a		
Part 801 Labeling	To ensure compliance to	Sec 801.15	Pass	
	Subpart A	All labels and IEUs must		
		meet Subpart A 801 15		
	Compare Subpart A			
	requirements against labels and			
	IFÛ			
Labeling	To ensure labeling meets	All product labeling	Pass	
Verification	Product Specifications	specifications must be		
		met		
	Compare labeling requirements			
	within the Product Specification			
	to IF Us and product labels			

Test	Objective	Standard/Acceptance	Results
Human Factors		· · · · · · · · · · · · · · · · · · ·	
Summative Usability Study	To assess labeling and Instructions for Use (IFU) effectiveness in a simulated clinical environment.	ANSI/AAMI/IEC 62366:2007/(R) 2013 ANSI/AAMI HE75:2009/(R) 2018	Pass
	Evaluate users as they perform use case scenarios of HeartSync electrodes. Participants' ability to properly use the device was evaluated against the intended use.	FDA – Applying Human Factors and Usability Engineering to Medical Devices, Guidance for Industry and Food and Drug Administration Staff, Feb 3, 2016	
		Each participant's accuracy and overall study accuracy must meet minimum percentage per protocol.	
Human Factors Comparative Study	To ensure human factors of the use of HeartSync devices do not introduce any additional mis-use risk when compared to the OEM device.	Comparative risk assessment shows no increased risk with HeartSync devices	Pass
	Comparison of IFU for OEM to HeartSync.		

X. <u>SUMMARY OF PRIMARY CLINICAL STUDY</u>

As per FDA Final Order, "Effective Date of Requirement for Premarket Approval for Automated External Defibrillator Systems," published February 3, 2015, existing published clinical literature may be leveraged as part of the PMA submission. Published clinical studies, performance testing, and post market surveillance data were utilized to support the clinical application of the subject devices.

A. Study Design

Oscilloscope Waveform Testing was completed to compare the waveform parameters of HeartSync Defibrillation Electrodes with the waveform of the OEM electrodes. Defibrillation shocks were applied through the electrodes with a range of impedance from 25Ω to 200Ω (when applicable) at maximum energy or highest energy protocol and at lowest energy or energy protocol. Energy was found to be within 15%, or 3J of the OEM, whichever is greater.

All other parameters (listed below) were nearly identical to the OEM and met predetermined acceptance criteria.

- Peak current of the leading edge of the first and second phase
- Peak voltage of the leading edge of the first and second phase
- First and second phase duration
- First and second phase tilt
- Selected energy and delivered energy

Nissha Medical performed an oscilloscope waveform comparison study, where the peak voltage, peak current, phase duration, and tilt were evaluated for both phases of the biphasic waveform as well as total energy delivered. The data generated for this study proved that the HeartSync electrodes performed the same as the OEM electrodes, thus indicating the same therapy is being delivered to the patient. The output waveforms were evaluated for clinical relevance in the respective SSED of the defibrillator as detailed by each manufacturer below.

B. Safety and Effectiveness Results

PUBLISHED CLINICAL DATA

The results of the Physio-Control OEM electrode to HeartSync comparisons demonstrated a comparable performance to the OEM electrodes. Therefore, as the Physio-Control electrodes have been approved by the FDA and considering¹ the clinical evaluation in the van Alem's 2003 study², the HeartSync electrodes for Physio-Control defibrillators can be inferred to be as effective as the OEM Physio-Control electrodes for electrophysiological treatments using the Physio-Control biphasic truncated exponential waveform.

The results of the Cardiac Science OEM electrode to HeartSync comparison demonstrated a comparable performance to the OEM electrodes. Therefore, as the OEM Cardiac Science electrodes have been demonstrated as clinically relevant^{3,4} with the clinical evaluation referenced in Kerber's 1997 study⁵, the HeartSync electrodes for Cardiac Science defibrillators can be inferred to be as effective as the OEM Cardiac Science strengthered for electrophysiological treatments using the Cardiac Science STAR® biphasic waveform.

The results of the Zoll OEM electrode to HeartSync comparison demonstrated a comparable performance to the OEM electrodes. Therefore, as the OEM Zoll electrodes have been demonstrated as clinically relevant⁶ with the clinical evaluation referenced in Hess's 2012 study⁷, the HeartSync electrodes for Zoll defibrillators can be inferred to be as effective as the OEM Zoll electrodes for electrophysiological treatments using the Zoll rectilinear biphasic waveform.

ELECTRICAL TEST RESULTS

In addition to the Oscilloscope Waveform analysis identified above that was conducted using the HeartSync electrodes as compared to OEM electrodes, and the supporting clinical data provided for the OEM electrodes; the HeartSync product was evaluated against the current international standard: ANSI/AAMI/IEC 60601-2-4:2010/A1:2018 for electrical performance testing in order to demonstrate safety and effectiveness:

- 1. Surface & Volume Resistivity testing was conducted to ensure sufficient resistivity to disperse the energy across the effective area of the electrode.
- 2. Defibrillator Recovery testing was conducted to ensure adequate energy dispersion in a timely manner.
- 3. Direct Current Offset testing was performed to ensure proper baselining of measurement equipment.
- 4. Alternating Current Signal Impedance test ensuring proper measurement of impedance of the patient for determination of shock delivery energy.
- 5. Energy Discharge Verification testing was conducted to ensure electrodes can withstand 50 defibrillation shocks.
- 6. Pacing Confirmation testing ensures electrodes can perform per published pacing charts.
- 7. Dielectric Strength Testing was conducted to ensure cable can withstand high-voltage shocks and protect the user.

HUMAN FACTORS TEST RESULTS

The safe and effective use of the HeartSync product was demonstrated with the successful completion of the following human factor studies:

- 1. Summative Usability Study, used to assess labeling and Instructions for Use (IFU) in a simulated clinical environment. Users were evaluated as they performed use case scenarios of HeartSync electrodes. Participants' ability to properly use the device was evaluated against the intended use.
- 2. Human Factors Comparative Study, conducted to ensure human factors of the HeartSync devices do not introduce any additional mis-use risk when compared to the OEM device.
- 1. <u>Pediatric Extrapolation</u>

In this premarket application, existing clinical data was not leveraged to support approval of a pediatric patient population.

C. Financial Disclosure

The Financial Disclosure by Clinical Investigators regulation (21 CFR 54) requires applicants who submit a marketing application to include certain information concerning the compensation to, and financial interests and arrangement of, any clinical investigator conducting clinical studies covered by the regulation. There was no pivotal clinical study. None of the clinical investigators had disclosable financial interests/arrangements as defined in sections 54.2(a), (b), (c), and (f). The information provided does not raise any questions about the reliability of the data.

XI. <u>SUMMARY OF SUPPLEMENTAL CLINICAL INFORMATION</u>

Since the Acquisition of HeartSync by Graphic Controls in May 2018 and through June 2023, there have been a total of 61 complaints related to product quality for HeartSync Defibrillation Electrodes. These complaints were evaluated and compared to the total amount of units sold in this period. The evaluation determined that the HeartSync Defibrillation Electrodes product family yielded a complaint rate of only .0013%. In addition, these complaints were reviewed regarding adverse event reporting. This review concluded that the complaint rate for reportable events was 0.0003%.

XII. PANEL MEETING RECOMMENDATION AND FDA'S POST-PANEL ACTION

In accordance with the provisions of section 515(c)(3) of the act as amended by the Safe Medical Devices Act of 1990, this PMA was not referred to the Cardiovascular Devices Panel, an FDA advisory committee, for review and recommendation because the information in the PMA substantially duplicates information previously reviewed by this panel.

XIII. CONCLUSIONS DRAWN FROM PRECLINICAL AND CLINICAL STUDIES

A. Effectiveness Conclusions

The outcomes of the Nonclinical Studies (Performance Testing) have demonstrated that the HeartSync Defibrillation Electrodes waveform delivery is comparable to the OEM defibrillation pads and have proven via design verification the reasonable effectiveness of the device. This includes the applicable specifications outlined in IEC 60601-2-4:2010/A1 :2018.

Conclusions of the Design Verification related to wearability of the device demonstrated that HeartSync Defibrillation Electrodes remain adhered to the patient for 24 hours after application. Additionally, the Human Factors Assessment performed compared the usability of HeartSync Defibrillation Electrodes to the OEM counterparts, and found that in all cases, the information on how to use the HeartSync device is both pictorially clear and the instruction as complete as the OEM devices. Since the Acquisition of HeartSync by Graphic Controls in May 2018 and through June 2023, there have been a total of 61 complaints related to product quality for HeartSync Defibrillation Electrodes. This represents a complaint rate of .0013%.

B. Safety Conclusions

The risks of the device are based on Nonclinical Studies (Performance Testing), existing clinical data, design verification, and post market surveillance conducted to support PMA approval as described above. Due to the longstanding commercial use of the HeartSync product line, the results of performance testing, and the low complaint rate and adverse event rate the risks and safety of the device have been found to be appropriately mitigated.

C. Benefit-Risk Determination

The probable benefits of the device are also based on data collected in a clinical study conducted to support PMA approval as described above. HeartSync Defibrillation Electrodes were originally reviewed and FDA-cleared via FDA 510(k) notifications and have been commercially available since 2008.

The established benefits and risks of the electrodes have been evaluated using post market surveillance information, existing clinical data, and design verification/performance testing.

The HeartSync Defibrillation Electrodes are accessories to AEDs. The benefit of AED use is increased chance of survival in a cardiac arrest by acutely treating ventricular fibrillation and pulseless ventricular tachycardia. Defibrillation in emergency events is life- saving and the ability to deliver defibrillation therapy quickly is imperative.

The probable risks of the device are also based on post market surveillance information and existing clinical data. The benefits of defibrillation to sustain life far outweigh the risks associated with it.

1. Patient Perspective

This submission either did not include specific information on patient perspectives or the information did not serve as part of the basis of the decision to approve or deny the PMA for this device.

In conclusion, given the available information above, the data supports that for external pacing, cardioversion, defibrillation, and monitoring applications the probable benefits outweigh the probable risks.

D. Overall Conclusions

The data in this application support the reasonable assurance of safety and effectiveness of this device when used in accordance with the indications for use. For patients in cardiopulmonary arrest who are unconscious, not breathing spontaneously, or without circulation (without a pulse), or in need of external pacing, cardioversion, or monitoring, the benefits of HeartSync Defibrillation Electrodes outweigh the risks.

XIV. CDRH DECISION

CDRH issued an approval order on 6/23/2023.

The applicant's manufacturing facilities have been inspected and found to be in compliance with the device Quality System (QS) regulation (21 CFR 820).

XV. <u>APPROVAL SPECIFICATIONS</u>

Directions for use: See device labeling.

Hazards to Health from Use of the Device: See Indications, Contraindications, Warnings, Precautions, and Adverse Events in the device labeling.

Post-approval Requirements and Restrictions: See approval order.

XVI. <u>REFERENCES</u>

- 1. Physio-Control LIKEPAK FDA Summary of Safety and Effectiveness Data (SSED), PMA P160026
- van Alem AP, Sanou BT, Koster RW. Interruption of cardiopulmonary resuscitation with the use of the automated external defibrillator in out-ofhospital cardiac arrest. Ann Emerg Med. 2003 Oct;42(4):449-57. doi: 10.1067/s0196-0644(03)00383-4. PMID: 14520315. <u>https://doi.org/10.1016/S0300-9572(03)00216-8</u>
- 3. Zoll Medical Corporation FDA Summary of Safety and Effectiveness Data (SSED), PMA P160033
- 4. Zoll Medical Corporation FDA Summary of Safety and Effectiveness Data (SSED), PMA P160034
- 5. RE Kerber, et al. "Automated External PMA P160034: FDA Summary of Safety and Effectiveness Data Page 16 Defibrillators for Public Access Defibrillation: Recommendations for Specifying and Reporting Arrhythmia Analysis Algorithm

Performance, Incorporation New Waveforms, and Enhancing Safety." Circulation (1997), Vol 95, No 6, 1677-1681. https://www.ahajournals.org/doi/10.1161/01.CIR.95.6.1677

- 6. ZOLL Medical Corporation FDA Summary of Safety and Effectiveness Data (SSED), PMA P160022
- Hess EP, White RD. Automated external defibrillation. Crit Care Clin. 2012 Apr;28(2):143-53. doi: 10.1016/j.ccc.2011.10.009. Epub 2011 Dec 1. PMID: 22433479. <u>https://doi.org/10.1016/j.ccc.2011.10.009</u>