

November 5, 2021

Bestsafe Glove CO., LTD. Piyawat Chirasakulkarun Chief Executive Officer 52/43 Wat Khot Hin-Khao Phai Road, T. Tubma Muang Rayong, Rayong 21000 Thailand

Re: K210253

Trade/Device Name: Best Glove-Latex Powder Free Examination Glove Regulation Number: 21 CFR 880.6250 Regulation Name: Non-powdered patient examination glove Regulatory Class: Class I, reserved Product Code: LYY Dated: May 25, 2021 Received: August 30, 2021

Dear Piyawat Chirasakulkarun:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray, III, Ph.D. Assistant Director DHT4B: Division of Infection Control and Plastic Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

**Indications for Use** 

Form Approved: OMB No. 0910-0120 Expiration Date: January 31, 2017 See PRA Statement below.

**510(k) Number (***if known***)** K210253

#### Device Name

BEST GLOVE -LATEX POWDER FREE EXAMINATION GLOVE

#### Indications for Use (Describe)

A powder-free patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

X Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRA Staff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



52/43 WATKODHIN-KAOPHAI RD. T.TUBMA, A.MUANG, RAYONG, RAYONG. THAILAND 21000 TEL OFFICE : 038-949877-79 FAX : (+66)038-694493, 038-949899 TAX ID : 0215563006105 TEL H/P : +66-0-86-3173770 , +66-0-81-4022099 ,+66-0-62-325-6558 EMAIL : SALES@BESTSAFE.CO.TH WEBSITE : WWW.BESTSAFE.CO.TH I WWW.BESTSAFEGLOVE.COM

## 510(k) SUMMARY K210253 Latex Powder Free Examination Glove

#### **1.0** Applicant Information

	Applicant:	BESTSAFE GLOVE CO., LTD
	Address	52/43 Watkodhin-Kaophai Rd. T. Tubma,
		A. Muang Rayong Rayong. THAILAND 21000
	Phone Number:	+66 (038)-949880
	Fax Number:	+66 (038)-949850
	Name of Contact Person:	Piyawat Chirasakulkarun
	Contact Number:	+66 81 4022099
	Contact Email:	sales@bestsafe.co.th
	Preparation date:	September 6, 2021
2.0	Device Identification:	
	Trade/Proprietary Name(s):	BEST GLOVE -LATEX POWDER FREE EXAMINATION GLOVE
	Common Name:	Latex Powder Free Examination Gloves
	Classification Name:	Patient Examination Gloves
	510(k):	K210253
	Device Class:	I
	Product Code:	LYY
	Registration Number	21 CFR 880.6250
	Review Panel	General Hospital
3.0	Predicate Device	
	Device Name:	MPXX <sup>™</sup> Powder Free Natural Rubber Latex Examination Gloves
	Manufacturer:	Total Glove Company Sdn. Bhd.
	510(k):	K110250
	Device Class:	I
	Product Code:	LYY



52/43 WATKODHIN-KAOPHAI RD. T.TUBMA, A.MUANG, RAYONG, RAYONG. THAILAND 21000 TEL OFFICE : 038-949877-79 FAX : (+66)038-694493, 038-949899 TAX ID : 0215563006105 TEL H/P : +66-0-86-3173770 , +66-0-81-4022099 ,+66-0-62-325-6558 EMAIL : SALES@BESTSAFE.CO.TH WWW.BESTSAFE.CO.TH I WWW.BESTSAFEGLOVE.COM

#### 4.0 Description of the Device:

Latex Examination Powder Free Gloves are manufactured to meet all the current specifications listed under the ASTM Specification D3578 - 19, Standard Specification for Rubber Examination Gloves. They are made from Natural Rubber Latex. These gloves are natural in color (no color is added) and are powder free.

#### 5.0 Indication for Use of the Device:

A powder-free patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner

#### 6.0 Comparison of Technological characteristics between predicate and subject devices

Characteristics	Poforoncos/Standards	Device per			
Characteristics	References/Standarus	Predicate	Current	Comparison	
510(k) Number	-	K110250	K210253		
Manufacturer(s)	-	Total Glove Company Sdn. Bhd.	BESTSAFE GLOVE CO., LTD		
Name of device	-	MPXX <sup>™</sup> Powder Free Natural Rubber Latex Examination Gloves	BEST GLOVE -LATEX POWDER FREE EXAMINATION GLOVE		
Indication for Use	Medical Gloves Guidance Manual, Issued on January 22, 2008	MPXX <sup>™</sup> Powder Free Natural Rubber Latex Examination Gloves are single use device intended for medical purposes that is worn on the hand of healthcare and similar personnel to prevent contamination between healthcare and the patient.	A powder-free patient Examination Glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner	Same	
Material	ASTM D3578-19	Natural Rubber Latex	Natural Rubber Latex	Identical	
Color		Natural Colo	Natural Colo	Same	
Size	Medical Glove Guidance Manual- Labeling- Issued on January 22, 2008	Extra Small Small Medium Large Extra Large	Small Medium Large Extra Large	Different	



 52/43
 WATKODHIN-KAOPHAI RD. T.TUBMA, A.MUANG, RAYONG, RAYONG. THAILAND 21000

 TEL OFFICE : 038-949877-79
 FAX : (+66)038-694493, 038-949899
 TAX ID : 0215563006105

 TEL H/P : +66-0-86-3173770 , +66-0-81-4022099 ,+66-0-62-325-6558
 TAX ID : 0215563006105

EMAIL : SALES@BESTSAFE.CO.TH WEBSITE : WWW.BESTSAFE.CO.TH I WWW.BESTSAFEGLOVE.COM

		Device performance				
Characteristics	References/Standards	Predicate		Current		Comparison
Single Use	Medical Gloves Guidance Manual - Issued on January 22, 2008	Single use		Single use	9	Same
Sterile/Non-sterile	-	Non-sterile		Non-steril	e	Same
Dimension	ASTM D3578-19	<u>Length</u> 230 mm minimum	<u>Length</u> 230 mm min		Similar	
			Size	•	Actual value	
			Sma	1	240	
			Mediu	m	240	
			Large	e	240	
			Extra La	arge	241	
		<u>Width</u> No data is available	Mediur (fo	<u>Width</u> m: 95 mm : r Medium s	± 10 mm size)	Different
			Size	)	Actual value	
			Sma	II	85	
			Mediu	ım	93	
			Larg	e	105	
			Extra La	arge	115	
Thickness	ASTM D3578-19	Finger = 0.08 mm. min Palm = 0.08 mm. Min	Finger = 0.08 Palm = 0.08	3 mm. min mm. Min		Similar
				Palm	Finger	
			Size	(Actual value)	(Actual value)	
			Small	0.09	0.12	
			Medium	0.09	0.12	
			Large	0.10	0.13	
			Extra Large	0.09	0.12	



52/43 WATKODHIN-KAOPHAI RD. T.TUBMA, A.MUANG, RAYONG, RAYONG. THAILAND 21000 TEL OFFICE : 038-949877-79 FAX : (+66)038-694493, 038-949899 TAX ID : 0215563006105 TEL H/P : +66-0-86-3173770 , +66-0-81-4022099 ,+66-0-62-325-6558 EMAIL : SALES@BESTSAFE.CO.TH WEBSITE : WWW.BESTSAFE.CO.TH I WWW.BESTSAFEGLOVE.COM

	References/Standard	Devic			
Characteristics	S	Predicate	Current		Comparison
Physical Properties	ASTM D3578-19	Before Aging Tensile Strength	Before Tensile Strength	<u>Aging</u> (18 MPa min)	Similar
		18 MPa min <u>After Aging</u>	Size	Actual value	
		Tensile Strength	Small	18.04	
		14 MPa min	Medium	19.13	
			Large	18.26	
			Extra Large	18.64	Similar
			Tensile Strength	(14 MPa min)	Similar
			Size	Actual value	
			Small	15.02	
			Medium	18.12	
			Large	16.01	
			Extra Large	15.06	Similar
		Before Aging Ultimate Elongation	<u>Before Aging</u> Ultimate Elongation : 650% min		Similar
		650% min	Size	Actual value	
		After Aging	Small	651	
		500% min (after aging)	Medium	654	
			Large	650	
			Extra Large	650	Similar
			<u>After Aging</u> Ultimate Elongation : 500% min		
			Size	Actual value	
			Small	501	
			Medium	601	
			Large	501	
			Extra Large	502	



Γ

## BESTSAFE GLOVE CO.,LTD. | บริษัท เบสท์เซฟ โกลฟ จำกัด

 52/43
 WATKODHIN-KAOPHAI RD. T.TUBMA, A.MUANG, RAYONG, RAYONG. THAILAND 21000

 TEL OFFICE
 : 038-949877-79
 FAX : (+66)038-694493, 038-949899
 TAX ID : 0215563006105

 TEL H/P : +66-0-86-3173770 , +66-0-81-4022099 ,+66-0-62-325-6558
 TAX ID : 0215563006105
 TAX ID : 0215563006105

Dovic	o porformanco	
EMAIL : SALES@BESTSAFE.CO.TH	WEBSITE : WWW.BESTSAFE.CO.TH I WWW.BESTS	SAFEGLOVE.COM

Characteristics	Poforoncos/Standards	Devic			
Characteristics	References/Standards	Predicate	Current		Comparison
Water Tight Test (1000 ml)	ASTM D5151 – 19	Pass AQL 1.5	Pass	AQL 1.5	Same
Powder Residual	ASTM D6124-06 (Reapproved 2017)	Meet <u>&lt;</u> 2.0 mg/glove	<u>&lt;</u> 2.0 mg/glove		Similar
			Size	Residual powder content (mg/glove)	
			Small	0.62	
			Medium	0.46	
			Large	0.61	
			Extra Large	0.61	
	Primary Skin Irritation – ISO 10993-10 Third Edition 2010-08-01	Not a primary skin irritant under the conditions of the study	Not a primary the conditions of	skin irritant under the study	Same
Biocompatibility	Dermal Sensitization – ISO 10993-10 Third Edition 2010-08-01	Not a contact sensitizer under the conditions of the study	Not a contact se conditions of the	nsitizer under the study	Same
	In vitro cytotoxicity ISO10993-5 :2009(E)	No data is available	Under the conditional cytotoxic for under the conditional transmission of the cytotoxic for under the study, 1:8, 1:16 and Moreover, under the study, non a toxic.	tions of the study liluted (neat) and noncytotoxic for d 1:32 dilutions. the conditions of cute systemic	Different
	Acute Systemic Toxicity ISO10993-11:2017(E)	No data is available	Under the condit did not induce a toxicity.	tions of the study, ny systemic	Different



52/43 WATKODHIN-KAOPHAI RD. T.TUBMA, A.MUANG, RAYONG, RAYONG. THAILAND 21000

 TEL OFFICE
 : 038-949877-79
 FAX: (+66)038-694493, 038-949899
 TAX ID: 0215563006105

 TEL H/P:
 +66-0-86-3173770, +66-0-81-4022099, +66-0-62-325-6558
 EMAIL:
 SALES@BESTSAFE.CO.TH
 WEBSITE:
 WWW.BESTSAFE.CO.TH
 WWWW.BESTSAFE.CO.TH
 WWW.BESTSAFE

7.0 Summary of Non-clinical performance Tests

Test Method	Standard	Purpose of testing	Acceptance Criteria	Result	Status
	ASTM D3578-19 Standard Specification for Rubber Examination Gloves	To determine the length of the gloves	Min 230 mm for all sizes	Small: 240 mm Medium: 240 mm Large: 240 mm Extra Large: 241 mm	Pass
Dimension	ASTM D3578-19 Standard Specification for Rubber Examination Gloves	To determine the width of the gloves	Small: 80 ± 10 mm Medium: 95+/-10mm Large: 111± 10 mm Extra Large: 120 ± 10 mm	Small: 85 mm Medium: 93 mm Large: 105 mm Extra Large: 115 mm	Pass
	ASTM D3578-19 Standard Specification for Rubber Examination Gloves	To determine the thickness of the gloves	Palm 0.08 mm min Finger 0.08 mm min for all sizes	Small: Palm 0.09 mm, Finger: 0.12 mm Medium: Palm: 0.09 mm, Finger: 0.12 mm Large: Palm 0.10 mm, Finger: 0.13 mm Extra Large: Palm 0.09 mm, Finger: 0.12 mm	Pass
Physical Properties	ASTM D3578-19 Standard Specification for Rubber Examination Gloves	To Determine the physical properties- Tensile strength	Before Ageing Tensile Strength 18Mpa Minimal for all sizes After Ageing Tensile Strength 14Mpa Minimal for all sizes	Before AgeingSmall: 18.04 MPaMedium: 19.13 MPaLarge: 18.26 MPaExtra Large: 18.64 MPaAfter ageing:Small: 15.02 MPaMedium: 18.12 MPaLarge: 16.01 MPaExtra Large: 15.06 MPa	Pass



# BESTSAFE GLOVE CO.,LTD. | บริษัท เบสท์เซฟ โกลฟ จำกัด 52/43 watkodhin-kaophai rd. t.tubma, a.muang, rayong, rayong. thailand 21000

TEL OFFICE : 038-949877-79 FAX : (+66)038-694493, 038-949899 TEL H/P: +66-0-86-3173770, +66-0-81-4022099,+66-0-62-325-6558 TAX ID: 0215563006105

EMAIL : SALES@BESTSAFE.CO.TH WEBSITE : WWW.BESTSAFE.CO.TH I WWW.BESTSAFEGLOVE.COM

Test Method	Standard	Purpose of testing	Acceptance Criteria	Result	Status
Physical Properties	ASTM D3578-19 Standard Specification for Rubber Examination Gloves	To Determine the physical properties- Ultimate Elongation	Before Ageing Ultimate Elongation 650% Min for all sizes After Ageing Ultimate Elongation 500% Min for all sizes	Before Ageing Small: 651%Medium: 654 % Large:650%Extra Large:650%After ageing: Small:501%Medium: 601% Large:501%Extra Large:502%	Pass
Watertight test	ASTM D5151-19 Standard Test Method for Detection of Holes in Medical Gloves	To determine the holes in the gloves	Sample size: 200 pcs Inspection level : GI AQL 1.5 Acceptance Number 7 Rejection Number 8	The batch size for this sampling is 35,001- 150,000. Hence, according to the single sampling plan GI, the sample to be drawn is under code L equivalent to 200 pcs with accept 7 and reject 8 to be accept under AQL 1.5. Small: 0 (Zero) Medium: 0 (Zero) Large:0 (Zero) Extra Large:0 (Zero)	Pass
Residual powder	ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves	To determine the residual powder in the gloves	2 mg per glove or less	Sample size : 5 pcs Requirement: 2 mg per glove or less Result : Small:0.62 mg/glove Medium: 0.46 mg/glove Large: 0.61 mg/glove Extra Large:0.61 mg/glove	Pass



52/43 WATKODHIN-KAOPHAI RD. T.TUBMA, A.MUANG, RAYONG, RAYONG. THAILAND 21000

 TEL OFFICE
 : 038-949877-79
 FAX:
 (+66)038-694493, 038-949899
 TAX ID:
 0215563006105

 TEL H/P:
 +66-0-86-3173770, +66-0-81-4022099, +66-0-62-325-6558
 EMAIL:
 SALES@BESTSAFE.CO.TH
 WEBSITE:
 WWW.BESTSAFE.CO.TH
 WWW.BESTSAFE.CO

Purpose of testing **Acceptance Criteria Test Method** Standard Result Status Protein content ASTM D5712 - 15, Standard Test To determine the Less than 200 µg/dm<sup>2</sup> Sample size : 3 pcs Pass Method for Analysis of Aqueous extractable protein in Requirement: Less than 200 µg/dm<sup>2</sup> Extractable Protein in Natural the gloves Result Small:124.36 µg /dm<sup>2</sup> Rubber. Medium: 140.78 µg /dm<sup>2</sup> Large: 134.26 µg /dm<sup>2</sup> Extra 159.46 µg /dm<sup>2</sup>

The performance test data of the non-clinical tests meet following standards:

ASTM D 3578 – 19 Standard Specification for Rubber Examination Gloves

ASTM D 5151-19: Standard Test Method for Detection of Holes in Medical Gloves

ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves

ASTM D5712 – 15, Standard Test Method for Analysis of Aqueous Extractable Protein in Natural Rubber.



 BESTSAFE GLOVE CO.,LTD. | บริษัท เบสท์เซฟ โกลฟ จำกัด

 52/43 WATKODHIN-KAOPHAI RD. T.TUBMA, A.MUANG, RAYONG, RAYONG. THAILAND 21000

 TEL OFFICE : 038-949877-79

 FAX : (+66)038-694493, 038-949899

 TAX ID : 0215563006105

 TEL H/P : +66-0-86-3173770, +66-0-81-4022099, +66-0-62-325-6558

 EMAIL : SALES@BESTSAFE.CO.TH

#### 8.0 Summary of Clinical Performance Tests:

Not Applicable

#### 9.0 Conclusion:

The conclusion drawn from the non-clinical tests demonstrates that the subject device Latex Examination Powder Free Gloves is as safe, as effective, and performs as well as or better than the legally marketed predicate device K110250

--END---